Stahl Bohrer • Steel Burs
Instrumentes en Acier • Instrumentos de Acero
Strumenti in Acciaio

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Correct application

- Ensure that only technically perfect, serviced and clean turbines, handpieces and contra-angles are used with the rotary instruments.
- The instruments must be inserted as far as possible.
- The instruments must be running before being placed on the surface.
- Avoid wedging or levering the instrument as this increases the risk of breakage.

Pressure

- With rotary grinding instruments excessive pressure may lead to breakouts or cause them to smear, which increases heat build-up (max. 0,3 - 2 N).

Cool adequately with water

- To prevent undesirable heat build-up during preparation, provide for adequate water irrigation (at least 50ml/min.).

Dispose of worn instruments

- Fractured and deformed blades cause vibration.
- Bent or non-concentric rotary instruments must also be disposed of.

Disinfection, cleaning and sterilization

- Rotary instruments must be disinfected, cleaned and sterilized prior to first use on patients and immediately after use.
- The method of use, reaction time and suitability of disinfectants and cleaning agents for certain types of instruments are covered by the manufacturers' instructions.

Recommended speeds

- To produce optimum results, run the rotary instruments at their recommended speeds.
- The recommended speeds and maximum permissible working speeds are included in the manufacturer's instructions.

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