

OsteoBiol[®]
by Tecnos

YEARS
20
ANNIVERSARY

Bone Grafting Materials


REGENERATION SCIENCE

INSPIRED BY NATURE


SEM image of an OsteoBiol® Gen-Os® granule colonised by osteoblasts
Author: Prof Ulf Nånmark, University of Göteborg, Sweden




4-months histology showing remnants of porcine dermal matrix surrounded by connective tissue
Experimental study on dogs. Htx-eosine
Author: Prof Ulf Nånmark, University of Göteborg, Sweden



LM image of an OsteoBiol® Lamina hydrated with blood: vascularisation enhanced by the presence of the original vascular canals
Author: Prof Ulf Nånmark, University of Göteborg, Sweden



Part of a biopsy showing newly formed bone after treatment with OsteoBiol® Putty. Htx-eosine. Magnification x20
Author: Prof Ulf Nånmark, University of Göteborg, Sweden



OUR MISSION

«To produce a xenogenic bone substitute as similar as possible to autogenous bone»

Giuseppe Oliva MD
R&D Director
Tecnoss S.r.l.




THE OSTEObIOL® DUAL-PHASE HETEROLOGOUS BONE MATRIX

OsteoBiol® is the family of biomaterials produced by TecnoSS® for the dental and maxillo-facial surgeons.

In each OsteoBiol® granule, besides its mineral phase, the TecnoSS® process retains the xenogenic collagen phase with its precious biological properties, making it biocompatible and ideal for grafting and augmentation purposes.

Avoiding high process temperatures, the OsteoBiol® bone matrix prevents ceramization, maintaining a chemical composition extremely similar to autogenous bone⁽¹⁾, and therefore gradually resorbable and replaceable by newly formed bone.

(1) Figueiredo M, Henriques J, Martins G, Guerra F, Judas F, Figueiredo H
Physicochemical characterization of biomaterials commonly used in dentistry as bone substitutes - comparison with human bone
Journal of Biomedical Materials Research Part B: Applied Biomaterials, 2010 Feb; 92(2):409-19




SEM image of an OsteoBiol® Gen-Os® granule: osteoblastic colonisation. Magnification x3000
Author: Prof Ulf Nånberg, University of Göteborg, Sweden

HIGH BIOCOMPATIBILITY

The chemical structure of each OsteoBiol® dual-phase granule, its ideal porosity and collagen content, make it a valid scaffold and substrate for osteoblasts anchorage, proliferation and new bone apposition⁽²⁾.

(2) Nannmark U, Sennerby L

The bone tissue responses to prehydrated and collagenated cortico-cancellous porcine bone grafts: a study in rabbit maxillary defects
Clinical Implant Dentistry and Related Research, 2008 Dec;10(4):264-70



Biopsy from upper jaw region retrieved at four months showing gradual osteoclastic resorption of an OsteoBio!® GTO® granule
Author: Prof Ulf Nannmark, University of Göteborg, Sweden

GRADUAL RESORPTION

Autogenous bone is gradually replaced by newly formed bone: similarly, the OsteoBiol® bone matrix allows progressive osteoclastic resorption, with simultaneous new bone apposition.

Cells receive nutrients from newly formed vessels, that are able to colonize adequately the grafted site.

New bone grows in and around the OsteoBiol® granules⁽³⁾, which are partially but significantly replaced by vital bone at re-entry time.

(3) Giuliani A, Iezzi G, Mazzoni S, Piattelli A, Perrotti V, Barone A

Regenerative properties of collagenated porcine bone grafts in human maxilla: demonstrative study of the kinetics by synchrotron radiation microtomography and light microscopy
Clinical Oral Investigations, 2018 Jan;22(1):505-513




VASCULARIZATION IS THE KEY FOR CLINICAL SUCCESS

Dual-phase biomaterials are progressively resorbed by osteoclasts and replaced by new vital bone produced by osteoblasts, similarly to autogenous bone grafts. Both types of cells live thanks to blood supply, which is critical and essential for the success of any bone regeneration procedure.

The progressive resorption of OsteoBiol[®] granules allows an adequate colonization of the grafting site by new vessels, and is therefore a positive and significant factor within the regenerative process⁽⁴⁾.

(4) Rombouts C, Jeanneau C, Camilleri J, Laurent P, About I
Characterization and angiogenic potential of xenogeneic bone grafting materials: Role of periodontal ligament cells
Dental Materials Journal, 2016 Dec 1;35(6):900-907



THE ROLE OF COLLAGEN

Collagen favours MSC differentiation and enhances osteoblasts proliferation^(5,6): it is considered as the ideal substrate for bone forming cells. OsteoBiol® dual-phase particulate bone substitutes contain approximately 22% collagen.

Furthermore, collagen gel mixed with dual-phase collagenated granules packed in syringes improves the handling and the stability of the graft, reducing also operatory time and risk of contamination.

(5) Brunelli G, Sollazzo V, Carinci F, Palmieri A, Girardi A, Monguzzi R
OsteoBiol® influences osteogenic differentiation of adipose derived stem cells
European Journal of Inflammation, 2011, Vol. 9, no. 3 (S), 103-107

(6) Jeanneau C, Le Fournis C, About I
Xenogeneic bone filling materials modulate mesenchymal stem cell recruitment: role of the complement c5a
Clinical oral investigations; 2019 oct 23




Fig. 1




Fig. 2




Fig. 3




Fig. 4



Fig. 5



Fig. 6



Fig. 7



Fig. 8



Fig. 9

A SPECIFIC PRODUCT FOR EVERY CLINICAL INDICATION

OsteoBiol® is not only a marvellous collagenated bone matrix: it is a complete family of biomaterials specifically designed for bone and soft tissue augmentation in dentistry. For every clinical indication a dedicated product has been developed, with the goal of providing the best handling, the ideal granulometry and consistency, and finally optimal regenerative results in adequate re-entry time.

Enjoy one of the widest and most complete product ranges, with the security and support of 15 years of clinical research: you will experience that today it is finally possible to achieve predictable clinical success⁽⁷⁾ without the availability limitations of autogenous bone.

(7) Checchi V, Felice P, Zucchelli G, Barausse C, Piattelli M, Pistilli R, Grandi G, Esposito M


Wide diameter immediate post-extractive implants vs delayed placement of normal-diameter implants in preserved sockets in the molar region: 1-year post-loading outcome of a randomised controlled trial
European Journal of Oral Implantology, 2017;10(3):263-278

Images authors pag 14: Fig.1, Dr Patrick Palacci, Marseille, France | Fig. 2-9, Dr Roberto Rossi, Genova, Italy

OsteoBiol® and natural human bone have the same density and very similar physico-chemical properties

Figueiredo et al. J Biomed Mater Res B: Appl Biomater, 2010 Feb; 92(2):409-19


Mineral content



Figueiredo et al. J Biomed Mater Res B: Appl Biomater, 2010 Feb; 92(2):409-19


Gen-Os® has a higher angiogenic potential compared to anorganic xenografts

Rombouts et al. Dent mater J, 2016 Dec 1;35(6):900-907



In ridge preservation collagenated biomaterials show significant smaller volume reduction and basal area shrinkage compared to slowly resorbable xenografts

Barone et al. Clin Oral Implants Res, 2016 Nov;27(11):E105-E115



OsteoBiol® bone matrix promotes osteoblast differentiation and bone regeneration

Brunelli et al. Eur J Inflamm, 2011, Vol. 9, no. 3 (S), 103-107



OsteoBiol® bone scaffolds absorb growth factors secreted by MSCs and improve bone tissue repair

Mijiritsky et al. Materials, 2017 Sep 8;10(9)

KEY SCIENTIFIC DATA

Over 190 articles have been published on peer-reviewed journals during the last 15 years, proving with in-vitro, experimental and clinical studies the outstanding biological properties and clinical performance of the OsteoBiol® collagenated biomaterials.

For further information see the complete literature on p. 92

Images author pag 16: Fig.1, 2, 3 - Prof Ulf Nannmark, University of Göteborg, Sweden



PATIENTS FIRST

Combining the best skills and the best materials, within the limits and guidelines provided by scientific evidence, is the key for clinical success: however let us all remember that the patients are and will always be the center of all our attentions.

Meeting their expectations, helping them to recover function and aesthetics with long term success^(8,9) is the greatest reward for any surgeon and fulfillment of our company mission.

(8) Barone A, Orlando B, Tonelli P, Covani U

Survival rate for implants placed in the posterior maxilla with and without sinus augmentation: a comparative cohort study
Journal of Periodontology, 2011 Feb; 82(2):219-26

(9) Scarano A, Piattelli A, Assenza B, Quaranta A, Perrotti V, Piattelli M, Iezzi G

Porcine bone used in sinus augmentation procedures: a 5-year retrospective clinical evaluation
Journal of Oral and Maxillofacial Surgery, 2010 Aug;68(8):1869-73

OsteoBiol® products clinical evidence*

Gen-Os®

Collagenated heterologous cortico-cancellous bone mix
Granulometry 250-1000 µm
For information on OsteoBiol® Gen-Os® see page 24

mp3®

Pre-hydrated collagenated heterologous cortico-cancellous bone mix
Granulometry 600-1000 µm
For information on OsteoBiol® mp3® see page 32

GTO®

Pre-hydrated collagenated heterologous cortico-cancellous bone mix
Granulometry 600-1000 µm
For information on OsteoBiol® GTO® see page 36

Putty

Pre-hydrated collagenated heterologous cortico-cancellous bone paste
Granulometry up to 300 µm
For information on OsteoBiol® Putty see page 40


Gel 40

Pre-hydrated collagenated heterologous cortico-cancellous bone gel
Granulometry up to 300 µm
For information on OsteoBiol® Gel 40 see page 44

ALVEOLAR REGENERATION



MAXILLARY SINUS LIFT




CRESTAL ACCESS ONLY



PERI-IMPLANT DEFECTS



HORIZONTAL AUGMENTATION



VERTICAL AUGMENTATION

INLAY TECHNIQUE



PERIODONTAL REGENERATION



3-WALL DEFECTS



SOFT TISSUE AUGMENTATION

* Based on published scientific literature and clinical experience of expert surgeons.

Apatos

Cortico-cancellous and cortical bone
Granulometry 600-1000 μm
For information on OsteoBio® Apatos
see page 48

Sp-Block

Collagenated heterologous
cancellous block
For information on OsteoBio® Sp-Block
see page 54

Evolution


Heterologous collagen membrane
For information on OsteoBio® Evolution
see page 62

Lamina

Collagenated heterologous cortical bone
For information on OsteoBio® Lamina
see page 70

Derma

Collagen dermal matrix
For information on OsteoBio® Derma
see page 66



WITH TSV GEL



IN ASSOCIATION



BONE LAYER TECHNIQUE

WITH LAMINA AND TSV GEL



BONE SUBSTITUTES



Fig. 1




Fig. 2



Fig. 3




Fig. 4




Fig. 5



Fig. 6





Fig. 7




Fig. 8

OsteoBiol® Dual-Phase bone substitutes


HETEROLOGOUS BONE MATRIX



Cortico-cancellous collagenated matrix					Apatos Cortical cortical bone
	TSV Gel	Collagen gel			
	Pre-hydrated	Pre-hydrated	Pre-hydrated	Pre-hydrated	
Gen-Os®	GTO®	mp3®	Putty	Gel 40	Apatos Mix
100% collagenated bone mix	~80% collagenated bone mix ~20% TSV Gel	~90% collagenated bone mix ~10% collagen gel	~80% collagenated bone mix ~20% collagen gel	~60% collagenated bone mix ~40% collagen gel	cortico-cancellous bone mix
					
Heterologous cortico-cancellous collagenated bone mix	Heterologous cortico-cancellous collagenated pre-hydrated bone mix	Heterologous cortico-cancellous collagenated pre-hydrated bone mix	Heterologous cortico-cancellous collagenated pre-hydrated bone paste	Heterologous cortico-cancellous collagenated pre-hydrated bone gel	Heterologous microcrystalline hydroxyapatite
For more information on OsteoBiol® Gen-Os® see page 24	For more information on OsteoBiol® GTO® see page 36	For more information on OsteoBiol® mp3® see page 32	For more information on OsteoBiol® Putty see page 40	For more information on OsteoBiol® Gel 40 see page 44	For more information on OsteoBiol® Apatos see page 48



Gen-Os[®]



The advantages of a dual-phase biomaterial

Collagenated heterologous cortico-cancellous bone mix

Characteristics and handling



Tissue of origin
Cortico-cancellous heterologous bone mix

Tissue collagen
Preserved

Physical form
Slightly radiopaque granules

Composition
100% granulated mix

Granulometry
250-1000 μm
1000-2000 μm

Re-entry time
4/5 months, depending on grafting site characteristics

Packaging
Vial: 0.25 g, 0.5 g, 1.0 g, 2.0 g

Product codes
250-1000 μm
M1052FS | 1 Vial | 0.25 g | Porcine
M1052FE | 1 Vial | 0.25 g | Equine
M1005FS | 1 Vial | 0.5 g | Porcine
M1005FE | 1 Vial | 0.5 g | Equine
M1010FS | 1 Vial | 1.0 g | Porcine
M1010FE | 1 Vial | 1.0 g | Equine
M1020FS | 1 Vial | 2.0 g | Porcine
M1020FE | 1 Vial | 2.0 g | Equine

1000-2000 μm
M0210FS | 1 Vial | 1.0 g | Porcine
M0220FS | 1 Vial | 2.0 g | Porcine

GMDN code **CND code**
46425 P900402


CHARACTERISTICS

A natural replicate of autologous bone, Gen-Os[®] conserves the same intimate structures⁽¹⁾ (matrix and porous form) and presents highly osteoconductive properties^(2,3). It is biocompatible and bioavailable, as recognized by tests made according to the ISO 10993 method conducted at Eurofins Biolab. Gen-Os[®] is gradually resorbable and provides support in bone neoformation helping to preserve the original graft shape and volume⁽⁴⁾.

Moreover, thanks to its collagen content, the product facilitates blood clotting and the subsequent invasion of repairing and regenerative cells, favouring *restitutio ad integrum* of missing bone. Because of its marked hydrophilia⁽⁵⁾, it can function as a carrier for selected medications and drugs⁽⁶⁾ and it is ideal to mix with GFs⁽⁷⁾.

HANDLING

Gen-Os[®] must always be hydrated and thoroughly mixed with either a few drops of sterile physiological solution (or patient's blood) to activate its collagen matrix and to enhance its adhesivity or with TSV Gel to increase graft stability in not self-contained defects. If necessary, it can as well be mixed with the drug selected for surgery.




SEM image of OsteoBio[®] Gen-Os[®] granules. Magnif. x50
Author: Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecnos[®] Dental Media Library

Gen-Os[®], a cortico-cancellous bone mix, has been the first product developed with the TecnoSS[®] innovative biotechnology and, due to its universal use, still is today the most demanded from the market. Gen-Os[®] has been successfully used and documented for alveolar ridge preservation⁽⁸⁾ in combination with *Evolution* membranes: the application of this biomaterial limits the alveolar ridge width reduction that would naturally occur with spontaneous healing, preserving thus the alveolar ridge volume and allowing a correct second stage implant placement⁽⁹⁾. Gen-Os[®] has been used for lateral access maxillary sinus lift⁽³⁾ and dehiscence regeneration⁽¹⁰⁾, always in association with *Evolution* membranes. Gen-Os[®] has been documented in periodontal regeneration of deep infrabony defects⁽¹¹⁾. Due to its collagen content, once hydrated Gen-Os[®] becomes very sticky and hydrophilic⁽⁵⁾: it combines therefore extremely well with blood and is very stable once applied into the grafting site.

Its cortico-cancellous composition allows a progressive resorption of osteoclastic type, with in parallel a similar rate of new bone formation⁽²⁾: these unique properties allow a very good graft volume preservation, a healthy and well vascularized new bony tissue and, ultimately, a successful implant rehabilitation. Gen-Os[®] is in fact able to boost vascularization: in vitro⁽¹²⁾ assays proved an increase in the secretion of VEGF by periodontal ligament cells (PDL) in the presence of Gen-Os[®], as well as an enhanced proliferation of endothelial cells.



LATERAL ACCESS SINUS LIFT
maxillary sinus floor augmentation




PERIODONTAL REGENERATION
intrabony defects



HORIZONTAL AUGMENTATION
two-wall defects



DEHISCENCES AND FENESTRATIONS
peri-implant lesions



CRESTAL ACCESS SINUS LIFT
osteotome technique




ALVEOLAR REGENERATION
socket preservation

BIBLIOGRAPHY

- (1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19
- (2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (3) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4
- (4) CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42
- (7) MIJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B
PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR
MATERIALS, 2017 SEP 8;10(9)
- (8) CHECCHI V, FELICE P, ZUCCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(3):263-278
- (9) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13
- (10) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (11) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRA-BONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (12) ROMBOUITS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907

For further information see the complete literature on p. 92



TSV Gel



The resorbable solution for ideal graft stability
Thermosensitive resorbable gel for graft stabilization

Characteristics and handling



Composition

Heterologous type I and III collagen gel
Thermogelling synthetic biocompatible copolymer

Physical form

LV phase at $+4^{\circ}\text{C}$
Gel viscosity at >math>+13^{\circ}\text{C}</math>

Packaging

Syringe: 0.5 cc, 1.0 cc

Available only in combination with OsteoBio!® Gen-Os® and Apatos
0.5 g, 1.0 g

Product codes

TSV005S | 1 Syringe | 0.5 cc | Porcine
TSV005E | 1 Syringe | 0.5 cc | Equine
TSV010S | 1 Syringe | 1.0 cc | Porcine
TSV010E | 1 Syringe | 1.0 cc | Equine

GMDN code

46425

CND code

P900402



CHARACTERISTICS

The purpose of *TSV Gel* is to provide mechanical stability to bone substitutes and barrier membranes.


TSV Gel is sterilized by Gamma irradiation and is radio-transparent. It contains heterologous type I and III collagen gel with polyunsaturated fatty acids diluted in aqueous solution containing a biocompatible synthetic copolymer that gives *TSV Gel* thermo-reversible and thermo-gelling properties. At low temperature (+4°C) the gel is relatively flowable and easy to mix and manipulate with the graft but it becomes more viscous when *in situ* and exposed to body temperature.

HANDLING

TSV Gel must be refrigerated for at least 20 minutes at +4°C before use, in order to reach the low viscosity (LV) phase, which makes it easier to mix with *Gen-Os®* or *Apatos*.

At room temperature, the product remains at LV phase for few minutes, whereas once *in situ* its viscosity quickly increases with body temperature. *TSV Gel* in LV phase can be used instead of saline for hydrating and mixing with *Gen-Os®* or *Apatos*. The result will be a sticky mixture easy to place and extremely stable once *in situ*.

TSV Gel can also be applied to the rough side of the *Evolution* membrane to stabilize it during graft covering and whilst suturing.



Part of a biopsy showing newly formed bone around a particle of OsteoBio!® Gen-Os® mixed with OsteoBio!® TSV Gel two weeks after grafting in rabbit. Htx-eosine.

Author: Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecnos® Dental Media Library



Source: Tecnos® Dental Media Library

TSV Gel can be used in GBR procedures together with OsteoBiol® bone substitutes and membranes to enhance graft stability. The viscosity reached by TSV Gel at body temperature improves the stability of Gen-Os® or Apatos granules and has proven particularly beneficial in cases where there is little bony support around the defect i.e. lateral augmentation, sockets with a compromised buccal wall, dehiscences and periodontal two and one wall defects.

Additionally, the viscosity of OsteoBiol® TSV Gel improves the stability and handling of Evolution membranes, particularly during the delicate phase of flap closure.

The above clinical information is based on the experience of expert surgeons



DEHISCENCES AND FENESTRATIONS
peri-implant lesions




ALVEOLAR REGENERATION
socket preservation



PERIODONTAL REGENERATION
intrabony defects




HORIZONTAL AUGMENTATION
two-wall defects



Peri-implant defect treated with OsteoBiol® Gen-Os® mixed with TSV Gel

Author: Dr Roberto Rossi, Genova, Italy

OsteoBiol® TSV Gel GELIFICATION KINETICS




Source: Politecnico di Torino, Italy

The graph shows the effect of temperature change on 3 TSV Gel samples.


As temperature increases from 0°C (1°C/min), the viscosity of the gel reaches its minimum at 4°C.

It then increases rapidly until it plateaus at 13°C. At room and body temperature TSV Gel is gel-like. It does not harden but keeps a soft consistency that allows the mixture with Gen-Os® or Apatos granules. Thanks to the hydrophilic properties of OsteoBiol® bone substitutes, the mixture becomes a sticky, stable conglomerate that can easily be placed in the defect site.

TSV Gel is biocompatible and rapidly resorbed.



mp3[®]



Ultimate performance and handling

Pre-hydrated collagenated heterologous cortico-cancellous bone mix



Tissue of origin

Cortico-cancellous heterologous bone mix

Tissue collagen

Preserved plus an additional 10% collagen gel

Physical form

Pre-hydrated granules and collagen gel

Composition

90% granulated mix, 10% collagen gel

Granulometry

600-1000 μm

1000-2000 μm

Re-entry time

About 5 months

Packaging

Syringe: 0.5 cc, 1.0 cc, 3x0.25 cc, 3x0.5 cc, 3x1.0 cc

Wide tip syringe: 2.0 cc

Product codes

600-1000 μm

A3095FS | 1 Syringe | 0.5 cc | Porcine

A3095FE | 1 Syringe | 0.5 cc | Equine

A3005FS | 1 Syringe | 1.0 cc | Porcine

A3005FE | 1 Syringe | 1.0 cc | Equine

A3075FS | 3 Syringes | 3x0.25 cc | Porcine

A3015FS | 3 Syringes | 3x0.5 cc | Porcine

A3015FE | 3 Syringes | 3x0.5 cc | Equine

A3030FS | 3 Syringes | 3x1.0 cc | Porcine

A3030FE | 3 Syringes | 3x1.0 cc | Equine

A3010FS | 1 Wide tip syringe | 2.0 cc | Porcine

A3010FE | 1 Wide tip syringe | 2.0 cc | Equine

1000-2000 μm

A3210FS | 1 Wide tip syringe | 2.0 cc | Porcine

A3210FE | 1 Wide tip syringe | 2.0 cc | Equine

GMDN code

46425

CND code

P900402

Characteristics and handling

CHARACTERISTICS

Heterologous origin biomaterial made of 600-1000 μm or 1000-2000 μm pre-hydrated collagenated cortico-cancellous granules, properly mixed with collagen gel. Thus, it is possible both skipping the hydration phase and decreasing the risk of accidental exposure of the material to pathogens during manipulation and grafting phases; furthermore, the syringe is flexible and ideal to simplify grafting in the receiving site.

The granules are endowed with characteristics very similar to human mineral bone, and can be used as an alternative to autologous bone.

Their natural micro-porous consistency facilitates new bone tissue formation⁽¹⁾ in defect sites and accelerates the regeneration process.


Gradually resorbable^(2,3), it preserves the original graft shape and volume (osteoconductive property)^(4,5).

Moreover, thanks to its collagen content, the product facilitates blood clotting and the subsequent invasion of repairing and regenerative cells.

HANDLING

mp3[®] is available in ready-to-use syringes and can be easily grafted avoiding the hydration and manipulation phases.

After adapting the material to the defect shape, it is necessary to remove non-stable residues before proceeding to soft tissue suture. It is recommended to always compact *mp3*[®] after grafting to achieve optimal stabilization.



Histology on maxillary sinus biopsy taken at 24 months. 48% new bone formation, 13% residual granules

Author: Biopsy by Dr Roberto Rossi, Genova, Italy. Histology by Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecnos[®] Dental Media Library

mp3® is a pre-hydrated cortico-cancellous bone mix with 10% collagen gel. It has been developed with this innovative biotechnology and is a "ready-to-use" product.

mp3® is commonly used for lateral access maxillary sinus lift^(1,6), always in association with *Evolution* membranes, to cover the antrotomy: the mp3® syringe can be directly applied into the bony window without having to mix the mp3® granules with saline.


Due to its collagen gel content, mp3® allows an excellent graft stability while its hydrophilia guarantees quick blood absorption and therefore the necessary graft vascularization. mp3® has also been successfully used in combination with *Evolution* membranes for alveolar ridge preservation^(3,7,8): the application of this biomaterial limits the alveolar ridge width and height reduction that would naturally occur with spontaneous healing, preserving thus the alveolar ridge volume and allowing a correct second stage implant placement.

mp3® has been documented for horizontal augmentation (two wall defects) in combination with autogenous bone blocks or with OsteoBioL® Lamina^(9,10): its cortico-cancellous composition allows a progressive resorption of osteoclastic type, and in parallel a similar rate of new bone formation⁽²⁾.

These unique properties allow a very good graft volume preservation⁽¹¹⁾, a healthy new bony tissue and ultimately, a successful implant rehabilitation.



Socket grafted with OsteoBioL® mp3®



Periodontal defect grafted with OsteoBioL® mp3®

Author: Dr Roberto Rossi, Genova, Italy

free animated videos
on OsteoBioL® APP



LATERAL ACCESS SINUS LIFT
maxillary sinus floor augmentation



ALVEOLAR REGENERATION
post-extractive sockets





HORIZONTAL AUGMENTATION
two-wall defects

BIBLIOGRAPHY

- (1) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAPTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30
- (2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (3) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLIN ORAL INVESTIG, 2017 2018 JAN;22(1):505-513
- (4) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP
BIOMED RES INT, 2016;2016:4086870
- (5) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZONI L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017, 18, 33-44
- (6) SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9
- (7) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13
- (8) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- (9) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- (10) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
MINERVA STOMATOL, 2016 OCT;65(5):257-68
- (11) BARONE A, TOTI P, MENCHINI-FABRIS GB, DERCHI G, MARCONCINI S, COVANI U
EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL
CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759

For further information see the complete literature on p. 92



GTO®



The new standard of excellence in biomaterials

Collagenated heterologous cortico-cancellous bone mix + TSV Gel

Characteristics and handling



Tissue of origin

Cortico-cancellous heterologous bone mix

Tissue collagen

Preserved

Physical form

Pre-hydrated granules and TSV Gel

Composition

80% granulated mix, 20% TSV Gel

Granulometry

600-1000 μm

Re-entry time

About 5 months

Packaging

Syringe: 0.5 cc

Wide tip syringe: 2.0 cc

Product codes

MU0005S | 1 Syringe | 0.5 cc | Porcine

MU0005E | 1 Syringe | 0.5 cc | Equine

MU0020S | 1 Wide tip syringe | 2.0 cc | Porcine

MU0020E | 1 Wide tip syringe | 2.0 cc | Equine

GMDN code

46425

CND code


P900402

CHARACTERISTICS

Heterologous bone grafting material made of a mix of collagenated cortico-cancellous granules of size ranging from 600 to 1000 μm , properly blended with TSV Gel, which is a mixture of heterologous type I and III collagen gel with polyunsaturated fatty acids and a biocompatible synthetic copolymer diluted in aqueous solution. GTO[®] is gradually resorbed and it is extremely osteoconductive. Moreover, the preserved collagen matrix characterizing the granules facilitates blood clotting and the subsequent invasion of repairing and regenerative cells. These unique properties guarantee an excellent rate of new bone formation, delivering adequate graft volume preservation, a healthy new bony tissue and, ultimately, a successful implant rehabilitation. The presence of the same kind of granules of its progenitor, mp3[®], which are very similar to human mineral bone, assures a similar biological response of the host tissue. GTO[®] can be used as alternative to autologous bone.

HANDLING

GTO[®] is available in two sizes (0.5 and 2.0 cc) as ready-to-use pre-hydrated biomaterial and can be easily grafted to the defect site. Thus, clinicians can skip the hydration step with saline or blood, saving time and decreasing the risk of accidental exposure to pathogens. The presence of TSV Gel ensures optimal stickiness of the material, which is also easily adaptable to the recipient site and extremely stable.



Part of a biopsy showing newly formed bone 4 months after treatment with OsteoBio[®] GTO[®]

Author: Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecno[®] Dental Media Library



Author: Dr Patrick Palacci, Marseille, France





GTO® has been conceived as a universal biomaterial, easily adaptable to any bone defect, in association with *Evolution* membranes or *Lamina* to protect the graft. Nonetheless, thanks to its stickiness, has proved particularly effective for horizontal augmentation procedures (e.g.: two-walls defects, when the crest is resorbed) and for socket preservation cases with compromised buccal plate.

GTO® can also be successfully used to treat peri-implant lesions.

In case of open defects, GTO® should be grafted in consecutive layers compacting each layer with a sterile gauze.

The above clinical information is based on the experience of expert surgeons



Peri-implant defect treated with OsteoBiol® GTO®
Author: Dr Patrick Palacci, Marseille, France



DEHISCENCES AND FENESTRATIONS
peri-implant grafting




ALVEOLAR REGENERATION
post-extractive sockets



HORIZONTAL AUGMENTATION
two-wall defects

OsteoBiol® TSV Gel GELIFICATION KINETICS




Source: Politecnico di Torino, Italy

The graph shows the effect of temperature change on 3 TSV Gel samples.

As temperature increases from 0°C (1°C/min), the viscosity of the gel reaches its minimum at 4°C.

It then increases rapidly until it plateaus at 13°C. At room and body temperature TSV Gel is gel-like. It does not harden but keeps a soft consistency that allows the mixture with Gen-Os® or Apatos granules. Thanks to the hydrophilic properties of OsteoBiol® bone substitutes, the mixture becomes a sticky, stable conglomerate that can easily be placed in the defect site.

TSV Gel is biocompatible and rapidly resorbed.



Putty



Engineered for peri-implant defects
Pre-hydrated collagenated heterologous cortico-cancellous bone paste

Characteristics and handling



Tissue of origin

Cortico-cancellous heterologous bone mix

Tissue collagen

Preserved plus an additional 20% collagen gel

Physical form

Plastic consistency composed of collagen gel loaded with 80% micronized bone mix

Composition

80% granulated mix, 20% collagen gel

Granulometry

Up to 300 μm

Re-entry time

About 4 months

Packaging

Syringe: 0.25 cc, 0.5 cc, 3x0.5 cc, 3x0.25 cc
Wide tip syringe: 1.0 cc

Product codes

HPT52S	1 Syringe	0.25 cc	Porcine
HPT09S	1 Syringe	0.5 cc	Porcine
HPT09E	1 Syringe	0.5 cc	Equine
HPT32S	3 Syringes	3x0.25 cc	Porcine
HPT32E	3 Syringes	3x0.25 cc	Equine
HPT35S	3 Syringes	3x0.5 cc	Porcine
HPT35E	3 Syringes	3x0.5 cc	Equine

HPT61S	1 Wide tip syringe	1.0 cc	Porcine
HPT61E	1 Wide tip syringe	1.0 cc	Equine

GMDN code

46425

CND code


P900402

CHARACTERISTICS

Putty is a bone paste with at least 80% micronized heterologous bone (granulometry up to 300 μm) and collagen gel. It is made with an exclusive process that provides the product with exceptional malleability and plasticity, making it easy to apply into peri-implant defects with walls. Thanks to its collagen component, the product facilitates blood clotting and the subsequent invasion of repairing and regenerative cells, showing an osteoconductive behaviour⁽¹⁾. Successful grafting needs complete stability of the biomaterial: for this reason *Putty* must be used only in cavities able to firmly contain it. Therefore, *Putty* must not be grafted in two wall defects or in lateral access sinus lift procedures.

HANDLING

Inject the product and adapt it to defect morphology without compression; any non-stable residue must be removed before soft tissue suture. An *Evolution* membrane is recommended to protect *Putty* grafted in peri-implant defects.



Part of a biopsy showing newly formed bone after treatment with OsteoBiol® Putty

Author: Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecnos® Dental Media Library



The extraordinary handling properties of *Putty* syringe make this product the ideal choice for self-contained peri-implant defects⁽²⁾ and all small defects that present a self-contained cavity. Furthermore, the TecnoSS® manufacturing process avoids granules ceramization, allowing a progressive resorption of the biomaterial and, at the same time, an adequate new-bone formation rate⁽³⁾. *Putty*'s "soft" consistency also guarantees an easy and healthy soft-tissues healing. Thanks to these unique characteristics, *Putty* has been effectively used for peri-implant defects regeneration: following immediate post-extractive implants placement, *Putty* can be injected between the defect walls and the implant, guaranteeing a perfect filling of the entire defect volume⁽⁴⁾.



SEM image of OsteoBio® Putty
Author: Prof Ulf Nannmark, University of Göteborg, Sweden


The product versatility also makes *Putty* the ideal solution when bone tissue has been lost due to peri-implantitis as long as the containing walls are present. In fact, the primary condition for gaining a successful regeneration is to achieve the biomaterial initial stability. Therefore, *Putty* must be used only in defects where the surrounding walls guarantee such condition: for example inside the bone crest when ridge-split technique is adopted⁽⁵⁾, or with horizontally resorbed crests, in association with OsteoBio® Lamina (Bone Layer technique)⁽⁶⁾.



Peri-implant lesion grafted with OsteoBio® Putty
Author: Dr Roberto Rossi, Genova, Italy




DEHISCENCES AND FENESTRATIONS
peri-implant defects



HORIZONTAL AUGMENTATION
ridge split

BIBLIOGRAPHY

- (1) ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
MINERVA STOMATOL, 2005 JUN;54(6):351-62
- (2) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106
- (3) NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (4) CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (5) SANTAGATA M, GUARINIELLO L, TARTARO G
A MODIFIED EDENTULOUS EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N.:114-9
- (6) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85



Gel 40



A unique heterologous bone gel
Collagenated heterologous cortico-cancellous bone mix



Characteristics and handling


CHARACTERISTICS

Gel 40 is made of a collagen matrix (type I and III) obtained using an exclusive Tecnos[®] process, loaded for 60% of its volume with micronized heterologous bone (granulometry up to 300 μm). Thanks to its collagen component, Gel 40 facilitates the formation of primary blood clot and the subsequent invasion of repairing and regenerative cells; moreover, the cortico-cancellous component provides the necessary scaffold function.

The collagen gel component contained in Gel 40 is rapidly and totally resorbed; it is also endowed with exceptional anti-inflammatory, eutrophic and cicatrizing properties. This lipophilia is due mainly to a percentage of polyunsaturated fatty acids of the oleic-linoleic series (to which Omega 3 also belongs) directly derived from the raw material. Such components possess a valuable antioxidant action on the free radicals and therefore aid tissue regeneration.

HANDLING

The distinctive characteristics of viscosity and density of Gel 40 facilitate the handling of the product by the operator, providing a glue-like support. If viscosity is excessive, add a few drops of sterile lukewarm saline and then re-mix thoroughly to obtain the desired density.



Part of a biopsy showing newly formed bone after treatment with OsteoBiol[®] Gel 40. Biopsies were taken 5 weeks after implantation in rabbit maxillae. Htx-eosine. Original magnification x20
Author: Prof Ulf Nannmark, University of Göteborg, Sweden



Source: Tecnos[®] Dental Media Library

Tissue of origin

Cortico-cancellous heterologous bone mix

Tissue collagen

Preserved plus an additional 40% collagen gel

Physical form

Collagen gel type I and III loaded with 60% bone mix

Composition

60% granulated mix, 40% collagen gel

Granulometry

Up to 300 μm

Re-entry time

About 4 months

Packaging

Syringe: 0.5 cc, 3x0.5 cc

Product codes

05GEL40S | 1 Syringe | 0.5 cc | Porcine
05GEL40E | 1 Syringe | 0.5 cc | Equine
15GEL40S | 3 Syringes | 3x0.5 cc | Porcine
15GEL40E | 3 Syringes | 3x0.5 cc | Equine

GMDN code

46425

CND code

P900402



The exclusive Tecnos® manufacturing process guarantees an exceptional malleability and plasticity: furthermore, the syringe packaging provides Gel 40 extraordinary handling properties making this product the ideal choice for crestal access sinus lift^(1,2), deep and narrow peri-implant defects⁽³⁾, three-wall intrabony defects and, in combination with Evolution membranes, for treating gingival recessions⁽⁴⁾.

Furthermore, the Tecnos® manufacturing process avoids granules ceramization, allowing a progressive resorption of the biomaterial and, at the same time, an adequate new-bone formation rate^(5,6).

Gel 40 "soft" consistency also guarantees an easy and healthy soft-tissues healing.



Crestal access sinus lift with OsteoBiol® Gel 40
Source: Tecnos® Dental Media Library



PERIODONTAL REGENERATION
intrabony defects and gingival recessions




CRESTAL ACCESS SINUS LIFT
crestal sinus floor augmentation


BIBLIOGRAPHY

- (1) BARONE A, CORNELINI R, CIAGLIA R, COVANI U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9
- (2) SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS
J ORAL IMPLANTOL, 2010; 36(6):485-489
- (3) COVANI U, CORNELINI R, BARONE A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6
- (4) CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESIONS USING A COLLAGEN MEMBRANE WITH A DEMINERALIZED XENOGRFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67
- (5) NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3
- (6) LORENZON G, BUTTARELLO GM, CHESSA G
CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION
DENTISTRY, 2015, 5:2

For further information see the complete literature on p. 92



Apatos



Microcrystalline hydroxyapatite

Heterologous cortico-cancellous and cortical bone

Characteristics and handling



Tissue of origin

Apatos Mix: cortico-cancellous heterologous bone mix
Apatos Cortical: heterologous cortical bone

Tissue collagen

Degraded

Physical form

Radiopaque granules of mineral hydroxyapatite

Composition

Apatos Mix: 100% cortico-cancellous mix
Apatos Cortical: 100% cortical bone

Granulometry

600-1000 μm
1000-2000 μm

Re-entry time

About 5 months

Packaging

Mix | Vial: 0.5 g, 1.0 g, 2.0 g
Cortical | Vial: 0.5 g, 1.0 g

Product codes

600-1000 μm
Mix | A1005FS | 1 Vial | 0.5 g | Porcine
Mix | A1005FE | 1 Vial | 0.5 g | Equine
Mix | A1010FS | 1 Vial | 1.0 g | Porcine
Mix | A1010FE | 1 Vial | 1.0 g | Equine
Mix | A1020FS | 1 Vial | 2.0 g | Porcine
Mix | A1020FE | 1 Vial | 2.0 g | Equine
Cortical | AC1005FS | 1 Vial | 0.5 g | Porcine
Cortical | AC1010FS | 1 Vial | 1.0 g | Porcine
1000-2000 μm
Mix | A0210FS | 1 Vial | 1.0 g | Porcine
Mix | A0210FE | 1 Vial | 1.0 g | Equine

GMDN code

46425

CND code


P900402

CHARACTERISTICS

Apatos is a biocompatible^(1,2), osteoconductive^(3,4) biomaterial of heterologous origin with characteristics similar to mineralized human bone⁽⁵⁾; it can therefore be used as an alternative to autologous bone. The natural microporous consistency of Apatos facilitates the formation of new bone tissue in bone defect area⁽⁶⁾, accelerating the process. Apatos microcrystalline hydroxyapatite is available in cortical and mixed granules.

HANDLING

Apatos must always be hydrated and thoroughly mixed with a few drops of sterile saline or with *TSV Gel* to increase graft stability in not self-contained defects; it can also be mixed with patient's blood. Finally it can be mixed if necessary with the drug selected for surgery; the mixture thus obtained should be positioned with a sterile spatula or syringe for biomaterials.



SEM image of OsteoBiol® Apatos, cancellous granules

Source: Nobil Bio Ricerche, Villafranca d'Asti, Italy



Source: Tecnos® Dental Media Library




Apatos is a universal filler, that can be used to treat peri-implant defects and two-wall defects^(7,8). Because of its granulometry, Apatos cannot be used in narrow defects, but it fits well in big sockets, e.g. after molar extractions⁽⁹⁾. Both sinus lift procedures (with crestal or lateral access)^(2,10) can be performed with Apatos as bone substitute, as well as surgeries for horizontal regenerations.

Apatos Cortical is characterized by a very long resorption time⁽¹¹⁾, guaranteeing adequate preservation of the grafted volume.

When needed, Apatos grafts can be protected with OsteoBiol® Evolution membrane⁽¹²⁾ or stabilized with Cortical Lamina.

The above clinical information is based on the experience of expert surgeons



LATERAL ACCESS SINUS LIFT
maxillary sinus floor augmentation




ALVEOLAR REGENERATION
socket preservation



HORIZONTAL AUGMENTATION
two-wall defects



DEHISCENCES AND FENESTRATIONS
peri-implant grafting






CRESTAL ACCESS SINUS LIFT
osteotome sinus floor augmentation

BIBLIOGRAPHY

- (1) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91
- (2)ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC;77(12):1984-90
- (3) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMM, 2011, VOL. 9, NO. 3(S), 103-107
- (4) CAKIR M, KARACA IR, AYSEGÜL F, KAYMAZ F, BOZKAYA S
EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION
CLIN ORAL IMPLANTS RES, 2015 MAR-APR;30(2):279-85
- (5) KOLMAS J, SZWAJA M, KOLODZIEJSKI W
SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENIC BIOMATERIALS USED AS BONE SUBSTITUTES
J PHARM BIOMED ANAL, 2012 MAR 5;61:136-41
- (6) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRI B, DI FELICE R, MARCHIONNI S, CALVO GUIRADO JL, COVANI U, NANNMARK U
CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAFTS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL
J CLIN PERIODONTOL, 2017 FEB;44(2):204-214
- (7) BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106
- (8) BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26
- (9) BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U
VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAFTS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL
CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26
- (10) IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V
COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76
- (11) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- 12 | MARCONCINI S, GIAMMARINARO E, DERCHI G, ALFONSI F, COVANI U, BARONE A
CLINICAL OUTCOMES OF IMPLANTS PLACED IN RIDGE-PRESERVED VERSUS NONPRESERVED SITES: A 4-YEAR RANDOMIZED CLINICAL TRIAL
CLIN IMPL DENT RELAT RES, 2018 Dec;20(6):906-914

For further information see the complete literature on p. 92


BLOCKS



OsteoBiol® bone blocks

Sp-Block

collagenated cancellous bone




SEM image of OsteoBiol® Sp-Block. Magnification 200x.
Source: Politecnico di Torino, Italy
For more information on OsteoBiol® Sp-Block see page 54

Dual-Block

collagenated cortico-cancellous bone



SEM image of OsteoBiol® Dual-Block. Magnification 20x.
Source: Politecnico di Torino, Italy
For more information on OsteoBiol® Dual-Block see page 54




Sp-Block

Cancellous block for the inlay technique in the mandible




Highly osteoconductive properties



Dual-Block

Cortico-cancellous scaffold for horizontal augmentation in the maxilla



Characteristics, handling and clinical information

free animated videos
on OsteoBiol® APP



CHARACTERISTICS


Sp-Block is a cancellous block of xenogenic bone produced with an exclusive TecnoSS® process which avoids ceramization of the hydroxyapatite crystals, thus accelerating physiological resorption. *Sp-Block* supports new bone formation⁽¹⁾: thanks to its rigid consistency it is able to maintain the original graft volume, which is particularly important in case of large regenerations. Moreover, its collagen content facilitates blood clotting and the subsequent invasion of regenerative and repairing cells, favoring the *restitutio ad integrum* of missing bone.

HANDLING

Sp-Block must be hydrated before use for 5/10 minutes with sterile lukewarm physiological solution or with antibiotics. Afterwards, it can be adapted to the receiving site; the block must always be fixed with osteosynthesis microscrews and should be protected with a resorbable membrane (*Evolution*).

CLINICAL INFORMATION

Sp-Block has been documented in cases where a vertical gain in the posterior mandible is required, to achieve an augmentation of maximum 5 mm, by means of the inlay technique⁽²⁻⁴⁾. The gaps around the block can be filled with a biomaterial in granules; the augmented area is stabilized with mini-plates and screws and covered with an *Evolution* membrane.



SEM image of OsteoBiol® cancellous block
Author: Prof Ulf Nannmark, University of Göteborg, Sweden

Tissue of origin

Cancellous bone

Tissue collagen

Preserved

Physical form

Rigid dried block

Composition

Collagenated cancellous bone

Re-entry time

About 8 months, variable depending on characteristics and irradiation grade of grafting site and on clinical conditions of the patient

Packaging

Sterile blister

Product codes

BNOE | 10x10x10 mm | Equine
BN1E | 10x10x20 mm | Equine
BN2E | 10x20x20 mm | Equine
BN8E | 35x10x5 mm | Equine

GMDN code

46425

CND code

P900402

BIBLIOGRAPHY

(1) SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP
BIOMED RES INT, 2016;2016:4086870

(2) ESPOSITO M, BARAUSSE C, PISTILLI R, PIATTELLI M, DI SIMONE S, IPPOLITO DR, FELICE P

POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
INT J OF ORAL IMPLANTOL, 2019;12(1):39-54

(3) FELICE P, BARAUSSE C, BARONE A, ZUCHELLI G, PIATTELLI M, PISTILLI R, IPPOLITO DR, SIMION M
INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP
INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480

(4) BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U
EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY
J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447




VERTICAL AUGMENTATION
inlay technique

CHARACTERISTICS

Dual-Block is a cortico-cancellous block of xenogenic bone with osteoconductive characteristics. It can be used when the regeneration of big volumes is needed: thanks to the collagen content that promotes blood clotting and migration of regenerative and repairing cells⁽¹⁾, the graft offers an adequate support for tissue reconstruction and is gradually resorbed, while new bone is produced by osteoblasts.

HANDLING

Dual-Block must be hydrated before use with sterile lukewarm physiological solution or with antibiotics (5/10 minutes for Soft version; up to 40 minutes for Norm version). Afterwards, the block can be adapted to the receiving site which must be accurately decorticated in order to guarantee maximum contact; the block should always be fixed with osteosynthesis microscrews and protected with *Evolution* membrane.



SEM image of OsteoBiol® *Dual-Block*
Source: Politecnico di Torino, Italy

CLINICAL INFORMATION

Dual-Block can be grafted with the onlay technique only to augment horizontally heavily resorbed maxilla. The gaps around the block can be filled with a biomaterial in granules to achieve the desired volume and contour of the augmented recipient site.


The above clinical information is based on the experience of expert surgeons

BIBLIOGRAPHY


(1) MANESCU A, GIULIANI A, MOHAMMADI S, TROMBA G, MAZZONI S, DIOMEDE F, ZINI N, PIATTELLI A, TRUBIANI O
OSTEOGENIC POTENTIAL OF DUAL-BLOCKS CULTURED WITH HUMAN PERIODONTAL LIGAMENT STEM CELLS: IN VITRO AND SYNCHROTRON
J PERIODONTAL RES, 2016 Feb;51(1):112-24



OsteoBiol® *Dual-Block*
Source: TecnoSS® Dental Media Library



HORIZONTAL AUGMENTATION
onlay technique



Tissue of origin

Cortico-cancellous bone

Tissue collagen

Preserved

Physical form

Rigid dried block

Composition

Collagenated cortico-cancellous bone

Re-entry time

About 8 months, variable depending on characteristics and irradiation grade of grafting site and on clinical conditions of the patient

Packaging

Sterile blister

Product codes

STS7S | 20x15x5 mm | Soft | Porcine curved
STN5S | 20x10x5 mm | Norm | Porcine curved


GMDN code

46425

CND code


P900402

Histological results in alveolar regeneration








- **A)** Barone A et al. - **Flap versus flapless procedure for ridge preservation in alveolar extraction sockets: a histological evaluation in a randomized clinical trial**
Clinical Oral Implants Research, 2015 Jul;26(7):806-13
- **B)** Giuliani A et al. - **Regenerative properties of collagenated porcine bone grafts in human maxilla: demonstrative study of the kinetics by synchrotron radiation microtomography and light microscopy**
Clinical Oral Investigations, 2018 Jan;22(1):505-513
- **C)** Barone A et al. - **Clinical and histological changes after ridge preservation with two xenografts: preliminary results from a multicenter randomized controlled clinical trial**
Journal of Clinical Periodontology, 2017 Feb;44(2):204-214
- **D)** Crespi R et al. - **Corticocancellous porcine bone in the healing of human extraction sockets: combining histomorphometry with osteoblast gene expression profiles in vivo**
Int Journal of Oral and Maxillofacial Implants, 2011 Jul - Aug; 26(4):866-72
- **E)** Crespi R et al. - **Comparison of magnesium-enriched hydroxyapatite and porcine bone in human extraction socket healing: a histologic and histomorphometric evaluation**
Int Journal of Oral and Maxillofacial Implants, 2011 Sep-Oct;26(5):1057-62
- **F)** Barone A et al. - **Xenograft versus extraction alone for ridge preservation after tooth removal: a clinical and histomorphometric study**
Journal of Periodontology, 2008 Aug; 79(8):1370-7

Histological results in sinus lift



- **A)** Cassetta M et al. - **Bone formation in sinus augmentation procedures using autologous bone, porcine bone, and a 50 : 50 mixture: a human clinical and histological evaluation at 2 months**
Clinical Oral Implants Research, 2015 Oct; 26(10):1180-4
- **B)** Scarano A et al. - **Maxillary sinus augmentation in humans using cortical porcine bone: a histological and histomorphometrical evaluation after 4 and 6 months**
Clinical Implant Dentistry and Related Research, 2011 Mar; 13(1):13-18
- **C)** Orsini G et al. - **Histologic and ultrastructural analysis of regenerated bone in maxillary sinus augmentation using a porcine bone-derived biomaterial**
Journal of Periodontology, 2006 Dec; 77(12):1984-90
- **D)** Iezzi G et al. - **Comparative histological results of different biomaterials used in sinus augmentation procedures: a human study at 6 months**
Clinical Oral Implants Research, 2012 Dec;23(12):1369-76
- **E)** Tanaka K et al. - **Sinus floor elevation and antrostomy healing: a histomorphometric clinical study in humans**
Implant dentistry, 2019 Dec; 28(6):537-542

MEMBRANES AND BARRIERS



OsteoBiol® membranes and barriers

MEMBRANES

BARRIERS

Evolution

*Heterologous
mesenchymal tissue*



**Dried membrane with one
smooth side and one
micro-rough side**

For more information on OsteoBiol® Evolution
see page 62

Derma

Porcine derma



**Dried
membrane**

For more information on OsteoBiol® Derma
see page 66

Special

*Heterologous
pericardium*



**Translucent
dried membrane**

For more information on OsteoBiol® Special
see page 74

Lamina


Cortical bone




**Semi-rigid and rigid
dried lamina**

For more information on OsteoBiol® Lamina
see page 70

SEM image showing collagenic matrix of OsteoBiol® membranes
Source: Courtesy of Nobil Bio Ricerche, Villafranca D'Asti, Italy



Evolution



The natural Evolution of collagen membranes
Heterologous mesenchymal tissue

Characteristics and handling



Tissue of origin

Heterologous mesenchymal tissue

Tissue collagen

Preserved

Physical form

Dried membrane with one smooth side and one micro-rough side

Thickness

X-Fine: 0.2 mm
Fine: 0.3 mm
Standard: 0.4 mm

Estimated resorption time

X-Fine: about 2 months
Fine: about 3 months
Standard: about 4 months

Size

20x20 mm, 30x30 mm, 25x35 mm (oval), 40x40 mm, 80x60 mm

Product codes

EM33XS | 3 pcs | 30x30 mm | X-Fine | Porcine
EV02LLE | 20x20 mm | Fine | Equine
EV03LLE | 30x30 mm | Fine | Equine
EVOLLE | 25x35 mm (oval) | Fine | Equine
EV04LLE | 40x40 mm | Fine | Equine
EV06LLE | 80x60 mm | Fine | Equine
EM02HS | 20x20 mm | Standard | Porcine
EV02HHE | 20x20 mm | Standard | Equine
EM03HS | 30x30 mm | Standard | Porcine
EV03HHE | 30x30 mm | Standard | Equine
EM00HS | 25x35 mm (oval) | Standard | Porcine

GMDN code

47184

CND code

P900402

CHARACTERISTICS


Obtained from heterologous mesenchymal tissue, the *Evolution* membrane is gradually resorbable⁽¹⁾. Its structure is made of dense collagen fibers of high consistency and of extraordinary resistance that offer the specialist surgeon:

- maximum adaptability to bone tissue and soft tissues
- easy and secure suturability to nearby tissues
- best membrane-bone and membrane-periosteum interface
- stability and prolonged protection of the underlying graft
- clot stabilization and isolation⁽²⁾

HANDLING

The membrane can be shaped with sterile scissors until the desired size is reached; unless the grafting site is already bleeding, the membrane should be rehydrated with lukewarm physiological solution. Once it acquires the desired plasticity, it must be adapted to the grafting site.

NB: in case of accidental exposure, the dense collagenic matrix of *Evolution* protects the graft from infection; the membrane itself will also not be infected, allowing second intention healing⁽³⁻⁵⁾.



SEM image of an OsteoBiol® Evolution standard membrane

Source: Politecnico di Torino, Italy



Source: Tecnos® Dental Media Library




Evolution is obtained from heterologous mesenchymal tissue and is completely resorbable. Experimental studies have shown histological evidence of the prolonged barrier effect of this membrane, which lasts at least eight weeks⁽¹⁾, protecting the graft from external agents.

This property is particularly important in case of flapless regeneration⁽³⁾ of large posterior sockets⁽⁵⁾: in these cases, the standard model has proved to be the most effective.

In lateral access sinus lift, *Evolution* membranes have been documented for antrostomy coverage (standard model)^(6,7) and for protection of the sinus membrane from cutting risk due to graft pressure (fine model)⁽⁸⁾.


Evolution can be used to protect peri-implant regenerations⁽⁹⁾ and periodontal grafts⁽¹⁰⁾. Furthermore, *Evolution* fine has been successfully used to protect *Sp-Block* in vertical augmentation with the inlay technique⁽¹¹⁾.




LATERAL ACCESS SINUS LIFT
maxillary sinus floor augmentation




PERIODONTAL REGENERATION
intra-bony defects



HORIZONTAL AUGMENTATION
two-wall defects



DEHISCENCES AND FENESTRATIONS
peri-implant lesions



ALVEOLAR REGENERATION
graft protection




VERTICAL AUGMENTATION
inlay technique


BIBLIOGRAPHY

- (1) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70
- (2) KILINC A, ATAOL M
HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY
BMC ORAL HEALTH, 2017 OCT 5;17(1):126
- (3) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13
- (4) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7
- (5) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLIN ORAL INVEST, 2017 2018 JAN;22(1):505-513
- (6) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6
- (7) SCARANO A, PIATTELLI A, PERROTTI V, MANZON L, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRICAL EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR; 13(1):13-18
- (8) CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8
- (9) BARONE A, MARCONCINI S, GIAMMARINARO E, MIJIRITSKY E, GELPI F, COVANI U
CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY
CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112
- (10) ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2015;8(3):233-244
- (11) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

For further information see the complete literature on p. 92



Derma



A xenogenic matrix for soft tissue augmentation
Collagen dermal membrane

Characteristics and handling



Tissue of origin

Porcine derma

Tissue collagen

Preserved

Physical form

Dried membrane

Composition

100% derma

Thickness

0.5 mm (± 0.1 mm)

0.9 mm (± 0.1 mm)

2.0 mm (± 0.2 mm)

Estimated resorption time

0.5 mm: about 1 month

0.9 mm: about 3 months

2.0 mm: about 5 months

Size

0.5 mm: 20x20 mm

0.9 mm: 25x25 mm, 12x8 mm, 50x50 mm

2.0 mm: 7x5 mm, 15x5 mm, 30x30 mm, 50x50 mm

Product codes

ED02LS | 20x20 mm | 0.5 mm | Porcine

ED21FS | 12x8 mm | 0.9 mm | Porcine

ED25FS | 25x25 mm | 0.9 mm | Porcine

ED05FS | 50x50 mm | 0.9 mm | Porcine

ED75SS | 7x5 mm | 2.0 mm | Porcine

ED15SS | 15x5 mm | 2.0 mm | Porcine

ED03SS | 30x30 mm | 2.0 mm | Porcine

ED05SS | 50x50 mm | 2.0 mm | Porcine

GMDN code

47184

CND code


P900402

CHARACTERISTICS

Obtained from derma of porcine origin, using an exclusive Tecnos[®] process that preserves the natural collagen fibers⁽¹⁾, *Derma* membranes are gradually integrated⁽²⁾ with the autologous soft tissues. Their strong consistency and resistance allow a perfect stabilization and a prolonged protection of the underlying graft⁽³⁾ in socket regeneration procedures, together with a strong barrier action to guide the growth of epithelium and preventing its invagination.


HANDLING

Derma membrane can be shaped with scissors until the desired size is reached; then it must be thoroughly hydrated in sterile lukewarm physiological solution until the desired consistency is obtained. Once it acquires the desired plasticity, it must be adapted to the grafting site. It is always recommendable to prepare a pocket with an elevator in order to stabilize the membrane in the site after stitching the flaps.




SEM image of OsteoBiol[®] Derma

Source: Politecnico di Torino, Italy



4-months histology showing remnants of porcine dermal matrix surrounded by connective tissue. Experimental study on dogs. Htx-eosine

Author: Prof Ulf Nannmark, University of Göteborg, Sweden




Source: Tecnos[®] Dental Media Library

Derma membrane is a collagen resorbable matrix useful to augment soft tissues and to protect and stabilize bone grafting materials; only in this specific indication it can be used also in open healing⁽³⁾ situations due to its perfect tissue integration characteristics.

If a residual band of keratinized tissue is still present around teeth or implants, *Derma* membrane can be used as an alternative to connective tissue graft⁽²⁾ to improve the quality of keratinized tissues⁽⁴⁾. *Derma* has been also documented for horizontal soft tissue augmentation around implants^(5,8).


Mild gingival recessions^(6,7) can be treated with *Derma* to avoid patient morbidity and discomfort due to connective tissue graft harvesting. To avoid membrane exposure, usually *Derma* is completely covered by the coronally advanced flap. A properly shaped *Derma* membrane with rounded edges has been also documented for the tunnel technique⁽⁶⁾.



Graft protection using OsteoBiol® *Derma*



Positioning of OsteoBiol® *Derma* with the tunneling technique



OsteoBiol® *Derma* shaped for a gingival recession treatment




Treatment of a gingival recession using OsteoBiol® *Derma*
Author: Dr Roberto Rossi, Genova, Italy



SOFT TISSUE AUGMENTATION
soft tissue improvement



PERIODONTAL REGENERATION
gingival recessions



ALVEOLAR REGENERATION
graft protection

BIBLIOGRAPHY

(1) DE MARCO P, ZARA S, DE COLLI M, RADUNOVIC M, LAZOVIC V, ETTORRE V, DI CRESCENZO A, PIATELLI A, CATALDI A, FONTANA A
GRAPHENE OXIDE IMPROVES THE BIOCOMPATIBILITY OF COLLAGEN MEMBRANES IN AN IN VITRO MODEL OF HUMAN PRIMARY GINGIVAL FIBROBLASTS
BIOMED MATER, 2017 SEP 13;12(5):055005

(2) FICKL S, NANNMARK U, SCHLAGENHAUF U, HÜRZELER M, KEBSCHULL M
PORCINE DERMAL MATRIX IN THE TREATMENT OF DEHISCENCE-TYPE DEFECTS – AN EXPERIMENTAL SPLIT-MOUTH ANIMAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):799-805

(3) TALLARICO M, XHANARI E, PISANO M, DE RIU G, TULLIO A, MELONI SM
SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2016;9(3):263-275

(4) FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N
STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS
QUINTESSENCE INT, 2014 NOV-DEC;45(10):853-60


(5) FISCHER K R, TESTORI T, WACHTEL H, MÜHLEMANN S, HAPPE A, DEL FABBRO M, DEL FABBRO M
SOFT TISSUE AUGMENTATION APPLYING A COLLAGENATED PORCINE DERMAL MATRIX DURING SECOND STAGE SURGERY: A PROSPECTIVE MULTICENTER CASE SERIES
CLIN IMPLANT DENT RELAT RES, 2019;1-8

(6) FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
PORCINE DERMAL MATRIX FOR COVERING OF RECESION TYPE DEFECTS: A CASE SERIES
QUINTESSENCE INT, 2013;44(3):243-6


(7) MATOH U, PETELIN M, GASPERIC R
SPLIT-MOUTH COMPARISON OF CORONALLY ADVANCED FLAP WITH CONNECTIVE TISSUE GRAFT OR COLLAGEN MATRIX FOR TREATMENT OF ISOLATED GINGIVAL RECESIONS
INT J PERIODONTICS RESTORATIVE DENT, 2019;39(3):439-446

(8) VERARDI S, ORSINI M, LOMBARDI T, AUSENDA F, TESTORI T, PULICI A, OREGLIA F, VALENTE NA, STACCHI C
COMPARISON BETWEEN TWO DIFFERENT TECHNIQUES FOR PERI-IMPLANT SOFT TISSUE AUGMENTATION: PORCINE DERMAL MATRIX GRAFT VS. TENTING SCREW
J PERIODONTOL. 2020; ACCEPTED, IN PUBLICATION

For further information see the complete literature on p. 92



Lamina



A unique cortical bone barrier
Heterologous collagenated cortical bone

Characteristics and handling



Tissue of origin

Cortical bone

Tissue collagen

Preserved

Physical form

Lamina soft: semi-rigid flexible dried lamina

Lamina: rigid dried lamina, flexible after re-hydration

Composition

100% cortical bone

Thickness

0.5 mm (± 0.1 mm)

0.7 mm (± 0.1 mm)

1.0 mm (± 0.1 mm)

3.0 mm (± 1 mm)

Estimated re-entry time

0.5 mm: about 5 months

0.7 mm: about 7 months

1.0 mm: about 6 months

3.0 mm: about 8 months

Size

0.5 mm: 25x25 mm, 25x35 mm (oval)

0.7 mm: 35x15 mm

1.0 mm: 35x35 mm (Curved), 20x40 mm

3.0 mm: 30x30 mm

Product codes

LS25FS | 25x25 mm | 0.5 mm | soft | Porcine

LS25FE | 25x25 mm | 0.5 mm | soft | Equine

LS23FS | 25x35 mm (Oval) | 0.5 mm | soft | Porcine

LS23FE | 25x35 mm (Oval) | 0.5 mm | soft | Equine

LS24LS | 20x40 mm | 1.0 mm | soft | Porcine

LS10HS | 35x35 mm (Curved) | 1.0 mm | soft | Porcine

LS10HE | 35x35 mm (Curved) | 1.0 mm | soft | Equine

LS03SS | 30x30 mm | 3.0 mm | soft | Porcine

LS03SE | 30x30 mm | 3.0 mm | soft | Equine

LS15LS | 35x15 mm | 0.7 mm | rigid | Porcine

LS35LS | 35x35 mm | 1.0 mm | rigid | Porcine

GMDN code


46425

CND code

P900402

CHARACTERISTICS

Lamina barriers are made of cortical bone of heterologous origin produced with an exclusive TecnoSS® process that avoids the ceramization of hydroxyapatite crystals, thus allowing gradual resorption. After a process of superficial decalcification, *Lamina soft* acquires an elastic consistency, nevertheless maintaining the typical compactness of the bone tissue from which it originates; the margins are soft in order not to cause micro-traumas to the surrounding tissues. *Curved soft Lamina* has a semi-rigid consistency and should be grafted without hydration, provided that it is previously shaped to fit the defect morphology. *Rigid Lamina* undergoes a process of superficial semi-decalcification (50% vs *Lamina soft*) therefore increasing its consistency, typical of the cortical bone tissue⁽¹⁻²⁾.



SEM image of OsteoBiol® Lamina

Source: Politecnico di Torino, Italy

HANDLING

Lamina soft can be shaped with sterile scissors until the desired size is reached, then it must be hydrated for 5/10 minutes in sterile physiological solution. Once it acquires the desired plasticity, it must be adapted to the grafting site; it should always be immobilized either with titanium microscrews or sutured (fine model) directly to the surrounding tissues with a triangular section non-traumatic needle. ***Curved soft Lamina should not be hydrated*** in order to maintain its tenting effect but can also be shaped with sterile scissors, and must be fixated with osteosynthesis screws. In case of exposure, *Lamina* should only be removed if there is a clear suprainfection, because its consistency allows to achieve a complete second intention healing of the wound.



Source: TecnoSS® Dental Media Library



Lamina soft becomes flexible after hydration and can be shaped⁽³⁾ and adapted to the defect morphology creating, once fixated with osteosynthesis screws, a semi-rigid covering to the underlying graft⁽⁴⁻⁶⁾. This property is particularly useful when it is necessary to maintain the graft volume in aesthetic areas, as well as in horizontal augmentation⁽⁶⁻⁸⁾ of two wall defects and in lateral access sinus lift procedures^(5,9,10). *Lamina* can also be used in regenerations with risks of exposure. *Curved soft Lamina* has a ~1.0 mm thickness and must be directly grafted without hydration⁽¹¹⁾; it can be particularly effective in association with GTO® for regeneration of ridges with compromised buccal plate.


Rigid Lamina (ref. LS35LS) has been documented for orbital floor and wall reconstruction⁽¹²⁾.

The new 0.7 mm thickness rigid *Lamina* (ref. LS15LS) represents a viable alternative to autogenous cortical bone plates in the reconstruction of three-dimensional crestal defects with the shell technique.

The above clinical information is based on the experience of expert surgeons



OsteoBiol® Lamina positioning
Source: TecnoSS® Dental Media Library



HORIZONTAL AUGMENTATION
two-wall defects



HORIZONTAL AUGMENTATION
bone-layer technique




ORBITAL FLOOR RESTORATION

BIBLIOGRAPHY

- (1) DI CARLO R, ZARA S, VENTRELLA A, SIANI G, DA ROS T, IEZZI G, CATALDI A, FONTANA A
COVALENT DECORATION OF CORTICAL MEMBRANES WITH GRAPHENE OXIDE AS A SUBSTRATE FOR DENTAL PULP STEM CELLS
NANOMATERIALS, 2019;9:604
- (2) CABALLÉ-SERRANO J, MUNAR-FRAUJA A, DELGADO L, PÉREZ R, HERNÁNDEZ-ALFARO F
PHYSICO-CHEMICAL CHARACTERIZATION OF BARRIER MEMBRANES FOR BONE REGENERATION
J MECH BEHAV BIOMED, 2019;97:13-20
- (3) ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C
THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY
MINERVA STOMATOL, 2016 OCT;65(5):257-68
- (4) PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58
- (5) FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT AND RELAT RES, 2013 OCT;15(5):707-13
- (6) WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7
- (7) LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D
THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES
J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85
- (8) POLIS-YANES C, CADENAS-SEBASTIÁN C, GUAL-VAQUÉS P, AYUSO-MONTERO R, MARÍ-ROIG A, LÓPEZ-LÓPEZ J
GUIDED BONE REGENERATION OF AN ATROPHIC MAXILLA USING HETEROLOGOUS CORTICAL LAMINA
CASE REP DENT, 2019; 5216362
- (9) HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W
ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY
ORAL CRANIOFAC TISSUE ENG, 2011;1:188-197
- (10) SCARANO A, MURMURA G, MASTRANGELO F, LORUSSO F, GRECO LUCCHINA A, CARINCI F
A NOVEL TECHNIQUE TO PREVENT SINUS MEMBRANE COLLAPSE DURING MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT: TECHNICAL NOTE
J BIOL REGUL HOMEOST AGENTS, 2018 NOV-DEC;32(6):1589-1592
- (11) ROSSI R, FOCE E, SCOLAVINO S
THE CORTICAL LAMINA TECHNIQUE: A NEW OPTION FOR ALVEOLAR RIDGE AUGMENTATION. PROCEDURE, PROTOCOL, AND CASE REPORT
J LEBANESE DENTAL ASS, 2017 JAN-JUN; 52(1):35-41
- (12) RINNA C, REALE G, FORESTA E, MUSTAZZA MC
MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX
J CRANIOFAC SURG, 2009 MAY; 20(3):881-4

For further information see the complete literature on p. 92




Special



***A translucent membrane to
separate bone and soft tissues***

Engineered to protect hard and soft tissue grafts

Characteristics, handling and clinical information



Tissue of origin

Heterologous pericardium

Tissue collagen

Preserved

Physical form

Translucent dried membrane

Composition

100% pericardium

Thickness

0.2 mm

Resorption time

About 40 days

Size

20x20 mm, 30x30 mm

Product codes

EM02LE | 20x20 mm | Equine

EM03LE | 30x30 mm | Equine

GMDN code

47184

CND code

P900402

CHARACTERISTICS

Obtained from pericardium of heterologous origin, using an exclusive TecnoSS® process, the dried *Special* membranes are completely resorbable. Once hydrated, they become translucent and flexible, guiding the growth of epithelium and preventing its invagination: their action favors therefore an optimal regeneration of the underlying bone tissue.

HANDLING


The membrane can be shaped with sterile scissors until the desired size is reached; it must then be rehydrated with lukewarm physiological solution. Once it acquires the desired plasticity, it must be adapted to the grafting site. It is recommended to prepare a pocket with an elevator in order to stabilize the membrane in the site after stitching the flaps. If this is not possible, the membrane can be stabilized with envelope sutures which bridle it with the gingival flaps.

CLINICAL INFORMATION

In periodontology, the *Special* membrane can be used to protect and stabilize the biomaterial in the treatment of intrabony defects.

Special can be used to protect the sinus membrane before the insertion of the grafting material, to close sinus membrane perforations. Grafts placed in post-extractive sockets with closed healing procedure can also be protected with this membrane.

The above clinical information is based on the experience of expert surgeons



SEM images of OsteoBio® Special
Source: Nobil Bio Ricerche, Villafranca d'Asti, Italy



PERIODONTAL REGENERATION
intrabony defects



LATERAL ACCESS SINUS LIFT
Schneider membrane protection

Bone, Biomaterials & Beyond

Prof Antonio Barone, Prof Ulf Nannmark

CONTENTS

The introduction of osseointegrated dental implants soon 50 years ago has indeed revolutionized dentistry.


The scientific evaluation of their use has shown good and increasingly successful treatment outcomes.

A prerequisite though is the availability of sufficient bone volumes to ensure integration and acceptable aesthetic results.

In this book, various surgical techniques using different augmentation materials are described and explained.

The aim is to highlight minimally invasive surgical techniques, which lead to less risk of morbidity and reduce treatment time.

Readers will enjoy a comprehensive atlas providing some practical advice for every day surgical practice based on solid scientific evidence.



EDRA Editions

Publication: March 2014

Pages: 200

Images: 786 color

Binding: hardcover

Size: 21x29,7 cm

Euro: 59,00

ISBN 978.88.214.3758.8

For information and orders
please contact: ordiniedra@lswr.it

Available in the
following languages:

English
Italian
German
French
Spanish
Russian
Korean

Success through innovation:

Products launch



Gen-Os®



Apatos



Gel 40



Special



Lamina



Sp-Block



Dual-Block



Derma

Year

2000

2001

2002

2003

2004

2005

2006

2007

2008

2009

2010

Worldwide distribution countries

1

2

8

12

19

28

41

44

Publications on international journals

1

4

9

12

19

26

32

history of the **OsteoBiol**[®] brand by Tecross



Putty



Evolution



mp3[®]




TSV Gel



GTO[®]



Lamina
for shell technique



INNOVATION

A close-up photograph of a scientist wearing a white lab coat, a white surgical cap, and a white face mask. The scientist is looking through the eyepieces of a white and black microscope. The background is a blurred laboratory setting with various pieces of equipment. The word "INNOVATION" is overlaid in large, bold, blue capital letters on the left side of the image.

Tecnoss® bone vs human bone

Studies and researches have demonstrated that gold standard in bone regeneration is autologous bone^(1,2).

It is also well known, though, what disadvantages are related to the harvesting and grafting of autogenous bone⁽²⁻⁴⁾.

The goal of bone regeneration is to heal bone deficits with newly-formed quality tissue, in order to achieve a functional recovery and esthetics. To obtain these results, hundreds of studies have been conducted about the clinical performance of biomaterials. The examination of clinical results and the commercial diffusion of various kinds of products developed by the biomedical industry show a

clear superiority of products of natural origin over those of synthetic derivation.

The structure of animal bone is morphologically more similar to human bone than any synthesized product, the latter presenting a morphological pattern and properties artificially created, which differ in various ways from the structure of natural bone⁽⁵⁾.

Over the last thirty years several processes have been developed to allow the grafting of heterologous

origin products in the human body without adverse reaction^(6,7).

The first products developed through these technologies have shown encouraging clinical results, even if made of bone mineral matrix only.

The OsteoBiol® new generation of biomaterials, thanks to a revolutionary technology, goes beyond the simple role of aiding natural bone regrowth by stimulating and accelerating contact osteogenesis, with a behaviour similar to that of autogenous bone⁽⁸⁻¹⁰⁾.



BIBLIOGRAPHY

- (1) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC;77(12):1984-90
- (2) BARONE A, ALFONSI F, BORGIA V, IEZZI G, PIATTELLI A, COVANI U, TONELLI P
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS
CURR PHARM BIOTECHNOL, 2017;18(1):64-75
- (3) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZONI L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS
CURR PHARM BIOTECHNOL, 2017, 18, 33-44
- (4) BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25
- (5) FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
MATER SCI ENG C MATER BIOL APPL, 2013 AUG 1;33(6):3506-13
- (6) SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73
- (7) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAPTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30
- (8) CASSETTA M, RICCI L, IEZZI G, DELLAQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9
- (9) CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4
- (10) JEANNEAU C, LE FOURNIS C, ABOUT I
XENOGENIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT C5A
CLINICAL ORAL INVESTIGATIONS, 2019 OCT 23

Why xenografts?



Xenografts are the most used biomaterials worldwide.

This is because:

- tissues of origin are extremely safe and available in unlimited quantities
- xenogenic bone surface and porosity are extremely similar to autogenous bone
- there is no need to harvest autogenous bone in extraoral sites, with the related risk of morbidity and post-operative complications
- sterile xenografts are completely biocompatible and safe
- no adverse reactions after grafting deriving from biomaterial degradation
- easy to handle, quick learning curve
- collagenated xenografts enhance osteoblasts and osteoclasts activity
- wide scientific documentation
- excellent clinical performance
- storage can be done at room temperature
- long shelf life (5 years from production date)
- excellent price/quality ratio

Characteristics of Tecnos® process

Tecnos® has developed manufacturing processes for the treatment of tissues from various animal species, allowing to obtain the biocompatibility of these tissues, preserving at the same time their collagen matrix⁽¹⁾.

The protein components of animal tissues are determinant to make every individual unique. They activate the cells of the immune system of the receiving organism by interacting with receptors of the Major Histocompatibility Complex (MHC).

Their neutralization/denaturation allows collagen mineral bone matrix to be transferred from animal to man without any dangerous adverse reaction outbreak.

Successful Guided Bone Regeneration (GBR) depends both on stimulation of tissues involved in new bone formation and on the characteristics of grafted biomaterials, which can determine the quality of bone/graft interface. The development of OsteoBio® product line has thus been driven by the ideal

biomaterial concept: a material with the highest affinity to the new endogenous bone.

To pursue this aim, Tecnos® developed a biotechnology able to preserve the structure of natural hydroxyapatite, avoiding the high temperature ceramization phase, therefore allowing a bone turnover time of the grafted site similar to the one of the physiologic natural process⁽²⁾.

Thanks to this innovative technology, the OsteoBio® line has the following important characteristics:

1. Cell growth support and differentiation⁽³⁾
2. Absence of a foreign body response^(4,5)
3. Gradual resorption over time^(2,6)
4. Stimulation of the physiological tissue regeneration process^(7,8)
5. Protection of the grafting site from infection (membranes)^(5,9)
6. Capability of carrying medication to the surgical site⁽¹⁰⁾
7. Absorption and release over time of growth factors⁽¹¹⁾
8. Enhancement of endothelial cells proliferation⁽⁷⁾



BIBLIOGRAPHY

(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

(3) TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20(1 SUPPL 1):87-91

(4) BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN;24(1):1-6

(5) BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

(6) BARONE A, RICCI M, COVANI U, NANNMARK U, AZARMEHR I, CALVO GUIRADO JL
MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICOCANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2012 JUN;14(3):373-9

(7) ROMBOULTS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS
DENT MATER J, 2016 DEC 1;35(6):900-907

(8) JEANNEAU C, LE FOURNIS C, ABOUT I
XENOGENIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT C5A
CLINICAL ORAL INVESTIGATIONS, 2019 OCT 23

(9) BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

(10) FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING – A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42

(11) MJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B
PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR MATERIALS, 2017 SEP 8;10(9)

Collagen: a key factor for clinical success

BIBLIOGRAPHY

(1) FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H

PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE

J BIOMED MATER RES B APPL BIOMATER, 2010 FEB; 92(2):409-19

(2) ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A

HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL

J PERIODONTOL, 2006 DEC;77(12):1984-90

(3) RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C

ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION

CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

(4) FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M

VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT

INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR-APR;33(2):159-66

(5) IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZON L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V

MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS

CURR PHARM BIOTECHNOL, 2017,18,33-44

(6) MIZUNO M, FUJISAWA R, KUBOKI Y

TYPE I COLLAGEN-INDUCED OSTEOBLASTIC DIFFERENTIATION OF BONE-MARROW CELLS MEDIATED BY COLLAGEN-A2B1 INTEGRIN INTERACTION

J CELL PHYSIOL. 2000 AUG;184(2):207-13

(7) HSU FY, CHUEH SC, WANG YJ

MICROSPHERES OF HYDROXYAPATITE/RECONSTITUTED COLLAGEN AS SUPPORTS FOR OSTEOBLAST CELL GROWTH

BIOMATERIALS 1999, 20:1931-1936

(8) ABDELGAWAD ME, SØE K, ANDERSEN TL, MERRILL DM, CHRISTIANSEN P, KJÆRGAARD-ANDERSEN P, DELAISSE JM

DOES COLLAGEN TRIGGER THE RECRUITMENT OF OSTEOBLASTS INTO VACATED BONE RESORPTION LACUNAE DURING BONE REMODELING?

BONE, 2014 OCT;67:181-8

(9) ROMBOUITS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I

CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENEIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS

DENT MATER J, 2016 DEC 1;35(6):900-907

(10) JEANNEAU C, LE FOURNIS C, ABOUT I

XENOGENEIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT C5A

CLINICAL ORAL INVESTIGATIONS, 2019 OCT 23

Tecnoss[®] exclusive manufacturing process is able to neutralize the antigenic components present in heterologous bone achieving biocompatibility and preserving the collagen matrix inside the biomaterial. Moreover, the molecular structure of natural hydroxyapatite is not significantly altered thanks to the mild process temperature⁽¹⁾.

These characteristics of OsteoBiol[®] products allow a consistent bone neo-formation and a close contact between mature neo-formed bone and biomaterial particles⁽²⁻⁵⁾.

Collagen has a key role in bone regeneration process in that:

- it acts as a valid substrate for platelet activation and aggregation
- it serves to attract and differentiate the mesenchymal stem cells present in the bone marrow⁽⁶⁾
- it increases the proliferation rate of the

osteoblasts up to 2/3 times⁽⁷⁾


• it stimulates the activation of the platelets, osteoblasts and osteoclasts in the bone healing process⁽⁸⁾.

The presence of collagen inside each granule makes OsteoBiol[®] Gen-Os[®] hydrophilic and facilitates further mixing with collagen gel and TSV Gel.


This technology has permitted the development of several versatile and innovative products: OsteoBiol[®] GTO[®], OsteoBiol[®] mp3[®], OsteoBiol[®] Putty and OsteoBiol[®] Gel 40. Their consistency allows an ideal filling of bone defects and guarantees simple handling and fast application.

The OsteoBiol[®] new generation of biomaterials, thanks to a revolutionary technology, goes beyond the simple role of aiding natural bone regrowth by stimulating and accelerating this vital physiological process^(9,10).

Composition of **OsteoBiol[®] Gen-Os[®]**



Source: University of Duisburg-Essen, Germany



Guided bone regeneration (GBR) is necessary to treat bone deficits due to lesions or bacterial infections.

Bone defect recovery occurs through the general mechanisms of tissue healing: complex dynamic mechanisms directed towards the repair of tissue function and anatomic integrity.

The discovery of the events pathway leading to tissue healing has helped to clearly identify the main actors in bone healing process; the concomitant presence of the following three components is necessary for the formation of “*de novo*” bone tissue:

- the platelets represent the principal actors during the first phase of the healing process, when, subsequent to a lesion, an initial deposition of fibrin and the formation of blood clot take place. This phase is characterized by significant activation of the chemical signals mediated by cytokines and growth factors.

In fact, the primary post-haemorrhagic clot formation process through platelet aggregation and lysis causes the release of both the coagulation cascade factors and growth factors, such as PDGF, IGF 1, IGF 2 and VEGF which are known for their activating effect on osteoblasts and osteoclasts, and TGF- β (Bone Morphogenetic Proteins belong to this superfamily) which starts bony callus formation.

- the osteoblastic precursors deriving from bone marrow mesenchymal stem cells are responsible, after cell differentiation in osteoblasts, for the second phase of the healing process (enchondral and/or intramembranous ossification) thanks to the synthesis of collagen and other components of the

extracellular matrix.

- an insoluble substrate, suitable carrier for osteoinductive signal and able to support and guide new bone tissue formation. Sampath and Reddi (1980) demonstrated crosslinked type I collagen to be the most appropriate carrier for promoting osteoinductive signal activity. The continuous progresses in comprehension of biological mechanisms regulating bone tissue morphogenesis can be exploited also for elaboration of natural or artificial products able to restore or maintain the function of damaged tissues and organs (tissue engineering)⁽¹⁻³⁾.

In vitro studies demonstrated that heterologous collagen is able to induce differentiation of mesenchymal osteoprogenitor stem cells into osteoblasts⁽⁴⁾, and that association of collagen type I with a scaffold of hydroxyapatite significantly enhances osteoblasts proliferation rate.

This important scientific evidence provides the rationale behind OsteoBioI[®] product line: a complete series of biomaterials with collagen base.

Collagen, in addition to its well-known structural action carried on connective tissues, is endowed with the following important properties, useful in tissue reparation processes:

1. Haemostasis

Collagen is able to activate the receptors present on cellular membranes of platelets, responsible for their aggregation and lysis process; moreover, during the first week, it reinforces the action of fibrin in the formation of the primary clot, and then, in the second week, it replaces the

function of fibrin.

2. Debridement

Collagen has a chemotactic action on monocyte/macrophage cell lines, from which osteoclasts derive; these cells, through their action on mineral component resorption of both bone tissue and OsteoBioI[®] biomaterials, can draw, activate and collaborate with osteoblasts in bone rearranging and remodeling.

3. Angiogenesis

The drawn monocytes/macrophages, in their turn, stimulate both osteoblastic activity and angiogenesis process in grafting site.

4. Osteoblastic activity

Collagen, binding to fibronectin, promotes the anchorage of mesenchymal stem progenitors, on which it exerts its chemotactic action, and induces differentiation into osteoblasts^(4,5).

5. Receiving site remodeling


Exogenous collagen grafting can contribute in decreasing remodeling times of immature bone tissue.

6. Osteoconduction and guided regeneration

Naturally integrated with mineral component, collagen is able to increase osteoblasts proliferation rate while as a resorbable membrane it is able to guide connective tissue regeneration.

BIBLIOGRAPHY

- (1) GRIFFITH LG, NAUGHTON G
TISSUE ENGINEERING-CURRENT CHALLENGES AND EXPANDING OPPORTUNITIES
SCIENCE 2002, 295:1009-14
- (2) REDDI AH
MORPHOGENESIS AND TISSUE ENGINEERING OF BONE AND CARTILAGE: INDUCTIVE SIGNALS, STEM CELLS, AND BIOMIMETIC BIOMATERIALS
TISSUE ENG 2000, 6(4):351-59
- (3) NAKASHIMA N, REDDI AH
THE APPLICATION OF BONE MORPHOGENETIC PROTEINS TO DENTAL TISSUE ENGINEERING
NAT BIOTECHNOL 2003, 9:1025-32
- (4) SALASZNYK RM, WILLIAMS WA, BOSKEY A, BATORSKY A, PLOPPER GE
ADHESION TO VITRONECTIN AND COLLAGEN I PROMOTES OSTEOGENIC DIFFERENTIATION OF HUMAN MESENCHYMAL STEM CELLS
J BIOMED BIOTECHNOL 2004, 1:24-34
- (5) BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL[®] INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
EUR J INFLAMMAT, 2011, VOL. 9, NO. 3(S), 103-107



From heterologous bone to biomaterial

RESULTS OF **INORGANIC** CHEMICAL ANALYSES PERFORMED ON OSTEObIOL® GEN-OS®

Chemical element	OsteoBiol® Gen-OS® (% in weight)	
Ca	25.7%	Mineral component 73.6%
PO ₄ ³⁻	35.2%	
C	13.6%	
H	2.2%	
N	2.9%	
O (not in PO ₄ ³⁻)	20.4%	
TOTAL	100.0%	Organic matrix 22.4%
Ca/P (n:n)	1.73	
		Water 4.0%

Inorganic chemical analyses results
Source: University of Duisburg-Essen, Germany

RESULTS OF **ORGANIC** CHEMICAL ANALYSES PERFORMED ON OSTEObIOL® GEN-OS®



“The separated proteins (one lane) were fractionated in ten portions and analysed with nano-LC-ESI MS/MS. In the fractions 1-5 in the range from 20-200kDa we found ONLY COLLAGEN. In the fractions 6-10 we identify NO PROTEIN”

Organic chemical analyses results
Source: Proteome Factory, Germany

A biomaterial for the reconstruction of bone defects must be biocompatible and have good handling and modeling properties; in specific clinical situations, it must also provide sufficient mechanical resistance. Tecnos® laboratories are specialized in processing heterologous bony and collagenic tissues. OsteoBiol® bone process, in particular, has been developed to modify while preserving the original collagen matrix of heterologous tissue, in order to maintain its positive biological functions, and complete biocompatibility⁽¹⁾. Most biomaterials are inert products that do not interfere, or rather, do not take

part in the physiology of bone remodeling: since they have been developed according to the sole concept of biocompatibility, their function is limited to the preservation of the graft volume (scaffold). The biocompatibility concept has an essential purpose in the implant of permanent prosthetic elements inside the human body, but it is extremely restrictive in case of materials used for bone reconstruction. OsteoBiol® biomaterials, being gradually resorbed and replaced by abundant newly formed bone over time, create the ideal conditions for the osseointegration of dental implants at re-entry⁽²⁾.

(1) NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70.

(2) CHECCHI V, FELICE P, ZUCCELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL
EUR J ORAL IMPLANTOL, 2017;10(3):263-278



OsteoBiol®: the most complete products range



The extensive OsteoBiol® range of products are engineered to help surgeons making the right decision when it comes to choose the perfect product for a specific clinical indication, both in dental and maxillofacial surgery.

Tecnoss® development of new products or improvement of existing ones, is focussed on supporting the technical capabilities of the practitioner to improve both intraoperative techniques and clinical results.

Specialists and researchers share their experience, blending clinical background and hands-on experience with the most advanced bio-technologies: the main goal is to obtain a specific solution to satisfy each clinical need.

OsteoBiol® collagenated grafting materials contribute to mineral deposition, vascular ingrowth and growth factor binding, thus providing a favourable environment for bone regeneration. The scientific literature has demonstrated that OsteoBiol® bone matrix is similar to human bone, and it has been reported in humans to be osteoconductive, well integrated in the host site and partially resorbed after 3-6 months, with no signs of adverse reaction⁽¹⁾.

(1) GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY
CLIN. ORAL INVESTIG., 2017 2018 JAN;22(1):505-513

CERTIFICATIONS





APPROVAL CERTIFICATE Bone Matrix
Source: TecnoSS® s.r.l.



APPROVAL CERTIFICATE Membranes
Source: TecnoSS® s.r.l.



UNI EN ISO 13485 KIWA CERMET quality certificate
Source: TecnoSS® s.r.l.



EC DESIGN-EXAMINATION CERTIFICATE Bone Matrix
Source: TecnoSS® s.r.l.



EC DESIGN-EXAMINATION CERTIFICATE Membranes
Source: TecnoSS® s.r.l.



FREE APP
for smartphone, iPhone, tablet
and iPad including:

8 animations to show your patients
the main GBR techniques

Information about the full range of
OsteoBiol® biomaterials

Direct access to the database of videos
and cases on osteobiol.com



This App may be too large to download over a mobile connection,
or may exceed data usage limits. Wi-Fi connection recommended.





• international distributors

Collagenated biomaterials

Distributed in 81 countries
200 international scientific publications
20 years of clinical success
Over 1.000.000 surgeries performed

OsteoBiol[®]
by Tecnos



REGENERATION SCIENCE

INSPIRED BY NATURE

LITERATURE



OsteoBiol®
by TecnoSS

SCIENTIFIC ABSTRACTS

REGENERATION SCIENCE

INSPIRED BY NATURE

Evolution:

LATERAL ACCESS SINUS LIFT

- Michela Cassella
- Laura Ricci
- Giovanna Izzzi
- Sabrina Calosso
- Adriano Pignelli
- Victoria Perrotti

ORIGINAL ARTICLE

International Journal of Periodontics & Endodontics
September 2012
2012; 2(4): 102-6

filled with
SUBSTITUTE
BioBiol® Gen-Ox

LASL

filled with
SUBSTITUTE
BioBiol® Gen-Ox

LAE

filled with
SUBSTITUTE
BioBiol® mp2®
MEMBRANE
OsteoBiol® Evolution

ALR

COVANI U, AMERI S, CRESPI R, BARONE A
PRESERVAZIONE DEL PROCESSO ALVEOLARE CON OSSO ETEROLOGO. CONSIDERAZIONI ISTOLOGICHE
ITALIAN ORAL SURGERY, 2004, VOL 3, 1: 17-23

CASSETTA M, CALASSO S, VOZZA I, DELL'AQUILA D
REHABILITATION OF ATROPHIC ALVEOLAR CRESTS WITH CYLINDRICAL SANDBLASTED AND ACID ETCHED IMPLANTS: A PILOT STUDY
EUR J IMPLANT PROSTHODONTICS, 2005;(3)1:133-144

ARCURI C, CECCHETTI F, GERMANO F, MOTTA A, SANTACROCE C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
MINERVA STOMATOL, 2005 JUN;54(6):351-62

BARONE A, CRESPI R, ALDINI NN, FINI M, GIARDINO R, COVANI U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
INT J ORAL MAXILLOFAC IMPLANTS, 2005 JUL-AUG; 20(4):519-25

RINNA C, UNGARI C, SALTARELLA A, CASSONI A, REALE G
ORBITAL FLOOR RESTORATION
J CRANIOFAC SURG, 2005 NOV; 16(6):968-72

BARONE A, AMERI S, COVANI U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
EUR J IMPLANT PROSTHODONTICS, 2006,2: 99-106

BARONE A, SANTINI S, SBORDONE L, CRESPI R, COVANI U
A CLINICAL STUDY OF THE OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION
INT J ORAL MAXILLOFAC IMPLANTS, 2006 JAN-FEB; 21(1):81-5

COVANI U, BARONE A, CORNELINI R, CRESPI R
CLINICAL OUTCOME OF IMPLANTS PLACED IMMEDIATELY AFTER IMPLANT REMOVAL
J PERIODONTOL, 2006 APR;77(4):722-7

ORSINI G, SCARANO A, PIATTELLI M, PICCIRILLI M, CAPUTI S, PIATTELLI A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
J PERIODONTOL, 2006 DEC;77(12):1984-90

TRUBIANI O, SCARANO A, ORSINI G, DI IORIO D, D'ARCANGELO C, PICCIRILLI M, SIGISMONDO M, CAPUTI S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
INT J IMMUNOPATHOL PHARMACOL, 2007 JAN-MAR; 20 (1 SUPPL 1):87-91

BARONE A, COVANI U
MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NON- VASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS
J ORAL MAXILLOFAC SURG, 2007 OCT;65(10):2039-46

DEL CORSO M
SOFT TISSUE RESPONSE TO PLATELET RICH FIBRIN: CLINICAL EVIDENCES
COSMETIC DENTISTRY, 2008, 3:16-20

BARONE A, SANTINI S, MARCONCINI S, GIACOMELLI L, GHERLONE E, COVANI U
OSTEOTOMY AND MEMBRANE ELEVATION DURING THE MAXILLARY SINUS AUGMENTATION PROCEDURE. A COMPARATIVE STUDY: PIEZOELECTRIC DEVICE VS. CONVENTIONAL ROTATIVE INSTRUMENTS
CLIN ORAL IMPLANTS RES, 2008 MAY;19(5):511-5

BARONE A, CORNELINI R, CIAGLIA R, COVANI U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
INT J PERIODONTICS RESTORATIVE DENT, 2008 JUN; 28(3):283-9

BARONE A, ALDINI NN, FINI M, GIARDINO R, CALVO GUIRADO JL, COVANI U
XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY
J PERIODONTOL, 2008 AUG;79(8):1370-7

COVANI U, CORNELINI R, BARONE A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
INT J ORAL MAXILLOFAC IMPLANTS, 2008 SEP-OCT; 23(5):841-6

CARDAROPOLI D, CARDAROPOLI G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
INT J PERIODONTICS RESTORATIVE DENT, 2008 OCT; 28(5):469-77

NANNMARK U, SENNERBY L
THE BONE TISSUE RESPONSE TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2008 DEC;10(4):264-70

SCARANO A, PIATTELLI M, CARINCI F, PERROTTI V
REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT
J OSSEOINTEGR, 2009;1(1):35-40

COVANI U, MARCONCINI S, CRESPI R, BARONE A
IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT
J ORAL IMPLANTOL, 2009; 35(4):189-95

FIGUEIREDO M, HENRIQUES J, MARTINS G, GUERRA F, JUDAS F, FIGUEIREDO H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
J BIOMED MATER RES B APPL BIOMATER, 2010FEB; 92(2):409-19

GRENGA PL, REALE G, COFONE C, MEDURI A, CERUTI P, GRENGA R
HESS AREA RATIO AND DIPLOPIA: EVALUATION OF 30 PATIENTS UNDERGOING SURGICAL REPAIR FOR ORBITAL BLOW-OUT FRACTURE
OPHTHAL PLAST RECONSTR SURG, 2009 MAR-APR; 25(2):123-5

CRESPI R, CAPPARÈ P, GHERLONE E
DENTAL IMPLANTS PLACED IN EXTRACTION SITES GRAFTED WITH DIFFERENT BONE SUBSTITUTES: RADIOGRAPHIC EVALUATION AT 24 MONTHS
J PERIODONTOL, 2009 OCT; 80(10):1616-1621

RINNA C, REALE G, FORESTA E, MUSTAZZA MC
MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX
J CRANIOFAC SURG, 2009 MAY; 20(3): 881-4

CARDAROPOLI D, CARDAROPOLI G
HEALING OF GINGIVAL RECESSIONS USING A COLLAGEN MEMBRANE WITH A THE MINERALIZED XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
INT J PERIODONTICS RESTORATIVE DENT, 2009 FEB; 29(1):59-67

NANNMARK U, AZARMEHR I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
CLIN IMPLANT DENT RELAT RES, 2010 JUN 1; 12(2):161-3

SCARANO A, PIATTELLI A, ASSENZA B, QUARANTA A, PERROTTI V, PIATTELLI M, IEZZI G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
J ORAL MAXILLOFAC SURG, 2010 AUG; 68(8):1869-73

ROSSI R, MORALES RS, FRASCARIA M, BENZI R, SQUADRITO N
PLANNING IMPLANTS IN THE ESTHETIC ZONE USING A NEW IMPLANT 3D NAVIGATION SYSTEM
EUR J ESTHETIC DENT, 2010 SUMMER; 5(2):172-88

SCARANO A, CARINCI F, ASSENZA B, PIATTELLI M, MURMURA G, PIATTELLI A
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE WITH A XENOGRAFT WITHOUT MINISCREWS AND MINIPLATES: CASE SERIES
CLIN ORAL IMPLANTS RES, 2011 OCT;22(10):1125-30

PAGLIANI L, ANDERSSON P, LANZA M, NAPPO A, VERROCCHI D, VOLPE S, SENNERBY L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
CLIN IMPLANT DENT RELAT RES, 2012 OCT;14(5):746-58

SANTAGATA M, GUARINIELLO L, TARTARO G
A MODIFIED EDENTULOUS RIDGE EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT
J ORAL IMPLANTOL, 2011 MAR;37 SPEC N°:114-9

SCARANO A, PIATTELLI A, PERROTTI V, MANZON I, IEZZI G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRIC EVALUATION AFTER 4 AND 6 MONTHS
CLIN IMPLANT DENT RELAT RES, 2011 MAR;13(1):13-18

CRESPI R, CAPPARÈ P, ROMANOS GE, MARIANI E, BENASCIIUTTI E, GHERLONE E
CORTICOCANCELLOUS PORCINE BONE IN THE HEALING OF HUMAN EXTRACTION SOCKETS: COMBINING HISTOMORPHOMETRY WITH OSTEOBLAST GENE EXPRESSION PROFILES IN VIVO
INT J ORAL MAXILLOFAC IMPLANTS, 2011 JUL-AUG; 26(4):866-72

HINZE M, VRIELINCK L, THALMAIR T, WACHTEL H, BOLZ W
ZYGOMATIC IMPLANT PLACEMENT IN CONJUNCTION WITH SINUS BONE GRAFTING: THE "EXTENDED SINUS ELEVATION TECHNIQUE". A CASE-COHORT STUDY
ORAL CRANIOFAC TISSUE ENG, 2011; 1:188-197

IEZZI G, DEGIDI M, PIATTELLI A, MANGANO C, SCARANO A, SHIBLI JA, PERROTTI V
COMPARATIVE HISTOLOGICAL RESULTS OF DIFFERENT BIOMATERIALS USED IN SINUS AUGMENTATION PROCEDURES: A HUMAN STUDY AT 6 MONTHS
CLIN ORAL IMPLANTS RES, 2012 DEC;23(12):1369-76

SLOTTE C, LINDFORS N, NANNMARK U
SURGICAL RECONSTRUCTION OF PERI-IMPLANT BONE DEFECTS WITH PREHYDRATED AND COLLAGENATED PORCINE BONE AND COLLAGEN BARRIERS: CASE PRESENTATIONS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):714-23

BARONE A, RICCI M, GRASSI RF, NANNMARK U, QUARANTA A, COVANI U
A 6-MONTH HISTOLOGICAL ANALYSIS ON MAXILLARY SINUS AUGMENTATION WITH AND WITHOUT USE OF COLLAGEN MEMBRANES OVER THE OSTEOTOMY WINDOW: RANDOMIZED CLINICAL TRIAL
CLIN ORAL IMPLANTS RES, 2013 JAN; 24(1):1-6

SANTAGATA M, GUARINIELLO L, RAUSO R, TARTARO G
IMMEDIATE LOADING OF DENTAL IMPLANT AFTER SINUS FLOOR ELEVATION WITH OSTEOTOME TECHNIQUE: A CLINICAL REPORT AND PRELIMINARY RADIOGRAPHIC RESULTS
J ORAL IMPLANTOL, 2010 DEC; 36(6):485-489

FESTA VM, ADDABBO F, LAINO L, FEMIANO F, RULLO R
PORCINE-DERIVED XENOGRAFT COMBINED WITH A SOFT CORTICAL MEMBRANE VERSUS EXTRACTION ALONE FOR IMPLANT SITE DEVELOPMENT: A CLINICAL STUDY IN HUMANS
CLIN IMPLANT DENT RELAT RES, 2013 OCT;15(5):707-13

RAMIREZ FERNANDEZ MP, CALVO GUIRADO JL, MATÉ SANCHEZ DE VAL JE, DELGADO RUIZ RA, NEGRI B, BARONA DORADO C
ULTRASTRUCTURAL STUDY BY BACKSCATTERED ELECTRON IMAGING AND ELEMENTAL MICROANALYSIS OF BONE-TO-BIOMATERIAL INTERFACE AND MINERAL DEGRADATION OF PORCINE XENOGRAFTS USED IN MAXILLARY SINUS FLOOR ELEVATION
CLIN ORAL IMPLANTS RES, 2013 MAY;24(5):523-30

CASSETTA M, RICCI L, IEZZI G, DELL'AQUILA D, PIATTELLI A, PERROTTI V
RESONANCE FREQUENCY ANALYSIS OF IMPLANTS INSERTED WITH A SIMULTANEOUS GRAFTING PROCEDURE: A 5-YEAR FOLLOW-UP STUDY IN MAN
INT J PERIODONTICS RESTORATIVE DENT, 2012 OCT;32(5):581-9

BARONE A, ORLANDO B, CINGANO L, MARCONCINI S, DERCHI G, COVANI U
A RANDOMIZED CLINICAL TRIAL TO EVALUATE AND COMPARE IMPLANTS PLACED IN AUGMENTED VS. NON-AUGMENTED EXTRACTION SOCKETS. A 3-YEAR EVALUATION
J PERIODONTOL, 2012 JUL;83(7):836-46

ESPOSITO M, CANNIZZARO G, SOARDI E, PISTILLI R, PIATTELLI M, CORVINO V, FELICE P
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM-LONG, 4 MM-WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2012 SPRING;5(1):19-33

FELICE P, PIANA L, CHECCHI L, PISTILLI R, PELLEGRINO G
VERTICAL RIDGE AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH A 2-STAGE INLAY TECHNIQUE: A CASE REPORT
 IMPLANT DENT, 2012 JUN;21(3):190-5

BARONE A, RICCI M, TONELLI P, SANTINI S, COVANI U
TISSUE CHANGES OF EXTRACTION SOCKETS IN HUMANS: A COMPARISON OF SPONTANEOUS HEALING VS. RIDGE PRESERVATION WITH SECONDARY SOFT TISSUE HEALING
 CLIN ORAL IMPLANTS RES, 2013 NOV;24(11):1231-7

CASSETTA M, RICCI L, IEZZI G, CALASSO S, PIATTELLI A, PERROTTI V
USE OF PIEZOSURGERY DURING MAXILLARY SINUS ELEVATION: CLINICAL RESULTS OF 40 CONSECUTIVE CASES
 INT J PERIODONTICS RESTORATIVE DENT, 2012 DEC;32(6):E182-8

BRUNELLI G, SOLLAZZO V, CARINCI F, PALMIERI A, GIRARDI A, MONGUZZI R
OSTEOBIOL® INFLUENCES OSTEOGENIC DIFFERENTIATION OF ADIPOSE DERIVED STEM CELLS
 EUR J INFLAMMAT, 2011, VOL. 9, NO. 3 (S), 103-107

FELICE P, PIANA L, CHECCHI L, CORVINO V, NANNMARK U, PIATTELLI M
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE WITH AN INLAY TECHNIQUE AND CANCELLOUS EQUINE BONE BLOCK: A CASE REPORT
 INT J PERIODONTICS RESTORATIVE DENT, 2013 MAR;33(2):159-66

FICKL S, JOCKEL-SCHNEIDER Y, LINCKE T, BECHTOLD M, FISCHER KR, SCHLAGENHAUF U
PORCINE DERMAL MATRIX FOR COVERING OF RECESSION TYPE DEFECTS: A CASE SERIES
 QUINTESSENCE INT, 2013;44(3):243-6

SILVESTRI M, MARTEGANI P, D'AVENIA F, FARNETI M, CAPRI D, PAOLANTONI G, LANDI L
SIMULTANEOUS SINUS AUGMENTATION WITH IMPLANT PLACEMENT: HISTOMORPHOMETRIC COMPARISON OF TWO DIFFERENT GRAFTING MATERIALS. A MULTICENTER DOUBLE-BLIND PROSPECTIVE RANDOMIZED CONTROLLED CLINICAL TRIAL
 INT J ORAL MAXILLOFAC IMPLANTS, 2013 MAR-APR;28(2):543-9

WACHTEL H, FICKL S, HINZE M, BOLZ W, THALMAIR T
THE BONE LAMINA TECHNIQUE: A NOVEL APPROACH FOR LATERAL RIDGE AUGMENTATION - A CASE SERIES
 INT J PERIODONTICS RESTORATIVE DENT, 2013 JUL-AUG;33(4):491-7

RODRIGUEZ JG, ELIDIBANY RM
VERTICAL SPLITTING OF THE MANDIBULAR BODY AS AN ALTERNATIVE TO INFERIOR ALVEOLAR NERVE LATERALIZATION
 INT J ORAL MAXILLOFAC SURG, 2013 SEP;42(9):1060-6

FIGUEIREDO A, COIMBRA P, CABRITA A, GUERRA F, FIGUEIREDO M
COMPARISON OF A XENOGENEIC AND AN ALLOPLASTIC MATERIAL USED IN DENTAL IMPLANTS IN TERMS OF PHYSICO-CHEMICAL CHARACTERISTICS AND IN VIVO INFLAMMATORY RESPONSE
 MATER SCI ENG C, MATER BIOL APP, 2013 AUG 1;33(6):3506-13

FELICE P, PISTILLI R, PIATTELLI M, SOARDI E, CORVINO V, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. PRELIMINARY RESULTS FROM A RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, SUMMER;5(2):149-61

TRAINI T, PIATTELLI A, CAPUTI S, DEGIDI M, MANGANO C, SCARANO A, PERROTTI V, IEZZI G
REGENERATION OF HUMAN BONE USING DIFFERENT BONE SUBSTITUTE BIOMATERIALS
 CLIN IMPLANT DENT RELAT RES, 2015 FEB;17(1):150-62

KOLMAS J, SZWAJA M, KOLODZIEJSKI W
SOLID-STATE NMR AND IR CHARACTERIZATION OF COMMERCIAL XENOGENEIC BIOMATERIALS USED AS BONE SUBSTITUTES
 J PHARM BIOMED ANAL, 2012 MAR 5;61:136-41

PISTILLI R, FELICE P, PIATTELLI M, GESSAROLI M, SOARDI E, BARAUSSE C, BUTI J, CORVINO V, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2013 WINTER;6(4):343-357

FICKL S, NANNMARK U, SCHLAGENHAUF U, HÜRZELER M, KEBSCHULL M
PORCINE DERMAL MATRIX IN THE TREATMENT OF DEHISCENCE-TYPE DEFECTS - AN EXPERIMENTAL SPLIT-MOUTH ANIMAL TRIAL
 CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):799-805

BARONE A, BORGIA V, COVANI U, RICCI M, PIATTELLI A, IEZZI G
FLAP VERSUS FLAPLESS PROCEDURE FOR RIDGE PRESERVATION IN ALVEOLAR EXTRACTION SOCKETS: A HISTOLOGICAL EVALUATION IN A RANDOMIZED CLINICAL TRIAL
 CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):806-13

BARONE A, RICCI M, ROMANOS GE, TONELLI P, ALFONSI F, COVANI U
BUCCAL BONE DEFICIENCY IN FRESH EXTRACTION SOCKETS: A PROSPECTIVE SINGLE COHORT STUDY
 CLIN ORAL IMPLANTS RES, 2015 JUL;26(7):823-30

BARONE A, TOTI P, QUARANTA A, DERCHI G, COVANI U
THE CLINICAL OUTCOMES OF IMMEDIATE VERSUS DELAYED RESTORATION PROCEDURES ON IMMEDIATE IMPLANTS: A COMPARATIVE COHORT STUDY FOR SINGLE-TOOTH REPLACEMENT
 CLIN IMPLANT DENT RELAT RES, 2015 DEC;17(6):1114-26

CASSETTA M, PERROTTI V, CALASSO S, PIATTELLI A, SINJARI B, IEZZI G
BONE FORMATION IN SINUS AUGMENTATION PROCEDURES USING AUTOLOGOUS BONE, PORCINE BONE, AND A 50 : 50 MIXTURE: A HUMAN CLINICAL AND HISTOLOGICAL EVALUATION AT 2 MONTHS
 CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1180-4

GHENO E, PALERMO A, BUFFOLI B, RODELLA LF
THE EFFECTIVENESS OF THE USE OF XENOGENEIC BONE BLOCKS MIXED WITH AUTOLOGOUS CONCENTRATED GROWTH FACTORS (CGF) IN BONE REGENERATION TECHNIQUES
 J OSSEOINTEGRATION 2014;6(2):37-42

FALISI G, GALLI M, VITTORINI-VELASQUEZ P, GALLEGOS-RIVERA JC, MINASI R, DE BIAZE A, DI PAOLO C
USE OF 3D CARTILAGE SCAFFOLDS FOR THE STABILIZATION OF IMPLANTS AND BONE REGENERATION WITH THE FIT-LOCK TECHNIQUE
 ACTA ODONTOL LATINOAM 2013;26(3):167-172

FISCHER KR, FICKL S, MARDAS N, BOZEC L, DONOS N
STAGE-TWO SURGERY USING COLLAGEN SOFT TISSUE GRAFTS: CLINICAL CASES AND ULTRASTRUCTURAL ANALYSIS
 QUINTESSENCE INT, 2014 NOV-DEC; 45(10):853-60

FISCHER KR, STAVROPOULOS A, CALVO GUIRADO JL, SCHNEIDER D, FICKL S
INFLUENCE OF LOCAL ADMINISTRATION OF PAMIDRONATE ON EXTRACTION SOCKET HEALING - A HISTOMORPHOMETRIC PROOF-OF-PRINCIPLE PRE-CLINICAL IN VIVO EVALUATION
 CLIN ORAL IMPLANTS RES, 2015 OCT;26(10):1135-42


SCARANO A, MURMURA G, SINJARI B, ASSENZA B, SOLLAZZO V, SPINELLI G, CARINCI F
EXPANSION OF THE ALVEOLAR BONE CREST WITH ULTRASONIC SURGERY DEVICE: CLINICAL STUDY IN MANDIBLE
 INT J IMMUNOPATHOL PHARMACOL, 2011 APR-JUN; 24(2 SUPPL):71-5

SCARANO A, PIATTELLI A, MURMURA G, IEZZI G, ASSENZA B, MANCINO C
DELAYED EXPANSION OF THE ATROPHIC MANDIBLE BY ULTRASONIC SURGERY: A CLINICAL AND HISTOLOGICAL CASE SERIES
 INT J ORAL MAXILLOFAC IMPLANTS, 2015 JAN-FEB;30(1):144-9

LORENZON G, BUTTARELLO GM, CHESSA G
CASE REPORT: IMPLANT PLACEMENT AND IMMEDIATE LOADING WITH SIMULTANEOUS BONE REGENERATION FOLLOWING JAW ODONTOGENIC CYST ENUCLEATION
 DENTISTRY, 2015, 5:2

THALMAIR T, FICKL S, SCHNEIDER D, HINZE M, WACHTEL H
DIMENSIONAL ALTERATIONS OF EXTRACTION SITES AFTER DIFFERENT ALVEOLAR RIDGE PRESERVATION TECHNIQUES - A VOLUMETRIC STUDY
 J CLIN PERIODONTOL, 2013 JUL;40(7):721-7

MANESCU A, GIULIANI A, MOHAMMADI S, TROMBA G, MAZZONI S, DIOMEDE F, ZINI N, PIATTELLI A, TRUBIANI O
OSTEOGENIC POTENTIAL OF DUAL-BLOCKS CULTURED WITH HUMAN PERIODONTAL LIGAMENT STEM CELLS: IN VITRO AND SYNCHROTRON
 J PERIODONTAL RES, 2016 FEB;51(1):112-24



SEM image of OsteoBioL® cancellous block
 Author: Prof Ulf Nannmark, University of Göteborg, Sweden



SCARANO A, PIATTELLI A, IEZZI G, VARVARA G
SPONTANEOUS BONE FORMATION ON THE MAXILLARY SINUS FLOOR IN ASSOCIATION WITH SURGERY TO REMOVE A MIGRATED DENTAL IMPLANT: A CASE REPORT

MINERVA STOMATOL, 2014 OCT;63(10):351-9

BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, CALVO GUIRADO JL, NEGRI B, DI FELICE R, COVANI U
VOLUMETRIC ANALYSIS OF REMODELLING PATTERN AFTER RIDGE PRESERVATION COMPARING USE OF TWO TYPES OF XENOGRAPHS. A MULTICENTRE RANDOMIZED CLINICAL TRIAL

CLIN ORAL IMPLANTS RES, 2016 NOV;27(11):E105-E115

ESPOSITO M, GRUSOVIN MG, LAMBERT F, MATOS S, PIETRUSKA M, ROSSI R, SALHI L, BUTI J
THE EFFECTIVENESS OF A RESORBABLE BONE SUBSTITUTE WITH A RESORBABLE MEMBRANE IN THE TREATMENT OF PERIODONTAL INFRABONY DEFECT - A MULTICENTER RANDOMISED CONTROLLED TRIAL

EUR J ORAL IMPLANTOL, 2015;8(3):233-244

OZEL B, FINDIKCIOGLU K, SEZGIN B, GUNAY K, BARUT I, OZMEN S

A NEW OPTION FOR THE RECONSTRUCTION OF ORBITAL FLOOR DEFECTS WITH HETEROLOGOUS CORTICAL BONE

J CRANIO-MAXILLOFAC SURG, 2015 OCT;43(8):1583-8

CORBELLA S, TASCHIERI S, WEINSTEIN R, DEL FABBRO M
HISTOMORPHOMETRIC OUTCOMES AFTER LATERAL SINUS FLOOR ELEVATION PROCEDURE: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS

CLIN ORAL IMPLANTS RES, 2016 SEP;27(9):1106-22

BARONE A, MARCONCINI S, GIAMMARINARO E, MIJIRITSKY E, GELPI F, COVANI U

CLINICAL OUTCOMES OF IMPLANTS PLACED IN EXTRACTION SOCKETS AND IMMEDIATELY RESTORED: A 7-YEAR SINGLE-COHORT PROSPECTIVE STUDY

CLIN IMPLANT DENT RELAT RES, 2016 DEC;18(6):1103-1112

CAKIR M, KARACA IR, AYŞEGÜL F, KAYMAZ F, BOZKAYA S
EXPERIMENTAL EVALUATION OF THE EFFECTS OF ANKAFERD BLOOD STOPPER AND COLLAGENATED HETEROLOGOUS BONE GRAFT ON BONE HEALING IN SINUS FLOOR AUGMENTATION

CLIN ORAL IMPLANTS RES, 2015 MAR-APR;30(2):279-85

LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F
REGENERATION OF ATROPHIC CRESTAL RIDGES WITH RESORBABLE LAMINA: TECHNICAL NOTE

J BIOL REGUL HOMEOST AGENTS 2015 JUL-SEP;29(3 SUPPL 1):97-100

ETTORRE V, DE MARCO P, ZARA S, PERROTTI V, SCARANO A, DI CRESCENZO A, PETRINI M, HADAD C, BOSCO D, ZAVAN B, VALBONETTI L, SPOTO G, IEZZI G, PIATTELLI A, CATALDI A, FONTANA A

IN VITRO AND IN VIVO CHARACTERIZATION OF GRAPHENE OXIDE COATED PORCINE BONE GRANULES CARBON, JULY 2016, VOLUME 103, PAGES 291-298

ROSSI R, RANCITELLI D, POLI PP, RASIA DAL POLO M, NANNMARK U, MAIORANA C

THE USE OF A COLLAGENATED PORCINE CORTICAL LAMINA IN THE RECONSTRUCTION OF ALVEOLAR RIDGE DEFECTS. A CLINICAL AND HISTOLOGICAL STUDY

MINERVA STOMATOL, 2016 OCT;65(5):257-68

SCARANO A, LORUSSO F, RAVERA L, MORTELLARO C, PIATTELLI A
BONE REGENERATION IN ILIAC CRESTAL DEFECTS: AN EXPERIMENTAL STUDY ON SHEEP

BIOMED RES INT, 2016;2016:4086870

FELICE P, ZUCCHELLI G, CANNIZZARO G, BARAUSSE C, DIAZZI M, TRULLENQUE-ERIKSSON A, ESPOSITO M
IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 4-MONTH POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL

EUR J ORAL IMPLANTOL, 2016;9(3):233-247

ROMBOUS C, JEANNEAU C, CAMILLERI J, LAURENT P, ABOUT I
CHARACTERIZATION AND ANGIOGENIC POTENTIAL OF XENOGENEIC BONE GRAFTING MATERIALS: ROLE OF PERIODONTAL LIGAMENT CELLS

DENT MATER J, 2016 DEC 1;35(6):900-907

BARONE A, TOTI P, MARCONCINI S, DERCHI G, MARCHIONNI S, COVANI U
ESTHETIC OUTCOME OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS BY CLINICIANS WITH OR WITHOUT EXPERIENCE: A MEDIUM-TERM PROSPECTIVE EVALUATION

INT J ORAL MAXILLOFAC IMPLANTS, 2016;31(6)

TALLARICO M, KHANARI E, PISANO M, DE RIU G, TULLIO A, MELONI SM

SINGLE POST-EXTRACTIVE ULTRA-WIDE 7 MM-DIAMETER IMPLANTS VERSUS IMPLANTS PLACED IN MOLAR HEALED SITES AFTER SOCKET PRESERVATION FOR MOLAR REPLACEMENT: 6-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL

EUR J ORAL IMPLANTOL, 2016;9(3):263-275

LOPEZ MA, MANZULLI N, CASALE M, ORMIANER Z, CARINCI F

THE USE OF RESORBABLE HETEROLOGOUS CORTICAL LAMINA AS A NEW SINUS LIFT FLOOR: A TECHNICAL NOTE

J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):75-79

LOPEZ MA, ANDREASI BASSI M, CONFALONE L, CARINCI F, ORMIANER Z, LAURITANO D

THE USE OF RESORBABLE CORTICAL LAMINA AND MICRONIZED COLLAGENATED BONE IN THE REGENERATION OF ATROPHIC CRESTAL RIDGES: A SURGICAL TECHNIQUE. CASE SERIES

J BIOL REGUL HOMEOST AGENTS, 2016 APR-JUN;30(2 SUPPL 1):81-85

ESPOSITO M, ZUCCHELLI G, BARAUSSE C, PISTILLI R, TRULLENQUE-ERIKSSON A, FELICE P

FOUR MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN POSTERIOR ATROPHIC JAWS: 4-MONTH POST-LOADING RESULTS FROM A RANDOMISED CONTROLLED TRIAL

EUR J ORAL IMPLANTOL, 2016;9(4):393-409

BARONE A, TOTI P, QUARANTA A, ALFONSI F, CUCCHI A, NEGRI B, DI FELICE R, MARCHIONNI S, CALVO GUIRADO JL, COVANI U, NANNMARK U

CLINICAL AND HISTOLOGICAL CHANGES AFTER RIDGE PRESERVATION WITH TWO XENOGRAPHS: PRELIMINARY RESULTS FROM A MULTICENTER RANDOMIZED CONTROLLED CLINICAL TRIAL

J CLIN PERIODONTOL, 2017 FEB;44(2):204-214

ALFONSI F, BORGIA V, IEZZI G, PIATTELLI A, COVANI U, TONELLI P, BARONE A
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING THE MANAGEMENT OF EXTRACTION SOCKETS

CURR PHARM BIOTECHNOL, 2017;18(1):64-75

BARONE A, TOTI P, MENCHINI FABRIS GB, MARCHIONNI S, COVANI U
EARLY VOLUMETRIC CHANGES AFTER VERTICAL AUGMENTATION OF THE ATROPHIC POSTERIOR MANDIBLE WITH INTERPOSITIONAL BLOCK GRAFT VERSUS ONLAY BONE GRAFT: A RETROSPECTIVE RADIOLOGICAL STUDY

J CRANIO-MAXILLOFAC, 2017 SEP;45(9):1438-1447

BARONE A, TOTI P, FUNEL N, CAMPANI D, COVANI U
EXPRESSION OF SP7, RUNX1, DLX5, AND CTNBN1 IN HUMAN MESENCHYMAL STEM CELLS CULTURED ON XENOGENEIC BONE SUBSTITUTE AS COMPARED WITH MACHINED TITANIUM

IMPLANT DENT, 2014 AUG;23(4):407-15

ESPOSITO M, ZUCCHELLI G, CANNIZZARO G, CHECCHI L, BARAUSSE C, TRULLENQUE-ERIKSSON, FELICE P
IMMEDIATE, IMMEDIATE-DELAYED (6 WEEKS) AND DELAYED (4 MONTHS) POST-EXTRACTIVE SINGLE IMPLANTS: 1-YEAR POST-LOADING DATA FROM A RANDOMISED CONTROLLED TRIAL

EUR J ORAL IMPLANTOL, 2017;10(1):11-26

SCARANO A, CRINCOLI V, DI BENEDETTO A, COZZOLINO V, LORUSSO F, PODALIRI VULPIANI M, GRANO M, KALEMAJ Z, MORI G, GRASSI FR

BONE REGENERATION INDUCED BY BONE PORCINE BLOCK WITH BONE MARROW STROMAL STEM CELLS IN A MINIPIG MODEL OF MANDIBULAR "CRITICAL SIZE" DEFECT STEM CELLS INT, 2017;2017:9082869

SCARANO A
TRADITIONAL POSTEXTRACTIVE IMPLANT SITE PREPARATION COMPARED WITH PRE-EXTRACTIVE INTERRADICULAR IMPLANT BED PREPARATION IN THE MANDIBULAR MOLAR REGION, USING AN ULTRASONIC DEVICE: A RANDOMIZED PILOT STUDY

INT J ORAL MAXILLOFAC IMPLANTS, 2017 MAY/JUN;32(3):655-660

BARONE A, TOTI P, MENCHINI-FABRIS GB, DERCHI G, MARCONCINI S, COVANI U
EXTRA ORAL DIGITAL SCANNING AND IMAGING SUPERIMPOSITION FOR VOLUME ANALYSIS OF BONE REMODELING AFTER TOOTH EXTRACTION WITH AND WITHOUT 2 TYPES OF PARTICULATE PORCINE MINERAL INSERTION: A RANDOMIZED CONTROLLED TRIAL

CLIN IMPLANT DENT RELAT RES, 2017 AUG;19(4):750-759

GIULIANI A, IEZZI G, MAZZONI S, PIATTELLI A, PERROTTI V, BARONE A
REGENERATIVE PROPERTIES OF COLLAGENATED PORCINE BONE GRAFTS IN HUMAN MAXILLA: DEMONSTRATIVE STUDY OF THE KINETICS BY SYNCHROTRON RADIATION MICROTOMOGRAPHY AND LIGHT MICROSCOPY

CLINICAL ORAL INVESTIGATIONS, 2018 JAN;22(1):505-513

IEZZI G, PIATTELLI A, GIULIANI A, MANGANO C, BARONE A, MANZONI L, DEGIDI M, SCARANO A, FILIPPONE A, PERROTTI V
MOLECULAR, CELLULAR AND PHARMACEUTICAL ASPECTS OF FILLING BIOMATERIALS DURING MAXILLARY SINUS-LIFT PROCEDURES. PART 2: DETAILED CHARACTERISTICS OF THE MATERIALS

CURR PHARM BIOTECHNOL, 2017, 18, 33-44

FELICE P, BARAUSSE C, BARONE A, ZUCCHELLI G, PIATTELLI M, PISTILLI R, IPPOLITO DR, SIMION M
INTERPOSITIONAL AUGMENTATION TECHNIQUE IN THE TREATMENT OF POSTERIOR MANDIBULAR ATROPHIES: A RETROSPECTIVE STUDY COMPARING 129 AUTOGENOUS AND HETEROLOGOUS BONE BLOCKS WITH 2 TO 7 YEARS FOLLOW-UP

INT J PERIODONTICS RESTORATIVE DENT, 2017 JUL/AUG;37(4):469-480

IIDA T, CARNEIRO MARTINS NETO E, BOTTICELLI D, APAZA ALCCAYHUAMAN KA, LANG NP, XAVIER SP
INFLUENCE OF A COLLAGEN MEMBRANE POSITIONED SUBJACENT THE SINUS MUCOSA FOLLOWING THE ELEVATION OF THE MAXILLARY SINUS. A HISTOMORPHOMETRIC STUDY IN RABBITS

CLIN IMPLANT DENT RELAT RES, 2017 JUN 7

DE MARCO P, ZARA S, DE COLLI M, RADUNOVIC M, LAZOVIC V, ETTORRE V, DI CRESCENZO A, PIATTELLI A, CATALDI A, FONTANA A

GRAPHENE OXIDE IMPROVES THE BIOCOMPATIBILITY OF COLLAGEN MEMBRANES IN AN IN VITRO MODEL OF HUMAN PRIMARY GINGIVAL FIBROBLASTS

BIOMED MATER, 2017 SEP 13;12(5):055005

MIJIRITSKY E, FERRONI L, GARDIN C, BRESSAN E, ZANETTE G, PIATTELLI A, ZAVAN B

PORCINE BONE SCAFFOLDS ADSORB GROWTH FACTORS SECRETED BY MSCS AND IMPROVE BONE TISSUE REPAIR

MATERIALS, 2017 SEP 8;10(9)

ROSSI R, FOCE E, SCOLAVINO S
THE CORTICAL LAMINA TECHNIQUE: A NEW OPTION FOR ALVEOLAR RIDGE AUGMENTATION. PROCEDURE, PROTOCOL, AND CASE REPORT

J LEBANESE DENTAL ASS, 2017 JAN-JUN; 52(1):35-41

CHECCHI V, FELICE P, ZUCCHELLI G, BARAUSSE C, PIATTELLI M, PISTILLI R, GRANDI G, ESPOSITO M
WIDE DIAMETER IMMEDIATE POST-EXTRACTIVE IMPLANTS VS DELAYED PLACEMENT OF NORMAL-DIAMETER IMPLANTS IN PRESERVED SOCKETS IN THE MOLAR REGION: 1-YEAR POST-LOADING OUTCOME OF A RANDOMISED CONTROLLED TRIAL

EUR J ORAL IMPLANTOL, 2017;10(3):263-278

CRESPI R, CAPPARÈ P, GHERLONE E
COMPARISON OF MAGNESIUM-ENRICHED HYDROXYAPATITE AND PORCINE BONE IN HUMAN EXTRACTION SOCKET HEALING: A HISTOLOGIC AND HISTOMORPHOMETRIC EVALUATION

INT J ORAL MAXILLOFAC IMPLANTS, 2011 SEP-OCT;26(5):1057-62

CORBELLA S, TASCHIERI S, FRANCETTI L, WEINSTEIN R, DEL FABBRO M

HISTOMORPHOMETRIC RESULTS AFTER POSTEXTRACTION SOCKET HEALING WITH DIFFERENT BIOMATERIALS: A SYSTEMATIC REVIEW OF THE LITERATURE AND META-ANALYSIS

INT J ORAL MAXILLOFAC IMPLANTS, 2017 SEP/OCT;32(5):1001-1017

RADUNOVIC M, DE COLLI M, DE MARCO P, DI NISIO C, FONTANA A, PIATTELLI A, CATALDI A, ZARA S
GRAPHENE OXIDE ENRICHMENT OF COLLAGEN MEMBRANES IMPROVES DPSCS DIFFERENTIATION AND CONTROLS INFLAMMATION OCCURRENCE

J BIOMED MATER RES A, 2017 AUG;105(8):2312-2320

KILINC A, ATAOL M
HOW EFFECTIVE IS COLLAGEN RESORBABLE MEMBRANE PLACEMENT AFTER PARTIALLY IMPACTED MANDIBULAR THIRD MOLAR SURGERY ON POSTOPERATIVE MORBIDITY? A PROSPECTIVE RANDOMIZED COMPARATIVE STUDY
 BMC ORAL HEALTH, 2017 OCT 5;17(1):126

TROIANO G, ZHURAKIVSKA K, LO MUZIO L, LAINO L, CICCÌU M, LO RUSSO L
COMBINATION OF BONE GRAFT AND RESORBABLE MEMBRANE FOR ALVEOLAR RIDGE PRESERVATION: A SYSTEMATIC REVIEW, META-ANALYSIS AND TRIAL SEQUENTIAL ANALYSIS
 J PERIODONTOL, 2017 SEP 12;117:1-17. Epub ahead of print

ROSSI R, LONGO E, MIJIRITSKY E
A NEW INTERPRETATION OF GUIDED IMPLANT SURGERY TO ACHIEVE AN OPTIMAL RESULT IN THE AESTHETIC ZONES
 MEDICAL RESEARCH ARCHIVES, 2017 APRIL, VOL. 5, ISSUE 4

BOLLE C, FELICE P, BARAUSSE C, PISTILLI V, TRULLENQUE-ERIKSSON A, ESPOSITO M
FOUR MM-LONG VERSUS LONGER IMPLANTS IN AUGMENTED BONE IN POSTERIOR ATROPHIC JAWS: 1-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2018;11(1):31-47

ESPOSITO M, DAVÒ R, MARTI PAGES C, FERRER FUENTES A, BARAUSSE C, PISTILLI R, IPPOLITO DR, FELICE P
IMMEDIATELY LOADED ZYGOMATIC IMPLANTS VS CONVENTIONAL DENTAL IMPLANTS IN AUGMENTED ATROPHIC MAXILLAE: 4 MONTHS POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2018;11(1):11-28

GASTALDI G, FELICE P, PISTILLI V, BARAUSSE C, IPPOLITO DR, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 × 5 MM IMPLANTS WITH A NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. 3-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2018;11(1):49-61

DIOMEDE F, D'AURORA M, GUGLIANDOLO A, MERCIARO I, ORSINI T, GATTA V, PIATTELLI A, TRUBIANI O, MAZZON E
BIOFUNCTIONALIZED SCAFFOLD IN BONE TISSUE REPAIR
 INT J OF MOLECULAR SCIENCES, 2018, 19, 1022

SCARANO A, LORUSSO F, SANTOS DE OLIVEIRA P, MURMURA G, CARINCI F
DENTAL IMPLANTS DISPLACED INTO THE MANDIBULAR CORPUS: CLINICAL NOTE
 JOURNAL OF DENTISTRY AND ORAL CARE, 2018, 4(1):8-12

FORABOSCO A, GHENO E, SPINATO S, GARUTI G, FORABOSCO E, CONSOLO U
CONCENTRATED GROWTH FACTORS IN MAXILLARY SINUS FLOOR AUGMENTATION: A PRELIMINARY CLINICAL COMPARATIVE EVALUATION
 INT J OF GROWTH FACTORS AND STEM CELLS IN DENTISTRY, 2018;1:2-7

PISTILLI R, FELICE P, CANNIZZARO G, PIATTELLI M, CORVINO V, BARAUSSE C, BUTI J, SOARDI E, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM LONG, 4 MM WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. ONE-YEAR POST-LOADING RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2013;6(4):359-372

FELICE P, BARAUSSE C, PISTILLI V, PIATTELLI M, IPPOLITO DR, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM LONG X 4 MM WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. 3-YEAR POST-LOADING RESULTS FROM A PILOT RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2018;11(2):175-187

DAVÒ R, FELICE P, PISTILLI R, BARAUSSE C, MARTI PAGES C, FERRER FUENTES A, IPPOLITO DR, ESPOSITO M
IMMEDIATELY LOADED ZYGOMATIC IMPLANTS VS CONVENTIONAL DENTAL IMPLANTS IN AUGMENTED ATROPHIC MAXILLAE: 1-YEAR POST-LOADING RESULTS FROM A MULTICENTRE RANDOMISED CONTROLLED TRIAL
 EUR J ORAL IMPLANTOL, 2018;11(1):145-161

OMORI Y, SILVA ER, BOTTICELLI D, APAZA ALCCAYHUAMAN KA, LANG NP, XAVIER SP
REPOSITION OF THE BONE PLATE OVER THE ANTROSTOMY IN MAXILLARY SINUS AUGMENTATION: A HISTOMORPHOMETRIC STUDY IN RABBITS
 CLIN ORAL IMPLANTS RES, 2018 AUG;29(8):821-834

CANULLO L, GENOVA T, NAENNIC N, NAKAJIMA Y, MASUDA K, MUSSANO F
PLASMA OF ARGON ENHANCES THE ADHESION OF MURINE OSTEOBLASTS ON DIFFERENT GRAFT MATERIALS
 ANN ANAT, 2018 JUL;218:265-270

RAMÍREZ FERNÁNDEZ MP, MAZÓN P, GEHRKE SA, CALVO GUIRADO JL, DE AZA PN
COMPARISON OF TWO XENOGRFT MATERIALS USED IN SINUS LIFT PROCEDURES: MATERIAL CHARACTERIZATION AND IN VIVO BEHAVIOR
 MATERIALS, 2017 JUN 7;10(6)

DIOMEDE F, ZINI N, GATTA V, FULLE S, MERCIARO I, D'AURORA M, LA ROVERE RM, TRAINI T, PIZZICANNELLA J, BALLERINI P, PIATTELLI A, TRUBIANI O
HUMAN PERIODONTAL LIGAMENT STEM CELLS CULTURED ONTO CORTICO-CANCELLOUS SCAFFOLD DRIVE BONE REGENERATIVE PROCESS
 EUR CELL MATER, 2016 SEP 16;32:181-201

DEVELIOGLU H, OZCAN G, GULTEKIN SE, SENGUVEN B, YILDIRIM A
THE SHORT-TERM EFFECTS OF VARIOUS XENOGRFTS ON BONE HEALING IN RATS CRANIAL DEFECTS
 BIOMEDICAL RESEARCH, 2018;29(8): 1598-1602

BRUNELLI G, CARINCI F, GIRARDI A, PALMIERI A, CACCIANIGA GL, SOLLAZZO V
OSTEOBIOL® EFFECT ON DENTAL PULP DERIVED STEM CELLS
 EUR J INFLAMMATION, VOL. 10, NO. 1 (S), 27-30

MAZZONI S, MOHAMMADI S, TROMBA G, DIOMEDE F, PIATTELLI A, TRUBIANI O, GIULIANI A
ROLE OF CORTICO-CANCELLOUS HETEROLOGOUS BONE IN HUMAN PERIODONTAL LIGAMENT STEM CELLS XENO-FREE CULTURE STUDIED BY SYNCHROTRON RADIATION PHASE-CONTRAST MICRO TOMOGRAPHY
 INT J MOL SCI, 2017;18:364

LAURITANO D, CARINCI F, ZOLLINO I, HASSANIPOUR A, SAGGESE V, PALMIERI A, GIRARDI A, CURA F, PIRAS A, ZAMBONI P, BRUNELLI G
OSTEOBIOL® ENHANCES OSTEOGENIC DIFFERENTIATION IN BONE MARROW DERIVED STEM CELLS
 EUR J INFLAMMATION, VOL. 10, NO. 1 (S3), 83-88

MATÉ SANCHEZ DE VAL J, MAZON P, PIATTELLI A, CALVO GUIRADO J, MAREQUE BUENO J, GRANERO MARIN J, DE AZA P
COMPARISON AMONG THE PHYSICAL PROPERTIES OF CALCIUM PHOSPHATE-BASED BONE SUBSTITUTES OF NATURAL OR SYNTHETIC ORIGIN
 INT J APPLIED CERAMIC TECHNOL, 2018;15(4):930-937

NAKAJIMA Y, PIATTELLI A, IEZZI G, FORTICH MESA N, FERRI M, BOTTICELLI D
INFLUENCE OF THE PRESENCE OF ALVEOLAR MUCOSA AT IMPLANTS: A HISTOLOGICAL STUDY IN HUMANS
 IMPLANT DENT, 2018;27(2):193-201

CASCONI P, VELLONE V, RAMIERI V, BASILE E, TARSITANO A, MARCHETTI C
RECONSTRUCTION OF THE ADULT HEMIFACIAL MICROSMIA PATIENT WITH TEMPOROMANDIBULAR JOINT TOTAL JOINT PROSTHESIS AND ORTHOGNATIC SURGERY
 CASE REP SURG, 2018;1-6

AMR AEH, ABDEL GHAFAR KA, ABUEL-ELA HA, ABD ELHAMID ES
XENOGENIC FLEXIBLE BONE LAMINA GRAFT: A SUCCESSFUL ALTERNATIVE TO THE AUTOGENOUS ONLY BONE BLOCK GRAFT IN ALVEOLAR RIDGE AUGMENTATION: A CLINICAL, RADIOGRAPHIC AND HISTOLOGICAL EVALUATION
 JOURNAL OF DENTAL TREATMENT AND ORAL CARE, 2017;1(1):104

ATTIA AM
CLINICAL AND RADIOGRAPHIC EVALUATION OF CORTICO-CANCELLOUS BONE MIX XENOGRFT (OSTEOBIOL GEN-OS) IN THE TREATMENT OF HUMAN PERIODONTAL INTRABONY DEFECTS
 EGYPTIAN DENTAL JOURNAL, 2017;63:1-10

BECHARA S, KUBILIUS R, VERONESI G, PIRES JT, SHIBLI JA, MANGANO FG
SHORT (6-MM) DENTAL IMPLANTS VERSUS SINUS FLOOR ELEVATION AND PLACEMENT OF LONGER (>10-MM) DENTAL IMPLANTS: A RANDOMIZED CONTROLLED TRIAL WITH A 3-YEAR FOLLOW-UP
 CLIN ORAL IMPLANTS RES, 2017;28:1097-1107

CHANDRASEKARAN B, SURESH N, MUTHUSAMY S
PLATELET-RICH FIBRIN WITH BONE GRAFTS FOR REGENERATION OF BONY DEFECT FOLLOWING EXTRACTION OF SUPERNUMERARY TEETH: A CASE REPORT
 CHIN J DENT RES, 2017;20(4):231-234

CHIRILA L, ROTARU C, FILIPOV I, SANDULESCU M
MANAGEMENT OF ACUTE MAXILLARY SINUSITIS AFTER SINUS BONE GRAFTING PROCEDURES WITH SIMULTANEOUS DENTAL IMPLANTS PLACEMENT - A RETROSPECTIVE STUDY
 BMC INFECT DIS, 2016 MAR 8;16 SUPPL 1:94

EKSTEIN J, TANDELICH M, NART J, CALVO GUIRADO JL, SHAPIRA L
MARGINAL BONE LEVEL AROUND CONICAL CONNECTION TAPERED IMPLANTS WITH PLATFORM SWITCHING: A MULTICENTER RETROSPECTIVE STUDY AT 14 MONTHS FOLLOW-UP
 J OSSEOINTEGRATION, 2016 JAN-APR;8(1)

DEL CORSO M, DOHAN EHRENFEST DM
IMMEDIATE IMPLANTATION AND PERI-IMPLANT NATURAL BONE REGENERATION (NBR) IN THE SEVERELY RESORBED POSTERIOR MANDIBLE USING LEUKOCYTE- AND PLATELET-RICH FIBRIN (L-PRF): A 4-YEAR FOLLOW-UP
 POSEIDO, 2013;1(2):85-92

NOAMI SA, ELMOSY K, ASKAR N
EVALUATION OF PRE-HYDRATED COLLAGENATED CORTICO-CANCELLOUS GRANULES (MP3®) IN AUGMENTATION OF THE MAXILLARY SINUS (PRELIMINARY STUDY)
 J OF DENT, ORAL DISORD AND THERAPY, 2014;2(3): 1-8

BARONE A, TOTI P, PIATTELLI A, IEZZI G, DERCHI G, COVANI U
EXTRACTION SOCKET HEALING IN HUMANS AFTER RIDGE PRESERVATION TECHNIQUES: COMPARISON BETWEEN FLAPLESS AND FLAPPED PROCEDURES IN A RANDOMIZED CLINICAL TRIAL
 J PERIODONTOL, 2014 JAN;85(1):14-23

LOPEZ MA, LICO S, CASALE M, ORMANIER Z, CARINCI F
THE USE OF VARIOUS BIOMATERIALS IN COMPUTER-GUIDED CRESTAL SINUS LIFT PROCEDURES. A REPORT ON TWO CASE STUDIES WITH VOLUME COMPARISON
 ORAL IMPLANTOL, 2016;APR-JUN 9(2)

MEHL C, BOSCH T
IMPLANTOLOGY AND PROTHODONTICS AT CROSSROADS - CLASSIC VS MODERN TREATMENT CONCEPTS
 JSM DENTAL SURGERY, 2016;1(1):1003


NEMTOI A, DANILA V, DRAGAN E, PASCA S, NEMTOI A, CONSTANTIN M, SAVA A, HABA D
THE EFFECTS OF INSULIN AND STRONTIUM RANELATE ON GUIDED BONE REGENERATION IN DIABETIC RATS
 REVISTA DE CHIMIE, 2017;68(4):693-697

RAMIREZ FERNANDEZ MP, GEHRKE SA, ALBACETE MARTINEZ CP, CALVO GUIRADO JL, DE AZA PN
SEM-EDX STUDY OF THE DEGRADATION PROCESS OF TWO XENOGRFT MATERIALS USED IN SINUS LIFT PROCEDURES
 MATERIALS, 2017;10:542

ZITA GOMES R, DE VASCONCELOS MR, LOPES GUERRA IM, DE ALMEIDA RAB, DE CAMPOS FELINO AC
IMPLANT STABILITY IN THE POSTERIOR MAXILLA: A CONTROLLED CLINICAL TRIAL
 BIOMED RES INT, 2017;6825213

COVANI U, CANULLO L, TOTI P, ALFONSI F, BARONE A
TISSUE STABILITY OF IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS: A 5-YEAR PROSPECTIVE SINGLE-COHORT STUDY
 J PERIODONTOL, 2014 SEP;85(9):E323-32

DIOMEDE F, D'AURORA M, GUGLIANDOLO A, MERCIARO I, ETORRE V, BRAMANTI A, PIATTELLI A, GATTA V, MAZZON E, FONTANA A, TRUBIANI O
A NOVEL ROLE IN SKELETAL SEGMENT REGENERATION OF EXTRACELLULAR VESICLES RELEASED FROM PERIODONTAL-LIGAMENT STEM CELLS
 INT J NANOMEDICINE, 2018 JUN 29;13:3805-3825



Histology showing a granule of OsteoBiol® mp3® partially resorbed and in contact with new bone tissue. Htx-eosine
Author: Prof Ulf Nanmark, University of Göteborg, Sweden

KAWAKAMI S, LANG NP, IIDA T, FERRI M, APAZA ALCCAYHUAMAN KA, BOTTICELLI D
INFLUENCE OF THE POSITION OF THE ANTROSTOMY IN SINUS FLOOR ELEVATION ASSESSED WITH CONE-BEAM COMPUTED TOMOGRAPHY: A RANDOMIZED CLINICAL TRIAL
JOURNAL INVESTIG CLIN DENT, 2018 NOV;9(4):E12362

SCARANO A, DE OLIVEIRA PS, TRAINI T, LORUSSO F
SINUS MEMBRANE ELEVATION WITH HETEROLOGOUS CORTICAL LAMINA: A RANDOMIZED STUDY OF A NEW SURGICAL TECHNIQUE FOR MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT
MATERIALS, 2018 AUG 17;11(8)

FIGLIUZZI MM, GIUDICE A, CRISTOFARO MG, PACIFICO D, BIAMONTE P, FORTUNATO L
POSTEXTRACTIVE IMPLANTS IN AESTHETIC AREAS: EVALUATION OF PERIMPLANT BONE REMODELING OVER TIME
ANN STOMATOL, 2015 MAY 18;6(1):29-34

KIVOVICS M, SZABÓ BT, NÉMETH O, TARI N, DORI F, NAGY P, DOBÓ-NAGY C, SZABÓ G
MICROARCHITECTURAL STUDY OF THE AUGMENTED BONE FOLLOWING RIDGE PRESERVATION WITH A PORCINE XENOGRAFT AND A COLLAGEN MEMBRANE: PRELIMINARY REPORT OF A PROSPECTIVE CLINICAL, HISTOLOGICAL, AND MICRO-COMPUTED TOMOGRAPHY ANALYSIS
INT J ORAL MAXILLOFAC SURG, 2017 FEB;46(2):250-260

MARCONCINI S, GIAMMARINARO E, DERCHI G, ALFONSI F, COVANI U, BARONE A
CLINICAL OUTCOMES OF IMPLANTS PLACED IN RIDGE-PRESERVED VERSUS NONPRESERVED SITES: A 4-YEAR RANDOMIZED CLINICAL TRIAL
CLIN IMPL DENT RELAT RES, 2018 DEC;20(6):906-914

MARCONCINI S, COVANI U, GIAMMARINARO E, VELASCO-ORTEGA E, DE SANTIS D, ALFONSI F, BARONE A
CLINICAL SUCCESS OF DENTAL IMPLANTS PLACED IN POSTERIOR MANDIBLE AUGMENTED WITH INTERPOSITIONAL BLOCK GRAFT: 3-YEAR RESULTS FROM A PROSPECTIVE COHORT CLINICAL STUDY
J ORAL MAXILLOFAC SURG, 2019 FEB;77(2):289-298

IIDA T, RICARDO SILVA E, LANG NP, APAZA ALCCAYHUAMAN KA, DANIELE BOTTICELLI D, XAVIER SP
HISTOLOGICAL AND MICRO-COMPUTED TOMOGRAPHY EVALUATIONS OF NEWLY FORMED BONE AFTER MAXILLARY SINUS AUGMENTATION USING A XENOGRAFT WITH SIMILAR DENSITY AND MINERAL CONTENT OF BONE: AN EXPERIMENTAL STUDY IN RABBITS
CLIN EXP DENT RES, 2018;1-7

BERNARDI S, GATTO R, SEVERINO M, BOTTICELLI G, CARUSO S, RASTELLI C, LUPI E, ROIAS AQ, IACOMINO E, FALISI G
SHORT VERSUS LONGER IMPLANTS IN MANDIBULAR ALVEOLAR RIDGE AUGMENTED USING OSTEOGENIC DISTRACTION: ONE-YEAR FOLLOW-UP OF A RANDOMIZED SPLIT-MOUTH TRIAL
J ORAL IMPLANTOL, 2018 JUN;44(3):184-191

DIKER N, SARICA H, CUMBUL A, KILIC E
EFFECTS OF SYSTEMIC ERYTHROPOIETIN TREATMENT AND HETEROGENEOUS XENOGRAFT IN COMBINATION ON BONE REGENERATION OF A CRITICAL-SIZE DEFECT IN AN EXPERIMENTAL MODEL
J CRANIOMAXILLOFAC SURG, 2018 NOV;46(11):1919-1923

SCARANO A, MURMURA G, MASTRANGELO F, LORUSSO F, GRECO LUCCHINA A, CARINCI F
A NOVEL TECHNIQUE TO PREVENT SINUS MEMBRANE COLLAPSE DURING MAXILLARY SINUS FLOOR AUGMENTATION WITHOUT BONE GRAFT: TECHNICAL NOTE
J BIOL REGUL HOMEOST AGENTS, 2018 NOV-DEC;32(6):1589-1592

SENESE O, BOUTREMAN E, GOSSIAUX C, LOEB I, DEQUANTER D
RETROSPECTIVE ANALYSIS OF 79 PATIENTS WITH ORBITAL FLOOR FRACTURE: OUTCOMES AND PATIENT-REPORTED SATISFACTION
ARCH CRANIOFAC SURG, 2018 JUN;19(2):108-113

GENOVA T, PESCE P, MUSSANO F, TANAKA K, CANULLO L
THE INFLUENCE OF BONE-GRAFT BIO-FUNCTIONALIZATION WITH PLASMA OF ARGON ON BACTERIAL CONTAMINATION
J BIOMED MATER RES A, 2019 JAN;107(1):67-70

ESPOSITO M, BARAUSSE C, PISTILLI R, PIATTELLI M, DI SIMONE S, IPPOLITO DR, FELICE P
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 5 X 5 MM IMPLANTS WITH A NOVEL NANOSTRUCTURED CALCIUM-INCORPORATED TITANIUM SURFACE OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR RESULTS FROM A RANDOMISED CONTROLLED TRIAL
INT J OF ORAL IMPLANTOL, 2019;12(1):39-54

FELICE P, PISTILLI R, BARAUSSE C, PIATTELLI M, BUTI J, ESPOSITO M
POSTERIOR ATROPHIC JAWS REHABILITATED WITH PROSTHESES SUPPORTED BY 6 MM LONG X 4 MM WIDE IMPLANTS OR BY LONGER IMPLANTS IN AUGMENTED BONE. FIVE-YEAR POST-LOADING RESULTS FROM A WITHIN-PERSON RANDOMISED CONTROLLED TRIAL
INT J OF ORAL IMPLANTOL, 2019;12(1):57-72

DI CARLO R, ZARA S, VENTRELLA A, SIANI G, DA ROS T, IEZZI G, CATALDI A, FONTANA A
COVALENT DECORATION OF CORTICAL MEMBRANES WITH GRAPHENE OXIDE AS A SUBSTRATE FOR DENTAL PULP STEM CELLS
NANOMATERIALS, 2019;9:604

CABALLÉ-SERRANO J, MUNAR-FRAUA A, DELGADO L, PÉREZ R, HERNÁNDEZ-ALFARO F
PHYSICO-CHEMICAL CHARACTERIZATION OF BARRIER MEMBRANES FOR BONE REGENERATION
J MECH BEHAV BIOMED, 2019;97:13-20

MATOH U, PETELIN M, GASPERIC R
SPLIT-MOUTH COMPARISON OF CORONALLY ADVANCED FLAP WITH CONNECTIVE TISSUE GRAFT OR COLLAGEN MATRIX FOR TREATMENT OF ISOLATED GINGIVAL RECESSIONS
INT J PERIODONTICS RESTORATIVE DENT, 2019;39(3):439-446

AMBROZEWICZ E, MUSZYNSKA M, TOKAJUK G, GRYNKIEWICZ G, ŻARKOVIC N, SKRZYDLEWSKA
BENEFICIAL EFFECTS OF VITAMINS K AND D3 ON REDOX BALANCE OF HUMAN OSTEOBLASTS CULTURED WITH HYDROXYAPATITE-BASED BIOMATERIALS CELLS
2019;8:325

POLIS-YANES C, CADENAS-SEBASTIÁN C, GUAL-VAQUÉS P, AYOUSO-MONTERO R, MARÍ-ROIG A, LÓPEZ-LÓPEZ J
GUIDED BONE REGENERATION OF AN ATROPHIC MAXILLA USING HETEROLOGOUS CORTICAL LAMINA
CASE REP DENT, 2019; 5216362

KAWAKAMI S, LANG NP, FERRI M, APAZA ALCCAYHUAMAN KA, BOTTICELLI D
INFLUENCE OF THE HEIGHT OF THE ANTROSTOMY IN SINUS FLOOR ELEVATION ASSESSED BY CONE BEAM COMPUTED TOMOGRAPHY: A RANDOMIZED CLINICAL TRIAL
INT J ORAL MAXILLOFAC IMPLANTS, 2019;34(1):223-232

CHECCHI V, GASPARRO R, PISTILLI R, CANULLO L, FELICE P
CLINICAL CLASSIFICATION OF BONE AUGMENTATION PROCEDURE FAILURES IN THE ATROPHIC ANTERIOR MAXILLAE: ESTHETIC CONSEQUENCES AND TREATMENT OPTIONS
BIOMED RES INT., 2019:4386709

FISCHER K R, TESTORI T, WACHTEL H, MÜHLEMANN S, HAPPE A, DEL FABBRO M
SOFT TISSUE AUGMENTATION APPLYING A COLLAGENATED PORCINE DERMAL MATRIX DURING SECOND STAGE SURGERY: A PROSPECTIVE MULTICENTER CASE SERIES
CLIN IMPLANT DENT RELAT RES.,2019;1-8

HIROTA A, LANG NP, FERRI M, FORTICH MESA N, APAZA ALCCAYHUAMAN KA, BOTTICELLI D
TOMOGRAPHIC EVALUATION OF THE INFLUENCE OF THE PLACEMENT OF A COLLAGEN MEMBRANE SUBJACENT TO THE SINUS MUCOSA DURING MAXILLARY SINUS FLOOR AUGMENTATION: A RANDOMIZED CLINICAL TRIAL
INT J IMPLANT DENT,2019; 5(1):31

TANAKA K, IEZZI G, PIATTELLI A, FERRI M, MESA NF, APAZA ALCCAYHUAMAN KA, BOTTICELLI D
SINUS FLOOR ELEVATION AND ANTROSTOMY HEALING: A HISTOMORPHOMETRIC CLINICAL STUDY IN HUMANS
IMPLANT DENT, 2019; AUG 14

JEANNEAU C, LE FOURNIS C, ABOUT I
XENOGENIC BONE FILLING MATERIALS MODULATE MESENCHYMAL STEM CELL RECRUITMENT: ROLE OF THE COMPLEMENT CSA
CLIN. ORAL INVESTIG.; 2019 OCT 23

RAMANAUSKAITE A, BORGES T, ALMEIDA BL, CORREIA A
DENTAL IMPLANT OUTCOMES IN GRAFTED SOCKETS: A SYSTEMATIC REVIEW AND META-ANALYSIS
J ORAL MAXILLOFAC RES., 2019 SEP 5;10(3):E8

FARIA-ALMEIDA R, ASTRAMSKAITE-JANUSEVICIENE I, PUISYS A, CORREIA F
EXTRACTION SOCKET PRESERVATION WITH OR WITHOUT MEMBRANES, SOFT TISSUE INFLUENCE ON POST EXTRACTION ALVEOLAR RIDGE PRESERVATION: A SYSTEMATIC REVIEW
J ORAL MAXILLOFAC RES., 2019 SEP 5;10(3):E5

ADIOGLU S, GIRAY CB, KULAC I, USUBUTUN A, AKTAS A
CLINICAL AND HISTOPATHOLOGICAL COMPARATIVE STUDY OF TWO EQUINE-DERIVED BONE GRAFT: A HUMAN STUDY
J PAK MED ASSOC., 2019; 69: 1617

ROSSI R, FOCE E
RECONSTRUCTION OF A HORIZONTAL AND VERTICAL BONE DEFECT USING THE CORTICAL LAMINA TECHNIQUE
MED. RES. ARCH, 2019; VOL 7 (11)


CANULLO L, GENOVA T, RAKIC M, SCULEAN A, MIRON R, MUZZI M, CAROSSA S, MUSSANO F
EFFECTS OF ARGON PLASMA TREATMENT ON THE OSTEOCONDUCTIVITY OF BONE GRAFTING MATERIALS
CLIN ORAL INVESTIG, 2019 NOV 20

OsteoBiol® product codes




PRODUCT	PACKAGING	TYPE	SIZE	PORCINE CODE	EQUINE CODE
BONE SUBSTITUTES					
Gen-Os®	1 Vial	DRIED GRANULES	0.25 g	M1052FS	M1052FE
Gen-Os®	1 Vial	DRIED GRANULES	0.5 g	M1005FS	M1005FE
Gen-Os®	1 Vial	DRIED GRANULES	1.0 g	M1010FS	M1010FE
Gen-Os®	1 Vial	DRIED GRANULES	2.0 g	M1020FS	M1020FE
Gen-Os® 1000-2000	1 Vial	DRIED GRANULES	1.0 g	M0210FS	
Gen-Os® 1000-2000	1 Vial	DRIED GRANULES	2.0 g	M0220FS	
TSV Gel	1 Syringe	TSV GEL	0.5 g	TSV005S <small>in kit with M1005FS or A1005FS</small>	TSV005E <small>in kit with M1005FE or A1005FE</small>
TSV Gel	1 Syringe	TSV GEL	1.0 g	TSV010S <small>in kit with M1010FS or A1010FS</small>	TSV010E <small>in kit with M1010FE or A1010FE</small>
mp3®	1 Syringe	BONE MIX	0.5 cc	A3095FS	A3095FE
mp3®	1 Syringe	BONE MIX	1.0 cc	A3005FS	A3005FE
mp3®	3 Syringes	BONE MIX	3x0.25 cc (0.75 cc)	A3075FS	
mp3®	3 Syringes	BONE MIX	3x0.5 cc (1.5 cc)	A3015FS	A3015FE
mp3®	3 Syringes	BONE MIX	3x1.0 cc (3.0 cc)	A3030FS	A3030FE
mp3®	1 Syringe (wide tip)	BONE MIX	2.0 cc	A3010FS	A3010FE
mp3® 1000-2000	1 Syringe (wide tip)	BONE MIX	2.0 cc	A3210FS	A3210FE
GTO®	1 Syringe	BONE MIX + TSV Gel	0.5 cc	MU0005S	MU0005E
GTO®	1 Syringe	BONE MIX + TSV Gel	2.0 cc	MU0020S	MU0020E
Putty	1 Syringe	BONE PASTE	0.25 cc	HPT52S	
Putty	1 Syringe	BONE PASTE	0.5 cc	HPT09S	HPT09E
Putty	3 Syringes	BONE PASTE	3x0.25 cc (0.75 cc)	HPT32S	HPT32E
Putty	3 Syringes	BONE PASTE	3x0.5 cc (1.5 cc)	HPT35S	HPT35E
Putty	1 Syringe (wide tip)	BONE PASTE	1.0 cc	HPT61S	HPT61E
Gel 40	1 Syringe	BONE GEL	0.5 cc	05GEL40S	05GEL40E
Gel 40	3 Syringes	BONE GEL	3x0.5 cc (1.5 cc)	15GEL40S	15GEL40E
Apatos Mix	1 Vial	DRIED GRANULES	0.5 g	A1005FS	A1005FE
Apatos Mix	1 Vial	DRIED GRANULES	1.0 g	A1010FS	A1010FE
Apatos Mix	1 Vial	DRIED GRANULES	2.0 g	A1020FS	A1020FE
Apatos Cortical	1 Vial	DRIED GRANULES	0.5 g	AC1005FS	
Apatos Cortical	1 Vial	DRIED GRANULES	1.0 g	AC1010FS	
Apatos Mix 1000-2000	1 Vial	DRIED GRANULES	1.0 g	A0210FS	A0210FE

PRODUCT	PACKAGING	TYPE	SIZE	PORCINE CODE	EQUINE CODE
BLOCKS					
Sp-Block	1 Blister	DRIED BLOCK / NORM	10x10x10 mm		BN0E
Sp-Block	1 Blister	DRIED BLOCK / NORM	10x10x20 mm		BN1E
Sp-Block	1 Blister	DRIED BLOCK / NORM	10x20x20 mm		BN2E
Sp-Block	1 Blister	DRIED BLOCK / NORM	35x10x5 mm		BN8E
Dual-Block CURVED	1 Blister	DRIED BLOCK / SOFT	20x15x5 mm	STS7S	
Dual-Block CURVED	1 Blister	DRIED BLOCK / NORM	20x10x5 mm	STN5S	
MEMBRANES AND BARRIERS					
Evolution	3 Blister	DRIED / X-FINE	30x30x (0.2) mm	EM33XS	
Evolution	1 Blister	DRIED / FINE	20x20x (0.3) mm		EV02LLE
Evolution	1 Blister	DRIED / FINE	30x30x (0.3) mm		EV03LLE
Evolution	1 Blister	DRIED / FINE	Oval 25x35x (0.3) mm		EVOLLE
Evolution	1 Blister	DRIED / FINE	40x40x (0.3) mm		EV04LLE
Evolution	1 Blister	DRIED / FINE	80x60x (0.3) mm		EV06LLE
Evolution	1 Blister	DRIED / STANDARD	20x20x (0.4) mm	EM02HS	EV02HHE
Evolution	1 Blister	DRIED / STANDARD	30x30x (0.4) mm	EM03HS	EV03HHE
Evolution	1 Blister	DRIED / STANDARD	Oval 25x35x (0.4) mm	EM00HS	
Derma	1 Blister	DRIED	20x20x (0.5) mm	ED02LS	
Derma	1 Blister	DRIED	Oval 12x8x (0.9) mm	ED21FS	
Derma	1 Blister	DRIED	25x25x (0.9) mm	ED25FS	
Derma	1 Blister	DRIED	50x50x (0.9) mm	ED05FS	
Derma	1 Blister	DRIED	7x50x (2.0) mm	ED75SS	
Derma	1 Blister	DRIED	15x5x (2.0) mm	ED15SS	
Derma	1 Blister	DRIED	30x30x (2.0) mm	ED03SS	
Derma	1 Blister	DRIED	50x50x (2.0) mm	ED05SS	
Soft Cortical Lamina	1 Blister	DRIED	25x25x (0.5) mm	LS25FS	LS25FE
Soft Cortical Lamina	1 Blister	DRIED	Oval 25x35x (0.5) mm	LS23FS	LS23FE
Soft Cortical Lamina	1 Blister	DRIED	20x40x (1.0) mm	LS24LS	
Curved Lamina	1 Blister	DRIED	35x35x (1.0) mm	LS10HS	LS10HE
Soft Cortical Lamina	1 Blister	DRIED	30x30x (3.0) mm	LS03SS	LS03SE
Cortical Lamina	1 Blister	DRIED	35x15x (0.7) mm	LS15LS	
Cortical Lamina	1 Blister	DRIED	35x35x (1.0) mm	LS35LS	
Special	1 Blister	DRIED	20x20x (0.2) mm		EM02LE
Special	1 Blister	DRIED	30x30x (0.2) mm		EM03LE





OsteoBiol® by Tecross




Tecross s.r.l. is an innovative, globally active company that develops, produces and documents premium-quality xenogenic biomaterials by the brands Tecross® and OsteoBiol®.

Its 20 years of research led to its patent-protected production process that ensures neutralization of antigenic components in order to achieve biocompatibility, while preserving the natural collagen matrix inside the biomaterial.

Tecross® products comply with highest quality standards such as ISO 13485 and European laws.

osteobiol.com

DISTRIBUTED BY



MK-CAT0120EN

Tecross® s.r.l.
info@tecross.com
www.tecross.com

BIOMATERIALS ENGINEERING

Tecross® Dental s.r.l.
info@tecross-dental.com
www.osteobiol.com

INTERNATIONAL SALES & MARKETING