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OsteoBio® Collagenated membranes and fetls – English version -

Class III Medical Device, collagenated barrier membranes of heterologous origin, resorbable, in guided tissue regeneration surgery.

Special: Made from extremely fine (pericardium), translucent, resorbable, mesenchymal, biocompatible tissue, thickness<1mm.

ESTIMATED RESCRIPTION TIME: about 40 days. Resorption of the material may take longer depending on the graft site. **CLINICAL INDICATIONS:** to protect the bone grafts, to separate the different tissue compartments and as a temporary barrier to protect the damaged tissues against external agents.

DIRECTIONS FOR USE: the grafting set must be properly washed and cleaned, and the product must be well stabilized and protected. Hydrate before use with sterile solutions and apply to the grafting site. Proceed with flap surgery.

PACKAGING: sterile disposable blister.

SIDE EFFECTS: none known. Although the grafting procedure adopted has demonstrated the inhibition of adverse antibody reactions, the remote possibility of graft rejection cannot be excluded.

CONTRAINDICATIONS: do not use in patients who have pathologies or conditions that contraindicate surgery or who have an acute or chronic infection (osteoemilias). Not indicated for patients with immunodeficiency or suffering from autoimmune diseases, particularly children and the elderly. Do not use directly on large vessel aneurysms or nerves. Do not use on poorly vascularized, neonate or infected receiving sites. Do not use on sites that do not guarantee the stability and protection of the product after grafting. The product must not be used if the packaging is damaged or expired.

WARNING: This product has not been tested on pregnant women. The product is for single use only and, as such, should not be re-sterilized. The graft site must be properly washed and cleaned, and the product must be well stabilized and protected.

Product reserved for surgeons with experience in biomaterial implants: strictly adhere to the instructions given in the scientific literature on heterologous grafts, especially the scientific literature on OsteoBio® devices available on the website www.osteobiol.com, as well as the instructions herein. For every single application, the specialist surgeon has the responsibility to ensure proper clinical management of the biomaterials and the patient.

EXPIRY DATE: the product is valid for five years from the date of sterilization and the expiry date is indicated on the label.

STORAGE: The product must be stored away from direct heat sources and sunlight, in a dry and clean place, under the conditions indicated on the label.

DISPOSAL: the product does not contain hazardous substances and it may be disassembled and disposed of as regular municipal waste.

PRODUCT ORIGIN: the last letter on the label of the product code REF indicates the source of the material (S = porcine, E = equine).

Symbols used on the label:

single use year and month of expiry oot number carefully read the instructions Sterile. Sterilization performed by irradiation with gamma rays. If the packaging is open or damaged the sterility of the product is no longer guaranteed limitation of the temperature

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