





control and comfort

multi-frequency apex locator technology

- · works in dry and wet canals
- no calibration, no zero adjustment are necessary

propex•pixi™ is designed to be placed on the patient's chest while working

- · less disturbance by tangling wires
- improved visual control of the file progression by placing the propex pixi™ in the line of vision
- clip-on measuring wire and no-slip pads keep the propex•pixi™ in place

progressive sound control with 4 volume levels

· dual control of the file progression

THINK BIG

"... the canal terminus must be detected accurately during canal preparation and precise control of working length during the process must be maintained."

Nekoofar M. H., Ghandi M. M., Hayes S. J. & Dummer P. M. H.

The fundamental operating principles of electronic root canal length measurement devices. International Endodontic Journal 2006; 39: 595-609

small size, big benefits

· takes up minimal space during treatment, in storage and transport

while

· preserving its accuracy and reliability




SUCCESS

"The establishment of a correct working length is one of the fundamental parameters for endodontic success."

Ebrahim A. K., Wadachi R., Suda H. Electronic Apex Locators — A Review. Journal of Medical and Dental Sciences 2007; 54: 125–136

propex•pixi™ operating features



The file enters the canal 1 double signal

The file is in the pre-apical zone «2.0» to «1.5» intermittent signal

Apical zone: «1.0» to «0.0» **APEX: «0.0»**

solid signal

PRECISION

"A precise working length prevents transportation of the apical foramen, possible overfilling... inadequate debridement, which, in turn, could lead to retained pulp tissue and lack of healing"

Jakobson S. J. M., Westphalen V. P. D., da Silva Neto U. X., Fariniuk L- F., Picoli F., Carneiro E. The Accuracy in the Control of the Apical Extent of Rotary Canal Instrumentation Using RootZX II and ProTaper Instruments: An In Vivo Study. Journal of Endodontics 2008; 34: 1342-1345

OVER

The file is through the apical foramen rapid intermittent signal

Attention

The scale «2.0» to «0.0» does not represent a distance in mm – it indicates the file progression towards the apex

propex•pixi™





propex•pixi™ including 1 measuring wire,

1 connection hook, 2 lip clips and 1 charger

propex•pixi™ measuring wire

REF A 1031-100

RFF A 1030

propex•pixi™ charger

REF A 1031-200

propex® accessory:

- lip clip, 5 pces REF A 1029-300



- connection hook, 2 pces REF A 1029-400

- fork, 2 pces

REF A 1029-600

propex•pixi™ operates on NiMH rechargeable battery

AAA 1.2 V 1000 mAh from GP Batteries or equivalent available in the shops