

OsteoBiol[®]

by TecnoSS

mp3

ULTIMATE PERFORMANCE AND HANDLING

Pre-hydrated collagenated heterologous cortico-cancellous bone mix

REGENERATION SCIENCE

INSPIRED BY NATURE



A unique biotechnology

TECNOSS®: A UNIQUE PROCESS THAT ACCELERATES AND GUIDES NATURAL BONE REGENERATION

Tecnoss® developed and patented a unique biotechnology that prevents the ceramization phase of natural bone and preserves the tissue collagen, allowing an osteoclastic-type remodelling of the biomaterial similar to physiological bone turnover and delivering a product endowed with characteristics very similar to human mineral bone⁽¹⁾.

The combination of these factors allows a consistent new bone formation and a close contact between neo-formed bone and biomaterial granules^(A).

COLLAGEN: A KEY FACTOR FOR BONE REGENERATION

Collagen has a key role in bone regeneration process in that:

- it acts as a valid substrate for platelet activation and aggregation
- it serves to attract and differentiate the mesenchymal stem cells present in the bone marrow⁽²⁾
- it increases the proliferation rate of the osteoblasts up to 2/3 times⁽³⁾
- it stimulates the activation of the platelets, osteoblasts and osteoclasts in the tissue healing process



OSTEOBIOL®: UNIQUE COLLAGENATED BIOMATERIALS

Thanks to the innovative Tecnoss® technology, the OsteoBiol® line has the following important characteristics:

- absence of a foreign body response
- gradual resorption over time⁽⁴⁾
- stimulation/acceleration of physiological tissue healing process
- protection of the grafting site from infection (membranes)

The Tecnoss® new generation of biomaterials, thanks to a revolutionary technology, goes beyond the simple role of aiding natural bone regrowth by stimulating and accelerating this vital physiological process.

OsteoBiol®
by Tecnoss



A | Histology at 6 months. Human sinus grafted with OsteoBiol® mp3. Biopsy courtesy of Dr P Palacci, Marseille, France. Histology courtesy of Dr U Nannmark, University of Göteborg, Sweden

B | SEM image of an OsteoBiol® bone matrix colonized by osteoblasts from a cell-line (MG63). Courtesy of Dr U Nannmark, University of Göteborg, Sweden

C | SEM image showing OsteoBiol® mp3 particles, granulometry 600-1000 microns. Courtesy of Dr U Nannmark, University of Göteborg, Sweden

1 | Figueiredo M et al. J Biomed Mater Res B Appl Biomater, 2010

2 | Salasznyk RM, et al. J Biomed Biotechnol, 2004

3 | Hsu FY, et al. Biomaterials, 1999

4 | Nannmark U et al. Clin Implant Dent Relat Res, 2008

5 | Barone A et al. Clin Implant Dent Relat Res, 2010

6 | Barone A et al. Clin Oral Implants Res, 2008

7 | Barone A et al. J Periodontol, 2008

8 | Barone A et al. J Oral Maxillofac Surg, 2007



Ultimate performance and handling




CHARACTERISTICS

Heterologous origin biomaterial made of 600-1000 microns collagenated cortico-cancellous bone mix properly mixed with collagen gel. Gradually resorbable⁽⁴⁾, it preserves the original graft shape and volume (osteoconductive property). Moreover, the preserved collagen matrix facilitates blood clotting and the subsequent invasion of repairing and regenerative cells. These unique properties allow an excellent rate of new bone formation⁽⁵⁾, delivering graft volume preservation, a healthy new bony tissue and ultimately, a successful implant rehabilitation.

HANDLING

Available in ready-to-use pre-hydrated syringes, mp3 can be easily grafted avoiding the hydration and manipulation phases decreasing the risk of accidental exposure to pathogens.



OsteoBiol® mp3 at 6 months from grafting with lateral access sinus lift procedure in human patients⁽⁵⁾

CLINICAL INDICATIONS

mp3 main clinical indication is lateral access maxillary sinus lift^(5,6), always in association with Evolution membranes: the mp3 syringe can be directly applied into the bony window without having to mix the product with saline or blood. Due to its collagen gel component, mp3 allows an excellent graft stability while its hydrophilia guarantees quick blood absorption and therefore the necessary graft vascularization.



Maxillary sinus grafted with OsteoBiol® mp3. Courtesy of Dr Antonio Barone, Lido di Camaiore, Italy

mp3 can also be used in combination with Evolution membranes for alveolar ridge preservation⁽⁷⁾: the application of this biomaterial limits significantly the alveolar ridge width and height reduction that would naturally occur with spontaneous healing, preserving thus the alveolar ridge volume and allowing a correct second stage implant placement.

Finally, mp3 is indicated for horizontal augmentation (two wall defects) in combination with autogenous bone blocks⁽⁸⁾ or with OsteoBiol® Cortical Lamina (curved model).

Tissue of origin

Heterologous cortico-cancellous bone mix

Tissue collagen

Preserved plus an additional 10% collagen gel (OsteoBiol® Gel 0)

Physical form

Pre-hydrated granules and collagen gel

Composition

90% granulated mix, 10% collagen gel

Granulometry

600-1000 microns

Re-entry time

About 5 months

Packaging

Syringe: 1.0cc, 3x 0.5cc, 3x 1.0cc



MAXILLARY SINUS FLOOR AUGMENTATION






POST-EXTRACTIVE SOCKETS



TWO-WALL DEFECTS



Excellent clinical performances



CASE REPORT

■ LATERAL ACCESS SINUS LIFT

Lateral access sinus lift with simultaneous implant and horizontal augmentation

Sex: **Female** | Age: **42**

Fig. 1 Initial x-ray showing a 3mm residual bone

Fig. 2 Flap opening, a substantial vestibular bone resorption can be determined

Fig. 3 Antrostomy performed with Piezo surgery technique

Fig. 4 A OsteoBio® Evolution membrane is inserted through the antrostomy to protect the Schneider membrane from grafting material

Fig. 5 Maxillary sinus grafted with OsteoBio® mp3

Fig. 6 Immediate implant placement

Fig. 7 A OsteoBio® Evolution membrane is stabilised with osteosynthesis screws above the antrostomy

Fig. 8 Cortical bone stimulation

Fig. 9 OsteoBio® mp3 is grafted on the vestibular side of the defect for horizontal augmentation

Fig. 10 The OsteoBio® Evolution membrane is stabilised into position with a transpalatal suture


Fig. 11 Final situation

Fig. 12 Post-operative x-ray

Documentation provided by
Dr **Rosario Sentineri**
Private practitioner in Genova, Italy
e-mail: rosario.sentineri@gmail.com

Bone substitute: **OsteoBio® mp3**
Membrane: **OsteoBio® Evolution**

- Covani U, Ameri S, Crespi R, Barone A
PRESERVAZIONE DEL PROCESSO ALVEOLARE CON OSSEO ETEROLOGO. CONSIDERAZIONI ISTOLOGICHE
 ITALIAN ORAL SURGERY, 2004
- Cassetta M, Calasso S, Voza I, Dell'Aquila D
REHABILITATION OF ATROPHIC ALVEOLAR CRESTS WITH CYLINDRICAL SANDBLASTED AND ACID ETCHED IMPLANTS: A PILOT STUDY
 EUROPEAN JOURNAL OF IMPLANT PROSTHODONTICS, 2005
- Arcuri C, Cecchetti F, Germano F, Motta A, Santacroce C
CLINICAL AND HISTOLOGICAL STUDY OF A XENOGENIC BONE SUBSTITUTE USED AS A FILLER IN POSTEXTRACTIVE ALVEOLUS
 MINERVA STOMATOLOGICA, 2005
- Barone A, Crespi R, Aldini NN, Fini M, Giardino R, Covani U
MAXILLARY SINUS AUGMENTATION: HISTOLOGIC AND HISTOMORPHOMETRIC ANALYSIS
 INTERNATIONAL JOURNAL OF ORAL AND MAXILLOFACIAL IMPLANTS, 2005
- Rinna C, Ungari C, Saltarel A, Cassoni A, Reale G
ORBITAL FLOOR RESTORATION
 JOURNAL OF CRANIOFACIAL SURGERY, 2005
- Barone A, Ameri S, Covani U
IMMEDIATE POSTEXTRACTION IMPLANTS: TREATMENT OF RESIDUAL PERI-IMPLANT DEFECTS. A RETROSPECTIVE ANALYSIS
 EUROPEAN JOURNAL OF IMPLANT PROSTHODONTICS, 2006
- Barone A, Santini S, Sbordone L, Crespi R, Covani U
A CLINICAL STUDY OF THE OUTCOMES AND COMPLICATIONS ASSOCIATED WITH MAXILLARY SINUS AUGMENTATION
 INTERNATIONAL JOURNAL OF ORAL AND MAXILLOFACIAL IMPLANTS, 2006
- Covani U, Barone A, Cornelini R, Crespi R
CLINICAL OUTCOME OF IMPLANTS PLACED IMMEDIATELY AFTER IMPLANT REMOVAL
 JOURNAL OF PERIODONTOLOGY, 2006
- Orsini G, Scarano A, Piattelli M, Piccirilli M, Caputi S, Piattelli A
HISTOLOGIC AND ULTRASTRUCTURAL ANALYSIS OF REGENERATED BONE IN MAXILLARY SINUS AUGMENTATION USING A PORCINE BONE-DERIVED BIOMATERIAL
 JOURNAL OF PERIODONTOLOGY, 2006
- Trubiani O, Scarano A, Orsini G, Di Iorio D, D'Arcangelo C, Piccirilli M, Sigismondo M, Caputi S
THE PERFORMANCE OF HUMAN PERIODONTAL LIGAMENT MESENCHYMAL STEM CELLS ON XENOGENIC BIOMATERIALS
 INTERNATIONAL JOURNAL OF IMMUNOPATHOLOGY AND PHARMACOLOGY, 2007
- Barone A, Covani U
MAXILLARY ALVEOLAR RIDGE RECONSTRUCTION WITH NONVASCULARIZED AUTOGENOUS BLOCK BONE: CLINICAL RESULTS
 JOURNAL OF ORAL AND MAXILLOFACIAL SURGERY, 2007
- Calvo Guirado JL, Pardo Zamora G, Saez Yuguero MR
RIDGE SPLITTING TECHNIQUE IN ATROPHIC ANTERIOR MAXILLA WITH IMMEDIATE IMPLANTS, BONE REGENERATION AND IMMEDIATE TEMPORISATION: A CASE REPORT
 JOURNAL OF IRISH DENTAL ASSOCIATION, 2007
- Del Corso M
SOFT TISSUE RESPONSE TO PLATELET RICH FIBRIN: CLINICAL EVIDENCES
 COSMETIC DENTISTRY, 2008
- Barone A, Santini S, Marconcini S, Giacomelli L, Gherlone E, Covani U
OSTEOTOMY AND MEMBRANE ELEVATION DURING THE MAXILLARY SINUS AUGMENTATION PROCEDURE. A COMPARATIVE STUDY: PIEZOELECTRIC DEVICE VS. CONVENTIONAL ROTATIVE INSTRUMENTS
 CLINICAL ORAL IMPLANTS RESEARCH, 2008
- Barone A, Cornelini R, Ciaglia R, Covani U
IMPLANT PLACEMENT IN FRESH EXTRACTION SOCKETS AND SIMULTANEOUS OSTEOTOME SINUS FLOOR ELEVATION: A CASE SERIES
 INTERNATIONAL JOURNAL OF PERIODONTICS AND RESTORATIVE DENTISTRY, 2008
- Barone A, Aldini NN, Fini M, Giardino R, Calvo Guirado JL, Covani U
XENOGRAFT VERSUS EXTRACTION ALONE FOR RIDGE PRESERVATION AFTER TOOTH REMOVAL: A CLINICAL AND HISTOMORPHOMETRIC STUDY
 JOURNAL OF PERIODONTOLOGY, 2008
- Covani U, Cornelini R, Barone A
BUCCAL BONE AUGMENTATION AROUND IMMEDIATE IMPLANTS WITH AND WITHOUT FLAP ELEVATION: A MODIFIED APPROACH
 INTERNATIONAL JOURNAL OF ORAL AND MAXILLOFACIAL IMPLANTS, 2008
- Cardaropoli D, Cardaropoli G
PRESERVATION OF THE POSTEXTRACTION ALVEOLAR RIDGE: A CLINICAL AND HISTOLOGIC STUDY
 INTERNATIONAL JOURNAL OF PERIODONTICS AND RESTORATIVE DENTISTRY, 2008
- Nannmark U, Sennerby L
THE BONE TISSUE RESPONSES TO PREHYDRATED AND COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS: A STUDY IN RABBIT MAXILLARY DEFECTS
 CLINICAL IMPLANT DENTISTRY AND RELATED RESEARCH, 2008
- Scarano A, Piattelli A, Perrotti V, Manzon L, Iezzi G
MAXILLARY SINUS AUGMENTATION IN HUMANS USING CORTICAL PORCINE BONE: A HISTOLOGICAL AND HISTOMORPHOMETRIC EVALUATION AFTER 4 AND 6 MTH
 CLINICAL IMPLANT DENTISTRY AND RELATED RESEARCH, 2009
- Perrotti V, Nicholls BM
RESORPTION PATTERN OF A PORCINE-DERIVED BONE SUBSTITUTE
 JOURNAL OF OSSEOINTEGRATION, 2009
- Calvo Guirado JL, Gomez Moreno G, Barone A, Cutando A, Alcaraz-Baños M, Chiva F, Lopez Mari L, Guardia J
MELATONIN PLUS PORCINE BONE ON DISCRETE CALCIUM DEPOSIT IMPLANT SURFACE STIMULATES OSTEOINTEGRATION IN DENTAL IMPLANTS
 JOURNAL OF PINEAL RESEARCH, 2009
- Scarano A, Piattelli M, Carinci F, Perrotti V
REMOVAL, AFTER 7 YEARS, OF AN IMPLANT DISPLACED INTO THE MAXILLARY SINUS. A CLINICAL AND HISTOLOGIC CASE REPORT
 JOURNAL OF OSSEOINTEGRATION, 2009
- Covani U, Marconcini S, Crespi R, Barone A
IMMEDIATE IMPLANT PLACEMENT AFTER REMOVAL OF A FAILED IMPLANT: A CLINICAL AND HISTOLOGICAL CASE REPORT
 JOURNAL OF ORAL IMPLANTOLOGY, 2009
- Calvo Guirado JL, Gomez Moreno G, Lopez Mari L, Ortiz Ruiz AJ, Guardia J
ATRAUMATIC MAXILLARY SINUS ELEVATION USING THREADED BONE DILATORS FOR IMMEDIATE IMPLANTS. A THREE-YEAR CLINICAL STUDY
 MEDICINA ORAL, PATOLOGIA ORAL Y CIRUGIA BUCAL, 2010
- Figueiredo M, Henriques J, Martins G, Guerra F, Judas F, Figueiredo H
PHYSICO-CHEMICAL CHARACTERIZATION OF BIOMATERIALS COMMONLY USED IN DENTISTRY AS BONE SUBSTITUTES - COMPARISON WITH HUMAN BONE
 JOURNAL OF BIOMEDICAL MATERIALS RESEARCH PART B: APPLIED BIOMATERIALS, 2010
- Grenga PL, Reale G, Cofone C, Meduri A, Ceruti P, Grenga R
HESS AREA RATIO AND DIPLOPIA: EVALUATION OF 30 PATIENTS UNDERGOING SURGICAL REPAIR FOR ORBITAL BLOW-OUT FRACTURE
 OPHTHALMIC PLASTIC AND RECONSTRUCTIVE SURGERY, 2009
- Crespi R, Cappare P, Gherlone E
DENTAL IMPLANTS PLACED IN EXTRACTION SITES GRAFTED WITH DIFFERENT BONE SUBSTITUTES: RADIOGRAPHIC EVALUATION AT 24 MONTHS
 JOURNAL OF PERIODONTOLOGY, 2009
- Rinna C, Reale G, Foresta E, Mustazza MC
MEDIAL ORBITAL WALL RECONSTRUCTION WITH SWINE BONE CORTEX
 THE JOURNAL OF CRANIOFACIAL SURGERY, 2009
- Cardaropoli D, Cardaropoli G
HEALING OF GINGIVAL RECESIONS USING A COLLAGEN MEMBRANE WITH A HEMINERALIZED XENOGRAFT: A RANDOMIZED CONTROLLED CLINICAL TRIAL
 INTERNATIONAL JOURNAL OF PERIODONTICS AND RESTORATIVE DENTISTRY, 2009
- Nannmark U, Azarmehr I
SHORT COMMUNICATION: COLLAGENATED CORTICO-CANCELLOUS PORCINE BONE GRAFTS. A STUDY IN RABBIT MAXILLARY DEFECTS
 CLINICAL IMPLANT DENTISTRY AND RELATED RESEARCH, 2010
- Barone A, Ricci M, Covani U, Nannmark U, Azarmehr I, Calvo-Guirado JL
MAXILLARY SINUS AUGMENTATION USING PREHYDRATED CORTICO-CANCELLOUS PORCINE BONE: HISTOMORPHOMETRIC EVALUATION AFTER 6 MTH
 CLINICAL IMPLANT DENTISTRY AND RELATED RESEARCH, 2010
- Calvo Guirado JL, Gomez Moreno G, Lopez Mari L, Guardia J, Marinez Gonzalez JM, Barone A, Tresguerras IF, Paredes SD, Fuentes Breto L
ACTIONS OF MELATONIN MIXED WITH COLLAGENIZED PORCINE BONE VERSUS PORCINE BONE ONLY ON OSTEOINTEGRATION OF DENTAL IMPLANTS
 JOURNAL OF PINEAL RESEARCH, 2010
- Scarano A, Piattelli A, Assenza B, Quaranta A, Perrotti V, Piattelli M, Iezzi G
PORCINE BONE USED IN SINUS AUGMENTATION PROCEDURES: A 5-YEAR RETROSPECTIVE CLINICAL EVALUATION
 JOURNAL OF ORAL AND MAXILLOFACIAL SURGERY, 2010
- Rossi R, Morales RS, Frascaria M, Benzi R, Squadrino N
PLANNING IMPLANTS IN THE ESTHETIC ZONE USING A NEW IMPLANT 3D NAVIGATION SYSTEM
 THE EUROPEAN JOURNAL OF ESTHETIC DENTISTRY, 2010
- Barone A, Orlando B, Tonelli P, Covani U
SURVIVAL RATE FOR IMPLANTS PLACED IN THE POSTERIOR MAXILLA WITH AND WITHOUT SINUS AUGMENTATION: A COMPARATIVE COHORT STUDY
 JOURNAL OF PERIODONTOLOGY, 2010
- Scarano A, Carinci F, Assenza B, Piattelli M, Murmura G, Piattelli A
VERTICAL RIDGE AUGMENTATION OF ATROPHIC POSTERIOR MANDIBLE USING AN INLAY TECHNIQUE WITH A XENOGRAFT WITHOUT MINISCREWS AND MINIPLATES: CASE SERIES
 In press 2010
- Pagliani L, Andersson P, Lanza M, Nappo A, Verrocchi D, Volpe S, Sennerby L
A COLLAGENATED PORCINE BONE SUBSTITUTE FOR AUGMENTATION AT NEOS IMPLANT SITES: A PROSPECTIVE 1-YEAR MULTICENTER CASE SERIES STUDY WITH HISTOLOGY
 CLINICAL IMPLANT DENTISTRY AND RELATED RESEARCH, 2010
- Santagata M, Guariniello L, Tartaro G
A MODIFIED EDENTULOUS RIDGE EXPANSION (MERE) TECHNIQUE FOR IMMEDIATE PLACEMENT OF IMPLANTS. A CASE REPORT
 THE INTERNATIONAL JOURNAL OF ORAL IMPLANTOLOGY, 2010
- Barone A, Ricci M, Calvo Guirado JL, Covani U
BONE REMODELLING AFTER REGENERATIVE PROCEDURES AROUND IMPLANTS PLACED IN FRESH EXTRACTION SOCKETS: AN EXPERIMENTAL STUDY IN BEAGLE DOGS
 CLINICAL ORAL IMPLANTS RESEARCH, 2010



Histology on maxillary sinus biopsy taken at 24 months. 48% new bone formation, 13% residual granules. Courtesy of Dr Roberto Rossi, Genova, Italy (biopsy) and Dr Ulf Nannmark, University of Göteborg, Sweden (histology).

mp3

ULTIMATE PERFORMANCE AND HANDLING

Pre-hydrated collagenated heterologous cortico-cancellous bone mix



Tecnoss s.r.l. is an innovative, globally active company that develops, produces and documents premium-quality xenogenic biomaterials by the brands Tecnoss® and OsteoBiol®.

Its 15 years of research led to its patent-protected production process that ensures neutralization of antigenic components in order to achieve biocompatibility, while preserving the natural collagen matrix inside the biomaterial.

Tecnoss® products comply with highest quality standards such as ISO13485 (notified body TÜV Rheinland), 93/42/EC (amended by 2007/47/EEC) and 03/32/EC (notified body CE 0373).

www.osteobiol.com


Tecnoss® Dental

Via Torino, 23
10044 Pianezza (TO) | Italy
Tel +39 011 9682823
Fax +39 011 9787577
info@tecnoss-dental.com

www.osteobiol.com

International Sales & Marketing

Authorized Distributor



MKT-MP3BRO031TEN