



ECO-XL
PLASTIC FREE WIPES
240 x 200mm

- Plastic free wipe
- Plant based fibres
- Renewable resources
- Reduced carbon footprint
- Alcohol free

V-WIPE ZERO

Surface Disinfectant Wipe BIODEGRADABLE[†]

EN14476
VIRUS TESTED

2 in 1
DISINFECTS
& CLEANS

STAY FRESH LID

CONTACT 60 SEC*

APPLE AROMA

PLANT BASED WIPE

KEY FEATURES

V-WIPE Zero ECO-XL wipes are the first alcohol free biodegradable[†] wipes in the Bossklein infection control range. Developed for use in clinical environments, on non invasive medical devices and general high traffic non-porous surfaces.

The brand new, exclusive wipe substrate material is constructed from plant based viscose fibres which are 100% renewable and sustainably sourced. Wipes are biodegradable[†] and totally plastic free.

The new flowpack style packaging helps to reduce the amount of plastic consumed by removing the dispenser tub. Each pack includes an easy open and sturdy stay fresh lid. This ensures wipes remain moist and ready for the next application.

The brand new, alcohol free Zero formulation is ideal to use on the majority of wipeable surfaces even those that are sensitive to alcohol, including dental chairs*

WHAT YOU NEED TO KNOW

- Biodegradable[†]. Plastic free wipes
- 100 extra large wipes
- Contact time from only 60 seconds
- Dual registration; Medical Device and Biocide
- EN14476 virus tested
- Safe to use on sensitive surfaces*
- 2 in 1; Cleans and disinfects
- Subtle apple aroma

Bossklein Zero contains a synergistic blend of two powerful biocides, is free from PHMB and contains detergent derived from renewable resources.

The new formula takes advantage of modern, safer and greener chemical components to create a more simple and effective product for the end user.

V-WIPE Zero ECO-XL wipes are ideal to use where time is critical, such as in between patients and also in high traffic areas.

COMPOSITION

Contains Aqua, Alkyl (C12-16) dimethylbenzyl ammonium chloride (ADBAC/BKC (C12-16)) (68424-85-1) 0.48% w/w (4.8g/kg), Didecylidimethylammonium Chloride (7173-51-5) 0.22% w/w (2.2g/kg), Perfumes.

INSTRUCTIONS FOR USE

Use sufficient wipes for each application. Wipe the clean surfaces to be disinfected so that they are completely moistened and let the disinfectant take action over 60 seconds (to achieve full efficacy leave for 2 minutes). Allow surface moisture to dry or wipe with clean tissue. For professional use only.

Use biocides safely. Always read the label and product information before use.

[†] Wipes are made from biodegradable plant fibres

* Always test an inconspicuous area first See website for full details. Do not use on copper or aluminium.



Complete hygiene solutions

DirectaDentalGroup
DIRECTA topdental ORSING parkell physics
RÖNIG TrollDental Sendoline CONTRACEZ

Topdental (Products) Ltd
12 Ryefield Way Silsden West Yorkshire BD20 0EF UK
Tel +44 (0)1535 652 750 www.topdental.co.uk

Surface Disinfectant Wipe BIODEGRADABLE†

ECO-XL

- Plastic free wipe material
- Made from renewable resources
- Extra large wipe 240 x 200mm
- Alcohol free formulation
- Pack of 100 wipes
- Safe to use on most surfaces**

BIODEGRADABLE† MATERIAL

Unlike most competitor wipes which are made from polypropylene, Bossklein ECO-XL wipes are constructed from 100% natural plant based fibres. After use the plastic free, biodegradable† wipes will naturally breakdown overtime.

The construction and large size of the 40gsm, plant based, viscose wipe material is extra durable and ideal for clinical disinfection. The plastic free substrate is designed to help transfer the alcohol free disinfectant solution on to the surface in question.



PLANT
BASED WIPE



ECO-XL
PLASTIC FREE WIPES
240 x 200mm

| MICROORGANISM | CONTACT TIME |
|-----------------------------|--------------|
| Bactericidal | |
| BS EN 16615:2015 | |
| BS EN 13727:2012 + A2:2015 | |
| Pseudomonas aeruginosa | 2 min |
| Staphylococcus aureus | 2 min |
| Enterococcus hirae | 2 min |
| Yeasticidal | |
| BS EN 16615:2015 | |
| BS EN 13624:2013 | |
| Candida albicans | 60 sec |
| Enveloped viruses | |
| BS EN 14476:2013 + A2:2019 | |
| Vaccinia virus | 60 sec |
| Hepatitis B virus (HBV) | 60 sec |
| Hepatitis C virus (HCV) | 60 sec |
| Hepatitis Delta virus (HDV) | 60 sec |
| HIV | 60 sec |
| Influenza virus | 60 sec |
| Coronavirus | 60 sec |
| Herpesviridae | 60 sec |
| Filoviridae (Ebola) | 60 sec |
| Rubella virus | 60 sec |
| Rabies virus | 60 sec |
| Measles virus | 60 sec |
| Poxviridae | 60 sec |

MICROBIOLOGICAL EFFICACY

Bactericidal, Yeasticidal and Virucidal against all enveloped viruses including Coronavirus, Hepatitis B and C and HIV. This highly effective formula has been tested to a variety of standards including EN 14476 for proven efficacy.

The new alcohol free Bossklein Zero surface disinfectant is effective against all listed (opposite) micro-organisms. This ensures complete peace of mind when it comes to clinical surface disinfection.

The following test standards have been achieved: EN 16615, EN 14885, EN 14476, EN 13727, EN 13624.

VARIANTS AND PACKAGING

BOS8001 Flowpack of 100 ECO-XL Wipes x 9 per box



BACTERICIDAL



YEASTICIDAL



VIRUCIDAL*

EN14476
VIRUS TESTED



† Wipes are made from biodegradable plant fibres

* Virucidal against all enveloped viruses

** Always test an inconspicuous area first. Do not use on copper or aluminium

