


Cavitron®

For hygienists who strive to deliver
efficiency, comfort, and performance,

**The Cavitron® family of inserts
offers the right tools for any
scaling environment.**


There's only one Cavitron®.



Cavitron® family of inserts work together to deliver the highest standard in ultrasonic scaling

Over the years, the Cavitron® brand has become synonymous with ultrasonic scaling. Dedicated to providing end-to-end solutions in dental hygiene, Cavitron® products have continued to evolve to best serve you, your patients, and your practice.

Today, Cavitron® products offer the latest technology in ultrasonic scaling units and a wide variety of scaling inserts to fit your needs. Each instrument in this comprehensive line of inserts is designed for the hygienist and the patient alike, to ensure efficiency, comfort, and performance during scaling procedures.



Ultrasonic scaling offers numerous benefits for practices, hygienists, and patients.

When selected and used appropriately, ultrasonic inserts are clinician- and patient-friendly, and offer ergonomic benefits over hand scaling.

Efficiency

Ultrasonic scaling can quickly and effectively remove even the most tenacious calculus.

Comfort


Both patients and dental professionals experience less pressure and greater comfort during ultrasonic scaling procedures compared to hand scaling. Without the repetitive motion of hand scaling, ultrasonic scaling can also lessen hand and wrist strain for dental professionals.

Performance

Cavitron® inserts are designed for ultimate longevity and reliability. Our inserts are made with durable materials and, when used properly, can withstand use and wear over time.

“Patients had a strong preference (74%) for ultrasonic scaling when compared to hand scaling. Particular preference for the ultrasonic scaling was registered for effective buildup removal, less irritating sound, clean feeling, less overall pain, more overall efficiency, and less mess.”

—Croft L. Int J Periodontics Restorative Dent. 2003.



In ultrasonic scaling, one insert is not enough

One insert can't do it all—that's why the Cavitron® family of inserts offers a broad selection of ultrasonic scaling inserts to fit the needs of your practice. Just like with hand scaling, some inserts are built for removing certain levels of deposit, while others are designed for working around specific types of tooth anatomy. Most patients require more than one type of insert during a single visit; following the guidelines for insert selection and use in this brochure can help with comfort and efficiency for your patients, and ultimately, your practice.

Use the chart to the right to help determine the best insert choices for certain procedures. Think about insert selection as a 2-step process based on the following assessments:

Type of deposit

Inserts are uniquely made to handle different levels of deposit. Consider the force needed to remove the deposits. A thicker insert with a larger tip is built for heavier levels of calculus, while thinner inserts with slim, adaptable tips are designed for the removal of light to moderate deposits or biofilm.










Anatomy of the treatment site

Look for inserts that are specially designed for the type of tooth anatomy you will be working on. Straight inserts are made for flat anatomy, while curved inserts will give you the best results in furcations and concavities.




Insert Use Guide*

Anatomy of Treatment Site

		For Flat Anatomy or Contoured Anatomy (With Heavier Deposits)		For Contoured Anatomy		For Flat Anatomy					
Calculus	Heavy	Dark Blue		Light Blue		Light Blue					
	Moderate	Dark Blue		Teal		Purple					
	Slight	Light Blue		Teal		Purple					
Biofilm	Heavy	Light Blue		Teal		Purple					
	Moderate	Light Blue		Teal		Purple					
	Slight	Light Blue		Teal		Purple					
Non-Visible Biofilm	Heavy	Light Blue		Teal		Purple					
	Moderate	Light Blue		Teal		Purple					
	Slight	Light Blue		Teal		Purple					
		 Powerline™-10 (FSI®-10)	 Powerline™-3 (FSI®-3)	 Slimline® L-Curve	 Slimline®-10	 THINsert®	 Powerline™-100 (FSI®-100)	 Powerline™-1000 (FSI®-1000)	 Slimline® R-Curve	 Slimline®-1000	

*These are only recommendations, based upon qualitative indications of calculus and biofilm levels.









All Cavitron® inserts are designed for optimal comfort and a range in style

Focused Spray® water delivery system

- Designed to deliver targeted lavage
- Enhances visibility of work area
- Improves patient comfort

Insert Design Guide

Cavitron® inserts are designed to give you options. Consider your needs in bend and cross-section when making your selection.

Bends	Cross-Section	Thicker	Thinner
 1000		Powerline™-1000 (FSI®-1000)	Slimline®-1000
 100		Powerline™-100 (FSI®-100)	n/a
 10		Powerline™-10 (FSI®-10)	Slimline®-10

Get your favorite insert in Bellissima® soft grip!

- Increased comfort and control with a soft, dimpled grip
- Available for Slimline®-10 and Curve inserts and all Powerline™ (FSI®) inserts



Standard-sized inserts are designed for removal of heavier deposits

Powerline™ (FSI®) series of inserts

- Specifically designed for efficiently removing heavier deposits
- The flat beaver tail tip on the Powerline™-3 (FSI®-3) insert is designed for use on calculus bridges and removal of heavy extrinsic staining



FOCUSED SPRAY®
PowerLINE™

“Ultrasonic scalers are often considered to be less strenuous for the operator and more comfortable for the patients than hand currettes.”

—Jepsen S. J Clin Periodontol. 2004.



Thinner Cavitron® inserts allow for improved access and adaptability

Slimline® series of inserts

- Slim design is 30% thinner than the Powerline™-1000 (FSI®-1000) insert
- Thin design provides improved subgingival access
- Slimline® Curve inserts allow adaptation to root anatomy and are designed for use in furcations and concavities
- Designed to remove light to moderate deposits



FOCUSED SPRAY®
slimLINE®

The newest insert in the Slimline® family

Slimline®-1000

Featuring a triple-bend design with beveled edges for better adaptation around line angles



THINsert® ultrasonic insert


- Delivers ultrasonic power with a tip that is 47% thinner than Slimline®-10
- Unique 9° backbend design allows for easier adaptation
- Provides access in difficult areas, including misalignments, interproximal surfaces, and areas of tight tissues, without loss of tactile sensation
- Designed for light to moderate calculus and deplaquing



THINsert®
FOCUSED SPRAY®

*“Subgingival debridement may be completed in less time
with ultrasonic/sonic than with hand instruments.”*

—Tunkel J. J Clin Periodontol. 2002.



The Cavitron® specialty line of inserts conform to your specific scaling needs

SofTip™ Ultrasonic Implant Insert for plaque removal around titanium implants and abutments

- Gently removes light to moderate plaque and calculus without damaging implant titanium abutment surfaces
- Provides benefits of ultrasonic lavage during implant maintenance procedures



DIAMONDCOAT® Ultrasonic Insert for surgical applications

- Unique diamond coating removes moderate to heavy calculus 3 times faster than traditional ultrasonic inserts†
- Provides excellent tactile sense, access, and control during periodontal surgical applications
- Effective for the removal of overhanging margins and re-contouring of dental restorations



Air Polishing available for stain removal

- **NEW** Cavitron® JET Plus units now feature prophylaxis mode auto cycles that automatically cycle between air polish and rinse



“It appears that air polishing is the least damaging and most efficient means of removing stain on enamel.”

—Gutmann ME. *J Dent Hygiene*. 1998.

Newly redesigned Cavitron® Plus and Cavitron® JET Plus ultrasonic scaling units including Tap-On™ Technology



†Data on file.

Monitoring the wear of your ultrasonic inserts is key to ensuring optimal performance

Using an insert that is excessively worn will cause it to operate at substandard efficiency. It can also cause you to use excessive pressure while scaling, which results in discomfort for both you and the patient.


Regularly checking the wear of your ultrasonic inserts and replacing them when necessary can help maximize efficiency and comfort for you and your patients.

Use your Efficiency Indicator to check the wear on your inserts.

Tip extends beyond blue line: **Optimum Efficiency**

Tip touches blue line: **Reorder** (25% efficiency loss)

Tip touches red line: **Discard** (50% efficiency loss)



Call your DENTSPLY Professional representative to request an Efficiency Indicator or to reorder inserts.

“Clinicians who retain their inserts for years may find their instrument tips become reduced in length through wear, potentially leading to a detrimental change in the performance of the scaling system.”

—Lea SC. *J Clin Periodontol.* 2006.

Recommended Inserts Guide

Use this guide to select inserts based on the types of procedures your practice generally performs on patients. The bolded insert in each category denotes the recommended insert.

	Prophylaxis	SPT (Supportive Periodontal Therapy)/Perio Maintenance	SRP (Scaling and Root Planing)/Debridement
Procedure Code	D1110	D4910	D4341
Scaling	Slimline®-1000	Slimline®-1000	Powerline™-100
	Powerline™-1000	Powerline™-1000	Powerline™-10
	Powerline™-100	Powerline™-100	Powerline™-1000
	Powerline™-10	Powerline™-10	Slimline®-1000
<i>Calculus Bridge</i>	Powerline™-3	N/A	Powerline™-3
Debridement	THINsert®	Slimline®-10	Slimline®-10
	Slimline®-10	THINsert®	THINsert®
	Slimline®-1000	Slimline®-1000	Slimline®-1000
Roots/Contours	Slimline® L/R Curve	Slimline® L/R Curve	Slimline® L/R Curve
Specialty Implants	SofTip™	SofTip™	SofTip™
<i>Flap</i>	N/A	N/A	DIAMONDCOAT®

Cavitron® products available through authorized distributors

PowerLINE™ Inserts

Item	SKU#
Powerline™-10	25K 80293
	30K 80294
	30K Bellissima® soft grip 81501
Powerline™-100	25K 80795
	30K 80798
	30K Bellissima® soft grip 81502
Powerline™-1000	25K 80796
	30K 80799
	30K Bellissima® soft grip 81503
Powerline™-3	25K 80794
	30K 80797
	30K Bellissima® soft grip 81500

slimLINE™ Inserts

Item	SKU#
Slimline®-10	25K 80392
	30K 80395
	30K Bellissima® soft grip 81510
Slimline®-1000	25K 81569
	30K 81570
Slimline® L-Curve	25K 80393
	30K 80396
	30K Bellissima® soft grip 81511
Slimline® R-Curve	25K 80394
	30K 80397
	30K Bellissima® soft grip 81512
Slimline®-10/ Slimline® L-Curve/ Slimline® R-Curve Combo Pack	25K 80398
	30K 80399
	30K Bellissima® soft grip 81513
Slimline® L-Curve/ Slimline® R-Curve Combo Pack	25K 80400
	30K 80401
	30K Bellissima® soft grip 81514

THINsert®

Item	SKU#
THINsert® 30K	81551

Specialty Inserts

Item	SKU#
SofTip™ Ultrasonic Implant Insert	30K 90411
DIAMONDCOAT®	30K 90243

For more information or to place an order, [contact your local DENTSPLY Professional distributor](#), call DENTSPLY Professional Customer Service at **1.800.989.8826**, or visit www.professional.dentsply.com/cavitron.