

Made in Sweden,
used worldwide.



Practice Box

The practice box is ideal for keeping interdental brushes within easy reach when trying out the correct size and instructing the patient. The box has eight removable cases that each holds 15 interdental brushes original or extra soft. Tip: Use the multi-pack to fill the practice box.



Multi-pack

TePe Interdental Brushes are available in packs of 25, all the same size and with plastic caps. Tip: Store the brushes in the practice box.



Mixed Pack

All models of TePe's interdental brushes are available in a mixed pack, convenient when trying out the correct brush size.



Information Material

TePe has developed a selection of instructive information material, useful when informing and educating patients.



Dental Model

The handcrafted dental model, intended for patient instruction, has colour coded interdental spaces matching the sizes of TePe's interdental brushes.



TePe Interdental Gel with Fluoride

TePe's fluoride gel offers extra caries protection between the teeth. The gel can easily be applied with an interdental brush.



Storage and Travel

The micro holder offers hygienic storage of interdental brushes in the bathroom.

The travel case fits six TePe Interdental Brushes original or extra soft.



© TePe Oral Hygiene Products AB. 180GB, 2011.



TePe Interdental Brushes
Choice of size, texture and design


TePe has been promoting healthy smiles since 1965. At our headquarters in Malmö, Sweden, we design and produce interdental brushes, toothbrushes and other oral hygiene products that are available in more than 50 countries worldwide. All TePe products are developed in collaboration with dental expertise.

TePe Oral Hygiene Products AB


Bronsåldersgatan 5, SE-213 76 Malmö, Sweden
Tel +46 (0)40-670 11 00 • Fax +46 (0)40-670 11 11
info@tepe.com • www.tepe.com

TePe Interdental Brushes


Use of interdental brushes has proven to be the most efficient way of removing plaque between the teeth. TePe's interdental brushes are developed in collaboration with Swedish and international dental professionals to certify the highest efficiency, comfort and quality in every detail. To minimise the effect on the environment, the handles are made of recyclable polypropylene.



It is important to select the correct size. The filaments should fill the interdental space.



All sizes and models have plastic coated wire for gentle cleaning.



TePe Interdental Brush – the original

The original range comprises nine colour coded sizes to fit narrow as well as wider interdental spaces. The user-friendly handle offers a stable and convenient grip facilitating cleaning with controlled movements.




The plastic cap can be used for hygienic protection of the brush or to extend the handle length.

TePe Interdental Brush **extra soft** – the softer choice

This unique range of interdental brushes has specially selected, extra soft filaments for gentle cleaning. These brushes are recommended after surgery or for patients with sore gingiva, sensitive teeth or dry mouth. They are also suitable for persons who prefer a softer brush. Available in six colour coded sizes.

Areas of use:



Curve the brush gently when cleaning between the posterior teeth. Do not straighten it again. Use another, straight brush for the anterior teeth.

Missing tooth.

Distal surface of the last tooth.

Implants.

Orthodontic appliances.

TePe Interdental Brushes comply with ISO Standard

TePe's interdental brushes are thoroughly tested in terms of quality and durability and fulfil all requirements of the international standard, ISO 16409:2006.

Original									
Extra soft									
Angle									
Colour	Pink	Orange	Red	Blue	Yellow	Green	Purple	Grey	Black
Brush size (0-7)	0	1	2	3	4	5	6	7	7
Wire size (mm)	0,4	0,45	0,5	0,6	0,7	0,8	1,1	1,3	1,5
PHD (mm)	0,6	0,7	0,9	1,1	1,3	1,6	2,0	2,6	3,4

The standard places the interdental brush in a brush size category (0-7) based on the measurement of the smallest hole which the brush can pass through; passage hole diameter, PHD.



TePe Angle™ – the angled option

An easy to use interdental brush, developed for excellent access to all interdental spaces. The slender brush head is angled for easy reach e.g., between the posterior teeth and from the palatal and lingual sides. The long, flat handle provides an ergonomic grip, allowing cleaning with steady and controlled movements.



Available in six colour coded sizes.

Areas of use:



Posterior cleaning from the buccal side.

Posterior cleaning from the lingual side.

Anterior cleaning from the lingual side.

Cleaning of implants from the buccal side.

Cleaning of implants from the palatal side.

Interdental brushes have a higher plaque reducing capacity on approximal surfaces than toothpicks or dental floss.

Särner B, et al.

Recommendations by dental staff and use of toothpicks, dental floss and interdental brushes for approximal cleaning in an adult Swedish population. *Oral Health Prev Dent 2010;8:185-194.*

Interdental brush use as an adjunct to toothbrushing removes more dental plaque than toothbrushing alone and is more effective than dental floss or wooden sticks in removing dental plaque. Studies have shown that the use of interdental brushes has a positive effect on plaque scores, bleeding scores and probing pocket depth.

Slot DE, et al.

The efficacy of interdental brushes on plaque and parameters of periodontal inflammation: a systematic review. *Int J Dent Hyg 2008;6:253-264.*