

Surface Disinfectant Wipe Comparison

All V-WIPE products have been developed for use inside clinical environments on both non invasive medical devices and general high traffic, non-porous surfaces.



**V-WIPE
Classic**



**V-WIPE
Ultra**



**V-WIPE ZERO
Classic**



**V-WIPE ZERO
Ultra**



**V-WIPE ZERO
Universal-XL**



**V-WIPE ZERO
Eco-XL**

	V-WIPE Classic	V-WIPE Ultra	V-WIPE ZERO Classic	V-WIPE ZERO Ultra	V-WIPE ZERO Universal-XL	V-WIPE ZERO Eco-XL
SKU(s)	BOS9783 Tub BOS9784 Refill	BOS1000 Tub BOS1001 Refill	BOS2000 Tub BOS2001 Refill	BOS3000 Tub BOS3001 Refill	BOS4001	BOS8001
Number of wipes per pack	200	160	200	160	200	100
Wipe Size	200x200mm	175x155mm	200x200mm	175x155mm	267x210mm	240x200mm
Wipe Thickness/Weight	18gsm	40gsm	18gsm	40gsm	34gsm	40gsm
Wipe Material	100% Polypropylene	100% Polypropylene	100% Polypropylene	100% Polypropylene	100% Polypropylene	100% Spunlace Viscose
Material Construction	SMS Triple Layer Spunbond/ Meltblown/Spunbond	Meltblown Diamond Pattern	SMS Triple Layer Spunbond/ Meltblown/Spunbond	Meltblown Diamond Pattern	SMS Triple Layer Spunbond/ Meltblown/Spunbond	Hydroentangled parallel lapped natural fibres
Total Disinfectant	620g (3.1g per wipe)	620g (3.9g per wipe)	600g (3g per wipe)	550g (3.4g per wipe)	1050g (5.25g per wipe)	543g (5.4g per wipe)
Area Coverage (per wipe)	Small to medium	Medium to large	Small to medium	Medium to large	Larger	Larger
Active Ingredient/s	Alcohol - Ethanol (63%)	Alcohol - Ethanol (63%)	Quaternary Ammonium Compounds	Quaternary Ammonium Compounds	Quaternary Ammonium Compounds	Quaternary Ammonium Compounds
Aroma	Fresh Lemon	Fresh Lemon	Subtle Apple	Subtle Apple	Subtle Apple	Subtle Apple
EN14476 Tested	Yes	Yes	Yes	Yes	Yes	Yes
EN16615 Tested	Yes	Yes	Yes	Yes	Yes	Yes
Contact Time	60 sec	60 sec	60 sec	60 sec	60 sec	2 mins
2 in 1: Cleans & Disinfects	No	No	Yes	Yes	Yes	Yes
Summary	<ul style="list-style-type: none"> Quick drying Streak free Most economical wipe Only for use on alcohol resistant surfaces 	<ul style="list-style-type: none"> Quick drying Streak free Premium feel wipe Only for use on alcohol resistant surfaces 	<ul style="list-style-type: none"> All non porous surfaces Cleans & disinfects Alcohol Free Low streaking Most economical wipe 	<ul style="list-style-type: none"> All non porous surfaces Cleans & disinfects Alcohol Free Low streaking Premium feel wipe 	<ul style="list-style-type: none"> All non porous surfaces Cleans & disinfects Alcohol Free Low streaking Largest wipe Strongest wipe 	<ul style="list-style-type: none"> Biodegradable wipe Renewable resources All non porous surfaces Cleans & disinfects Alcohol Free Low streaking