

Kerr News UK



1st April - 30th June 2024

Follow Us




@kerrdentaluk

Kerr[™]

Shop Now on store.kerrdental.com

- Choose from award winning product lines
- View exclusive online promotions
- Receive regular rewards



**Register today and receive
a fabulous welcome gift**



Kerr[™] Innovation Symposium — 2024


THE BELFRY
HOTEL & RESORT

Sutton Coldfield,
20th - 21st September 2024

Shaping the Future of Oral Health

Join us at the Kerr Innovation Symposium for a full day of hands-on sessions and lectures designed to elevate your dentistry innovation journey. With three tracks led by renowned speakers, including restorative and endodontic workshops, modern hygiene methods, and inclusive amenities, this event is your pathway to the next level in your career.



Prof. Philippe Sleiman




Dr. Jordi Manauta

STYLE ITALIANO[™]



Dr. Tracey Lennemann



Scan here and secure
your spot at the forefront.

Demi™ Plus

Light Emitting Diode (LED) Curing Light

Thanks to the PLS technology, Demi Plus intermittently shifts the output intensity from a base of 1,100 mW/cm² to a peak of 1,330 mW/cm². This enables the Demi Plus to deliver a 5-second cure for A3 or lighter shades.

Enhanced Body Design

This curing light's exterior is made from Valox resin, a high-performance plastic that will protect and extend the life of your unit. In addition, the upgraded closed housing allows for easier cleaning and disinfection. Improved Cooling System and fan free for quiet operation.

True 5-second cure

Consistent output of Periodic Level Shifting (PLS) delivers a true 5-second cure for composite shades A3 and lighter.

Never shuts down

Whisper-quiet fan and unique thermal control system keeps the unit cool for uninterrupted operation. Extends the overall life of the unit.

Advanced battery technology

The lithium-ion battery technology can yield 300 10-second cures with a single charge.

Lightweight, ergonomic design

Multiple holding options reduce hand and arm stress.



Buy 30
**refills of any
composite,
cement or bond
and receive
Demi Plus
FREE***

910860-2 - Demi Plus

OptiBond™ FL

Filled, Light-Cure Total-Etch Dental Adhesive



26684 E
OptiBond FL Kit

£ 153

Buy 1 get Gel Etchant
FREE*



OptiBond™ Universal

Single-Component Universal Dental Adhesive



36519
OptiBond Universal
Refill 5ml

£ 89

Buy 2
get NX3 Dual
Cure Syringe
FREE*



OptiBond™ Solo Plus

Single Component Total-Etch Dental Adhesive



29692
OptiBond Solo Plus Refill 5ml

£ 74


Buy 1 get Gel Etchant
FREE*



OptiDisc™

High Performing Polishing and Contouring Disc

OptiDisc™ finishing and polishing discs allow you to obtain a final polish equal to that of natural dentition on all kinds of restorations. Made of durable polyester impregnated with aluminium oxide particles, the four grit system achieves a long-lasting result. A colour-coded system enables immediate identification of the grit type.



AVAILABLE IN THREE
DIFFERENT DIAMETER SIZES



Buy 2
get 1
FREE*


OptiDisc Refills - 100 pcs (except Extra-Course, 80 pcs)

- | | |
|---|---|
| 4197 - 80x 9.6mm - Extra-Coarse | 4185 - 100x 12.6mm - Fine |
| 4181 - 100x 9.6mm - Coarse/Medium | 4186 - 100x 12.6mm - Extra-Fine |
| 4182 - 100x 9.6mm - Fine | 4199 - 80x 15.9mm - Extra-Coarse |
| 4183 - 100x 9.6mm - Extra-Fine | 4191 - 100x 15.9mm - Coarse/Medium |
| 4198 - 80x 12.6mm - Extra-Coarse | 4192 - 100x 15.9mm - Fine |
| 4184 - 100x 12.6mm - Coarse/Medium | 4193 - 100x 15.9mm - Extra-Fine |

£ **26**

Herculite™ XRV Ultra

Universal Nanohybrid Dental Composite



33860
Herculite XRV Ultra Mini Kit Syringe

3x 4g syringes in shades enamel A2, A3, dentine A2, 1x 5ml OptiBond™ Solo Plus bottle, 1x 3g gel etchant, 25x mixing wells, 50x applicator tips.

£ **94**



34071
Herculite XRV Ultra Mini Kit Unidose

30x 0.2g unidoses in shades enamel A2, A3, dentine A2, 1x 5ml OptiBond™ Solo Plus bottle, 1x 3g gel etchant, 25x mixing wells, 50x applicator tips.

£ **87**



Buy 1 get OptiDisc Mini Kit **FREE***

Adapt SuperCap™ Matrices

Universal Matrix Tensioning System



£ 50



Buy 4 SuperCap refills
get SuperMat Handle FREE*

Adapt SuperCap Pack of 50

2161 - Steel - 0.030mm/5.0mm
2162 - Steel - 0.030mm/ 6.3mm
2181 - Steel - 0.038mm/5.0mm

2182 - Steel - 0.038mm/6.3mm
2171 - Trans. - 0.050mm/5.0mm
2172 - Trans. - 0.063mm/5.0mm



SuperMat Intro Kit

Universal Matrix Tensioning System

2150
SuperMat Intro Kit

£ 163

Buy 1 get OptiDam
Anterior Intro Kit FREE*

TempBond™

Temporary Cement



61087 - TempBond Standard Package
60256 - TempBond NE Standard Package

£ 29

Buy Any 3 TempBond get Maxcem
Elite Mini Kit
FREE*



SimpliShade™

Simplified Universal Composite System

Simplify Shade-matching

Shade selection made easier and faster.
SimpliShade matches all 16 VITA classical shades with just 3 shades: Light, Medium, Dark.

Adaptive Response Technology (ART)


The ART filler network, a proven technology found in Harmonize, helps you achieve lifelike restorations with exceptional strength and unmatched esthetics.

Streamline Inventory Management

Save time and money with minimal shades to meet almost all your composite needs. Declutter your shelves, avoid the burden of re-ordering too many different part numbers, and reduce or eliminate expired products.

Universal For All Cavity Classes

Use SimpliShade with confidence in all cavity classes without the need for a blocker, unlike other competitors.




Excellent Radiopacity

The more radiopaque, the easier it is to distinguish the restoration from the tooth structure under a radiograph. SimpliShade at 300% AI is almost twice as radiopaque as competition.

Optimal Translucency

SimpliShade provides ideal opacity to mimic tooth structure, to better mask underlying tooth stains, and blend well with surrounding dentition, resulting in a more lifelike restoration.



Unidoses 20x 0,25g
37006 - Light
37007 - Medium

£ **59**

Unidoses 10x 0,25g
37008 - Dark
37009 - Bleach White
37010 - Universal Opaque

£ **39**

Syringe 4g
37001 - Light
37002 - Medium
37003 - Dark
37004 - Bleach White
37005 - Universal Opaque

£ **49**

Buy 2
**get OptiDisc
Mini Kit FREE***



Maxcem™ Elite Universal

Self-Etch, Self-Adhesive Resin Dental Cement

2x Automix Syringe 5g

34058 - Clear

34060 - White Opaque

34059 - White

34061 - Yellow

37003 - Dark

Buy 2
get 1
FREE*

£ 100




Maxcem™ Elite Chroma Universal

Self-Etch, Self-Adhesive Resin Dental Cement with Cleanup Indicator

Buy 2
get 1
FREE*

£ 128



2x Automix Syringe 5g

36372 - Clear

36374 - Yellow

36373 - White

NX3 Nexus™ Third Generation

Universal Adhesive Resin Dental Cement

Automix Syringe 1.8g

33646 - Bleach

33643 - Clear

33647 - Opaque

33644 - White

33645 - Yellow

£ 85

Light-Cure Syringe 5g

33651 - Bleach

33648 - Clear

33652 - Opaque

33649 - White

33650 - Yellow

£ 42

Buy 3 get
OptiBond
Universal
FREE*



SimpliShade™ BULK FILL & BULK FILL FLOW

ONE SHADE IS ALL IT TAKES



Single-shade Simplicity meets Superior Adaptability

Beautiful restorations just got simpler, with SimpliShade™ Bulk Fill and SimpliShade Bulk Fill Flow. Skip the capping layer with a one-step placement and reach every area of the tooth's anatomy — with exceptional durability and depth of cure. Requiring just one shade to match all 16 VITA® classical shades, SimpliShade Bulk Fill and Bulk Fill Flow reduce inventory requirements while saving you time for your patients.



ART³: A NEW WAY TO BULK FILL



SimpliShade™ BULK FILL



One Step

Eliminate the need for a capping layer, complete restorations with a SingleFill™ composite. With its enhanced curing mechanism, SimpliShade Bulk Fill makes it easy to achieve polymerization depths up to 5mm in a single step.



Simplicity with Durability

SimpliShade Bulk Fill offers low wear rate, high flexural strength, and comparable shrinkage stress relative to other bulk fill brands on the market, keeping your restorations stronger and long-lasting.



Chameleon Effect

Special nanofiller architecture mimics the natural enamel structure and interacts with light for virtually indistinguishable blending.



SimpliShade™ BULK FILL FLOW



Save Time

Simplify your workflow with just one step to cure up to 5 mm, before adding a capping layer.



Go with the Flow

Fill every area of the tooth's anatomy with a self-leveling bulk fill flowable that easily adapts to the cavity floor.




Low Shrinkage

The low polymerization-related shrinkage and shrinkage stress of SimpliShade Bulk fill Flow ensure long-lasting restorative results.



Add Confidence

Restore confidently with a flowable bulk fill composite that offers comparable or higher mechanical strength than the competitors.



COMING SOON


Contact Your Local Kerr Product Specialist to arrange a free sample.

elements™ Connect

– perfection in every detail



Dr. Maciej Mikołajczyk DDS, PhD is a senior Assistant in English Dental Division of Medical University in Lodz, Poland. Author of 5-time book series “Endodontics in dental practice” and author and co-author of multiple scientific papers and lectures concerning minimally invasive adhesive dentistry, use of lasers and microscopy in dentistry. He works in his own private practice in Lodz, specializes in esthetic reconstructions of teeth and lectures in many European countries.




Elements Connect is a cordless Endomotor with a set of quite unique features to make every day work of the endodontist on root canal preparation easier.

One notable feature to highlight initially is the excellent ergonomics of the device. The motor is lightweight and very well balanced. It fits the operator's hand stably, without causing any muscle strain, joint tension, or fatigue.

The grip has been shaped to make maneuvering the handpiece in the patient's mouth as easy as possible. The unique profile of the handpiece's surface offers excellent support for the fingers during preparation and vertical movements.

The large, contrasting and ergonomically located screen enables continuous observation of the operating parameters of the endodontic file. It also allows for precise control of the torque and position of the file within the canal, while simultaneously offering the option to activate auto-reverse or auto-reverse/stop programs after reaching the limit parameters for a given endodontic file system.

In my many years of practice, I have never worked with such an ergonomically well-designed endodontic motor that puts so little strain on my hand during the preparation of root canals.



The size of the head of the device, with an inserted endodontic file, is among the smallest in its class of devices, yet the design of the head itself provides a solid grip for any standard endodontic file. The handpiece has been designed in a way that allows for maximum control over the torque, RPM and the range of rotation angles. This level of control is particularly important when working with reciprocal and adaptive motion, both possible with elements Connect. Changes in angular value and rotation direction occur multiple times per second during file movement and must be executed in a repeatable and fully controlled manner.

I use two lengths of Zenflex files in my practice. When working with them and the elements Connect micromotor, I encountered no difficulties accessing any root canal, regardless of the tooth's position in the patient's mouth.



The built-in ability to work with an external Apex Connect endometer enhances the root canal preparation process by providing continuous control of the endodontic file's position in relation to the apex of the root canal. Additionally, Apex Connect is a fully functional endometer that enables precise determination of the working length using manual canal tools during canal negotiations before rotational preparation begins.

When connected to elements Connect, it allows for continuous monitoring of the file's position in the canal during canal preparation and activation of the auto-reverse program upon reaching the working length.

I like the versatility of combining an endometer and an endodontic micromotor. This combination offers the usability and ergonomics of a free-standing endometer while providing continuous control of the working length during root canal preparation when connected to the motor.

The built-in programming module has saved operating parameters for reciprocal movement, adaptive movement, and all Kerr endodontic file systems. All settings are quickly accessible with a single press of a button, and they can be easily modified in any way as needed. Additionally, regardless of the saved file systems, it's possible to program operating parameters for any endodontic files available on the market, whether rotary and reciprocal. This makes elements Connect a versatile tool for every dentist performing endodontic treatment.

Easy access to each program and the ability to quickly change operating parameters greatly simplify my work with this motor and various endodontic files, whether preparing the entrance to the canal, obtaining a glide path, or final preparation of the root canal.



The unique system of two batteries and a properly designed charger/stand guarantees the constant readiness of the motor for my work. During root canal preparation, one of the batteries can always be charged on the base, ensuring uninterrupted usage. Thanks to the quick contact charging of the spare battery, the micromotor does not have to be stored on the base at all times, it can remain available on the assistant's or on the doctor's table, ready for immediate use.

Kerr has already successfully used similar solutions in their curing lights, and the elements Connect motor continues to work perfectly in my practice. The micromotor's ability to remain readily available ensures seamless workflow. I no longer have to remember to place the motor on its stand; instead, it's always within reach whenever I need it.



elements™ Connect & Apex Connect

Cordless Endodontic Motor and Electronic Apex Locator

Your safe path

**Elevate your confidence
while enhancing patient safety.**

Greater visibility into canal depth and exceptional reliability during every shaping procedure.

Add confidence to your shaping, with elements Connect and Apex Connect. Designed to work seamlessly together to support patient safety while simplifying your procedures. The intuitive display provides real-time working length feedback to minimize error, while the motor delivers reliable performance with virtually any file system.

**All you need to support the success
of your shaping procedures from start to finish.**



scan the QR Code
for more information
on elements™ Connect
& Apex Connect



Versatile Motions

With Adaptive, Reciprocation and Rotary motions in one motor, elements Connect allows offices to maintain their existing techniques and protocols.

Preset

9 Kerr preset file settings,
1 reciprocating setting and
4 customizable settings.

Continuous Use



Two removable batteries are included to ensure one is always charged so you can complete your procedures safely without interruption.

Built-in Accuracy

Ensure measurement accuracy with Apex Connect's automatic self-calibration feature that activates every time the device is turned on.

Shape with Confidence

View real-time working length information during procedures with a smart light indicator that changes color, as well as audible indicators, as you reach the apex.

Buy **60 packs of**
Kerr NiTi Files**
(TF, TFA, K3 or Zenflex)
and receive
elements Connect
FREE*

Buy **100 packs of**
Kerr NiTi Files**
(TF, TFA, K3 or Zenflex)
and receive
elements Connect
& Apex Connect
FREE*



Buy **Apex Connect**
and receive
20 packs of K-Files
FREE*

K-Files™

Stainless steel files

**CONTACT
YOUR KERR REP
FOR SIZES**

£ **10**



*Buy 3
get 1
FREE**

K-Flex™

Stainless steel files

**CONTACT
YOUR KERR REP
FOR SIZES**

£ **10**




*Buy 3
get 1
FREE**

K3™ Files

NiTi Rotary Files

**CONTACT
YOUR KERR REP
FOR SIZES**

£ **47**



*Buy 3
get 1
FREE**

TF Adaptive™

NiTi Twisted Files

**CONTACT
YOUR KERR REP
FOR SIZES**

£ **45**



*Buy 3
get 1
FREE**

ZenFlex™

NiTi Rotary Shaping Files

Maximum Strength. Minimally Invasive

Featuring an ideal balance of strength and flexibility, ZenFlex lets you reach the full working length with reduced risk of breakage while preserving more tooth integrity.



High Cutting Efficiency

Designed with triangular cross section with proprietary heat treatment, ZenFlex exhibits high cutting efficiency that lets you reach the full working length of even complex root geometries.



Minimally Invasive

A maximum flute diameter of 1mm and a non-cutting tip work together to minimize transportation and preserve the tooth's structural integrity.



Exceptional Strength

Resists torsional stress and cyclic fatigue resulting in minimal breakage and reduced transportation.



Excellent Flexibility

Proprietary heat-treatment design delivers excellent controlled memory and flexibility for up to 90° curves.

*Buy 3
get 1
FREE**

£ **38**

**CONTACT
YOUR KERR REP
FOR SIZES**



SHARING KNOWLEDGE ADDING VALUE TO YOUR PRACTICE.

WEBINARS | HANDS-ON COURSES | LECTURES



KERR IS COMMITTED TO OFFERING DENTISTS CONTINUOUS LEARNING OPPORTUNITIES

Featuring renowned international speakers; register and join our events on both endodontic & restorative dentistry.

APPLY NOW! Each educational program is designed for dentists at all levels, whether you're still in your early years of training and want some practice with innovative equipment, or you're a dental specialist looking to brush up on specific approaches or widen your subspecialty knowledge.

UPCOMING HANDS-ON COURSES


Tips & Tricks for a Predictable Root Canal Treatment



May 17th | 9:00 - 17:00 | Edinburgh

Presented by Dr. Arvind Sharma

**MORE DETAILS
AND REGISTRATION** ▶



May 24th | 9:00 - 17:00 | Newport Wales

Presented by Dr. Adi Moran

**MORE DETAILS
AND REGISTRATION** ▶



June 7th | 9:00 - 17:00 | Manchester

Presented by Dr. Arvind Sharma

**MORE DETAILS
AND REGISTRATION** ▶



SPEAKERS ARE LOOKING FORWARD TO MEETING YOU!
FOR MORE INFORMATION, PLEASE CONTACT YOUR LOCAL REPRESENTATIVE

elements™ IC System

Cordless Obturation System

Worth
over
£650

Buy elements IC *get* 5 boxes
of Gutta Percha cartridges,
2 Heat Pluggers, 3 Hand Pluggers
and 1 TubliSeal EWT **FREE***



973-0600-TYPEF - elementsIC Obturation System

5x Gutta Percha cartridges

972-2500 - 23GA Gutta Percha cartridge (medium body)

972-2501 - 25GA Gutta Percha cartridge (light body)

972-2502 - 23GA Gutta Percha cartridge (heavy body)

3x Hand Pluggers

974-0060 - Hand plugger (size # 0)

974-0058 - Hand plugger (size # 1)

974-0059 - Hand plugger (size # 2)

1x Heat Plugger

952-0005 - Buchanan Heat Plugger .08 (FM)

25903 - TubliSeal EWT

£ **3,050**

SealTight™

Disposable Air/Water Syringe Tips

Seal-Tight™ is a disposable syringe tip with a built-in rubber O-ring seal making syringes work like new every time. The tips are fast to change thanks to the easy-to-install push-button converter. The interlock system eliminates wear and tear at the insertion point of the air/water syringe.



**CONTACT YOUR KERR REP
FOR YOUR FREE ADAPTOR**

Pack of 200

PE77200 - Seal-Tight White

PE72200 - Seal-Tight Spectrum

£ 36

OptiView™

Lip and Cheek Retractor

The innovative OptiView™ provides retraction of lips and cheeks for improved buccal and gingival access at the same time. It facilitates both lip and cheek retraction with one device without the need for additional hands or instruments.



£ 35

5500 - OptiView Standard Kit



Buy 1 get **Cleanic Mint FREE***

X-rays holders

Film, Phosphor Plate and Sensor Holders

Kwik-Bite™

X-Ray Holders for Horizontal Bite Wings

Endo-Bite™

X-Ray Holders for Diagnostic X-Rays

Super-Bite™

X-Ray Holders for Periapical Radiographs

Paro-Bite™

X-Ray for Perio Patients and Children

Kwik-Bite™ Senso


Universal Sensor Holder for Horizontal Bite Wings

Endo-Bite™ Senso

Universal Sensor Holder for Endodontic Radiographs

Super-Bite™ Senso

Universal Sensor Holder for Periapical Radiographs



*Buy 3 get 1 FREE**
on **ALL X-Ray Holders**

**CONTACT YOUR
KERR REP FOR MODELS**

from **£ 36**

PENTRON

FibreKleer™ 4x
Fiber Post

Everything you want to see
and nothing you don't



£ 143

Tapered Fiber Post Kit

N83B
FibreKleer 4x - Tapered Fiber Post Kit

Buy 1 get
10 Posts FREE*

Parallel Fiber Post Kit

N83A
FibreKleer 4x - Parallel Fiber Post Kit

Buy 1 get
10 Posts FREE*

Original Fiber Post Kit

N83C
FibreKleer 4x - Original Fiber Post Kit

Buy 1 get
10 Posts FREE*



The Power of Choice



Request a demo to discover more about
managing each step of your digital workflow

DEXIS.com



EVERYDAY DENTAL TRACK



9 - 11 MAY 2024

CCIB - INTERNATIONAL CONVENTION CENTER

Barcelona

Limited places left!

SCAN QR CODE OR CONTACT YOUR KERR SALES REP
TO GET MORE INFORMATION



UPCOMING **FREE** | **LIVE** WEBINARS

WEBINAR | ENDODONTICS



April 19th | 13.00 – 14.00

Minimally invasive endodontic treatment
- biomechanical approach to root canal preparation
Dr. Maciej Mikołajczyk DDS, PhD

WEBINAR | ENDODONTICS



June 7th | 13.00 – 14.00

Digital tools to manage complex endodontic cases
Dr. Mohamed Jamal

**FOR MORE
INFO ABOUT
OUR HANDS-ON
COURSES AND
WEBINARS**



[kerrdental.com/en-uk/
dental-courses-and-events](https://kerrdental.com/en-uk/dental-courses-and-events)

Can't watch live?

Don't worry once registered you will have access
to the recording for one week, and still receive receive your CPD!

Kerr UK Sales Team



Julie Todd
DSO Manager UK/Ireland
Mobile: 07933 418738
julie.todd@envistaco.com



Chris McAnena
Mobile: 07974 620612
christopher.mcanena@envistaco.com



Tommy Maguire
Mobile: 07711 750624
tommy.maguire@envistaco.com



Heidi Campbell
Mobile: 07711 750617
heidi.campbell@envistaco.com



Sonia Khela
Mobile: 07929 187102
sonia.khela@envistaco.com



Joanne Shannon
Mobile: 07711 750622
joanne.shannon@envistaco.com




Matthew Green
Mobile: 07855 161659
matthew.green@envistaco.com



Kuljeet Dhak
Mobile: 07711 750619
kuljeet.dhak@envistaco.com



Sarah Green
Mobile: 07855 178315
sarah.green@envistaco.com



Ian Brill
Regional Leader UK/Ireland
Mobile: 07711 750623
ian.brill@envistaco.com



Paul Byfleet
Marketing Manager UK/Ireland
Mobile: 07711 750618
paul.byleet@envistaco.com



Skye Hammond
Sales and Marketing Executive UK/Ireland
Mobile: 1733 892292
skye.hammond@envistaco.com

London

Office: 01733 892292

* Terms & Conditions. Promotion valid 1st April through 30th June 2024 via participating dealers in UK and Ireland only. Prices displayed are Kerr recommend promotional price per item and may vary with your participating dealer. All free stock will be shipped by participating dealers. Promotions subject to change without notice. E & OE.

** Not to be used in combination with any other offers or discount schemes. Free equipment shipped from dealer.

Kerr UK Ltd

c/o Orega Stockley Park | 4 Longwalk road | Stockley Park | Uxbridge | UB11 1FE | UK

Phone: 01733 892292 - www.kerrdental.com

MKT-23-1647 Rev0

Kerr™