

Kerr News UK

WINTER PROMOTIONS

1st January - 31st March 2024

NEW

elements™ Connect & Apex Connect

Cordless Endodontic Motor and Electronic Apex Locator

Seamless Connection. Greater Confidence.

Designed to work seamlessly together, interconnected system supports patient safety while simplifying your procedures.



Follow Us



@kerrdentaluk

Kerr™


Shop Now on **store.kerrdental.com**

- Choose from award winning product lines
- View exclusive online promotions
- All Kerr promotions are automatically applied, the free stock is shipped to your practice without the need to submit a copy invoice
- Receive regular rewards



**Register today and receive
a fabulous welcome gift**

Don't miss
the **Exclusive
January Sale!**



Demi™ Plus

Light Emitting Diode (LED) Curing Light

Thanks to the PLS technology, Demi Plus intermittently shifts the output intensity from a base of 1,100 mW/cm² to a peak of 1,330 mW/cm². This enables the Demi Plus to deliver a 5-second cure for A3 or lighter shades.

Enhanced Body Design

This curing light's exterior is made from Valox resin, a high-performance plastic that will protect and extend the life of your unit. In addition, the upgraded closed housing allows for easier cleaning and disinfection. Improved Cooling System and fan free for quiet operation.

True 5-second cure

Consistent output of Periodic Level Shifting (PLS) delivers a true 5-second cure for composite shades A3 and lighter.

Never shuts down

Whisper-quiet fan and unique thermal control system keeps the unit cool for uninterrupted operation. Extends the overall life of the unit.

Advanced battery technology

The lithium-ion battery technology can yield 300 10-second cures with a single charge.

Lightweight, ergonomic design

Multiple holding options reduce hand and arm stress.



*Save 40%**

910860-2 - Demi Plus

Only £ **895**

Adapt SuperCap™ Matrices

Universal Matrix Tensioning System



£ **50**



*Buy 4 SuperCap refills
get SuperMat Handle FREE**

Adapt SuperCap Pack of 50

2161 - Steel - 0.030mm/5.0mm
2162 - Steel - 0.030mm/ 6.3mm
2181 - Steel - 0.038mm/5.0mm

2182 - Steel - 0.038mm/6.3mm
2171 - Trans. - 0.050mm/5.0mm
2172 - Trans. - 0.063mm/5.0mm

OptiBond™ FL

Filled, Light-Cure Total-Etch Dental Adhesive



26684 E
OptiBond FL Kit

£ 153

*Buy 1 get Gel Etchant
FREE**



OptiBond™ eXTRa Universal

Dual-Component Universal Dental Adhesive



36658
OptiBond eXTRa Universal
Bottle Kit

£ 149

*Buy 1 get NX3 Clear
Refill FREE**



OptiBond™ Universal

Single-Component Universal Dental Adhesive



*Buy 2
get 1
*FREE**

36519
OptiBond Universal Refill 5ml

£ 89

OptiBond™ Solo Plus

Single Component Total-Etch Dental Adhesive



29692
OptiBond Solo Plus
Refill 5ml

£ 74

*Buy 1 get Gel Etchant
FREE**



Vertise™ Flow

Self-Adhering Flowable Composite

Vertise Flow is the first self-adhering flowable dental composite powered by OptiBond™. It greatly simplifies your direct restorative procedures by incorporating a bonding agent into a flowable.

Universal use

Unique self-adhering property offers new opportunities for use in a wide selection of indications and situations.

Reduced post-op sensitivity

Easy handling

Non-slumping property creates an ideal viscosity for different indications.

Low microleakage

Excellent marginal integrity and reduced voids.

Self-adhering.

No separate bonding step.

High bond strengths to dentin and enamel.

Incorporates clinically proven OptiBond adhesive technology for excellent bonds.



Syringe 2 x 2g

- 34401 - A1
- 34402 - A2
- 34403 - A3
- 34404 - A3.5
- 34405 - B1
- 34406 - B2
- 34407 - Universal Opaque
- 34408 - XL
- 34409 - Translucent

Buy 3
get 1
FREE*

£ **73**

*Save 20%**

34398
Vertise Flow Test-Me Kit
1x 2g syringe A2 shade

£ **31**

Herculite™ XRV Ultra

Universal Nanohybrid Dental Composite



Save
35%*

33860
Herculite XRV Ultra Mini Kit Syringe

3x 4g syringes in shades enamel A2, A3, dentine A2, 1x 5ml OptiBond™ Solo Plus bottle, 1x 3g gel etchant, 25x mixing wells, 50x applicator tips.

£ **61**



Save
35%*

34071
Herculite XRV Ultra Mini Kit Unidose

30x 0.2g unidoses in shades enamel A2, A3, dentine A2, 1x 5ml OptiBond™ Solo Plus bottle, 1x 3g gel etchant, 25x mixing wells, 50x applicator tips.


£ **57**

AMAZING VALUE!


Both kits contain **Gel Etchant**
& **5ml OptiBond Solo Plus!**

OptiDisc™

High Performing Polishing and Contouring Disc




AVAILABLE IN THREE
DIFFERENT DIAMETER SIZES



OptiDisc™ finishing and polishing discs allow you to obtain a final polish equal to that of natural dentition on all kinds of restorations. Made of durable polyester impregnated with aluminium oxide particles, the four-grit system achieves a long-lasting result. A colour-coded system enables immediate identification of the grit type. The flush-mounted mandrel ensures correct positioning of the disc and protects against scarring of the restoration. OptiDisc also features the choice of a short-size mandrel to facilitate access to posterior and other restricted areas.

Ready to use and long-lasting performance

OptiDisc is made of polyester coated with aluminium oxide particles. Its abrasive layer has no extra adhesive coating, meaning the moment you fix the disc to the mandrel you are ready to go!



Precise



Flexible



Transparent

Flexible and precise

The shorter mandrel and the disc's bendability allow you to reach posterior and other restricted areas with ease. The transparent material gives you great visibility while working, for extra precision in your finishing and polishing.

Safe


OptiDisc's active side is easy to recognise and the patented mandrel locks below the surface of the disc, protecting your restoration against scarring and scratching. It also guarantees optimal torque transmission, so the disc does not slide and is not rpm sensitive, giving you better and safer handling.



Active Side



Inactive Side



OptiDisc



Competitor Disc

OptiDisc Refills - 100 pcs (except Extra-Coarse, 80 pcs)

- 4197 - 80x 9.6mm - Extra-Coarse
- 4181 - 100x 9.6mm - Coarse/Medium
- 4182 - 100x 9.6mm - Fine
- 4183 - 100x 9.6mm - Extra-Fine
- 4198 - 80x 12.6mm - Extra-Coarse
- 4184 - 100x 12.6mm - Coarse/Medium

- 4185 - 100x 12.6mm - Fine
- 4186 - 100x 12.6mm - Extra-Fine
- 4199 - 80x 15.9mm - Extra-Coarse
- 4191 - 100x 15.9mm - Coarse/Medium
- 4192 - 100x 15.9mm - Fine
- 4193 - 100x 15.9mm - Extra-Fine

£ 26

Buy 2
get 1
FREE*

JET Diamonds™

Dental Burs

Streamline your workflow with only 23 burs

Improve your workflow and chair time with a trusted solution. Jet's diamond burs have outstanding product quality and a proven high satisfaction rate from clinicians all around the world.



Multi-Grit Technology



Outstanding cutting efficiency.



Unique diamond plating process ensures excellent diamond particle retention.



Shanks are designed to enable concentric running with **extremely low vibration** and **less wear and tear** on handpiece turbine.

**CONTACT
YOUR KERR REP
FOR BUR TYPES**

Special **PRICE**

£ 11

Pack of 5

TempBond™ NE

Non-Eugenol Temporary Cement



61087 - TempBond Standard Package

60256 - TempBond NE Standard Package

£ 29

*Buy Any 2 TempBond
get SimpliShade
Refill FREE**



1st January - 31st March 2024 | 7

SimpliShade™

Simplified Universal Composite System

Simplify Shade-matching

Shade selection made easier and faster.
SimpliShade matches all 16 VITA classical shades with just 3 shades: Light, Medium, Dark.

Adaptive Response Technology (ART)

The ART filler network, a proven technology found in Harmonize, helps you achieve lifelike restorations with exceptional strength and unmatched esthetics.

Streamline Inventory Management

Save time and money with minimal shades to meet almost all your composite needs. Declutter your shelves, avoid the burden of re-ordering too many different part numbers, and reduce or eliminate expired products.

Universal For All Cavity Classes

Use SimpliShade with confidence in all cavity classes without the need for a blocker, unlike other competitors.



SimpliShade™

Competitor 1

Competitor 2



Excellent Radiopacity

The more radiopaque, the easier it is to distinguish the restoration from the tooth structure under a radiograph. SimpliShade at 300% AI is almost twice as radiopaque as competition.

Optimal Translucency

SimpliShade provides ideal opacity to mimic tooth structure, to better mask underlying tooth stains, and blend well with surrounding dentition, resulting in a more lifelike restoration.

SimpliShade™

Competitor 1

Competitor 2



Unidoses 20x 0,25g

37006 - Light

37007 - Medium

£ **59**

Unidoses 10x 0,25g

37008 - Dark

37009 - Bleach White

37010 - Universal Opaque

£ **39**

Syringe 4g

37001 - Light

37002 - Medium

37003 - Dark

37004 - Bleach White

37005 - Universal Opaque

£ **49**

*Buy 2
get 1
FREE**

Maxcem™ Elite Universal

Self-Etch, Self-Adhesive Resin Dental Cement

Buy 2
get 1
FREE*



2x Automix Syringe 5g

34058 - Clear

34060 - White Opaque

34059 - White

34061 - Yellow

37003 - Dark

£ 100

Maxcem™ Elite Chroma Universal

Self-Etch, Self-Adhesive Resin Dental Cement with Cleanup Indicator

Buy 2
get 1
FREE*



2x Automix Syringe 5g

36372 - Clear

36373 - White

36374 - Yellow

£ 128

NX3 Nexus™ Third Generation

Universal Adhesive Resin Dental Cement

Buy 2
get 1
FREE*

Automix Syringe 1.8g

33646 - Bleach

33643 - Clear

33647 - Opaque

33644 - White

33645 - Yellow

£ 85

Light-Cure Syringe 5g

33651 - Bleach

33648 - Clear

33652 - Opaque

33649 - White

33650 - Yellow

£ 42



elements™ Connect & Apex Connect

Cordless Endodontic Motor and Electronic Apex Locator

NEW

An Interconnected System.

Versatile Motions

With Adaptive, Reciprocation and Rotary motions in one motor, elements Connect allows offices to maintain their existing techniques and protocols.

Preset

9 Kerr preset file settings, 1 reciprocating setting and 4 customizable settings.

Continuous Use

Two removable batteries are included to ensure one is always charged so you can complete your procedures safely without interruption.



Built-in Accuracy

Ensure measurement accuracy with Apex Connect's automatic self-calibration feature that activates every time the device is turned on.



Shape with Confidence

View real-time working length information during procedures with a smart light indicator that changes color, as well as audible indicators, as you reach the apex.



Launch Promo

815-1902-UK
Apex Connect

£ **625**

Save £ 265

815-1901-UK
elements Connect

£ **1,350**

Save £ 500

Endo Ring Kit

Endodontic Instruments Holder



Special PRICE

Save
20%*

£ **100**

825-0049 - Endo Ring Kit

K3™ Files

Nickel Titanium Files

**CONTACT
YOUR KERR REP
FOR SIZES**

£ **10**



*Buy 12
get 5
FREE**

K-Flex™

Stainless Steel Files

**CONTACT
YOUR KERR REP
FOR SIZES**

£ **10**



*Buy 12
get 5
FREE**

Tubli-Seal™ EWT

Root Canal Sealer



*Buy 3
get 1
FREE**

£ **53**


25903 - TubliSeal EWT

K3™ Files

Nickel Titanium Files

**CONTACT
YOUR KERR REP
FOR SIZES**

£ **47**



*Buy 12
get 5
FREE**

K3XF™ Files

Rotary Nickel Titanium Files

**CONTACT
YOUR KERR REP
FOR SIZES**

£ **48**



*Buy 12
get 5
FREE**

TF Adaptive™

NiTi Twisted Files

**CONTACT
YOUR KERR REP
FOR SIZES**

£ **45**



*Buy 12
get 5
FREE**

Irrigation Probe

Irrigation Needles

Special PRICE

*Save
20%**



£ **32**

Traverse™

Rotary Glide Path Files

Special **PRICE**

Save
20%*

**CONTACT
YOUR KERR REP
FOR SIZES**



£ **46**



ZenFlex™

NiTi Rotary Shaping Files

**CONTACT
YOUR KERR REP
FOR SIZES**



£ **38**

*Buy 3
get 1
FREE**

elements Gutta Percha Cartridges

Innovative Gutta Percha Cartridges

972-2500

23GA Gutta Percha cartridge (medium body)

972-2501

25GA Gutta Percha cartridge (light body)

972-2502

23GA Gutta Percha cartridge (heavy body)



Save
20%*

£ **57**

elements™ IC System

Cordless Obturation System

Buy elements IC *get* 5 boxes of Gutta Percha cartridges, 1 Heat Plugger, 3 Hand Pluggers and 1 TubliSeal EWT FREE*

Worth
over
£650

£ **3,050**

973-0600-TYPEF - elementsIC Obturation System

5x Gutta Percha cartridges

972-2500 - 23GA Gutta Percha cartridge (medium body)

972-2501 - 25GA Gutta Percha cartridge (light body)

972-2502 - 23GA Gutta Percha cartridge (heavy body)

3x Hand Pluggers

974-0060 - Hand plugger (size # 0)

974-0058 - Hand plugger (size # 1)

974-0059 - Hand plugger (size # 2)

1x Heat Plugger

952-0005 - Buchanan Heat Plugger .08 (FM)

25903 - TubliSeal EWT



OptiView™

Lip and Cheek Retractor

5500 - OptiView Standard Kit



BIG SAVING!

*Save 40%**

£ **21**

SealTight™

Disposable Air/Water Syringe Tips

*Buy 2
get 1
FREE**

Pack of 200

PE77200 - Seal-Tight White

PE72200 - Seal-Tight Spectrum

**CONTACT YOUR KERR REP
FOR YOUR FREE ADAPTOR**



£36

CaviWipes™

Wipe Surface Disinfectants

CaviWipes™ are broad-spectrum medical device surface disinfectant, with clinically relevant efficacy claims: bactericidal, virucidal, fungicidal, tuberculocidal and yeasticidal. These ready to use multipurpose intermediate level disinfectants are made for use on hard, non-porous surfaces of non-invasive medical devices within dental practices, including patient care areas and laboratories.

Fast, Easy and Multi-surface Compatible

Kill Most Pathogens in 3 minutes or Less

Faster change over time and helps lower risks of cross contamination.

Ready-To-Use Products

No dilution or mixing required.

Saves Time & Inventory

A cleaner and disinfectant in one.

Low Alcohol Content

Less than 20%, which offers excellent material compatibility.

Fragrance and Colorants-free.

Does not contain phenol, aldehyde, bleach or other toxic chemicals



*Buy 3
get 1
FREE**

4731160
CaviWipes Canisters
Pack of 160

£8.40

4731245
CaviWipes Flatpacks
Pack of 45

£4.70

A professional solution that you can trust at home

Perlite. Works professionally, works at home

Whether you are a movie star or an everyday action hero, a pop singer or a shower performer you deserve a beautiful, healthy and white smile.

Cleanicdent™ with Perlite technology is an effective solution from the market leader* in professional prophylaxis pastes. A gentle, effective stain removal toothpaste that helps to extend the benefits of professional tooth cleaning.

* Source: GfK DDM Jahresbericht 2017, February 2018 – Pastes segment
Key-Stone Fast Track Sell in data 2017



Hard on stains, gentle on teeth.

Cleanic™

Cleaning & Polishing Paste

*Buy 2 get 2
Cleanicdent
FREE**



- 3380** - Mint
- 3382** - Green Apple
- 3383** - Mint w/o Fluoride
- 3384** - Light
- 3386** - Berry Burst

£ 18

X-rays holders

Film, Phosphor Plate and Sensor Holders

Kwik-Bite™

X-Ray Holders for Horizontal Bite Wings

Endo-Bite™

X-Ray Holders for Diagnostic X-Rays

Super-Bite™

X-Ray Holders for Periapical Radiographs

Paro-Bite™

X-Ray for Perio Patients and Children

Kwik-Bite™ Senso

Universal Sensor Holder for Horizontal Bite Wings

Endo-Bite™ Senso

Universal Sensor Holder for Endodontic Radiographs

Super-Bite™ Senso

Universal Sensor Holder for Periapical Radiographs



*Buy 3 get 1 FREE**
on **ALL X-Ray Holders**

**CONTACT YOUR
KERR REP FOR MODELS**

from **£ 36**

Bib-Eze™

Single-use Disposable Bib Holders

Bib-Eze™ Disposable Bib Holders are lightweight stretchable paper holders that attach to the dental napkin with their easy-to-peel, self-adhesive ends. The bib holder can easily be removed after the procedure and thrown away alongside the dirty napkin, thus reducing the risk of cross-contamination.



Special **PRICE**

Save
20%*

35053 - Bib-Eze - Pack of 250

£ 11

PENTRON

FibreKleer™ 4x

Fiber Post



**Everything you want to see
and nothing you don't**

£ 143

Tapered Fiber Post Kit

N83B

FibreKleer 4x - Tapered Fiber Post Kit

**Buy 1 get
10 Posts FREE***

Parallel Fiber Post Kit

N83A

FibreKleer 4x - Parallel Fiber Post Kit

**Buy 1 get
10 Posts FREE***

Original Fiber Post Kit

N83C

FibreKleer 4x - Original Fiber Post Kit

**Buy 1 get
10 Posts FREE***



EVERYDAY DENTAL TRACK



9 - 11 MAY 2024

CCIB - INTERNATIONAL CONVENTION CENTER

Barcelona

SCAN QR CODE OR CONTACT YOUR KERR SALES REP
TO GET MORE INFORMATION



UPCOMING **FREE** | **LIVE** WEBINARS

WEBINAR

ENDODONTICS



February 23rd | 15.00 – 16.00

**How to minimize errors
in endodontic treatment**

MUDr. Roman Blahuta, PhD.

WEBINAR

RESTORATIVES



March 22nd | 10.00 – 11.00

**Posterior composite solutions
for everyday dentistry**

Dr. Ahmed Tadfi BDS, BSc

WEBINAR

ENDODONTICS




April 19th | 14.00 – 15.00 CET

**Minimally invasive endodontic treatment
- biomechanical approach to root canal preparation**

Dr. Maciej Mikolajczyk DDS, PhD

**FOR MORE
INFO ABOUT
OUR HANDS-ON
COURSES AND
WEBINARS**



[kerrendental.com/en-uk/
dental-courses-and-events](https://kerrendental.com/en-uk/dental-courses-and-events)

Can't watch live?

Don't worry once registered you will have access
to the recording for one week, and still receive receive your CPD!

Kerr UK Sales Team



Chris McAnena

Mobile: 07974 620612
christopher.mcanena@envistaco.com



Tommy Maguire

Mobile: 07711 750624
tommy.maguire@envistaco.com



Heidi Campbell

Mobile: 07711 750617
heidi.campbell@envistaco.com



Joanne Shannon

Mobile: 07711 750622
joanne.shannon@envistaco.com



Sonia Khela

Mobile: 07929 187102
sonia.khela@envistaco.com



Matthew Green

Mobile: 07855 16159
matthew.green@envistaco.com




Kuljeet Dhak

Mobile: 07711 750619
kuljeet.dhak@envistaco.com



Sarah Green

Mobile: 07855 178315
sarah.green@envistaco.com



Ian Brill

Regional Leader UK/Ireland

Mobile: 07711 750623
ian.brill@envistaco.com



Paul Byfleet

Marketing Manager UK/Ireland

Mobile: 07711 750618
paul.byleet@envistaco.com



Julie Todd

DSO Manager UK/Ireland

Mobile: 07933 418738
julie.todd@envistaco.com

London

Office: 01733 892292

* Terms & Conditions. Promotion valid 1st January through 31st March via participating dealers in UK and Ireland only. Prices displayed are Kerr recommend promotional price per item and may vary with your participating dealer. All free stock will be shipped by participating dealers. Promotions subject to change without notice. E & OE.

** Not to be used in combination with any other offers or discount schemes. Free equipment shipped from dealer.

Kerr UK Ltd

c/o Orega Stockley Park | 4 Longwalk road | Stockley Park | Uxbridge | UB11 1FE | UK

Phone: 01733 892292 - www.kerrdental.com

MKT-23-1647 Rev0

Kerr™