

DENTSPLY SIRONA SPECIAL DEAL



Buy an **X-Smart® Pro+**

Our latest endodontic motor
with integrated apex locator

£2,095



Buy an X-Smart Pro+ Endo Motor,
Get a Contra-Angle with Sheath (worth £367) FREE!

**OFFER ENDS
30TH JUNE
2025**

Why Choose X-Smart® Pro+?

1. Compared to contra-angles with larger head height
2. The Bluetooth® word mark is a registered trademark owned by Bluetooth SIG, Inc.

✓ Specifically designed to drive Dentsply Sirona's endodontic file systems

The X-Smart® Pro+ endo motor has been specifically designed to drive Dentsply Sirona's endodontic file systems in both continuous rotation and reciprocating modes.

✓ Mini-Contra Angle head with integrated 10-Lumen LED

Great visibility and access. Locate and enter root canal with confidence. Easy treatment of patients with limited mouth opening.¹ 360° adjustable position.

✓ User-friendly touchscreen interface

Large, easy-to-clean screen. Easy navigation with gloved hand. Up to 30 customisable file sequences with upgradable firmware. Optional Bluetooth® foot switch.²

✓ Integrated Apex Locator with Dynamic Accuracy™ Technology

Our Dynamic Accuracy™ technology reads the measurement signal accurately while shaping, with auto-reverse function to keep your file on your desired target.

✓ Removable/autoclavable handpiece metal sheath

Supports compliance with infection control guidelines. Reduces waste, making handpiece barrier sleeves obsolete.

✓ Portable and ergonomic design

Battery power to support a full day of work. Portable convenience with the power and functions of a tabletop unit. Drag and drop magnetic holder for enhanced comfort, reduced fatigue.