

Validity April - June 2025 | United Kingdom

promotion

 **HuFriedyGroup**
The Best In Practice

COME
TO VISIT US @

**ADI
Congress
Brighton**

**1 - 3 MAY
2025**

NEW!




Sterify Gel

**Expanding the
Boundaries of Care
and Protection**



- **DentaPure™**
Independent
Water Bottle
Cartridge..... 19



- **Needle
Holders**..... 10

- **Sutures**
Absorbable
and Non-absorbable.... 8

NEW!

Sterify Gel* HU-FRIEDY

Expanding the Boundaries of Care and Protection

Meet the reduced bacterial contamination solution to support the healing process in the treatment of periodontal disease

Periodontitis is a chronic multifactorial disease that can lead to the destruction of the structures supporting the teeth.

Affects about



70%



Leads to



of the adult population¹


the destruction of the supporting apparatus of the tooth¹

When treating periodontal disease, **Scaling & Root Planing (SRP)** might currently require the addition of local antibiotic or local antibacterial agents¹



Sterify Gel is a medical device that offers a new approach to periodontal care.

Mucoadhesive polymeric hydrogel in a pre-filled syringe for the treatment of periodontal disease²



Sterify Gel

THE MAIN BENEFITS:



Reduction of bacterial contamination thanks to its occluding action



Support to tissue healing with significant and proven constant improvement of pocket depth



No antibiotics, no antiseptics and no chlorhexidine but simultaneously acts as broad-spectrum antibacterial barrier, without generating resistance



No animal-derived components




Stable at room temperature with no need to manage the cold chain



Easy to use


Sterify Gel provides a unique solution in the treatment of periodontal disease by preventing bacterial contamination and supporting tissue healing.



Sterify Gel easily penetrates into the deepest and hard-to-reach areas of periodontal and peri-implant pockets, adhering to the gingival tissue, alveolar bone or root and dental implant surfaces, **providing complete coverage of the pockets**²

The occluding action at the level of the periodontal and peri-implant gingival pockets **creates an environment unsuitable for bacterial growth and supporting the tissue healing process**²

Changes in pocket depth, gingival recession, and clinical attachment level after treatment with Sterify Gel compared to srp alone¹



Sterify Gel in conjunction with SRP consistently demonstrated a **noteworthy improvement in PD, surpassing twofold enhancement compared to the control group** throughout all follow-up intervals¹



Start using Sterify Gel today and see the difference!
HuFriedyGroup.eu/Sterify-Gel

SPECIAL PROMO!

4 x PCPUNC156 Perio Probes
+ 1 x Sterify Gel

£ 99

Sterify Gel 0,3 ml,
1 syringe
STY-GEL0301

£ 49





NEW

Harmony™ Scalars and Curettes HU-FRIEDY

EVEREDGE™ 2.0

WORKING ENDS

COLOR CODED

SECURE AND NIMBLE GRASP

The silicone grip has been extended by 30% which can provide a secure and nimble grasp while using the instrument.³ The geometric pattern is specifically designed to maximize grip in all directions.



REDUCED PINCH FORCE

REDUCES PINCH FORCE UP TO 65%⁴

This may alleviate hand fatigue and injury risk due to repetitive motions.^{5,6,7,8}

LESS PRESSURE

REDUCES PRESSURE ON TOOTH BY 37%⁴

Harmony™ Ergonomic Scalars and Curettes may increase patient comfort, while reducing clinician fatigue because it requires less pressure to the tooth when scaling.

BALANCED

NOT TOO HEAVY OR TOO LIGHT


The Harmony™ Ergonomic Scaler and Curette's perfectly balanced handle is similar in weight to our best in class EverEdge™ 2.0 #9 Metal Handle.



SPECIAL PROMO!

4 x Harmony Scalars
+ 2 x Sterify Gel

£ 182



Learn How to Scale in Perfect Harmony at
[HuFriedyGroup.eu/Harmony-Scaler](https://www.HuFriedyGroup.eu/Harmony-Scaler)

3) Compared to the Hu-Friedy #9 metal handle, these are the nominal values. Data on file. Available on request. 4) Data on file. Available upon request. 5) Int J Dent Hygiene 7, 2009; 159-165 DOI: 10.1111/j.1601-5037.2009.00395.x, Hayes MJ, Cockrell D, Smith DR. A systematic review of musculoskeletal disorders among dental professionals. 6) Rempel, David, et al. "The Effects of Periodontal Curette Handle Weight and Diameter on Arm Pain." The Journal of the American Dental Association, vol. 143, no. 10, 2012, pp. 1105-1113, doi:10.14219/jada.archive.2012.0041. 7) Lalumandier, James A, and Scott D McPhee. "Prevalence and Risk Factors of Hand Problems and Carpal Tunnel Syndrome among Dental Hygienists." Journal of Dental Hygiene, vol. 75, no. II, 2001, pp. 130-134. 8) Mulmani P, Hoe VCW, Hayes MJ, Idiculla JJ, Abbas ABL, Karanth L. Ergonomic interventions for preventing musculoskeletal disorders in dental care practitioners. Cochrane Database of Systematic Reviews 2018, Issue 10. Art. No.: CD011261. DOI: 10.1002/14651858.CD011261.pub2. 9) See www. HuFriedyGroup.eu/Harmony-Scaler for terms & conditions



EverEdge™ 2.0 Instruments HU-FRIEDY

NEW



SG1/2XCE2



SG7/8XCE2



SG11/12XCE2



SG13/14XCE2



SAS1/2XCE2



SAS7/8XCE2



SAS11/12XCE2



SAS13/14XCE2



SCNEVI1XE2



SCNEVI2XE2




SCNEVI3XE2



SCNEVI4XE2

COMPARISON CHART



STANDARD GRACEY CURETTE

Description	Handle Harmony™
1/2 Standard Gracey	● SG1/2XCE2
7/8 Standard Gracey	● SG7/8XCE2
11/12 Standard Gracey	● SG11/12XCE2
13/14 Standard Gracey	● SG13/14XCE2

MINI FIVE™ GRACEY CURETTES

Description	Handle Harmony™
1/2 Mini Five™ Gracey	● SAS1/2XCE2
7/8 Mini Five™ Gracey	● SAS7/8XCE2
11/12 Mini Five™ Gracey	● SAS11/12XCE2
13/14 Mini Five™ Gracey	● SAS13/14XCE2

UNIVERSAL CURETTES

Description	Handle Harmony™
1/2 Langer	● SL1/2XE2
3/4 Langer	● SL3/4XE2
5/6 Langer	● SL5/6XE2
17/18 Langer	● SL17/18XE2

SICKLE SCALERS

Description	Handle Harmony™
H6/H7 Scaler	● SH6/7XE2
204S Scaler	● S204SXE2
204SD Scaler	● S204SDXE2
1 Nevi Anterior	● SCNEVI1XE2
2 Nevi Posterior	● SCNEVI2XE2
3 Nevi Scaler Posterior	● SCNEVI3XE2
4 Nevi Scaler Posterior	● SCNEVI4XE2

Curettes and Scalers
EverEdge™ 2.0 with
Harmony™ Handle

£ 48
£ 68

Save even
more when
bundled with
Sterify Gel



LC Implant Suite OMNIA

This multifunctional surgical trolley is developed to improve the organization of the dental office

- Properly stores electronic equipment
- Avoids cable and wiring clutter
- Can improve office hygiene
- Can save time and money in workflow organization

Essential for professionals using physiodispensers or piezosurgical equipment.


Multifunctional Surgical Trolley

50 x 54 x 94 cm, White

30.E0067.00

£ 1,250

£ 1,520



Oleodynamic Instrument Trolleys OMNIA



Discover
more

Equipped with single tray with raised edges.

Manually adjustable height from 86 to 120 cm

(tray: 73 x 47 x 2,5 cm) - Weight 31,5 kg

30.E0022.00

£ 1,200

£ 1,505



Surgical Over Patient Table OMNIA

Manually adjustable height from 98 to 150 cm

(tray: 60 x 90 cm)

Weight 15 kg

30.E0047.00

£ 950

£ 1,082






Irrigation Sets OMNIA

Physiodispenser and Piezoelectric Surgery

- Proper irrigation during the creation of the implant site is a fundamental element for successful bone integration of the implant and the overall outcome of the operation.
- Wide range of mechanical irrigation lines which are compatible with most implant motor models. Irrigation systems include a vented spike, a roller or clamp for flow control, fittings, and an integrated peristaltic pump section.
- Convenient terminal cap makes the connection of the irrigation line to the hand piece easy to achieve. The cap allows you to connect tubes to the wide variety of cooling needle models available.



For more information on all Irrigation Sets available, [download the dedicated brochure](#)



[Watch](#) our how-to video on YouTube



Code	Compatible with:		QTY	Promo	Sugg. Retail
32.F0134.00	NSK	Surgic XT/Plus, Vario Surg, SurgicPro	10 pcs	£ 47	£ 69
32.F0139.00	W&H	Implantmed, Piezomed	10 pcs	£ 47	£ 69
32.F0170.00	Bien Air	Chiropro 980	10 pcs	£ 47	£ 69
32.F0171.00	EMS	Piezon Master	10 pcs	£ 61	£ 89
32.F0194.00	Acteon	Implant Centre, Piezotome	10 pcs	£ 67	£ 99
32.F0113.00	Kavo / Mectron	Kavo Intrasept 905 - Mectron Piezosurgery II, Piezosurgery 3	10 pcs	£ 26	£ 38



Sterile Saline Solution OMNIA

The sterile solution is a product designed to cool rotating surgical instruments during implant surgery or oscillating instruments during piezoelectric surgery.

Sterile Saline Solution (sodium chloride 0.9 %) with connector for vented spike

Code	Description	Promo	Sugg. Retail
32.E0000.00	1000 ml bag, 10/Box (non medicinal)	£ 25	£ 36
32.E0001.00	500 ml bag, 15/Box (non medicinal)	£ 26	£ 39
32.E0002.00	250 ml bag, 20/Box (non medicinal)	£ 29	£ 43

Suction Tips OMNIA


Code	Description	Promo	Sugg. Retail
32.F6565.00	Omniasurg model sterile (each pack contains 1 adapter Adat. 08 - tip length: 20 cm) in PVC, 20/Box	£ 16	£ 24

Code	Description	Promo	Sugg. Retail
32.F6363.00	Yankauer sterile (tip length: 24,5 cm) in PVC, 10/Box	£ 13	£ 19

Code	Description	Promo	Sugg. Retail
32.F7070.00	Omniasurg model sterile with disposable bone trap (each pack contains 1 adapter Adat. 04 - tip length: 20 cm) in PVC, 5/Box	£ 23	£ 34




Absorbable Sutures




PGA HU-FRIEDY
Perma Sharp™, Violet thread
12 PER BOX

Code	Description	Promo	Sugg. Retail
PSN393V	3/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 40	£ 64
PSN389V	4/0, 45cm, needle 16 mm, 3/8 circle ReverseCutting	£ 46	£ 64
PSN388V	5/0, 45cm, needle 16 mm, 3/8 circle ReverseCutting	£ 46	£ 64
PSN387V	6/0, 45cm, needle 16 mm, 3/8 circle ReverseCutting	£ 46	£ 64




PGA FA HU-FRIEDY
Perma Sharp™ Fast Absorption, Undyed White thread
12 PER BOX

Code	Description	Promo	Sugg. Retail
PSN540FA	3/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 67	£ 126
PSN535FA	4/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 61	£ 126
PSN530FA	5/0, 70cm, needle 19 mm, 3/8 circle ExtraReverseCutting	£ 90	£ 185



PGLA RAPID OMNIA
Undyed White thread
12 PER BOX

Code	Description	Promo	Sugg. Retail
32.Z3450.00	4/0, 75cm, needle 16 mm, 3/8 circle ReverseCutting	£ 58	£ 86
32.Z3451.00	5/0, 75cm, needle 16 mm, 3/8 circle ReverseCutting	£ 58	£ 86
32.Z3453.00	6/0, 45cm, needle 12 mm, 3/8 circle ExtraReverseCutting	£ 70	£ 103



PGCL OMNIA
Violet thread
12 PER BOX

Code	Description	Promo	Sugg. Retail
32.Z3300.00	4/0, 75cm, needle 16 mm, 3/8 circle ReverseCutting	£ 58	£ 86
32.Z3302.00	5/0, 75cm, needle 16 mm, 3/8 circle ReverseCutting	£ 58	£ 86
32.Z3301.00	5/0, 45cm, needle 12 mm, 3/8 circle ExtraReverseCutting	£ 70	£ 103
32.Z3304.00	6/0, 45cm, needle 12 mm, 3/8 circle ExtraReverseCutting	£ 70	£ 103

Discover more
on the following
links:




Sutures
Catalogue



Sutures
Comparison




Non-absorbable Sutures




POLYPROPYLENE HU-FRIEDY Perma Sharp™, Blue thread 12 PER BOX

Code	Description	Promo	Sugg. Retail
PSN8683P	4/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 50	£ 71
PSN8661P	5/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 50	£ 71
PSN8614P	5/0, 45cm, needle 16 mm, 3/8 circle ReverseCutting	£ 41	£ 71
PSN8660P	6/0, 45cm, needle 16 mm, 3/8 circle ReverseCutting	£ 50	£ 71



SILK HU-FRIEDY Perma Sharp™, Black thread 12 PER BOX

Code	Description	Promo	Sugg. Retail
PSN7772S	3/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 27	£ 47
PSN683S	4/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 29	£ 47
PSN682S	5/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 34	£ 47




PTFE HU-FRIEDY PERMA SHARP™ 12 PER BOX

Code	Description	Promo	Sugg. Retail
PSNPTFE1	3/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 119	£ 135
PSNPTFE2	3/0, 45cm, needle 16 mm, 3/8 circle ReverseCutting	£ 119	£ 135
PSNPTFE3	4/0, 45cm, needle 19 mm, 3/8 circle ReverseCutting	£ 119	£ 134
PSNPTFE4	4/0, 45cm, needle 13 mm, 3/8 circle ExtraReverseCutting	£ 135	£ 152
PSNPTFE5	4/0, 45cm, needle 12 mm, 3/8 circle ReverseCutting	£ 119	£ 135
PSNPTFE6	4/0, 45cm, needle 16 mm, 3/8 circle ReverseCutting	£ 119	£ 135
PSNPTFE7	5/0, 45cm, needle 11 mm, 3/8 circle ExtraReverseCutting	£ 135	£ 153
PSNPTFE8	5/0, 45cm, needle 13 mm, 3/8 circle ExtraReverseCutting	£ 135	£ 152
PSNPTFE9	5/0, 45cm, needle 12 mm, 3/8 circle ReverseCutting	£ 119	£ 135
PSNPTFE10	6/0, 45cm, needle 11 mm, 3/8 circle ExtraReverseCutting	£ 135	£ 153
PSNPTFE11	6/0, 45cm, needle 13 mm, 3/8 circle ExtraReverseCutting	£ 135	£ 153
PSNPTFE12	6/0, 45cm, needle 12 mm, 3/8 circle ReverseCutting	£ 119	£ 135



The beauty of HuFriedyGroup Sutures is that you can count on the high quality of the suture material, the precise engineering of the needles and of course good value for money.

Amit Patel,
BDS MSc MClint Dent FDS RCSEd MRD
RCSEng Specialist in Periodontics





Needle Holders

HU-FRIEDY

Castroviejo, Perma Sharp™
Round Handle, 18 cm

5-0, 6-0, 7-0, 8-0, 9-0, 10-0

NH5024R

£ 299
£417



Crile-Wood, Perma Sharp™, 15 cm
8-0, 9-0, 10-0

NH5038


£ 99
£138




Tissue Pliers

HU-FRIEDY

PLAIN



1X2



Adson Plain
Perma Sharp™,
12 cm
TP5041

£ 96
£132



Adson 1x2
Perma Sharp™,
12 cm
TP5042

£ 118
£167



General Plain
Perma Sharp™,
15,5 cm
TP5070

£ 96
£128



General 1x2
Perma Sharp™,
15,5 cm
TP5071

£ 110
£156



Gerald Plain
Perma Sharp™,
18 cm
TP5050

£ 110
£150



Gerald 1x2
Perma Sharp™,
18 cm
TP5051


£ 124
£176



Scissors HU-FRIEDY

**#17 Iris, Straight, Super-Cut,
11,5 cm, 1 serrated blade**
S17SC

£ 58
£ 81



**#18 Iris, Curved, Super-Cut
11,5 cm, 1 serrated blade**
S18SC

£ 58
£ 81



**Cortellini, Super-Cut, Curved
12 cm, 1 serrated blade**
SCORTSC

£ 72
£ 101



**#14 LaGrange, Double Angled,
Super-Cut, 11,5 cm ,
1 serrated blade**
S14SC

£ 74
£ 104



**Iris, Perma Sharp™, Straight
11,5 cm**
S5082

£ 119
£ 168



**Iris, Perma Sharp™, Curved
11,5 cm**
S5083

£ 131
£ 185





Scalpel Blades & Handles HU-FRIEDY

SURGICAL

Round Scalpel Handle #5E,
European Style, Matte Finish
10-130-05E

£ 54
£73



Round Angled Scalpel Handle
#5E, European Style, Matte Finish
10-130-05AE

£ 52
£70



Scalpel
Blades #12,
Box of 100,
Expiry Date
10-255-12

£ 30
£40



Scalpel
Blades #15,
Box of 100,
Expiry Date
10-255-15

£ 30
£40



Sterile
Blade #15C,
Box of 100,
Expiry Date
10-256-15

£ 36
£48

MICROSURGICAL

Scalpel Blade/Mirror Handle,
SinusLine, Stainless Steel
SHDPV

£ 131
£176



Scalpel Blade/Mirror Handle,
SinusLine, Titanium
SHTMBH

£ 140
£189



Single Edge
Sterile
Scalpel
Blade,
box of 12
MB64

£ 57
£77



Mini Curved
Edge Scalpel
Blade,
box of 12
MB67

£ 57
£77



Full Radius
Round
Scalpel
Blade,
box of 12
MB69

£ 71
£77



Single Edge
Scalpel
Blade,
box of 12
MIM64


£ 49
£45



Urban Sinus Lift Kit HU-FRIEDY

Urban
Sinus
Lift Kit
URBANSLDIN

£ 530
£ 670



We can customise your kit by lasermarking your name, practice logo or kit designation.
Just ask when ordering!



Istvan Urban
Dmd, Md, PhD



Focus on:
Guided Bone Regeneration



Country:
Hungary



Favourite Instrument:
Mini Me PFIWDS1MK



Website:
implant.hu/en



Inspired by:
Henry H. Takei

KIT COMPOSITION:

Code	Description
IME12DIN8_	IMS™ Infinity Series™ Cassette, 10 Instruments, small
IMPURBA1S6	Sinuslift Instrument Urban #1 hdl #6
IMPURBA2S6	Sinuslift Instrument Urban #2 hdl #6
IMPURBAN36	Sinuslift Instrument Urban #3 hdl #6
PPR3S6	Periosteal Prichard #3 short hdl #6
PLGJO/46	Plugger Jovanovic #4 short hdl #6
CRURBAN	Cheek Retractor Urban


COLOUR CODE OF CASSETTES AND TUBS

Replace the _ column with the number of your desired colour in the code.






Luxating Elevator Kit HU-FRIEDY



Luxation Kit Handle #510,
includes: EL3S, EL5S, EL3C, EL5C, SSLUX
ELKIT

£ 250
£-350

Elevators HU-FRIEDY



Cryer #A30,
Handle #510,
right angle
EA30
left angle
EA31

£ 66 each
£-93




Coupland
Gouge,
Handle #510:
small
EA1C
medium
EA2C
large
EA3C

£ 66 each
£-93



Warwick-
James,
Handle #504:
straight
EWJS
left
EWJL
right
EWJR

£ 63 each
£-89



ORGANIZE YOUR ELEVATORS WITH IMS™ CASSETTES

IMS™ Infinity Series™ Cassette 1/2 DIN,
for 4 elevators, available in all colors.
IME12D8_E4

£ 125

COLOUR CODE OF CASSETTES AND TUBS

Replace the _ column with the number of your desired colour in the code.





NEW

Bio-Emulation™ Group HU-FRIEDY

The_alpha kit

Clinical tasks requiring high levels of manual and cognitive skill require precise fine motor control, so instruments designed to feel like a natural extension of the hand are an essential addition to your armamentarium.


The Bio-Emulation™ alphabet kit has been designed with a focus on proprioception and fine dexterity, to deliver outstanding direct composite restorations. **Designed by Dentists, for Dentists.**





1
Bio-Emulation™
Composite
Instrument,
Alpha
CIBIOA

Bio-
Emulation™
Single
Instrument


£ 48




2
Bio-Emulation™
Composite
Instrument,
Beta
CIBIOB




3
Bio-Emulation™
Composite
Instrument,
Gamma
CIBIOY



4
Bio-Emulation™
Composite
Instrument,
Delta
CIBIOD



5
Bio-Emulation™
Composite
Instrument,
Epsilon
CIBIOE



6
Bio-Emulation™
Composite
Instrument,
Zeta
CIBIOZ




The Bio-Emulation™ Philosophy can be seen as the natural next step to the foundations of biomimetic dentistry that our mentors established in the past. There is an apt quote from Albert Einstein that represents the Bio-Emulation™ philosophy very well: “Look deep into nature, and then you will understand everything better.”

Bio-Emulation™ Group



Discover
more



Bio Emulation Kit
BIODIN

£ 282

£382




Composcul Kit HU-FRIEDY

CompoSculp instruments were designed by Dr. Didier Dietschi with the intention of aiding clinicians in achieving superior esthetic restorations. The CompoSculp kit is comprised of 5 complimentary composite instruments for both anterior and posterior restorations.

CompoSculp Instrument Kit, Smooth Handles **PFIDDDIN**

£ 260
£ 322



Discover
more

Code	Description
IME14DIN54	IMS™ Infinity Series™ Cassette 1/4 DIN 5 instruments, orange
PFIDD1/2	Dietschi Composite #1/2, Smooth Handle
PFIDD3/4	Dietschi Composite #3/4, Smooth Handle
PFIDD5/6	Dietschi Composite #5/6, Smooth Handle
PFIDD7/8	Dietschi Composite #7/8, Smooth Handle
PFIDD9/10	Dietschi Composite #9/10, Smooth Handle



Success usually
results from an
ideal combination
of proper clinical
procedures,
good product and
well designed
instruments.



Dr. Didier Dietschi



PFIDD1/2



PFIDD3/4



PFIDD5/6



PFIDD7/8



PFIDD9/10



Black Line Rubber Dam Clamps

HU-FRIEDY

Anteriors





Premolars



Upper Molars



Lower Molars



General Molars



£ 16 £20

Single price valid across the entire range



Discover
more



Rubber Dam Instruments

HU-FRIEDY



Rubber Dam Punch,
Ainsworth
17,5cm

RDPA

£ 169


£ 213

Rubber Dam Forceps
17,5cm

RDF4E

£ 97

£ 122



Rubber Dam Frame
Adult, 15 cm

RDAF6

£ 27


£ 34

Rubber Dam Frame
Child, 13 cm

RDCF5

£ 27

£ 34



Rubber Dam Clamp Organizing Board
11,7 x 6,1 x 1 cm

Store, sterilize and organize with eight pegs that securely hold your clamps. Clamps not included.

RDCOB

£ 57

£ 72



DentaPure™ Cartridge CROSSTEX



I highly recommend this for maintaining clean water through DUWL's. DentaPure™ Cartridge, originally developed by NASA, uses Iodine to control harmful bacteria and is extremely cost effective compared to water treatment additives.



Pete Gibbons,
Decontamination
Specialist,
www.deconpete.co.uk



Watch now!
Easy to install

DentaPure™
365 Day Independent
Water Bottle Cartridge,
1 pcs

DPI365B10

£ 140

£215

or buy 2 for only
£ 125 each!

Stages of microbial growth in untreated waterlines



This picture is representative and not sold from Hu-Friedy Mfg. Co., LLC, Crosstex International, Inc., or Omnia Srl.



Compliant[†]

No harsh chemicals. no special disposal requirements



Effective

EPA registered to provide ≤ 10 CFU/mL for one-year or 240L of water^{†††}



Easy

No routine shocking, can be quickly and easily installed by office staff



Reliable

Reduces the possibility of human error. Technology developed for NASA



365 DAYS

One DentaPure™ Cartridge = 365 days^{††} compliant[†] dental unit water


For more information on how the multi-award-winning DentaPure™ Cartridge can help simplify your DUWL maintenance protocols, visit:
hufriedygroup.eu/en/dental-unit-waterlines



Comfort Plus™ Premium Saliva Ejectors

CROSSTEX

Comfort Plus™ Saliva Ejectors are the kinder, gentler saliva ejector that takes in fluid, not tissue. With a smooth, rounded tip designed to enhance patient comfort, it will not grab mucosal tissue. Facilitates uninterrupted fluid and light debris removal. Optimizes fluid removal as tip design prevents tissue from blocking suction.



Code	Description	Promo	Sugg. Retail
ZWGCPM	White with Green Tip Mint Scented, 100/Bag	£ 6	£10
ZWPCPBG	White with Pink Tip Bubble Gum Scented, 100/Bag	£ 6	£10
ZWWCP	White with White Tip, 100/Bag	£ 6	£10

MaxVac Plus Combo-Tip™ Evacuators

CROSSTEX

Manufactured from high density polyethylene. MaxVac Plus Combo-Tip™ Evacuators are two tips in one – vented on one end/non-vented on the other. Unique “S” tip design eliminates sharp edges enhancing patient safety and comfort. Stiff, firm construction allows for cheek retraction.



Code	Description	Promo	Sugg. Retail
ZETMVBL	Blue - 100 pc/pack	£ 3.50	£6
ZETMVGR	Green - 100 pc/pack	£ 3.50	£6
ZETMVLV	Lavender - 100 pc/pack	£ 3.50	£6
ZETMVPK	Pink - 100 pc/pack	£ 3.50	£6
ZETMVYE	Yellow - 100 pc/pack	£ 3.50	£6

Air & Water Syringe Tips

CROSSTEX

Air/Water Syringe Tips feature a Quad Core™ air channel to ensure complete air and water separation. The five translucent colors allow for visualization of water flow/blockage, and no adapter is required for easy and secure installation. A locking groove ensures tip is securely installed and a beveled edge provides increased patient comfort and safety. Single-use eliminates risk of cross-contamination and cleaning of syringe between patients.

Code	Description	Promo	Sugg. Retail
BCSAWSCL	Sparkle™ Air water syringe tips - clear - 8,9 cm - 250 pcs/pack	£ 21	£39
BCSAWSWH	Sparkle™ Air water syringe tips - white - 8,9 cm - 250 pcs/pack	£ 21	£39
BCSAWS	Sparkle™ Air water syringe tips - assorted colors - 8,9 cm - 250 pcs/pack	£ 21	£39
BCSAWSB	Sparkle™ Air water syringe tips - assorted colors - 8,9 cm - 1.500 pcs/pack	£ 112	£238





Lilac Utility Gloves HU-FRIEDY

Lilac gloves can be used in the sterilization area and during disinfection processes in the operatory.

- Can be autoclaved up to 121°C up to five times



Lilac Utility
Gloves,
3 pairs/
Package

£ 22

£ 35

Code	Description
40-060	Size S
40-062	Size M
40-064	Size L
40-066	Size XL

Lead-Free Steam Indicator Tape HU-FRIEDY

Made with adhesive and ink that is free of natural rubber latex, lead or solvents; safer for your patients, staff and the environment.

Not made with natural rubber latex, 18 mm x 55 m

IMS-1461

£ 6

£ 10



Chemical Monitoring

Type 5 IMS™ Steam Sterilization Integrators, Lead free and not made with natural rubber latex - **100 per box**

HU-FRIEDY

IMS-1241

£ 37

£ 58



Instrument Care HU-FRIEDY

Shine reNEW™ Instrument Wipes

Effectively removes stain and rust, restoring the luster to stainless steel and tungsten carbide instruments

20 wipes per canister

IMS-1455

£ 10

£ 15



Instrument Lubricant Spray

Keeps hinged instruments moving freely

237 ml

ILS

£ 18

£ 29



Pre-Cleaner HU-FRIEDY

Dual enzyme gel with corrosion inhibitor acts as a ready-to-use instrument pre-cleaner. Maintains moisture to help prevent debris from drying on instruments before cleaning and helps dissolve blood, tissue, and other debris. Can also be used as a linen spot remover.

Enzymax™ Spray Gel Ready-to-Use

710 ml

IMS-1229

£ 15


£ 24





IMS™ Cassettes are manufactured from high quality, electroplished stainless steel shells to protect against corrosion while providing a long-lasting, professional appearance. IMS™ Cassettes feature medical-grade silicone rails, robust and ergonomic latches with positive locking mechanism, and are available in multiple sizes and colors that can be configured to best fit the needs of your facility.

IMS™ Infinity Series™ Cassettes HU-FRIEDY




1/4 DIN CASSETTE
5 Instruments
72 x 184 x 34 mm
IME14DIN5_

£ 88
£124



1/2 DIN CASSETTE
10 instruments
152 x 184 x 34 mm
IME12DIN1_

£ 98
£139



DIN CASSETTE
290 x 184 x 34 mm

COLOUR CODE OF CASSETTES AND TUBS


Replace the _ column with the number of your desired colour in the code.

- 0
- 1
- 2
- 3
- 4
- 5
- 6
- 7
- 8
- 9
- 11

Code	Description	Promo	Sugg. Retail
IMEDIN13_	13 Instruments - with accessory space	£ 133	£188
IMEDIN16_	16 instruments - with accessory space	£ 133	£188
IMEDIN18_	18 instruments - with accessory space	£ 133	£188



Infinity Series™ Specialty Cassettes HU-FRIEDY



**Ultrasonic Insert/Handpiece
Cassette - Small**
38 x 203 x 32 mm
Holds 2 Ultrasonic Inserts
or 1 handpiece

IMN30HU_

£ 77
£109

COLOUR CODE OF CASSETTES AND TUBS

Replace the _ column
with the number of your
desired colour in
the code.

0

1

2

3

4

5


6

7

8

9

11



**Ultrasonic Insert/Handpiece
Cassette - Large**
76 x 203 x 32 mm
Holds up to 3 high-speed handpieces
or 5 ultrasonic inserts

IMN60H_

£ 86
£122

IMS™ Tubs HU-FRIEDY

**Complete Tub with lid and
accessories**

IMS-140_

£ 57
£79



Who We Are

We are a combination of companies that have come together under one roof to deliver the products, services, and solutions you need to be the best in practice. With our passion for infection prevention solutions, premium instrumentation, and service and support programs, HuFriedyGroup will help you perfect your craft and continue to deliver the highest level of care to each patient.



Our Contacts



HuFriedyGroupEurope



HuFriedyGroup.eu



Free Call 00800 48 37 43 39



For more information please contact:

Chris Mason (Country Manager UK & Ireland)

Tel.: 07703 186 474 | cmason@hu-friedy.com



@HuFriedyGroupUK

Hu-Friedy Mfg. Co., LLC • European Headquarters • Lyoner Str. 9 • 60528 Frankfurt am Main, Germany • HuFriedyGroup.eu

As some of the products included in this catalogue might not be available in some countries, we invite you to verify the availability of the product with your local authorised dealer or to contact our customer care at the free phone number 00800 48 37 43 39.

Available while stocks last. All prices indicated are the recommended prices in UK and don't include VAT and other taxes.

All company and product names are trademarks of Hu-Friedy Mfg. Co. LLC, its affiliates or related companies, unless otherwise noted. Sontara® is trademark of Sontara® South Asia and Tyvek® is trademark of DuPont Safety & Construction, Inc.

©2025 Hu-Friedy Mfg. Co., LLC All rights reserved. HFL-058/0325