

Special offers

Winter / Spring 2025

EXTENDED UNTIL 30/04/2025

~~VALID UNTIL 31/03/2025~~

Decon and IPC

New iClave 24



+ 12 months additional NSK
handpiece warranty (p12)

Implantology and Oral Surgery

New VarioSurg 4



Minimally invasive ultrasonic bone surgery system
Plus option to connect to Surgic Pro2 with one foot control.

Perio and Prevention



Restorative



Restorative - Turbines

BDIA DENTAL SHOWCASE

14 - 15 MARCH 2025

JOIN NSK ON STAND J50

Ti-Max Z HIGH PERFORMANCE TURBINES

2 + 1 **FREE*** + 1 **FREE Z micro**



NEW

ULTIMATE POWER for quicker preps and cutting Zirconia

• Lightweight titanium body with Duragrip coating to reduce hand fatigue



4 years for iClave users



- Optic
- **ULTIMATE Power 30-44w****
- Std head: Ø 12.1 x 12.7 mm

MODEL Z990L
CODE PA23730001
RRP £1,173

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®

**At 3 - 4.2 bar



- Optic
- **Power 28-41w****
- Mini head: Ø 10.6 x 10.8 mm

MODEL Z890L
CODE PA23720001
RRP £1,173

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®



- Optic
- **Power 20w**
- **Micro head: Ø 9 x 10.8 mm**

MODEL Z micro (ultra-mini head)
CODE PA23570001
RRP £1,173

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®

for regular or short shank burs

Scan/Click
For full info on models/order codes



S-Max M



3 years for iClave users



OPTIC 3 + 1 **FREE** + 1 **FREE Pico***



• Stainless steel • Enhanced power • 4 hole spray • DIY cartridge change • Clean head system



- Optic
- **Power 26w**
- **Std head: Ø 12.1 x 13.3 mm**

MODEL M900L
CODE P1254001
RRP £773

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®



- Optic
- **Power 23w**
- **Mini head: Ø 10.6 x 12.4 mm**

MODEL M800L
CODE P1253001
RRP £773

Also available in fittings for KaVo®



- Optic
- **Power 9w**
- **Ultra-mini head: Ø 8.6 x 9 mm**

MODEL Pico
CODE P1140001
RRP £1,121

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®

Used only with short shank burs.



- **Non-Optic**
- **Power 26w / 23w**
- **Std and mini heads**

MODEL M900
CODE P1256001
RRP £720

MODEL M800
CODE P1255001
RRP £720

M900 also available in fittings for KaVo®, W&H®

Scan/Click
For full info on models/order codes



Pana-Max2



2 years for iClave users



1 + 1 **FREE***

• Non-Optic turbine with single hole spray • Save time/money with DIY cartridge change • Clean head system • **12 month warranty**



- **Non-Optic**
- **Head: Ø 11.2 x 13.4 mm**
- **Coupling connection**

MODEL Pana-Max2
CODE P1211001
RRP £468

for NSK PTL coupling

- **Non-Optic**
- **Head: Ø 11.2 x 13.4 mm**
- **Coupling connection**

MODEL Pana-Max2 KV
CODE P1212001
RRP £468

for KaVo® MULTiflex® coupling



- **Non-Optic**
- **Head: Ø 11.2 x 13.4 mm**
- **Direct connection**

MODEL Pana-Max2 M4
CODE P1154001
RRP £425

For Midwest M4 hose

- **Non-Optic**
- **Head: Ø 11.2 x 13.4 mm**
- **Direct connection**

MODEL Pana-Max2 B2
CODE P1155001
RRP £425

for Borden hose

Scan/Click
For full info on models/order codes



*Terms and conditions apply. Offers only available through participating dealers.

Restorative - Contra Angles

BDIA DENTAL SHOWCASE

14 - 15 MARCH 2025
JOIN NSK ON STAND J50

Ti-Max Z

2 + 1 FREE + 1 FREE Z micro*



- Lightweight external titanium body • Optic with internal water • Small head/narrow neck
- Unique spray switch on Z95L, Z45L, Z25L



- Optic
 - 1:5 Increasing+Spray switch
 - Head size $\varnothing 8.9 \times 13.4\text{mm}$
- MODEL Z95L
CODE C1034001-SW
RRP £1,610



- Optic
 - 1:5 Increasing
 - 10% shorter & lighter than Z95L
 - Head size $\varnothing 8.9 \times 13.4\text{mm}$
- MODEL nano95LS
CODE C1099001
RRP £1,533



- Optic
 - 1:4.2 Increasing+Spray switch
 - 45° head
 - Head size $\varnothing 9.6 \times 14.5\text{mm}$
 - Ideal for minor oral surgery too
 - For regular or extra long (25mm) FG burs
- MODEL Z45L
CODE C1064001
RRP £1,610



- Optic
 - 1:1 Direct Drive +Spray switch
- MODEL Z25L
CODE C1038001-SW
RRP £1,276



Scan/Click

For a full range of NSK Contra-angles



3 years for iClave users



S-Max M

3 + 1 FREE + 1 FREE Pico*



- Stainless steel • All with internal water • Choice of optic & non-optic



- Optic
 - 1:5 Increasing
 - Head size $10 \times \text{H}14.6\text{mm}$
- MODEL M95L
CODE C1023001
RRP £1,040

- Non-Optic
 - 1:5 Increasing
 - Head size $10 \times \text{H}14.6\text{mm}$
- MODEL M95
CODE C1026001
RRP £863



- Optic
 - 1:1 Direct Drive
- MODEL M25L
CODE C1024001
RRP £842

- Non-Optic
 - 1:1 Direct Drive
- MODEL M25
CODE C1027001
RRP £553



- Optic
 - 4:1 Reducing
- MODEL M15L
CODE C1025001
RRP £985

- Non-Optic
 - 4:1 Reducing
- MODEL M15
CODE C1028001
RRP £658



- Non-Optic
 - 1:1 Straight
- MODEL M65
CODE H1008001
RRP £541

FXm

1 + 1 FREE*



2 years for iClave users



FX

3 FOR £195, 10 FOR £550*

New modern version of the classic NAC-EC



- Non-Optic Straight
 - 1:1 Direct Drive
 - Push button chuck
 - No internal water
- MODEL FX65m
CODE H1014001
RRP £247



- Non-Optic Contra-angle
 - 1:1 Direct Drive
 - Push button chuck
 - No internal water
- MODEL FX25m
CODE C1052001
RRP £273



- Non-Optic Propphy Contra-angle
 - 4:1 speed reducing
 - Max speed 1,000 rpm
 - Brushes/Cups: Screw in
- MODEL FX57m
CODE C1054001
RRP £307



- Non-Optic Contra-angle
 - 1:1 Direct Drive
 - Latch-type chuck
 - No internal water
- MODEL FX22
CODE C1129001
RRP £160

NEW

Scan/Click

For a full range of NSK Contra-angles



*Terms and conditions apply. Offers only available through participating dealers.

Restorative - Electric

BDIA **DENTAL SHOWCASE**

14 - 15 MARCH 2025

JOIN NSK ON STAND J50

NLX

RETRO FITTED MICROMOTORS

Electric micromotors are the driving force behind speed increasing handpieces and can be fully integrated into the delivery system of most new and existing dental chairs. They can be retro fitted by an engineer or come as an option when buying new chairs. Controlling the specific speed and torque of micromotors within existing platforms is made easy thanks to built-in control pads.

Speed increasing contra-angles not only do everything a turbine does using the same FG burs, they provide a constant torque from the micromotor that does not decrease with resistance and higher load with no loss of power as they come in contact with the tooth.

Scan/Click



For more info



To discuss the offers & options contact NSK:

OFFER 1

BUY 3 x



AND GET **OPTION A FREE**



OFFER 2

BUY 4 x



AND GET **OPTION B FREE**



OFFER 3

BUY BUILT-IN OPTION A or B AND GET FREE Ti-Max Z95L SPEED INCREASING CONTRA-ANGLE (1:5)



OPTION A

NLX nano micromotor & BUILT-IN KIT * (including MultiPad) RRP £2,307

OPTION B

NLX plus micromotor (includes Rotary Endo functions) & **BUILT-IN KIT *** (including MultiPad) RRP £2,774

OPTION C

Varios170 Lux ULTRASONIC HANDPIECE (tip options for scaling, perio, endo & restorative procedures) + **BUILT-IN KIT *** (including MultiPad) RRP £1,501 (Y1004145)

OPTION D

Option A + C RRP £3,654

OPTION E

Option B + C RRP £4,095

*Excludes install. Through participating dealers only.

NLX nano

AMAZING CLEARANCE OFFER

AVAILABLE WHILE STOCKS LAST!

PORTABLE MICROMOTOR

- Portable electric micromotor for connecting slow, straight and speed increasing handpieces
- Puts you in control - set your own speeds
- Connects to your existing cart

MODEL NLX nano Portable
CODE Y1001362
RRP £2,721



OFFER



+



£1,175 ~~£3,761~~
SAVE £2,586

NLZ E **NEW**

PORTABLE MICROMOTOR

Including reciprocating and rotary Endo features

NEW compact portable electric micromotor, for general & endodontic applications with contra angles (ie speed increasing/red band, 1:1 direct drive/blue band or speed reducing/green band).

- Lighter motor with enhanced power/torque
- Ideal for general (restorative, prophylaxis etc) and endo applications (rotary & reciprocating* modes)
- Enhanced patient safety function indicates when the handpiece needs repair to minimise it overheating
- Versatile installation options – portable or semi built-in

MODEL NLZ E
CODE Y1004004
RRP £3,250



*The 'Reciprocating Endo Mode' can be operated only with the NLZ Endo Contra-angle and WAVE ONE® file.

The NLZ requires access to a plug socket, and the use of an existing Midwest hose from the cart to gain access to coolant water/air and the existing foot control.

INCLUDES NLZ Endo 6:1 CONTRA ANGLE (for rotary & reciprocating)



Endo-Mate TC2

CORDLESS ENDO MOTOR/HANDPIECE

- Use with rotary files
- MODEL Endo-Mate TC2
CODE Y1001027
RRP £1,769



Endo-Mate DT2

CORDED ENDO MOTOR/HANDPIECE

- Use with rotary files



MODEL Endo-Mate DT2
CODE Y1004230
RRP £1,672

iPex II

APEX LOCATOR



MODEL iPex II
CODE Y1002208
RRP £1,050



OFFER BUY Endo-Mate DT2 AND GET iPex II FREE WORTH £1,050*

Varios Tips

ULTRASONIC TIPS FOR ENDO

Designed to aid agitation and irrigation. Compatible with stand alone Varios 970 / Varios Combi Pro, Varios 170 (built-in).



- E5 & E6 for lateral condensation
- E4D (diamond coated) for cleaning/smoothing or widening canals
- E7 / E8 for removing filling material from the canal
- E15D (diamond coated) for locating root canals & decalcification

MODEL Varios Endo tip kit
(tip stand + 6 tips above)
CODE Z221086/6

OFFER £255.10

Equivalent of 3+3 free tips + free tip stand

Additional tips RRP £83 **Buy 1 get 1 free**
(diamond tips £89.10 Buy 1 get 1 free)



Scan/Click
For info article



Surgic Pro and Surgic Pro2

MICROMOTORS FOR IMPLANTS / ORAL SURGERY

Scan/Click
For more about
Surgic Pro2



Including the new stylish Surgic Pro2 with a cordless foot control, offering wireless connectivity to VarioSurg4 and Osseo 100+.



MODEL Surgic Pro (optic)
CODE Y1003586
RRP £6,212

MODEL Surgic Pro (non-optic)
CODE Y1003587
RRP £4,637

Scan/Click
For more about
Surgic Pro



- Complete Kits contain:
- Control unit
 - Motor/Cord
 - 20:1 Implant contra angle
 - Foot control
 - 2 year warranty



MODEL Surgic Pro2 (optic)
CODE Y1004195
RRP £6,212

MODEL Surgic Pro2 (non-optic)
CODE Y1004196
RRP £4,637

FREE



Osseo 100+

OR



iSD900

OR 2 x



20:1 implant
contra-angles

ADDITIONAL MOTOR/CORD



• For Surgic Pro

MODEL SGL70M (optic)
CODE E1023001
RRP £2,758

MODEL SG70M (non-optic)
CODE E1025001
RRP £1,917



• For Surgic Pro2

MODEL SGL80M (optic)
CODE Y1004211
RRP £2,758

MODEL SG80M (non-optic)
CODE Y1004212
RRP £1,917

FREE 20:1 IMPLANT CONTRA-ANGLE OR SURGICAL STRAIGHT

WHEN YOU BUY AN ADDITIONAL MOTOR/CORD



SURGICAL STRAIGHTS 1:1

1 + 1 FREE*



- Used with standard HP burs
- **2 year warranty**
(3yrs for iClave users)

• Optic

MODEL Ti-Max X-SG65L
CODE H1009001
RRP £965

• Non-Optic

MODEL Ti-Max X-SG65
CODE H1038001
RRP £658

IMPLANT CONTRA-ANGLES 20:1

1 + 1 FREE*



- Max speed: 2,000min⁻¹
- **2 year warranty**
(3yrs for iClave users)

• Optic

MODEL Ti-Max X-SG20L
CODE C1003001
RRP £1,286

• Non-Optic

MODEL S-Max SG20
CODE C1010001
RRP £1,085

ZYGOMATIC CONTRA-ANGLE

1 + 1 FREE*


• 20:1 Reduction

- Twist chuck
- **12 Month warranty**
(2yrs for iClave users)

Handpiece for Zygomatic & Pterygoid implants. Its unique angulated shape offers optimal access and visibility.

• Non-Optic

MODEL SGX-E20R
CODE HA1200001
RRP £1,582



*Terms and conditions apply. Offers only available through participating dealers.

VarioSurg 4

ULTRASONIC BONE SURGERY SYSTEM

NEW

FREE 12 EXTRA TIPS OR AN OSSEO 100+

With its stylish design and wireless foot control, VarioSurg 4 builds on the reliability and control of previous models to deliver a modern ultrasonic bone surgery system, with Piezo technology to minimise damage to surrounding tissue and bone.

Contains

- Slimline lightweight LED handpiece
- Stylishly designed control box
- Wireless foot control**
- Tip Stand with 6 tips - SG11, SG6D and SG7D for sinus lifts. H-SG1 for cutting bone with SG3 and SG5 for scraping.



Scan/Click
Extensive range of tips



Scan/Click
For the core range of tips

MODEL VarioSurg 4
CODE Y1500706
RRP £8,369
(option with corded foot pedal Y1500712 + Z1102003)

VarioSurg tips 1 + 1 FREE*



IRRIGATION TUBING

- Used with Surgic Pro / Pro2 / XT and VarioSurg 3 / 4



Irrigation tubing x5 (disposable)
CODE Y900079
RRP £64.30



NSK Eco Aware!

Autoclavable Irrigation Tubing x1
CODE S900055
RRP £53.80

LINK OPTIONS

Already have a Surgic Pro2?

Link it to a VarioSurg 4 and control both with 1 foot control**

VarioSurg 4 (WITHOUT FOOT CONTROL)



50% OFF FOR Surgic Pro2 USERS

+ FREE LinkStand 3 OR iCart Duo



Link Stand for VarioSurg 4 / Surgic Pro2

MODEL Link Stand3
CODE ZA16230001
RRP £311



iCart

MODEL iCart Duo
CODE S9090
RRP £1,041



*Terms and conditions apply. Offers only via participating dealers. **Only the Surgic Pro2 foot control (FC-81) is designed to link and control both units.

Scan/Click
Watch VarioSurg 4 and
Surgic Pro2 link together



Implantology & oral surgery

BDIA **DENTAL SHOWCASE**

14 - 15 MARCH 2025

JOIN NSK ON STAND J50

Osseo 100+

OSSEOINTEGRATION MONITORING DEVICE

Osseo 100+ measures implant stability and osseointegration to enhance decisions about when to load the implant. This is invaluable when using protocols with shorter treatment times and treating higher risk patients.

The simple steps to measure the ISQ value allow the implant loading period to be planned in advance.

MODEL Osseo 100+
CODE Y1004176
RRP £2,660

£1,862



Reusable pegs!

Model: Multipeg
RRP £43.20

Scan/Click
More info &
Multipeg List



“The Osseo device allows us to easily measure implant stability by assessing the increasing ISQ readings during the healing period. We use it routinely on sinus graft cases, when we place simultaneously, and also in cases with low stability at placement.”

The Campbell Clinic, Nottingham

Mobile dentistry

VIVA ace

MOBILE DENTISTRY SYSTEM

- Ideal as a surgery back up or to visit shielded patients
- Lightweight micromotor with additional endo functions
- Multipurpose ultrasonic handpiece
- Quick setup and storage
- 4-level adjustable height and angle



Width 462 mm, Height 284 mm, Depth 292 mm



MODEL VIVA ace
CODE Y1003773
RRP £12,955

Weight 8.6kg

iProphy

CORDLESS PROPHY HANDPIECE

The complete kit comes with a 4:1 reduction head for screw-in prophy cups/brushes. Alternatively a PR-F head (optional extra) can be used that connects to disposable Prophy heads.



MODEL iProphy Complete Kit
CODE Y1500015
RRP £1,122

FREE* 1 x  M95L  M25L  M65

Monthly **RENTAL** £433* Fully maintained for 3 years

PACKAGE OFFER with Prophy-Mate neo

£1,375 ~~£2,285~~

SAVE £910



Perio & prevention

Scan/Click
For more info on
NSK Oral Hygiene



BDIA **DENTAL SHOWCASE**

14 - 15 MARCH 2025

JOIN NSK ON STAND J50

Prophy-Mate neo SUPRAGINGIVAL POWDER THERAPY

- Compact, lightweight & portable
- Includes 60° and 80° nozzles
- 360° rotation
- Includes 6 sachets of FLASH pearl® powder

CODE Y135029 for NSK PTL coupling*
CODE Y135030 for KaVo® MULTIFlex® LUX coupling
CODE Y135031 for Sirona® quick coupling
CODE Y135032 for W&H® Roto Quick® coupling
CODE Y135033 for Bien-Air® Unifix® coupling

RRP £1,163



Scan/Click
For more info



+  **FREE SECOND HANDPIECE WORTH £441***

Prophy-Mate neo nozzles



CODE T797012 60° nozzle
RRP £124.20

CODE T797013 80° nozzle
RRP £124.20

FLASH pearl



- Calcium carbonate (54µm) based powder for stain removal

4 x 300 g bottles
CODE Y900698
RRP £142.80

100 x 15 g sachets
CODE Y900693
RRP £195



SOFT pearl



- Glycine based powder (65µm), ideal for moderate staining. Water soluble, making clean up easier.

4 x 250 g bottles
CODE Y1500721
RRP £137.60

NEW

1 + 1 FREE*

**BUY ANY 2 PACKS OF FLASH OR SOFT pearl GET A
+ FREE*** 

COMBO OFFER Prophy-Mate neo + Perio-Mate £1,928* £2,488

SAVE £560

Perio-Mate

SUBGINGIVAL POWDER THERAPY / BIOFILM REMOVAL

- Subgingival version of ProphyMate
- Eradicates biofilm in just 20 seconds
- Powder flow adjuster for delicate pockets
- Includes flexible plastic nozzle tips for deeper pockets
- Packed in a handy case with 160g powder
- Also used supragingivally on moderate staining and around composite restorations/orthodontic brackets



Scan/Click
For more info



CODE Y1002654 for NSK coupling
CODE Y1002655 for KaVo® MULTIFlex® LUX coupling
CODE Y1002656 for Sirona® quick coupling
CODE Y1002657 for W&H® Roto Quick® coupling
CODE Y1002658 for Bien-Air® Unifix® coupling

RRP £1,325



+  **FREE SECOND HANDPIECE WORTH £441***

FREE



- Glycine based powder (25µm)
 - 2x 160 g bottles
- Perio-Mate Powder
CODE Y900938
RRP £99.30



- Pack of 40
Perio-Mate Nozzle Tips
CODE Y1002741
RRP £99.40

POWDER THERAPY PATIENT LITERATURE



Available on request with Prophy-Mate neo, Perio-Mate and Varios Combi Pro purchases.

Scan/Click



*Terms and conditions apply. Offers only available through participating dealers.

Varios Combi Pro iPiezo engine®

TWO ORAL CARE SYSTEMS IN ONE!



Ultrasonic Treatment

Powder Therapy

Combines ultrasonic treatments and powder therapy in one unit.
Comes with a detachable Prophy chamber.
Detachable Perio chamber available as an optional extra.

MODEL Varios Combi Pro
CODE Y1002843
RRP £5,151



Portable cart for Varios Combi Pro



MODEL iCart Prophy
CODE S9095
RRP £1,041

Note: Connection to external 6mm air & water line sockets are needed to run this unit.

FREE Perio Set WORTH £1,231*



Scan/Click



For more info

Varios 970 iPiezo engine®

ULTRASONIC SCALER

Ultrasonic scalers not only offer patients a more comfortable experience as less water is needed, but offer the user a wide range of 70 tips, not only for scaling/perio/implant maintenance, but for endo, extractions, restorative & MI applications.

Units come with:

- 1 handpiece
- G4, G8 & G16 tips
- 3 tip wrenches
- sterilisation case

• Optic

MODEL Varios 970 LUX
CODE Y1004161
RRP £2,999

• Non-Optic

MODEL Varios 970
CODE Y1004160
RRP £2,581



Scan/Click



For more info



OFFER

FREE

2 x



2 x



Varios2 handpiece

• Optic

MODEL VA2-LUX-HP
CODE E351050
RRP £549



Varios2 handpiece

• Non-Optic

MODEL VA2-HP
CODE E350050
RRP £307

1 + 1 FREE*

1 + 1 FREE*

SCALER TIPS 1 + 1 FREE*



Scan/Click
Full range
of tips



Scan/Click
To order



*Terms and conditions apply. Offers only available through participating dealers.

Scaler Tips

1 + 1 FREE*

TIPS FOR VARIOS
ULTRASONIC SCALER




Scan/Click
Full range
of tips

TIPS FOR AIR SCALERS



Scan/Click
Tips for
Air Scaler



AIR SCALERS



Direct connection
air scaler

- Oscillation frequency 5,600-5,850Hz
- Kit includes 3 scaler tips

MODEL AS2000
CODE T698001 (Midwest M4 hose)
CODE T697001 (Borden B2/B3 hose)
RRP £493 each



Optic air scaler with
coupling connection

- Oscillation frequency 6,200-5,850Hz
- Scratch resistant with titanium body
- Kit includes 3 scaler tips

MODEL Ti-Max S970L
CODE T1005001**
RRP £1,385

**also available in KaVo® fitting T1013001
and Sirona® fitting T1015001

FREE FX22



FREE 2 x FX22



IKIGAI

ORAL HYGIENE COMMUNITY

For all the latest in Perio &
Prevention visit IKIGAI on



FACEBOOK



INSTAGRAM
@ikigaioralhygiene




PROPHY CONTRA-ANGLES

1 + 1 FREE*



- Max speed 1,000 rpm
 - Brushes/Cups: **Screw in**
- MODEL FX57m
CODE C1054001
RRP £307

- Max speed 40,000 rpm
 - Brushes/Cups: **Disposable Prophy heads**
- MODEL FX65m
CODE H1014001
RRP £247



Paste defence system

3 FOR £195, 10 FOR £550*



AR-EC(k)

AR-EC(s)

FX22

- Max speed 5,000 rpm
 - Brushes/Cups: **Snap on cups**
- MODEL AR-EC(k)
CODE Y110165
RRP £152

- Max speed 5,000 rpm
 - Brushes/Cups: **Screw in**
- MODEL AR-EC(s)
CODE Y110164
RRP £152

NEW

- New modern version of
the classic NAC-EC
- Max speed 30,000 rpm
- Brushes/Cups: **Latch grip chuck**

MODEL FX22
CODE C1129001
RRP £160

*Terms and conditions apply. Offers only available through participating dealers.

Decontamination & IPC

BDIA **DENTAL SHOWCASE**

14 - 15 MARCH 2025

JOIN NSK ON STAND J50

AFFORDABLE, RELIABLE & RAPID STERILISATION

iClave 24

NEW

Large B-type vacuum autoclave 24 litre



- Perfectly combines performance, user experience, & technology, ensuring an effective, comfortable, and worry-free experience.
- 24L chamber with 5 trays of 2 different dimensions maximises capacity.
- High-efficiency steam generator allows for superior performance
- B-cycle 40 mins, S-cycle 31 mins
- User-friendly colour touchscreen
- Footprint: W451 x H492 x D669mm

CODE 92412B0100
RRP £5,725

FROM ONLY £4,995

+ install & delivery

RENT for just

£195pm[†]

iClave plus

B-type vacuum autoclave 18 litre



- Easy to use + USB data storage
- Fast 25-minute cycle for unwrapped loads
- 35 minutes for pouched loads (+econ cycle 20 minutes – one large tray; 15 minutes – two small trays)
- 18L cylinder-shaped chamber offers 20% more space than traditional chambers & is ideal for implant kits
- Includes Helix Test + Strips
- Footprint: W445mm x H428 x D532mm

CODE DXY1003077
RRP £6,413

FROM ONLY £3,995

+ install & delivery

RENT for just £155pm[†]

iMax S

S-type autoclave 18 litre



- Rapid 20-minute cycles including drying for unwrapped loads
- 25-minute handpiece cycle with vacuum drying
- 18L chamber for improved performance and space
- Footprint: W450 x H430 x D600mm

CODE DX91803S0200
RRP £3,910

ONLY £3,070

+ install & delivery

RENT for just £115pm[†]

12 MONTHS ADDITIONAL WARRANTY on new NSK handpieces or repairs for iClave users**

PROTECT IT

*Visit mynsk.co.uk/rental for more information. **Additional warranty applies to NSK handpieces (new or repaired by NSK), that are still within their warranty period.

All-Inclusive Maintenance and Breakdown Contracts for NSK autoclaves from just **£64pm**.



One annual maintenance visit or every 1500 cycles



Additional maintenance visit for units exceeding 1500 cycles



All service parts and service sundries included



All breakdown calls*



All breakdown parts including labour*



All engineer travel in the UK*



Free on-site operator training



Free technical software upgrades



More info on the NSK decon range and the Smart3 System of working

iCare

Handpiece lubrication unit



Precise dosing - less waste, better care

MODEL C2
CODE Y1002787
RRP £2,659

MODEL C3
CODE Y1002788
RRP £2,795

Connects to 6mm airline, 5-6 bar

NSK Eco Aware!

£2,000 OF FREE NSK HANDPIECES

2 × Optic Turbine 1 × Straight 1 × 1:1 Contra-angle

OR HALF PRICE £1,330

RENT for just £73.95pm[†]

CONTACT NSK

National Sales Manager
Mark Beckwith
mbeckwith@nsk-uk.com 07900 246529

N. Ireland & ROI
Stephanie Gilliland
sgilliland@nsk-uk.com 07435 840097

Scotland - Sales Manager
Angela Glasgow
aglasgow@nsk-uk.com 07525 911006

North West / North East / North Wales
Clare Morris
cmorris@nsk-uk.com 07901 117669

Derbys/Notts/Yorkshire/Lincs
Jo Dronfield
jdrnfield@nsk-uk.com 07545 786619

West Midlands
Andy Hargreaves
ahargreaves@nsk-uk.com 07872 824175

Leics / Northants / CV & PE postcodes
Graham Brown
gbrown@nsk-uk.com 07541 864641

South Wales
Robert Musgrove
rmusgrove@nsk-uk.com 07798 918741

South West
Ben Guest
bguest@nsk-uk.com 07350 361543

South Central
Steve Barnes
sbarnes@nsk-uk.com 07889 181099

Central London
Teruyuki Nakanishi
tnakanishi@nsk-uk.com 07570 293317

South East & East Anglia
Francois Faro
ffaro@nsk-uk.com 07900 246822

Field Service Engineer - North
Jamie Duncombe
jduncombe@nsk-uk.com 07780 105457

Field Service Engineer - South
Lee Morris
lmorris@nsk-uk.com 07435 783707

NSK United Kingdom LTD www.mynsk.co.uk

Office 4, Gateway 1000 Arlington Business Park, Whittle Way Stevenage, SG1 2FP, UK - info@nsk-uk.com - 0800 6341909



*Offers cannot be combined. NSK reserve the right to change prices and offers without notice. Published prices are RRP's and exclude VAT and installation. Offers only available through participating dealers. Additional warranties available from NSK UK on new products and NSK UK repairs. To qualify for specific offers please send a copy invoice to NSK UK at Office 4, Gateway 1000, Arlington Business Park, Whittle Way, Stevenage, Hertfordshire SG1 2FP or email sales@nsk-uk.com. Offers end 31/03/25.