

WOODPECKER Ultrasonic Scalers **U-600 Scaler**

- Intelligent power output for patient comfort
- Smart touch waterproof interactive panel for ease of use
- Includes: 2 handpieces, 2 boxes for autoclaving and 10 tips
- Features General, Endo and Perio functions

£595.00

UDS-E + LED Scalers

- Auto-water supply function
- Automatic frequency testing
- EMS compatible
- LED and Non-LED handpiece options
- PC controlled for easy operation and more efficient scaling

LED

£350.00

Non-LED

£325.00

UDS-L Scaler

- Complete with 7 tips and torque wrench
- Comes with self contained system water bottle
- Autoclavable detachable handpiece

£239.00

WOODPECKER EMS Compatible Tips

- G1, G2, G3, G4
EMS compatible scaler tips

OFFER £12.00



**Scalers
for all
budgets!**



* DTE/ Satelec compatible
scaler tips also available

WOODPECKER **iLED Curing Light**

- One second cure for 2mm of resin
- 360° fully rotatable head
- Ergonomic design
- Rechargeable lithium battery
- Two programs power output 1300 mW/cm² up to 2300 mW/cm²
- Available in White or Purple

£125.00

LED-B Curing Light

- Four different time settings 5, 10, 15 & 20 seconds
- Automatic sleep mode
- Lithium battery
- 70 minutes of use per charge
- Lightweight

£119.00

LED-E Curing Light

- Three different modes
- Automatic sleep mode
- 400 minutes per charge
- Features a replaceable/rechargeable battery
- Built in light meter

£155.00

LM-I Light Meter

- More than 95% accurate
- Measuring range: 0 - 3000 mW/cm²
- Wavelength range: 400 - 500 nm
- Light, compact design
- Measures LED and Halogen Curing Lights

£26.95



**WOODPECKER**

Woodpex III Apex Locator

- Automatic calibration ensures measurements are accurate
- LED indicator and clear bright 4.5" LCD display
- Autoclavable accessories: liphook, file clips and probe
- Calibration tester
- Rechargeable battery

£239.00**DENTAL ELECTRONIC**

Pulptester DP2000

- Cordless battery operated pulptester with digital readout
- Sends low current pulses through the tooth with automatically increasing amplitude
- Strength of the pulse is simultaneously indicated by a flashing LED
- A ground lead is available for use with gloved hands

**£235.00**

TRAUS Endo Cordless Motor System

- Simple interface design for easy control
- Ergonomic and compact body
- 5 memory programs
- RPM: 100 - 600 and auto reverse
- Torque: 1 - 4 N.cm

£440.00

TRAUS E-Cube Endo Motor System

- Complete with 16:1 contra angle
- 9 Programs & up to 5-hour battery life
- Auto reverse function
- Torque control from 0.6 to 5.2 N.cm
- Foot on/off switch
- Compact designed controller

£480.00

DELDENT

Miniblaster Clinical Sandblaster

- For left or right hand operators
- Index finger controlled operating valve
- In-line high pressure stopcock provided for ease of connection



Midwest



Sirona



Kavo



W&H

All Connectors
£44.00

£189.00



Dust-Inn 2000 Dust Collector

- Compact and lightweight chamber
- Fitted with a quiet fan extractor system
- Abrasive collects beneath the perforated work table for simple disposal
- Easily changeable filters

£189.00



Jetpolisher 2000 Prophy

- Unique advanced air polishing system
- Guaranteed non-clog nozzle
- Harmless to titanium implant surfaces
- Connects directly to existing high speed handpiece tubing
- Lightweight and easily transportable
- Uses the regular foot control of the dental unit
- Fully-detachable autoclavable handpiece

£749.00



All prices exclude VAT E & OE Valid 1st October - 31st December 2017

Trycare, Rialto Buildings, Station Road, Clayton, Bradford BD14 6JA
01274 88 55 44 • www.trycare.co.uk • dental@trycare.co.uk