

elexxion **pico**

The bare essentials:

der elexxion pico

Finally – a dental laser “Made in Germany” that combines mobility **and** professionalism.



Diode	max. 5 watt pulse	20.000 Hz pulse rate	808 nm diode
-------	-----------------------------	--------------------------------	------------------------



Image subject to model

Trycare

01274 885544 www.trycare.co.uk

elexxion
dental laser

Scope of delivery: Everything you need

The basic system comprises the pico laser unit, plus:

- 3 claros protect laser safety glasses
- 1 elexxion pico step foot switch
- 1 Ergoflex plus handpiece
- 1 ergo T handpiece
- 1 claros ergo T8 glass rod
- 1 bending tool
- 8 soft metal tips (two different sizes)
- 1 starter set of fibres (replaceable, multiple-use fibres: 2 x 200 µm, 2 x 300 µm, 3 x 400 µm, 2 x 600 µm)
- 1 table stand
- 1 charging cable
- 1 case
- 1 set of operating instructions, 1 document folder
- "Laser" wall sign

Technical data:

- Class 4 dental laser
- 230 V / 50Hz/60Hz power supply
- Wavelength: 808 nm +/- 10 nm
- Max. pulse output: 5 W
- Max. average output: 2.5 W
- Max. pulse rate: CW – 20,000 Hz
- Pulse duration: 26 µs / CW
- Dimensions: H 11.5 cm, W 9 cm, D 2.8 cm
- Weight: 590 g



Competitor's laser



elexxion laser

Five buttons.

One complete treatment programme.

So brilliantly simple – and professional, too! With just five "buttons", you can navigate your way through a wealth of scientifically recognised treatment options. It is well known that high-quality dental lasers enable you to achieve results more quickly and to work more effectively. That also spares the nerves of your patients. A positive side effect is that word gets around in no time that your practice is using a state-of-the-art dental laser for treatment, giving you a clear advantage among potential patients, too!

Essential: mobility

The pico allows you to work completely independently from the mains electricity supply. This means that you can get by with just one dental laser because the pico is mobile and can be used in any treatment room.

Practical: pre-programmed

By simply pressing on the touch-sensitive surface of the pico, you can choose from five different programs matched precisely to the needs of your practice and start treatment immediately. The pico is ideal for all practitioners who value straightforward installation and a short learning curve.

Ergoflex plus handpiece with 360° rotation

The ergonomically shaped handpiece is sterilisable and enables safe, controlled operation. Ergoflex plus means a single handpiece for all applications, with glass fibres of four different strengths and individually adjustable lengths. With its sterilisable and re-usable fibres, the complete application system corresponds to current RKI guidelines.

Plus: more than 25 soft laser indications

The pre-programmed indications for endodontics, periodontology and surgery alone reduce the amount of time spent on treatment and improve prognoses. In addition, the elexxion pico offers more than 25 "soft laser" indications (LLLT).

Bleaching

and PDT with perio green®: feel free to delegate!

Bleaching is a service that more and more patients are requesting – and, of course, one for which they pay extra, too. With the elexxion pico, bleaching becomes laser power bleaching. This means lasting results achieved quickly – and happy patients. perio green® avoids the need for antibiotics and can also be used preventatively. perio green® is a fantastic alternative to conventional periodontitis therapy. We will be happy to provide further information. www.periogreen.de



article number	description	price
----------------	-------------	-------

13512 claros pico 2.0



system contains:
 1* base unit
 1* table stand
 1* footswitch
 3* safety goggles
 1* battery charger
 1* starterset fibertips longlife
 1* handpiece Ergoflex plus
 1 bending tool
 8 soft metal tips

1* handpiece "ergo T"
 1* fiber rod T8
 1* primary adapter Europe pico
 1* primary adapter UK pico
 1* warning sign "Laser"

Promotional Price £5215.00

11199 claros nano system



system contains:
 1* base unit
 1* transmission system nano
 1* footswitch
 1* handpiece Ergoflex plus
 1* bending tool for handpiece Ergoflex plus
 8* soft metal tips (2 different sizes)
 1* handpiece "ergo T"
 1* fiber rod T8
 1* starterkit application fibres
 3* safety goggles claros protect (glass holder)
 1* quickcard holder

1* power cord
 1* warning sign laser
 1* set of documents
 1* user manual

Promotional Price £9990.00

13231 claros pico lite



system contains:
 1* base unit
 1* footswitch
 2* safety goggles claros protect (glass holder)
 1* handpiece Ergoflex plus
 8* soft metal tips (2 sizes)
 1* bending tool for handpiece Ergoflex plus
 1* handpiece "ergo T"
 1* fiber rod T8
 1* starterkit elexxion fibers 300 µ longlife
 1* battery charger
 1* primary adapter Europe pico
 1* primary adapter UK pico
 1* warning sign "Laser"
 1* set of documents
 1* user manual

Promotional Price £3990.00

article number	description	price
----------------	-------------	-------

14227 delos 4.0 System with SNORE3 software



system contains:

- | | |
|---|---|
| 1* base unit | 1* handpiece Ergoflex plus |
| 1* system tray "elexxion comfort" | 1* bending tool for handpiece Ergoflex plus |
| 1* footswitch "duros/delos step" | 8* soft metal tips (2 different sizes) |
| 1* handpiece duros 3.0 | 1* handpiece "ergo T" |
| 1* starterkit duros tips 3.0 | 1* fiber rod T8 |
| 3* safety goggles elexxion protect (glass holder) | 1* starterkit application fibres |
| 1* key-card "therapy" | 1* fiber rod SNORE3 |
| 1* power cord | 1* green card SNORE3 |
| 1* key-card "all applications" | |
| 1* warning sign "Laser" | |
| 1* set of documents | |
| 1* user manual | |

This item also includes flights and hotel accommodation on one of our snoring workshops, either in Germany or Dubai.

Promotional Price £77500.00

12259 duros 3.0 System



system contains:

- 1* base unit
- 1* system tray "elexxion comfort"
- 1* footswitch "duros/delos step"
- 1* handpiece duros 3.0
- 1* starterkit duros tips 3.0
- 3* safety goggles elexxion protect (glass holder)
- 1* key-card "all applications"
- 1* power cord
- 1* elexxion lite (warning light laser)
- 1* warning sign "Laser"
- 1* set of documents
- 1* user manual

Promotional Price £66500.00

14225 claros system EN 50W silver with SNORE3 software



system contains:

- | | |
|---|-------------------------|
| 1* base unit | 1* key-card "therapy" |
| 1* system tray "elexxion comfort" | 1* power cord |
| 1* footswitch "claros step" | 1* warning sign "Laser" |
| 3* safety goggles claros protect (glass holder) | 1* set of documents |
| 1* handpiece Ergoflex plus | 1* user manual |
| 1* bending tool for handpiece Ergoflex plus | 1* fiber rod SNORE3 |
| 8* soft metal tips (2 different sizes) | 1* green card SNORE3 |
| 1* handpiece "ergo T" | |
| 1* fiber rod T8 | |
| 1* starterkit application fibres | |
| 1* key-card "all applications" | |

This item also includes flights and hotel accommodation on one of our snoring workshops, either in Germany or Dubai.

Promotional Price £25500.00

article number	description	price
----------------	-------------	-------

13144 perio green (pack of 2 pcs) Starter Kit



photodynamic therapy

- 1 Perio Green pill, 0,2mg
- 1 bottle sterile water 2ml
- 1 blending bottle 10 ml (sterile)
- 1 disposable syringe 3ml (sterile)
- 1 aid tip 0,9 mm (sterile)
- 1 application tip 0,3 mm (sterile)

Promotional Price £105.95

13929	perio green tablets (10 pcs) refill pack	£105.95
13924	perio green sterile water 2ml refill pack	£19.50
13925	perio green blending bottle 10 ml sterile refill pack	£19.50
13926	perio green syringe 3ml sterile refill pack	£19.50
13927	perio green cannula 2er pink and green refill pack	£19.50

13499 odobleach



elexxion Laser-Powerbleaching

- 1 double syringe odobleach® + 3 mixing tips
- 1 Gingiva Protector + 3 application tips
- 1 After-Bleaching-Care

Promotional Price £59.95

14082 green card SNORE3



240 minutes of patient treatment e.g. 12x 20 minutes


Promotional Price £1350.00

P = pico P2= pico 2 PL= pico lite
N = nano C= claros D= delos/duros

to user for	article number	discription	price
-------------	----------------	-------------	-------

consumable

D	12639	starterkit saphir tips 3.0	£486.00
D	12825	saphir tip 3.0 400µm pack of 4 pcs	£365.00
D	12826	saphir tip 3.0 800µm pack of 4 pcs	£365.00
D	12827	saphir tip 3.0 1200µm pack of 4 pcs	£365.00
D	12828	saphir tip 3.0 paro 10mm pack of 4 pcs	£365.00
D	12829	saphir tip 3.0 paro 15mm pack of 4 pcs	£365.00
D	12830	saphir tip 3.0 nonkontakt pack of 4 pcs	£365.00



P	14232	Starter kit Duros Saphir tip 3.0	£299.00
P	14273	1200µ duros saphir tip 3.0-DP pack of 4	£299.00
P	14274	800µ duros saphir tip 3.0-DP pack of 4	£299.00
P	14275	400µ duros saphir tip 3.0-DP pack of 4	£299.00
P	14276	window duros saphir tip 3.0 Dp pack of 4	£235.00

P2, N, C, D	11056	starterkit fibres claros box of 9: 2x 200µ, 2x 300µ, 3x 400µ, 2x 600µ	£279.00
-------------	-------	--	---------



P2, N, C, D	10063	elexxion 200µ fiber longlife pack of 10 pcs	£299.00
PL, P2, N, C, D,	10391	elexxion 300µ fiber longlife pack of 10 pcs	£299.00
P2, N, C, D	10112	elexxion 400µ fiber longlife pack of 10 pcs	£299.00
P2, N, C, D	10120	elexxion 600µ fiber longlife pack of 10 pcs	£347.00

PL, P2, N, C, D,	11256	soft metal tip 18 pack of 20 pcs (für 400µ/600µ Fasern)	£15.70
PL, P2, N, C, D,	11257	soft metal tip 20 pack of 20 pcs (für 200µ/300µ Fasern)	£15.70



P	11855	starterset fibertips pico box of 7 : 1x S200µ, 1x S300µ, 1x S400µ, 1x S600µ, 1x E200µ, 1x P300µ, 1x P400µ	£369.00
---	-------	--	---------



P	11645	pico surgery 200µ pack of 5 pcs	£295.00
P	11646	pico surgery 300µ pack of 5 pcs	£295.00
P	11647	pico surgery 400µ pack of 5 pcs	£295.00
P	11648	pico surgery 600µ pack of 5 pcs	£295.00
P	11707	pico endo 200µ pack of 5 pcs	£295.00
P	11708	pico paro 300µ pack of 5 pcs	£295.00
P	11709	pico paro 400µ pack of 5 pcs	£295.00