

Build-It[®] Core Forms

Physical Properties:

Composition: Polyether carbonate urethane

Color: Transparent

Sizes:

- | | |
|----|--|
| #1 | Anterior Small (7.5mm H x 6.5mm L x 4mm W) |
| #2 | Anterior medium (7.7mm H x 7.5mm L x 4.5mm W) |
| #3 | Premolar (7mm H x 7.5mm L x 5mm W) |
| #4 | Large Canine (8mm H x *.5mm L x 6mm W) |
| #5 | Small/medium molar (6.5mm H x 8mm L x 6.5mm W) |
| #6 | Large molar (7mm H x 10.5mm L x 8.5mm W) |