



Directa Crowns Selecta Mould Guide



Opaque / Opak / Opaca / Opache
Translucent / Translucide / Translucida / Translucenti

Centrals/Centraler/Centrales/ÖK-Front/Centrales/Centrali

mm	↓	14.0	13.0	12.4	12.0	11.5	11.0	10.6	10.6	11.0	11.5	12.0	12.4	13.0	14.0
mm	↔	10.0	9.2	8.7	8.4	8.2	7.8	7.5	7.5	7.8	8.2	8.4	8.7	9.2	10.0



Upper bicuspids/ÖK-premolarer/Prémolaires supérieures/ ÖK Prämolaren/Premolares superiores/Bicuspidi superiori

mm	↓	9.8	9.4	8.9	8.6	8.4	7.5	8.1	8.5	9.0	9.5
mm	↔	7.4	6.9	6.7	6.3	6.0	6.0	6.3	6.7	7.2	7.4



Laterals/Lateraler/Latérales/OK-Front/Lateralest/Laterali

mm	↓	11.0	10.4	10.1	9.5	8.8	8.2	8.2	8.8	9.5	10.1	10.4	11.0
mm	↔	7.3	6.8	6.5	6.2	5.9	5.5	5.5	5.9	6.2	6.5	6.8	7.3



Lower bicuspids/UK-premolarer/Prémolaires inférieures/ UK Prämolaren/Premolares inferiores/Bicuspidi inferiori

mm	↓	9.8	9.5	9.2	8.8	8.2	7.1	7.3	8.0	8.1	8.4
mm	↔	7.4	7.2	6.8	6.6	6.3	6.4	6.8	7.1	7.5	8.0



Cuspids/Hörntänder/Canines/Eckzähne/Caninos/Cuspidi

mm	↓	13.0	12.4	11.6	11.4	11.1	10.8	10.5	10.5	10.8	11.1	11.4	11.6	12.4	13.0
mm	↔	9.0	8.5	8.1	7.9	7.7	7.5	7.3	7.3	7.5	7.7	7.9	8.1	8.5	9.0



Lower incisors/UK-incisiver/Incisives inférieures/ UK Front/Incisivos inferiores/Incisivi inferiori

mm	↓	12.0	11.5	11.0	10.5	10.0	8.0	8.5	9.1	9.6	10.3
mm	↔	6.2	5.6	5.3	5.1	4.9	4.5	4.9	5.2	5.5	5.8

