

www.wh.com



Reliable, durable,
cost-effective.

ALEGRA HS turbines
from W&H




DENTAL INSTRUMENTATION

ALEGRA HS

Reliable, durable, cost-effective.

ALEGRA HS instruments from W&H

The ALEGRA class with its sensational price-performance ratio has been developed using W&H know-how and the precision of Austrian manufacturing.

- High performance
- Low noise level
- Robust design
- Hygiene head system
- Triple spray (TE-98)
- Mini-turbine head (TE-96)
- Press-button chuck system
-  Can be sterilized up to 135 °C
- One-year warranty
- Manufactured in Austria/Europe



"ALEGRA TE-98 with triple spray"

Particularly fine spray provides for maximum cooling of the treatment site.

Hygiene head system

More security for dentist and patient and longer life span of instrument.

Mini-turbine head

The TE-96 BC/RM also offers a mini-turbine head with press-button chuck system.



Connections


All models are optionally equipped with fixed connections for 2(3) Borden hole systems or standard 4-hole systems.



ALEGRA

Design

Seamless hygienic design for the best grip - especially with gloves.



ALEGRA^{HS}



Technical data	TE-95 BC/RM	TE-96 BC/RM	TE-98 BC/RM
Power (W)	12	11	12
Free-running speed (rpm)	330,000	390,000	330,000
Operating pressure (bar)	2.2 - 2.5	2.2 - 2.5	2.2 - 2.5
Spray water quantity (ml/min)	>50	>50	>50
Air consumption (l/min)	42	42	42
Spray	single	single	triple
Head diameter (mm)	12.2	10.4	12.2
FG bur EN ISO 1797 (Ø/L mm)	1.6/25	1.6/21	1.6/25

Hygienic maintenance



F1 Service Oil, MD-400
 Spray can (400 ml) for quick and easy lubrication and for cleaning the internal lumens of instruments and air motors.



Assistina 301 plus
 Cleaning and care equipment for instruments, turbines and air motors. Fast, easy and economical.
 Only compressed air connection required.



DENTALWERK



ISO 13485
 Annex II
 93/42/EWG
 EC 0297

W&H Dentalwerk Bürmoos GmbH
Austria, A-5111 Bürmoos

Postfach 1, Ignaz-Glaser-Straße 53
 Telefon: +43/62 74/62 36-0 E-Mail: office@wnh.co.at
 Telefax: +43/62 74/62 36-55 www.wh.com