

Stay-Soft S.R.M.[™] Soft Reline Material

Physical Properties

StaySoft S.R.M. Liner

End Hardness	~ 30 Shore A
Tensile Strength	> 3.0 N/mm ² (squared)
Strain in Compression	6.5 - 7.5%
Working Time	1 - 1:20 Min.
Setting Time	3:30 Min.

StaySoft S.R.M. Varnish

End Hardness	~ 25 Shore A
Tensile Strength	>0.5 N/mm ² (squared)
Strain in Compression	10-13%
Working Time	1 - 1:20 Min.
Setting Time	3:30 Min.