

For all your  
**Aligner Solutions...**

Essentially **Trycare**





Dentalastics Separators

Blue for Posteriors. For quick and easy separation of teeth. Radiopaque. Inner diameter 2.1mm.

Blue	Pack 1000	Code: UWDS000
------	-----------	---------------



Separating Plier

Ideal for the placement of elastic separators.

Each	Code: UHFPOSP
------	---------------



Elastic Separators With Stops

Safe-Seps Elastomeric Separators, Regular Blue.

Blue	Pack 1000	Code: UPRO009
------	-----------	---------------



IPR Gauge

Includes 6 Stainless Steel gauges of 0.1 to 0.5mm.

Code: UFSIPRG
---------------



ApproxOpener

Designed for opening of contact points with a serrated strip edge and 8µm diamond coating for initial IPR. Safe grip thumb & finger holder.

Thickness			
Double-Sided	Pack 3	0.09mm	Code: UITV001
Single-Sided	Pack 3	0.10mm	Code: UITV002



ProxoContour

Intensiv ProxoContour metal strip with single sided diamond-coated grit sizes: 80µm (Blue) for contouring and 40µm (Red) for finishing.

Pack 6	Code: UITV046
--------	---------------



ProxoPolish

Single sided diamond coated metal strip with 2 grit sizes 15µm (yellow) for polishing and 8µm (orange) for opening.

Pack 6	Code: UITV047
--------	---------------



Dr. Francesco Garino, Turin, Italy




IPR Distance Control Set

- Stainless Steel and autoclavable.
- A complete and practical IPR Gauge.
- Safe grip thumb & finger holder.
- 8 sizes correlate to the Intensiv IPR range.
- 8 Sizes: 0.10mm to 1.0mm.

Code: UITV006
---------------



Separating Strips



Width	Single Sided	Thickness	Code
3mm	F303	0.07mm	UPRO004
4mm	C304	0.09mm	UPRO005
6mm	C306	0.09mm	UPRO006
6mm	416 Serrated	0.05mm	UPRO007




Corundum-bound, electro-deposited Stahlcarbo® strips. Can be cut to desired length. Packet of 12 strips.





Diamond Finishing Strips - Perforated - Single Sided

Gain space by gently contouring the sides of the overlapping teeth. Made of high grade Stainless Steel to resist breaking and stretching. Autoclavable. Available in a variety of widths and grain.



Perforated Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0,08mm	UEDFS18
Fine	3.75mm	30µm	0,10mm	UEDFS17
Medium	3.75mm	45µm	0,13mm	UEDFS16


Perforated Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	2.5mm	15µm	0,08mm	UEDFS12
Fine	2.5mm	30µm	0,10mm	UEDFS13
Medium	2.5mm	45µm	0,13mm	UEDFS11



Diamond Finishing Strips - Perforated Serrated - Single Sided

A serrated saw edge and diamond coating combined in one strip, allowing for simpler working procedures.



Perforated Pack 10


Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0,08mm	UEDFS09
Fine	3.75mm	30µm	0,10mm	UEDFS08
Medium	3.75mm	45µm	0,13mm	UEDFS07

Perforated Serrated Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	2.5mm	15µm	0,08mm	UEDFS15
Fine	2.5mm	30µm	0,10mm	UEDFS14
Medium	2.5mm	45µm	0,13mm	UEDFS10



Diamond Finishing Strips - Plain Serrated - Single Sided



Perforated Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0,08mm	EDNT0265


Perforated Serrated Pack 10

Colour	Width	Grain	Thickness	Code
Fine	3.75mm	30µm	0,10mm	UEDFS04



Diamond Finishing Strips - Plain - Single Sided

Diamond Interproximal Strips ideal for enamel reduction as recommended for aligner treatment. High grade Stainless Steel to resist breaking and stretching. Autoclavable.



Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0,08mm	UEDFS03
Fine	3.75mm	30µm	0,10mm	UEDFS05
Medium	3.75mm	45µm	0,13mm	UEDFS06

Pack 10

Colour	Width	Grain	Thickness	Code
Extra Fine	2.5mm	15µm	0,08mm	UEDFS02
Fine	2.5mm	30µm	0,10mm	UEDFS01
Medium	2.5mm	45µm	0,13mm	UEDFS00



Diamond Finishing Strips - Plain Wide - Single Sided




Pack 5

Colour	Width	Grain	Thickness	Code
Extra Fine	6mm	15µm	0,08mm	UEDCDS6
Fine	6mm	30µm	0,10mm	UEDFDS6
Medium	6mm	45µm	0,13mm	UEDFS19

Essix

Flex-View™  
Single Sided

These perforated Stainless Steel strips allow optimal visibility and minimise clogging. Flex-View™ Interproximal Strips also feature a gateway area for easy access for interproximal stripping and finishing. All strips are single sided and autoclavable.



2.5 mm Narrow Strips - Pack 10

X-Fine Grit Yellow	Code: UORC095
Fine Grit Red	Code: UORC093
Med Grit Blue	Code: UORC094

3.75 mm Wide Strips - Pack 10

X-Fine Grit Yellow	Code: UORC099
Fine Grit Red	Code: UORC097
Med Grit Blue	Code: UORC098

Pack 15

5 x Wide Extra Fine, 5 x Fine, 5 x Medium	
Assorted	<b>Code:</b> UORC092

Pack 15

5 x Wide Extra Fine, 5 x Fine, 5 x Medium	
Assorted	<b>Code:</b> UORC096

Essix  
Diamond Interproximal Strips

Fully autoclavable, surgical steel plain strips. These very thin strips (0.08mm/0.003") are made of electrolytically bonded diamonds on a vibro cushion base.



Single Sided Pack 10				
Colour	Width	Grain	Thickness	Code
Fine	2.5mm	30µm	0.08mm	UORC107
Fine	3.75mm	30µm	0.08mm	UORC110
Medium	2.5mm	45µm	0.08mm	UORC106
Medium	3.75mm	45µm	0.08mm	UORC109

Double Sided Pack 10				
Colour	Width	Grain	Thickness	Code
Medium	2.5mm	45µm	0.08mm	UORC111
Medium	3.75mm	45µm	0.08mm	UORC108

Essix  
Interproximal Serrated Strips 6mm

Stainless steel ultra thin interproximal serrated strips, used prior to reduction with diamond strips. Width: 6mm, Thickness: 0.05mm.

Pack 10	Code: UORC105
---------	---------------

Komet  
Serrated Saw Diamond Strips

Diamond Strips. 0.05mm Serrated Separating Saw. Width: 6mm, Thickness: 0.05mm.

Pack 10	Code: TD0703
---------	--------------

Komet  
Vision Flex Diamond Strips

These Strips have a more dead-soft quality than the typical springy feel of conventional diamond strips.

Pack 10 - Length 147mm				
Colour	Width	Grain	Thickness	Code
Extra Fine	3.75mm	15µm	0.08mm	UVFW37E
Fine	3.75mm	30µm	0.10mm	UVFW37F
Medium	3.75mm	45µm	0.13mm	UVFWS37

Assorted Pack 15	
Contains: 5 x medium (45µm), 5 x fine (30µm) and 5 x extra-fine (15µm)	
Assorted 3.75mm	Code: UVFW37A

Pack 10 - Length 147mm				
Colour	Width	Grain	Thickness	Code
Extra Fine	2.5mm	15µm	0.08mm	UVFW25E
Fine	2.5mm	30µm	0.10mm	UVFW25F
Medium	2.5mm	45µm	0.13mm	UVFWS25

Assorted Pack 15	
Contains: 5 x medium (45µm), 5 x fine (30µm) and 5 x extra-fine (15µm)	
Assorted 2.5mm	Code: UVFW25A



Solid Diamond Strips



Pack 10 - Length 148mm				
Colour	Width	Grain	Thickness	Code
Extra Fine	3.7mm	15µm	0.08mm	UVFD37E
Fine	3.7mm	30µm	0.10mm	UVFD37F
Medium	3.7mm	45µm	0.13mm	UVFDS37
Medium	3.7mm	45µm	0.21mm	UVFD37D
Coarse	3.7mm	90µm	0.18mm	UVFD37C

Pack 10 - Length 148mm				
Colour	Width	Grain	Thickness	Code
Extra Fine	2.5mm	15µm	0.08mm	UVFD25E
Fine	2.5mm	30µm	0.10mm	UVFD25F
Medium	2.5mm	45µm	0.13mm	UVFDS25



Code: UEDS04K  
PDL-4 KIT 4mm



Code: UEDS25K  
PDL-2.5 KIT 2.5mm

EVE Diastrips

Diastrips (Porcelain, Composite, Metal, Enamel).

- For manual surface processing
- Interdental application
- Perforation removes the cut material.

PACK 9

						Dimensions (mm)
0.11	0.11	0.10	0.10	0.07	0.07	Thickness (mm)
UEDS454	UEDS452	UEDS354	UEDS352	UEDS104	UEDS102	Code:

Kits Contain: 3 x of each grit (fine, medium and coarse) in the selected size of 2.5mm or 4mm widths.







**Garrison FitStrip Starter Kit** (FPSK01)

- Starter Kit Includes:**

  - 2 x Serrated Strip (0.05mm)

**Single Sided**

  - 1 x Super Fine Grit Strip (0.08mm)
  - 1 x Fine Grit Strip (0.10mm)
  - 1 x Medium Grit Strip (0.13mm)
  - 1 x Coarse Grit Strip (0.18mm)
- Double Sided**

  - 1 x Super Fine Grit Strip (0.11mm)
  - 1 x Fine Grit Strip (0.15mm)
  - 1 x Medium Grit Strip (0.21mm)
  - 1 x Coarse Grit Strip (0.30mm)
  - 2 x FitStrip Handle

Starter Kit    Code: UFPSK01



**Garrison FitStrip**

**FitStrips are your complete line of interproximal trimmers and finishers!**

Whether you're looking for safe and accurate IPR for your aligner cases, interproximal cement clean-up, or smoothing out cervical margins, FitStrip™ is a fast and convenient solution.

FitStrip™ is indicated for composite restorations (contouring/finishing all interproximal surfaces), orthodontics (IPR/tooth slenderizing), and crown and bridge (cement removal and clean-up).



**Garrison FitStrip Serrated**


FitStrips' simple, detachable handle makes a world of difference for both clinician ease of use and patient comfort.

Your hand stays outside the mouth, granting a clear field of view and easier operation. Twist the colour-coded barrel to adjust FitStrips' curvature and it's automatically locked in place. It doesn't get any easier than that. FitStrips are very durable and long lasting.

Serrated Grey 0.05mm    Pack of 4    Code: UFPSSAW



**Garrison FitStrip**



Double Sided	Pack Size	Code:	Single Sided	Pack Size	Code:
Superfine Yellow 0.11mm	Pack of 4	UFDSF11	Superfine Yellow 0.08mm	Pack of 4	UFPSFS4
Fine Red 0.15mm	Pack of 4	UFDF015	Fine Red 0.10mm	Pack of 4	UFPFSR4
Medium Blue 0.21mm	Pack of 4	UFPSMDS	Medium Blue 0.13mm	Pack of 4	UFPSMSS
Coarse Green 0.30mm	Pack of 4	UFDC030	Coarse Green 0.18mm	Pack of 4	UFPSCS4
Coarse White 0.40mm	Pack of 4	UFDMG04			
X-Coarse Black 0.50mm	Pack of 4	UFDMB24			



**Garrison FitStrip Universal Kit** (FPSK04)

- Universal Kit Includes:**

  - 1 x IPR gauge & 2 x FitStrip™ Handle.

**Single Sided**

  - 1 x Super Fine Grit Strip (0.08mm)
  - 1 x Fine Grit Strip (0.10mm)
  - 1 x Medium Grit Strip (0.13mm)
  - 1 x Coarse Grit Strip (0.18mm)
  - 1 x Super Fine Grit Strip (0.11mm)
- Double Sided**

  - 1 x Fine Grit Strip (0.15mm)
  - 1 x Medium Grit Strip (0.21mm)
  - 1 x Coarse Grit Strip (0.30mm)
  - 1 x Coarse Grit Strip (0.40mm)
  - 1 x X-Coarse Grit Strip (0.50mm)

Universal Kit    Code: UFSUNKT



**Garrison FitStrip IPR Kit** (FPSK05)

- IPR Kit Includes:**

  - 1 x IPR Gauge
  - 2 x FitStrip Handle

**Double Sided**

  - 2 x Super Fine Grit Strip (0.11 mm)
- 2 x Medium Grit Strip (0.21mm)
  - 2 x Coarse Grit Strip (0.30mm)
  - 2 x Coarse Grit Strip (0.40mm)
  - 2 x X-Coarse Grit Strip (0.50mm)

Universal Kit    Code: UFSIPRK



**QwikStrip Ultimate Assorted Kit**

- Industry Leading Qwik Strips.
- Easy to follow IPR sequence
  - Interproximal cement removal
  - Use for any aligner regime
  - Double or Single Sided
  - Autoclavable



Scan the QR Code:  
to watch a QwikStrip  
IPR Demonstration.

- Kit Contains:
- 5 x Single-sided (.07mm - .18mm)
  - 5 x Double-sided (.18mm - .40mm)

Code: UQWKS01

Ideal for breaking contacts and for interproximal cement removal.



**QwikStrip Single-Sided**

.07mm Yellow	Pack 10	Code: UQWKS03
.10mm Red	Pack 10	Code: UQWKS04
.13mm Blue	Pack 10	Code: UQWKS05
.16mm Green	Pack 10	Code: UQWKS06
.18mm Black	Pack 10	Code: UQWKS07



**QwikStrip Curved**

.07mm Yellow	Pack 10	Code: UQWKS08
.10mm Red	Pack 10	Code: UQWKS09
.13mm Blue	Pack 10	Code: UQWKS10
.16mm Green	Pack 10	Code: UQWKS11
.18mm Black	Pack 10	Code: UQWKS12




**QwikStrip Double-Sided**

.18mm Yellow	Pack 10	Code: UQWKS13
.20mm Red	Pack 10	Code: UQWKS14
.25mm Blue	Pack 10	Code: UQWKS15
.30mm Green	Pack 10	Code: UQWKS16
.40mm Black	Pack 10	Code: UQWKS17



**QwikStrip Serrated**

.05mm White	Pack 10	Code: UQWKS02
-------------	---------	---------------



**StripFix Finishing Strip System**

StripFIX finishing strip system Manual Diamond Strip Holders. Use with or without the handle. Hold between middle finger with your index finger on the diamond strip plate. Strip holder is swivel mounted. May be locked in 4 positions for optimal access to the oral cavity. Simple quick-plug system to hold the strip. All components autoclavable.

- Pack contains:
- 1 x .04mm Serrated Edge
  - 1 x .10mm Fine
  - 1 x .08mm Super Fine
  - 1 x .13mm medium

Pack of 4 Files and Handle      Code: USFX000

**Refills - Single-Sided**

White Serrated Edge .04mm	Pack 10	Code: USFX001
Yellow Super Fine .08mm	Pack 10	Code: USFX002
Red Fine .10mm	Pack 10	Code: USFX003
Blue Medium .13mm	Pack 10	Code: USFX004
2 x Handles .13mm	Pack 2	Code: USFX005







IPR Set Professional

Calibrated and precise stripping (IPR) for orthodontics. Mechanical and manual strips in a stainless steel tray. The professional kit contains both manual and mechanical IPR options all in an autoclavable cassette.

REF WG-69A * 1 Piece	Intensiv Swingle 4:1, reciprocating contra-angle (without light)
REF 053 1 Piece	Intensiv Ejector for safe file removal between IPR stages
REF OS08OP-DS 1 Piece	Intensiv Ortho-Strips, 8µm opener, serrated strip edge, for opening of contact points
REF OS40M-DS 1 Piece	Intensiv Ortho-Strips, 40µm medium, for step one creating additional space
REF OS60C-DS 1 Piece	Intensiv Ortho-Strips, 60µm coarse, for step two in creating additional space
REF OS80XC-DS 1 Piece	Intensiv Ortho-Strips, 80µm extra coarse, perforated, for high enamel reduction
REF OS40M-DS 1 Piece	Intensiv Ortho-Strips, 40µm medium, for pre-polishing of the enamel surface
REF OS15POL-DS 1 Piece	Intensiv Ortho-Strips, 15µm polishing, for fine polishing of treated surfaces

All in a single autoclavable box. Code: UITV005



Swingle Kit With Light & Strips



Swingle Kit Without Light & Strips




Swingle Kit

W&H Synea 4:1 Reciprocating Handpiece with a 6 file set.


Contains:

- Contra Angle Ref WG-69 LT
- Ejector Ref 053
- Spray Canal Cleaner Needle Ref 054
- 6 Intensiv Ortho Strips.

H/P With Light Plus Strips	Code: UITV012
H/P Without Light Plus Strips	Code: UITV011



- OS08OP-DS, Opener, 8µm
- OS80XC-DS, Extra coarse, 80µm
- OS60C-CEN-DS, Coarse, 60µm









- OS40M-CEN-DS, Medium, 40µm
- OS25F-CEN-DS, Fine, 25µm
- OS15POL-CEN-DS, Polishing, 15µm







Colour Code:	Micron	Description
494	8µm	Opener
504	15µm	Polishing
514	25µm	Fine
514	40µm	Medium
524	60µm	Coarse
524	80µm	Extra coarse









Ortho Strips

Fast, efficient and precise, double and single sided. Metal frames for reduction and polishing. Pack 3

Double-Sided diamond						
Code:						
	UITV028	UITV033	UITV034	UITV037	UITV040	UITV045
	Opening	Polishing	Finishing	Contouring	Reduction	Space creation

One-Sided diamond						
	Left	Right	Left	Right	Left	Right
One-Sided diamond						
	Code:	UITV029	UITV030	UITV031	UITV032	UITV035
	Opening	Opening	Polishing	Polishing	Finishing	Finishing

One-Sided diamond							
	Left	Right	Left	Right	Left	Right	
One-Sided diamond							
	Code:	UITV038	UITV039	UITV041	UITV042	UITV043	UITV044
		Contouring	Contouring	Reduction	Reduction	Space creation	Space creation



Ortho Strips

Contains an assortment of 6 oscillating files:

- Opening 08µm
- Space Creating 80µm
- Polishing 15µm
- Reduction 60µm
- Finishing 25µm
- Contouring 40µm

Double-Sided Set	Each	Code: UITV007
Single-Sided Set - Left	Each	Code: UITV008
Single-Sided Set - Right	Each	Code: UITV009





New Central System



Intensiv OST400Set01/6

Thanks to the diamond coating in the central operating area, the Intensiv Ortho-Strips System Central allows for high precision in approximal reduction, finishing and polishing of enamel, respecting the tooth morphology and dental tissue.

During the process of enamel reduction, trimming and polishing in orthodontics, parts of enamel or dentine below the tooth margin may be removed unintentionally, which may cause subsequent hypersensitivity of the tooth in the patient.

Colour Code:	Micron	Description
494	8µm	Opener
504	15µm	Polishing
514	25µm	Fine
514	40µm	Medium
524	60µm	Coarse

Contains:  
• 1 x Ortho-Strips Tray  
• 1 x Ortho-Strips System Opener.

Ortho-Strips system central double sided:  
• 2 x Contouring  
• 1 x Reduction  
• 1 x Finishing  
• 1 x Polishing.

Kit

Code: UITV010

Double-Sided diamond

Opening

Contouring

Reduction

Finishing

Polishing



DiaKleen

Pack of 3.

Code: UITV003



Ejector for Swingle

Intensiv File Ejector

Code: UITV004



Ortho-Strip Central

Pack of 3.

Colour Code:	Micron	Description
524	60µm	Coarse
514	40µm	Medium
514	25µm	Fine
504	15µm	Polishing

Double-Sided diamond

Reduction

Contouring

Finishing

Polishing

Code: UITV025   UITV022   UITV019   UITV016

Single-Sided diamond

Reduction

Contouring

Finishing

Polishing

Left Code: UITV026   UITV023   UITV020   UITV017

Right Code: UITV027   UITV024   UITV021   UITV018



Pop On Diamond Finishing Discs Double Sided

Double sided 0.10mm thick. Diamond Coated. Each.

10mm	Each	Code: UKMT001
14mm	Each	Code: UKMT002
18mm	Each	Code: UKMT003
Mandrel	Each	Code: UKMT004

Komet DISCtance IPR system kit

Rotary diamond discs for interproximal enamel reduction.


The Komet DISCtance IPR system features a range of new, sharp stripping discs that ensure the amount of enamel is easily, efficiently and precisely removed. All discs are mounted on a straight handpiece shank. Fully autoclavable bur holder.

Code: UKMT005



Comes with 5 new discs:

Ref:	Grit	Thickness	Coated on
RS15FV	Fine	0.15mm	Upper side
RS15FH	Fine	0.15mm	Bottom side
RS15F	Fine	0.15mm	Both sides
RS20M	Medium	0.20mm	Both sides
RS30M	Medium	0.30mm	Both sides



Superflex Diamond Open Mesh Disc

Superflex Diamond Open Meshed Disc Double Sided 22mm Diameter, 0.15mm Thickness.

22mm Diameter	Each	Code: USFDO22
---------------	------	---------------



IPR Diamond Bur

IPR Diamond Bur. Pack 5.

392 Series

Super Fine	Code: UEI392S
Fine	Code: UEI392F
Medium	Code: UEI392M



Komet

Orthodontic IPR Kit

The Komet IPR & Disc kit contains:

- Vision Flex Strips
- 3x WS37 .13mm 45µm Blue Medium
- 3x WS37F .10 mm 30µm Red Fine
- 3x WS37EF .08 mm 15µm Yellow Extra Fine
- 3x DS37D .21mm 45µm Blue Medium Double Sided
- 3x DS37C .18mm 90µm Green Coarse Single Sided
- 1x 934.104.180 Diamond Disc .18mm 18mm Diameter
- 1x 8934A.900.180 .15mm 14mm Diameter
- 1x Contra Angle Mandrel
- 1x De-Bonding Bur
- 1x H379AGK Egg De-Bonding Bur
- Non-autoclavable storage box

Code: UKMT006





# Precise Adjustments For Better Results

Clear Collection

Clear Collection



**Hu-Friedy**

## Tear Drop

Creates a notch at the gingival margin of clear aligners anywhere elastic hooks are required within the arch. Provides a reservoir to hold the elastics in place as the aligner is seated, permitting easier manipulation of elastics by the patient.

**Code:** UHFAPTD

**Ref:** 678-800



**Hu-Friedy**

## Hole Punch

Creates half-moon cut-outs for relief of clear aligners around bonded buttons used for elastics. Nips away plastic to provide soft tissue clearance to prevent aligner impingement of tissue.

**Code:** UHFAPPN

**Ref:** 678-803



**Hu-Friedy**

## Vertical

Emphasizes over-correction of rotations by producing an an indentation on the facial or lingual aspect of a specific tooth.

**Code:** UHFAPVT

**Ref:** 678-801



**Hu-Friedy**

## Petite Punch

The Petite Punch cuts a 4.5mm half-moon cut-out in an aligner for relief of the plastic around bonded buttons used for elastics.

**Code:** UHFAPP

**Ref:** 678-805



**Hu-Friedy**

## Horizontal

Enhances retention of either clear aligners or retainers by placing an indent into a tooth's undercut or below a bonded attachment.

**Code:** UHFAPHZ

**Ref:** 678-802



**Hu-Friedy**


## Spot

The Spot creates a tiny rounded dimple for spot contacts to accent small individual tooth movements or added retention.

**Code:** UHFAPSP

**Ref:** 678-804





**Hu-Friedy**

### Counter Clockwise Wedge

The Counter Clockwise Wedge is used to create a counter clockwise rotation on a single tooth within an aligner.

**Code:** UHFAPCC

**Ref:** 678-807




**Hu-Friedy**

### Clockwise Wedge

The Clockwise Wedge is used to create a clockwise rotation on a single tooth within an aligner.

**Code:** UHFAPCW

**Ref:** 678-808



**Hu-Friedy**

### Eraser

The Eraser is designed to reverse or remove unwanted indents placed by other instruments in the Clear Collection.

**Code:** UHFAPER

**Ref:** 678-811

**Essix**

### Aligner Chewies

Same size as standard Dental rolls - but better. Made of soft plastic, they are hygienic and easy to use. Two seaters per patient pack.

**10 Packs of 2 chewies.**

Bubblegum **Code:** UACMB00

Grape **Code:** UACMGR0

Mint Green **Code:** UACMG00

Pineapple **Code:** UACMPA0

Plain **Code:** UACMP00



**ReliaBites in classic roll and Y-shaped design.**



**Introducing ReliaBites® Aligner Chewies**



ReliaBites improve clear aligner treatment to help patients properly seat aligners, remove air bubbles, and reduce overall treatment time. Available in classic cylindrical rolls or Y-shaped design.

- Food grade elastomer material
- Applies uniform force to properly seat aligners
- Patient controlled forces – bite and seat
- Eliminates “halo” air bubbles in sockets
- Y-shaped available Soft or Firm, in various flavours.



**CLASSIC ROLL DESIGN**  
10 Pack's of 2



**Y-SHAPED DESIGN**  
Pack of 10

ReliaBites®

Flavour	Rolls Code:	Y-Shaped Code:	Firmness
Apple	-	UREL008	Soft
Grape	UREL001	UREL005	Firm
Mint	UREL002	UREL006	Firm
Orange	UREL003	UREL009	Soft
Peach	-	UREL007	Firm
Strawberry	UREL004	UREL010	Soft

\* Firmness is only applicable for Y-Shaped ReliaBites.





**Essix  
Retainer Brite**

Cleans and brightens retainers, night guards, mouth guards, splints and most other removable dental appliances.

- Helps remove plaque & tartar
- Kills germs & bacteria
- Keeps appliance minty-fresh.

Retainer Brite 36 - 1 month supply	<b>Code:</b> URBT000
Retainer Brite 92 - 3 month supply	<b>Code:</b> URBT001



**Retainer Boxes**

Assorted colours.

Small	Pack 24	<b>Code:</b> UPRO002
* Small Colours: Red, Green, Pistachio, Blue, Fuchsia, Orange. 29mmx76mmx72mm		
Glitter Retainer Boxes	Pack 24	<b>Code:</b> UPRO003
* Glitter Colours: Purple, Green, Yellow, Blue, Pink, Orange. 29mmx76mmx72mm		
Deep Pack	Pack 24	<b>Code:</b> UPRO001
* Deep Colours: Red, Green, Pistachio, Blue, Fuchsia, Orange. 42mmx82mmx72mm		

**ReliaPry**

**Introducing ReliaPry  
Clear Aligner Removal Tool**

Designed to offer an effortless solution to the daily task of aligner removal. With ReliaPry, your patients will have a reliable and gentle tool that allows removal of clear aligners easily and without stress.

Resealable and tamper-proof patient packs.

- Individually wrapped patient packs
- Each patient pack contains 1 ReliaPry tool
- Contains 100 patient packs - Assorted Colours.

Pack 100	<b>Code:</b> UREL011
----------	----------------------

**Essix  
Art Appliance  
Remover Tool**



ART makes it easy for patients to remove their clear appliances! Removes clear appliances from the cheek-side and/or from the tongue side. Increases patient comfort and compliance. Latex free.

Pack 25	<b>Code:</b> UORC100
---------	----------------------

**Retainer Box**

Bleaching/Retainer boxes in a variety of colours.

Shallow	Pack 10	<b>Code:</b> UBRBSAC
* 88x73x30mm		
Deep	Pack 10	<b>Code:</b> UBRBDSC
* 80x70x40mm		



**GC  
Tooth Mousse**

To reduce white spot lesions and a topical application after IPR.

Intro Kit	Pack 5	<b>Code:</b> UTM1000
Assorted	Pack 10	<b>Code:</b> UTM0006
Mint	Pack 10	<b>Code:</b> UTM0003
Strawberry	Pack 10	<b>Code:</b> UTM0001



**Fluorigard AF Daily Rinse**

A daily mouth rinse that prevents dental caries. Suitable for children over 6 years of age.

Listed on the British National Formulary (BNF) as Sodium Fluoride mouth wash 0.05% (225ppm F) Mint flavour. Alcohol free.

Blue 400ml	Pack 4	<b>Code:</b> UFAFDRB
Green 400ml	Pack 4	<b>Code:</b> UFGFMRG

**Fresh & Go Mouthwash**

Fresh & Go ® Mouthwash is specifically formulated for orthodontic appliance wearers as it is alcohol free, contains fluoride and has a refreshing mild after taste.

- Alcohol free
- 0.05% Sodium Fluoride 225ppm F
- Made in the UK

500ml	Pack 12	<b>Code:</b> UORC104
-------	---------	----------------------

**Clinisept+  
Dental**

**Outstanding clinical  
results for mucosa.**

**Clinisept**

Clinisept+ Mouthwash is an oral cleanser and hygienic rinse that delivers professional levels of antimicrobial protection against the harmful bacteria that cause gum disease and tooth decay, but without stinging or causing irritation.

400ml	Pack 12	<b>Code:</b> UCSPMW4
600ml	Pack 6	<b>Code:</b> UCSPMW6

- **Effective cleansing and antimicrobial protection** (EN1040:2005, EN1276:2019, EN14476:2013 +A1:2015)
- **Contains no chlorhexidine, alcohol, colouring or flavouring**
- **Non-irritant to teeth or gums**
- **Causes no staining**
- **Neutral pH 7.0**
- **Non-cytotoxic and hypo-allergenic**

Clinisept+ Dental Mouthwash has an 18-month shelf life and should be used within 3 months of opening.





# BioMin®

## Armour for Teeth

Protect your teeth from white spot lesions during Orthodontic treatment with BioMin® F toothpaste



BioMinF Toothpaste 75ml	Pack 8	Code: UBIOMINF
BioMinC Toothpaste 75ml	Pack 8	Code: UBIOMINC



- Helps reduce the risk of white spot lesions during Orthodontic treatment
- Active remineralisation of white spots post treatment
- Strengthens teeth and protects against decay
- Makes teeth feel cleaner and fresher
- Developed from research at Queen Mary's University of London
- First Oral Health Foundation accredited toothpaste for sensitivity relief and remineralisation

[www.biomin.co.uk](http://www.biomin.co.uk)

@BioMinTech /BioMinTechnologiesLimited @biomintechnologieslimited

Trycare Tel: 01274 88 10 44 | Web: [www.trycare.co.uk](http://www.trycare.co.uk) | Email: [Orthodontics@trycare.co.uk](mailto:Orthodontics@trycare.co.uk)

## Home Care Hygiene Range



### Clear Aligner Home Care Kit

This pack contains the essential items needed to help patients clean and care for their aligners as well as maintaining good Oral Hygiene.

Code: UORC090

#### Pack contents:

- Sonic cleaner
- Retainer Brite - aligner appliance cleanser (1 box 36 tablets)
- Interdental brushes (assorted pack)
- Aligner Chewies (help seat aligners correctly)
- Aligner remover
- Fresh & Go Mouthwash
- Aligner case
- Mirror
- All supplied in a convenient, large, heavy duty zip bag.



### Mini Aligner Home Care Kit

This pack contains the essential items needed to help patients clean and care for their aligners as well as maintaining good Oral Hygiene.

#### Pack contents:

- 1 x Aligner Remover
- 2 x Retainer Brite Tablets
- 2 x Aligner Chewies
- 1 x Storage Case

Code: UORC091

Tel: 01274 88 10 44 | Web: [www.trycare.co.uk](http://www.trycare.co.uk) | Email: [Orthodontics@trycare.co.uk](mailto:Orthodontics@trycare.co.uk) Trycare





## EverStick Ortho

everStick ORTHO is indicated for the aesthetic retention phase after active orthodontic treatment. The resin matrix holds the individual glass fibres in a bundle, which facilitates the handling and adaptation to the arch.

2 x 10 **Code:** UESO8CM



### Perfect A Smile - Pontic Paint

Quick and economical. Light Cure, apply smoothly to clear aligner sockets. Translucent to match corresponding enamel shade. High elasticity, bonds to any clear aligner.



4g Jar Bleach	<b>Code:</b> UREL021
4g Jar A1 Shade	<b>Code:</b> UREL012
4g Jar A2 Shade	<b>Code:</b> UREL013
4g Jar A3.5 Dark Shade	<b>Code:</b> UREL014
4g Jar B1 Shade	<b>Code:</b> UREL015
4g Jar C1 Shade	<b>Code:</b> UREL016
4g Jar C4 Shade	<b>Code:</b> UREL017
4g Jar D3 Shade	<b>Code:</b> UREL018
4g Jar Light Shade	<b>Code:</b> UREL019
4g Jar Medium Shade	<b>Code:</b> UREL020




### Ortho Flextech

- Use post aligner treatment.
- Flexible lingual retainer.
- Secure and stable long term passive retention.
- Low profile for patient comfort.
- Stainless Steel offers increased tensile and torsional strength.
- Chairside application.
- No Lab Fee or delay in delivery to the patient.

Popular Choice



Stainless Steel Etched	<b>Code:</b> UREL030
Stainless Steel Non Etched	<b>Code:</b> UREL048
14kt White Gold Etched	<b>Code:</b> UREL031
14kt White Gold Non Etched	<b>Code:</b> UREL032
14kt White Gold Nickel Free	<b>Code:</b> UREL033



Scan the QR Code: to watch the Reliance University clinic procedures.



### Perfect A Smile - Pontic Putty

Finger tip pontic fabrication. Light cure composite paste with putty like consistency. Facilitates tooth shape formation. Eliminates pontics, wax & PVS.

5g Syringe **Code:** UREL022



### Retainium XL

Pack 10 **Code:** UREL050



- Flat .027" x .016" ribbon arch
- Low profile – comfortable and functional
- Increased strength
- MRI compatible
- Titanium Molybdenum Alloy
- Nickel Free.



### Retainium

Pack 10 **Code:** UREL049



- Flexes Interproximally
- Lowest profile
- Flat .027" x .011" ribbon arch
- Adaptability
- MRI compatible
- Titanium Molybdenum Alloy
- Nickel Free.



Scan the QR Code: Video on Maxillary Retention Using Retainium



### Bond-A-Braid Lingual Retainer Wire

Pack 10 **Code:** UREL034



The dead soft wire easily adapts and is excellent for semi-permanent orthodontic splinting of the upper and lower incisors, diastema maintenance and other tooth splinting.

- Lingual retainer
- .027" x .011" ribbon arch
- 8-strand braided wire
- Each tube contains 10 (6") wires





- Kit contains:**
- 12 x Habit Control
  - 12 x Cuspid Ramp
  - 12 x Button Hook
  - 12 x Molar Tube
  - 12 x Bite Ramp 3 & 5mm
  - 12 x Wire Bonder 3 & 5mm
  - 2 x Brush handles.



**Reliamold Starter Kit**

Moulds are translucent and flexible to facilitate composite loading, curing and disengagement.

- A mould is filled with a light cure composite (GoTo®, T.K.O.®, or Flowtain™ recommended) then delivered to the tooth and cured.
- The desired auxiliary is delivered, formed, and functional immediately.

**Code:** UREL023



**Button Attachment**

Moulds are translucent and flexible to facilitate composite loading, curing and disengagement.

25 Moulds per Package **Code:** UREL024



**Reliamold Molar Pad Attachment**

Open the bite with posterior occlusal build-ups. T.K.O.® or Ultra Band-Lok® recommended.

25 Moulds per Package **Code:** UREL025

**Attachments**



**Lingual Button**

Pack 10.

- |                          |                      |
|--------------------------|----------------------|
| Flat Lingual             | <b>Code:</b> UPRO012 |
| Curved Lingual           | <b>Code:</b> UPRO010 |
| Flat Lingual Nickel Free | <b>Code:</b> UPRO011 |



**Square Composite Button**

Pack 10.

**Code:** UPRO008



**Freddy Ceramic Button**

Pack 10.

**Code:** UGSC001



**Freddy Hook**

Pack 10.

**Code:** UGSC002

MAKE ME  
SMILE.

**EMS**+



**GBT PREVENTS DEMINERALIZATION, WHITE SPOTS, CARIES, GINGIVITIS, CALCULUS FORMATION AND BAD BREATH FOR ORTHODONTIC PATIENTS**

- Patients who wear orthodontic appliances are more likely to be subject to biofilm-related oral pathologies, such as caries or periodontal disease, as some areas in the mouth are more difficult to reach when brushing at home.
- Professional oral hygiene on a regular basis is therefore highly recommended to avoid complications.
- Not only does regular professional oral hygiene reduce the build-up of debris on and around brackets and attachments, but it also helps minimize wire friction.

**MINIMALLY INVASIVE PROPHYLAXIS ON ORTHODONTIC PATIENTS**

**METAL BRACKETS**



BEFORE DISCLOSING



AFTER DISCLOSING



AFTER GBT

**CLEAR ALIGNERS**



BEFORE DISCLOSING



AFTER DISCLOSING



AFTER GBT

**SAFE FOR ALL ORTHODONTIC APPLIANCES**



Adhesives

Adhesives



**Starflow Flowable Composite (Starflow)**

Starflow is a light-cure, highly filled flowable composite for anterior and posterior restorations, low wear applications, cavity liner and conservative restorations.

- Radiopaque
- Sub-micron particle size gives superior polishability and resistance to staining
- True Vita shade match
- Non-sticky
- Fluoride slow release and natural fluorescence.

A1	<b>Code:</b> UEVI003
A2	<b>Code:</b> UEVI001
A3	<b>Code:</b> UEVI004
A3.5	<b>Code:</b> UEVI002
A4	<b>Code:</b> UEVI005
Replacement Tips	<b>Code:</b> UEVI006



**VOCO Aligner Flow**

Fluorescent composite for aligner attachments.

- Secure grinding and removal under UV-A light
- Precise application with perfect flow
- Light-curing through the transparent template
- Stable and abrasion resistant throughout the entire treatment period
- All Popular Vita Shades.

A1	2 x 2g	<b>Code:</b> UALFSA1
A2	2 x 2g	<b>Code:</b> UALFSA2
A3	2 x 2g	<b>Code:</b> UALFSA3
44 Gauge Tips		<b>Code:</b> UVDC044
46 Gauge Tips		<b>Code:</b> UVDC046

**GoTo Adhesive Paste in Push Syringe**



Now Fluoresces Under Black Light for Complete Clean-up! GoTo is a light cure adhesive with superior handling properties.

It does not allow bracket flotation once placed on the tooth or during flash removal. A fraction of the retail cost of leading, competitive light cure pastes.

1 x 4gm Syringe	<b>Code:</b> UREL036
10 x 0.04gm Tips	<b>Code:</b> UREL037



**Flow Tain**

- Needle Tip Dispenser
- Using tips provided, these materials are easily injected onto custom pads
- Low modulus of elasticity to resist breakage
- Highly filled to bridge gaps between custom pads and tooth surface
- Fluoresces under UV light

Flow Tain + 10 Tips	<b>Code:</b> UREL042
Flow Tain LV + 10 Tips	<b>Code:</b> UREL043
Flowtain 19 gauge Tips x 50	<b>Code:</b> UREL051



**Assure Plus**

All surface primer, eliminates hydrofluoric acid application. Enhances bond strength on all surfaces, tolerates moisture contamination. Primes and conditions any surface.



6ml Bottle	<b>Code:</b> UREL035
------------	----------------------



**L.C.R. Paste in Push Syringe**

L.C.R. is a flowable, light cure orthodontic adhesive intended to be used within an orthodontic, dental or paediatric dental practice for the bonding of lingual retainers, creation of occlusal buildups and for retention of thermoplastic aligners.

5g L.C.R. Paste in Push Syringe Micro+ nano fill for flowable consistency with durability of highly filled composite. Smooth viscosity. Cures to a glassy finish. Patient comfort – reduced wear.

1 x 5gm Syringe	<b>Code:</b> UREL039
1 x 5gm Syringe + 5 Tips	<b>Code:</b> UREL040
11 x 0.45gm Tips	<b>Code:</b> UREL041



**Bond Aligner Adhesive with Tips**

Bonds Attachments To Clear Aligners

- Secure aesthetic bonds
- Bond to any clear aligner material
- No primers or conditioners
- Clear light cure paste
- High modulus of elasticity
- Conforms to aligner distortions
- Secure bonds throughout
- Bond ceramic, metal or plastic attachments.

1 x 1.4gm Syringe +10 Tips	<b>Code:</b> UREL038
----------------------------	----------------------





**Essix  
C+ Retainer Sheets**

**Round 5"**

<b>Pack 100</b>				
.040"	1mm	120mm	1C120-50	<b>Code:</b> UORC064
<b>Pack 50</b>				
.040"	1mm	125mm	1C125-50	<b>Code:</b> UORC065

**Square 5" x 5"**

<b>Pack 50</b>				
.040"	1mm	127mm	1C500-50	<b>Code:</b> UORC063

**Essix  
Plus Retainer Sheets**

**Round 5"**

<b>Pack 80</b>				
.040"	120mm	1APL40C120		<b>Code:</b> UORC081
.040"	125mm	1APL40C125		<b>Code:</b> UORC082

**Square 5" x 5"**

<b>Pack 80</b>				
.040"	5"	127mm	1APL40S5	<b>Code:</b> UORC080

**Essix  
A+ Retainer Sheets**

**Round 5"**

<b>Pack 50</b>				
.80"	2mm	120mm	11081	<b>Code:</b> UORC041
.080"	2mm	125mm	11080	<b>Code:</b> UORC042

**Square 5" x 5"**

<b>Pack 50</b>				
.030"	.75mm	10035-50		<b>Code:</b> UORC030
.040"	1mm	10045-50		<b>Code:</b> UORC034
.060"	1.5mm	10065		<b>Code:</b> UORC038
.080"	2mm	10085		<b>Code:</b> UORC040

**Essix  
Ace Retainer Sheets**

**Round 5"**

<b>Pack 50</b>				
.030"	120mm .75mm	1A320-50		<b>Code:</b> UORC046
.030"	125mm .75mm	1A325-50		<b>Code:</b> UORC047
.040"	120mm 1mm	1A420-50		<b>Code:</b> UORC050
.040"	125mm 1mm	1A425-50		<b>Code:</b> UORC049

**Square 5" x 5"**

<b>Pack 50</b>				
.030"	.75mm	1A350-50		<b>Code:</b> UORC060
.040"	1mm	1A450-50		<b>Code:</b> UORC048



**3DAlign**  
Unmatched Simplicity  
A Complete Aligner System

- Full Treatment Control**  
Whatever the Malocclusion 3DAlign will deliver
- 3D Flex Polymer Plastic**  
Allows controlled and predictable tooth movements
- UK Based Production Team**  
No UK Customs import delays or duties
- Technicians On Hand For Full Clinical Support**  
Stage production in lots of 5 to enable patient tracking.
- Environmentally Friendly**  
Reduced plastic waste for aligner case modification. No need to discard sets of aligners if the patient is not tracking correctly
- Progressive Treatment Options**  
Staged level costs dependant on the case complexity, which is monitored on case submission.
- Average Treatment Plan**  
10 to 20 Aligners
- Patient Case Visualisation**  
Patient Empowerment.



Contact: 01142 721392

Find out more!



**Perfection Plus Crown Scissors Curved**

Stainless steel Bee Bee Crown Scissors Curved, high quality, excellent functionality. Fully autoclavable.

**Code:** UPPSCCC



# Tokuyama omniCHROMA

## Ortho Application

- Ideal for use with lingual retention.
- Perfect for creating buttons and attachments during aligner treatment.



## Just One Does It All

Colourless **Omnicroma** is the world's only composite that shades by structural colour. It's unique spherical filler (Smart Chromatic Technology), enables **infinite colour matching** through A1 to D4.

The colour of the surrounding tooth naturally appears in a semi-translucent final shade, creating a **perfect match every time!**

Tokuyama's spherical fillers ensure Omnicroma is unrivalled in **polishability, abrasion resistance & natural aesthetics.**

With **faultless shade matching**, Omnicroma is undoubtedly the most advanced composite available.

Use the F blocker to stop any intrinsic tooth shading to transfer through to the flowable composite.

## Colourless Omnicroma...



### Flow

Code: TK0326

### F Blocker

Code: TK0328



**WORKING TIME: 90 SECONDS**

**CURING TIME: 20 SECONDS**



## The teeth whitening for you

**BlancOne® HOME+ EVER Is a night treatment especially suited for use in aligners.**

Featuring 7% carbamide peroxide, it can be used for extended periods. The addition of vitamins D and K as catalysts, and nano-hydroxyapatite (nHAP), ensure that the gel is effective, whilst providing conditions for remineralisation of the teeth.



Type	Specification	Active Ingredients	Treatment cycle	Strong points
HOME+ EVER	Orthodontic aligners Sensitive teeth Maintenance	CP 7% (HP 2.4%) nHAP 1.0% Vitamin Booster	10-15 applications at night	Effectiveness Stability

### BlancOne Home+ Ever Single Kit

Code: UBOHESK

### BlancOne Home+ Ever Multi Kit

Code: UBOHEMK

