


A More Innovative & Cost Effective Alternative to Locators™

Kerator!

Kerator Denture Abutments are self locating with dual retention allowing easy placement and removal of dentures for your patients.

Denture Abutments For All Implants


Adin, Nobel Biocare)Branemark, Sterioss, Replace, Replace Select), 3i Astra, ITI. Ankylos, Biohorison, Biolok, Bicon, Camlog. Endopore, Frialit-2 & Xive, Lifecore, Simmer (Calcitek, Paragon, Swissplus), Uniti, Pitt-Easy, Osstem, Dentium, Dio, Dentis, Megagen, Cowell Medi, Warantec...



Excellent Retention & Resilience

Dual Retention(outside and inside) ensures the longer lasting performance while still allowing mobility for patient comfort. The self-aligning design allows easy insertion.

Collar Heights from Zero to 6.0mm




Low Vertical Height

The 15° angled abutment enables denture attachment to non-parallel implants and can correct more than 40° angles.



Versatility

Kerator abutments have much lower vertical height than regular ball attachments



Easy Handling

The plastic carrier supplied with each abutment makes it easier, safer to handle and check alignment in the mouth.



£83.00 +vat*

All abutments are supplied with denture cap and 3 inserts.

*£89.00 for angled abutments

