

OFFERS 2024

Valid until : 30.06.2024

Bien Air⁺
Dental



AT AN
UNBEATABLE
PRICE

Discover more on page 10

CLASSIC
1:1

SWISS  MADE

Bien-Air UK Ltd

Unit 19 Crawley Business Centre Stephenson Way Three Bridges Crawley RH10 1TN Tel. +44 (0)1293 550 200 ba-uk@bienair.com www.bienair.com

EDITO

Dear Practitioners,

It is with great satisfaction that we send you this editorial, bearing good news and strong commitments from BIEN-AIR for the first half of 2024.

At BIEN-AIR, we consider your comfort as oral health professionals and the safety of your patients to be top priorities. Our commitment to constantly improving your practice is reflected in significant investments in our industrial facilities.

By investing in sustainable energy sources such as photovoltaics and geothermal energy, we ensure the environmentally-friendly operation of our machines. The heat emitted by our equipment is judiciously used to heat our buildings, thus eliminating the use of fossil fuels. Bien-Air is actively committed to minimizing carbon emissions, thus contributing to the preservation of our planet.

These commitments strengthen our operational capacity, which translates into the ability to provide you with top-quality instruments at competitive prices, without compromising on quality. A case in point is the iconic 1:1 blue contra-angle ring on page 10, Bien-Air quality at an unbeatable price.

The New Year also marks an important step in our European expansion, with the transformation of our French subsidiary into Bien-Air Europe. This development marks our determination to achieve excellence and consolidate our market position, while continuing to strengthen our commitment to you and our business partners.

Thanks to these investments and this new organization, we are ready to take on new challenges and continue to meet your needs.

We thank you for your continued confidence. Together, we are building the future of oral health.

Your Bien-Air team.

CONTACT



CARL WOOD

AREA SALES MANAGER SOUTH
07967 593 655
carl.wood@bienair.com



PAUL TAYLOR

AREA SALES MANAGER NORTH
07553 106 005
paul.taylor@bienair.com




FIND US ON :



/bienairdental /bienairdental /bien-airdental310 /bien-air

 By telephone : +44 (0)1293 550 200

 Email us : ba-uk@bienair.com

 Our website : www.bienair.com

Throughout the catalogue, you'll find QR codes that will take you directly to the product pages on our website - just scan them! You'll find our maintenance guides, user manuals and regulatory information.



CLASSIC

PREMIUM



£ **445**

BORA 2


BORA 2 LED FOR BIEN-AIR UNIFIX COUPLING
REF 1601152-001
BORA 2 LK FOR KAVO® MULTIFLEX® COUPLING
REF 1601153-001

24 month warranty



MORE POWERFUL

With power and torque improved by 20%, the Bora 2 allows you to cut **twice as quickly**.



EXTENDED LIFESPAN

Wear-resistant **ceramic ball bearings** and a redesigned rotor extend the service life.



INCREASED PRECISION

The **Accu-Chuck PreciPlus** chuck design ensures the most secure clamping of the bur.



£ **725**

TORNADO

TORNADO LED FOR BIEN-AIR UNIFIX COUPLING
REF 1600908-001
TORNADO LK FOR KAVO® MULTIFLEX® COUPLING
REF 1600907-001
TORNADO S LED FOR BIEN-AIR UNIFIX COUPLING
REF 1600933-001
TORNADO S LK FOR KAVO® MULTIFLEX® COUPLING
REF 1600932-001

24 month warranty




WHEN POWER MEANS CONTROL

Thanks to **Steady Torque** technology, the Tornado turbine is the most powerful in the Bien-Air range



EXTENDED LIFESPAN

Wear-resistant **ceramic ball bearings** and a redesigned rotor extend the service life.



INCREASED PRECISION

The **Accu-Chuck PreciPlus** chuck design ensures the most secure clamping of the bur.



SMALLER IS BETTER !

The Tornado[®] with its **small head** allows better visibility in the mouth.



£ **850** ~~£ 890~~

OR £425 PER INSTRUMENT

DUOPACK

BORA 2

BORA 2 LED FOR BIEN-AIR UNIFIX COUPLING
REF 1601152-001
BORA 2 LK FOR KAVO® MULTIFLEX® COUPLING
REF 1601153-001

2x TURBINES



£ **175**

UNIFIX

OPTIC COUPLING
REF 1600363-001

12 month warranty



£ **210**

MULTIFLEX LK

LED OPTIC
REF 1600902-001

12 month warranty



£ **230**

UNIFIX REG

OPTIC WITH WATER SPRAY ADJUSTER
REF 1600086-001

12 month warranty



£ **235**

MULTIFLEX LK REG

LED OPTIC WITH WATER SPRAY ADJUSTER
REF. 1600866-001

12 month warranty



£ **1 295** ~~£ 1 450~~

OR £647.50 PER INSTRUMENT

DUOPACK

TORNADO

TORNADO LED TO FIT BIEN-AIR UNIFIX
REF 1600908-001
TORNADO LK TO FIT KAVO® MULTIFLEX®
REF 1600907-001
TORNADO S LED TO FIT BIEN-AIR UNIFIX
REF 1600933-001
TORNADO S LK TO FIT KAVO® MULTIFLEX®
REF 1600932-001

CONTRA - ANGLES

CLASSIC

PREMIUM



£ **580**

CA 1:5 L

CA 1:5 L STANDARD REF 1600386-001
CA 1:5 L MICRO-SERIES REF 1600690-001
24 month warranty



ENHANCED SAFETY
With the patented CoolTouch heat-arresting technology, heat accumulation is prevented to reduce the risk of burns.



ACCU-SPRAY WITH ILLUMINATION
This ultra-precise nebulization composed of three separate air/water sprays ensures complete cooling of the bur and treatment site.



SCAN ME



£ **320**

CA 1:1 L

CA 1:1 L STANDARD REF 1600384-001
CA 1:1 L MICRO-SERIES REF 1600691-001
24 month warranty



SOFT PUSH
Clamping mechanism with push button exclusive to Bien-Air. Light pressure is all that is needed.



PRÉCISION ACCRUE
The Accu-Chuck PreciPlus chuck design ensures the most secure clamping of the bur.



SCAN ME



£ **685**

NOVA 1:5 L

CA NOVA 1:5 L REF 1601139-001
CA NOVA 1:5 L MICRO-SERIES REF 1601138-001
36 month warranty



ENHANCED SAFETY
With the patented CoolTouch+ heat-arresting technology, heat accumulation is prevented to reduce the risk of burns.



EXTENDED LIFESPAN
Thanks to the new Cylcro gears designed to limit friction and reduce their wear, durability is increased.



4X MORE DURABLE HEAD
Made from stainless steel, it is more resistant to shocks and has an increased service life of 30% compared to a classic contra-angle.



SCAN ME



£ **545**

NOVA 1:1 L

CA NOVA 1:1 L REF 1601137-001
CA NOVA 1:1 L MICRO-SERIES REF 1601136-001
36 month warranty



EXTENDED LIFESPAN
Thanks to the new Cylcro gears designed to limit friction and reduce their wear, durability is increased.



4X MORE DURABLE HEAD
Made from stainless steel, it is more resistant to shocks and has an increased service life of 30% compared to a classic contra-angle.



£ **265**

PM 1:1

PM 1:1 STANDARD REF 1600383-001
PM 1:1 MICRO-SERIES REF 1600693-001
24 month warranty



SCAN ME

£ **1 095** ~~£1 160~~

OR £547.50 PER INSTRUMENT

DUOPACK

CA 1:5 L

CA 1:5 L STANDARD REF 1600386-001
CA 1:5 L MICRO-SERIES REF 1600690-001

2x CONTRA-ANGLES RED



£ **860** ~~£ 900~~

DUOPACK

CA MIXED

CA 1:1 L STANDARD REF 1600384-001
CA 1:1 L MICRO-SERIES REF 1600691-001
CA 1:5 L STANDARD REF 1600386-001
CA 1:5 L MICRO-SERIES REF 1600690-001

1x CONTRA-ANGLE RED
+ 1x CONTRA-ANGLE BLUE



£ **1 310** ~~£1 370~~

OR £655 PER INSTRUMENT

DUOPACK

NOVA 1:5 L

CA NOVA 1:5 L REF 1601139-001
OR CA NOVA 1:5 L MICRO-SERIES REF 1601138-001

2x CONTRA-ANGLES NOVA RED



£ **1 175** ~~£1 230~~

DUOPACK

NOVA MIXED

CA NOVA 1:5 L REF 1601139-001
CA NOVA 1:5 L MICRO-SERIES REF 1601138-001
CA NOVA 1:1 L REF 1601137-001
CA NOVA 1:1 L MICRO-SERIES REF 1601136-001

1x CONTRA-ANGLE NOVA RED
+ 1x CONTRA-ANGLE NOVA BLUE



PROPHY ENDO



£ **470**
PROPHY 15:1
 CAP 15:1
 REF 1600290-001
 24 month warranty



£ **630**
ENDO 4:1
 CA ENDO 4:1
 REF 1600955-001
 12 month warranty



£ **730**
MICROMOTOR MC2
 MC2 LED REF 1600681-001
 MC2 ISOLITE 300 REF 1600078-001
 18 month warranty

£ **755**
MICROMOTOR MC3
 MC3 LED REF 1600680-001
 MC3 LK REF 1600077-001
 18 month warranty

£ **275**
PROPHY 10:1
 PM 10:1
 REF 1600289-001
 24 month warranty



Prices are quoted excluding VAT.

LABORATORY



£ **975**
STATION STS-TRIMMER + TDS 890
 STATION STS-TRIMMER
 REF 1700296-001
 12 month warranty



£ **425**
STATION S001 + TD 783
 STATION S001
 REF 1700126-001
 12 month warranty



£ **320**
TURBINE TD 783
 TD 783
 REF 1600161-001
 12 month warranty



£ **399**
STATION SF 811 + TD 783
 STATION SF 811
 REF 1700083-001
 12 month warranty



£ **585**
TURBINE TDS 890
 TD S890
 REF 1600163-001
 12 month warranty



Prices are quoted excluding VAT.

FROM
£ 290

CLASSIC
CA 1:1 / **LEGENDARY RELIABILITY**

Rediscover a timeless handpiece. With close to **30 years of existence**, the durability of the CA 1:1 is beyond compare. Designed with cutting-edge innovation, it offers unrivalled precision guaranteed to last for years. Thanks to the CA 1:1 contra-angle, Bien-Air quality is now at your fingertips.




micro
SERIES

AT AN UNBEATABLE PRICE



DUALOOK™ lighting paves the way for exceptional visibility of the working area.



ENHANCED LIGHT INTENSITY thanks to multi-strand optical glass conductors.



ACCU-CHUCK™ technology provides outstanding axis stability, eliminating noise and vibration.



ACCU-SPRAY™ system ensures better irrigation of the treated area, and safer treatment for your patients.



SWISS  MADE

CLASSIC **1:1**



£ 320
CA 1:1
REF 1600424-001



£ 320
CA 1:1 L
REF 1600384-001



£ 290
CA 1:1 L Micro-Series
REF 1600691-001



SCAN TO HEAD
TO BIENAIR.COM

INSTRUMENT CARE



£ **1795**
LUBRICARE 2
 LUBRICARE 2 REF 1601156-001
 + DUOPACK MAINTENANCE
 24 month warranty



SCAN ME

LUBRICARE 2 COUPLINGS

BUY ONE COUPLING = GET ONE FREE
 WHEN PURCHASING LUBRICARE 2
 (UP TO 2 IDENTICAL CONNECTIONS)



£ **130**
BIEN-AIR
 REF 1600082-001

£ **150**
SIRONA
 REF 1308700-001



£ **150**
W&H
 REF 1308701-001

£ **150**
KAVO
 REF. 1308699-001

Prices are quoted excluding VAT.

ACCESSORIES



£ **138**
LUBRIFLUID
 LUBRICATION OF YOUR INSTRUMENTS
 BOX OF 6 BOTTLES
 LUBRIFLUID REF 1600064-006



SCAN ME



£ **138**
SPRAYNET
 CLEANING OF YOUR INSTRUMENTS
 BOX OF 6 BOTTLES
 SPRAYNET REF 1600036-006



SCAN ME



£ **95**
AQUACARE
 ELIMINATES SALINE FROM YOUR INSTRUMENTS
 BOX OF 6 BOTTLES
 AQUACARE REF 1600617-006



SCAN ME



£ **19**
LUBRIMED
 LUBRIMED REF 1000003-001



£ **19**
LUBRIMED CARTRIDGES
 PACK OF 6
 CARTOUCHES LUBRIMED REF 1600037-006



£ **13.50**
OIL ABSORBING PAD
 5 PIECES
 PAD REF 1306003-005



£ **17**
OIL ABSORBING SHEET
 5 PIECES
 FEUILLES REF 1308452-005



£ **26**
DOOR OIL ABSORBING PAD
 5 PIECES
 DOOR PAD REF 1308446-005

Prices are quoted excluding VAT.

CONTRA-ANGLE

1:2.5

SWISS  MADE

Bien Air
Dental

ULTIMATE RELIABILITY

ONLY

£ 955

INTRODUCTORY OFFER
INCLUDES 4 FREE BURS

LATCH TYPE

HIGH SPEED  SYSTEM

SETTING NEW STANDARDS

Discover a new world when you combine your high-speed contra-angle with latch-type burs. The result is an instrument whose durability is multiplied, demonstrating unrivaled versatility and reliability.

micro
SERIES


APICECTOMY
HEMI SECTION
CORONAL ELONGATION
**WISDOM TOOTH
EXTRACTION**

Offering you a high-tech contra-angle with unrivaled performance, the new version of the CA 1:2.5 is now compatible with extra-long 34 mm latch type surgical burs. Making this instrument the new standard in surgery and of periodontology.

Its internal irrigation line and smaller head allow ergonomic handling of the instrument and provide better visibility of the surgical area. The use of stainless-steel guarantees greater reliability and superior resistance to corrosion.



- 1** IMPROVED DURABILITY THROUGH THE DESIGN OF UNIQUE LATCH TYPE BURS
- 2** ENHANCED STABILITY AND REDUCED VIBRATIONS THANKS TO THE UNIQUE CHUCK DESIGN
- 3** SUPERIOR AND OPTIMIZED CUTTING SPEED FOR EXCEPTIONAL EFFICIENCY



NEW Type Zekrya

SET DE 5x RA BURS ZEKRYA L 28mm
REF. 1101384-005
SET DE 5x RA BURS ZEKRYA XL 34mm
REF. 1101385-005



NEW Type Lindemann

SET DE 5x RA BURS LINDEMANN L 28mm
REF. 1101386-005
SET DE 5x RA BURS LINDEMANN XL 34mm
REF. 1101387-005

SWISS  MADE

CA 1:2.5 L + STARTER KIT OF 4 BURS REF. 1601163-001
CA 1:2.5 L MS + STARTER KIT OF 4 BURS REF. 1601164-001

SCAN TO HEAD
TO BIENAIR.COM




CHIROPRO

CHIROPRO PLUS



Simplified handling and visibility of the surgical area thanks to the unique internal irrigation line.

Multistrand fibre optic technology.



Perfect balance and ease of use with unrivalled ergonomics when it is used with the MX-i PLUS micromotor.

£ **2 250**
CHIROPRO

CHIROPRO 3RD GEN REF 1700708-001
2 years warranty

£ **2 850** £3 940
CHIROPRO
+ 1 CONTRA-ANGLE 20:1 L

SET CHIROPRO 3RD GEN REF 1700707-001
2 years warranty

£ **3 390** £4 840
CHIROPRO
+ 2 CONTRA-ANGLES 20:1 L

SET CHIROPRO 3RD GEN REF 1700767-001
2 years warranty

£ **2 675**
CHIROPRO PLUS

CHIROPRO PLUS 3RD GEN REF 1700710-001
2 years warranty

£ **3 290** £4 311
CHIROPRO PLUS
+ 1 CONTRA-ANGLE 20:1 L

SET CHIROPRO PLUS 3RD GEN REF 1700709-001
2 years warranty

£ **3 830** £5 672
CHIROPRO PLUS
+ 2 CONTRA-ANGLES 20:1 L

SET CHIROPRO PLUS 3RD GEN REF 1700709-001
+ CA 20:1 L MICRO-SERIES REF 1600692-001
2 years warranty

iCHIROPRO



COMPATIBLE WITH
CO-DIAGNOSTIX
SOFTWARE


PRE PROGRAMMED
IMPLANT PROTOCOLS

MULTI IMPLANT
OPERATIONS

MULTI USER
INTERFACE

* iPad non-inclus

Simplified handling and visibility
of the surgical area thanks to the
unique internal irrigation line.



Perfect balance and ease of use with
unrivalled ergonomics when it is
used with the MX-i PLUS micromotor.

Multistrand fibre
optic technology.



£ **3 520**

iCHIROPRO

iCHIROPRO REF 1700439-001
2 years warranty

£ **4 060** ~~£4 968~~

iCHIROPRO
+ 1 CONTRA-ANGLE 20:1 L

SET iCHIROPRO REF 1700440-001
2 years warranty

£ **4 600** ~~£6 329~~

iCHIROPRO
+ 2 CONTRA-ANGLES 20:1 L

SET iCHIROPRO REF 1700440-001
+ CA 20:1 L MICRO-SERIES REF 1600692-001
2 years warranty

Prices are quoted excluding VAT.



£ **985**

CA 20:1 L

INTERNAL IRRIGATION
CA 20:1 L REF 1600598-001
CA 20:1 L MICRO-SERIES REF 1600692-001
24 month warranty



£ **295**

PM 1:1

EXTERNAL IRRIGATION
PM 1:1 REF 1600052-001
24 month warranty



£ **1 182**

SET MX-i PLUS
+ CABLE

(FOR iCHIROPRO OR CHIROPRO PLUS)
SET MX-i PLUS REF 1700758-001



£ **1 120**

SET MX-i LED
+ CABLE

(FOR CHIROPRO)
SET MX-i LED REF 1700744-001



£ **66**

IRRIGATION
LINES

BOX OF 10 LINES
LIGNES D'IRRIGATION REF 1500984-010



£ **44**

STERILE
PROTECTIVE
SCREENS

FOR IPAD iCHIROPRO REF 1501746-010
FOR CHIROPRO 3RD GEN REF 1502329-010



£ **70**

MAGNETIC
CONTROL
KNOB

KNOB REF 1307031-001



Prices are quoted excluding VAT.

ELECTRIC



£ **2850**

SET iOPTIMA^{INT}

SET iOPTIMA^{INT} REF 1700704-001
+ REF 6200217-001
(INCLUDES IPAD MINI)



SCAN ME

£ **2150**

MX2 OPTIMA^{INT} PRO

MX2 OPTIMA^{INT} PRO
REF 6200186-001

£ **1955**

MX2 OPTIMA^{INT} PRO W/O TRANSFORMER
REF 6200164-001



SCAN ME

£ **1675**

MCX OPTIMA^{INT}

MCX OPTIMA^{INT}
REF 6200222-001

£ **1425**

MCX OPTIMA^{INT} W/O TRANSFORMER
REF 6200183-001



SCAN ME

£ **1170**

MCX LED

MCX LED
REF 6200188-001

£ **955**

MCX MOTOR
REF 6200189-001



£ **1215**

OPTIMA MCX

SEMI-INTEGRATED
MICROMOTOR SOLUTION

OPTIMA MCX REF 1700595-001
OPTIMA MCX BLUE REF 1700596-001
OPTIMA MCX ORANGE REF 1700597-001
OPTIMA MCX GREEN REF 1700598-001
OPTIMA MCX PINK REF 1700599-001



SCAN ME