

# Kerr™ SELECT PROMOTIONS Q2

1<sup>st</sup> April - 30<sup>th</sup> June 2025

## OptiBond™ Universal

Single-Component Universal Dental Adhesive

OptiBond™ Universal is the one-bottle dental bonding agent that combines the reliability of patented OptiBond GPDM technology with the ease of use of a universal solution. It offers reliable bond strength to all surfaces and substrates, in any etching technique.

- One bottle. One coat
- Proven strength
- Complete versatility
- Protection against post-operative sensitivity



36519 - OptiBond Universal Refill 5ml

*Buy 3 get 1 FREE\**



## SimpliShade™ BULK FILL & BULK FILL FLOW

Using one shade that blends seamlessly with all 16 VITA® classical shades, SimpliShade Bulk Fill and Bulk Fill Flow saves you time selecting shades and managing inventory.

**SimpliShade Bulk Fill**  
37211 - Syringe - 4g  
37212 - Unidose 20x 0.25g

**SimpliShade Bulk Fill Flow**  
37217 - Syringe - 2x 2g  
37218 - Unidose 20x 0.25g

*Buy 3 get 1 FREE\**

# SimpliShade™

Simplified Universal Composite System

- Syringe 4g**  
**37001** - Light  
**37002** - Medium  
**37003** - Dark  
**37004** - Bleach White  
**37005** - Universal Opaque



- Unidoses 10x 0,25g**  
**37008** - Dark  
**37009** - Bleach White  
**37010** - Universal Opaque

- Unidoses 20x 0,25g**  
**37006** - Light  
**37007** - Medium



*Buy 3 get 1 FREE\**

## Herculite XRV Ultra

Nanohybrid Dental Composite



## Herculite™ XRV

Microhybrid Dental Composite



*Buy 3  
get 1  
FREE\**

**CONTACT  
YOUR KERR  
FOR SHADES**

## Herculite™ XRV Ultra Flow

Nanohybrid  
Flowable Composite



**CONTACT  
YOUR KERR REP  
FOR SHADES**

*Buy 3  
get 1 FREE\**

## SonicFill™ 3

Single Step Bulk Fill Composite

- Unidoses 20x 0,25g**  
**36711** - Refill A1  
**36712** - Refill A2  
**36713** - Refill A3  
**36714** - Refill B1



*Buy 5 get 1 FREE\**

## OptiDisc™ Refill

Polishing and Finishing Discs



*Buy 1 get 1 FREE\**

# Vertise™ Flow

Self-Adhering  
Flowable Composite



**CONTACT  
YOUR KERR REP  
FOR SHADES**

**Buy 2 get 1 FREE\***

# Metafix™

All-In-One Matrix System



Refill 50pcs  
**3601** - MetaFix small  
**3602** - MetaFix medium  
**3603** - MetaFix large

**Buy 2 get 1 FREE\***

# NX3 Nexus™ Third Generation

Universal Adhesive Resin Dental Cement

- Dual-cure automix syringe eliminates hand mixing
- For indications requiring unlimited work time
- Excellent adhesion to dentine, enamel, CAD/CAM blocks, ceramic, porcelain, resin and metal
- Self-etch or total-etch
- Superior colour stability for long-term aesthetics whether dual-cure or light-cure
- Easy clean-up in gel state



**Automix Syringe 1.8g**

- 33646** - Bleach
- 33643** - Clear
- 33647** - Opaque
- 33644** - White
- 33645** - Yellow

**Light-Cure Syringe 5g**

- 33651** - Bleach
- 33648** - Clear
- 33652** - Opaque
- 33649** - White
- 33650** - Yellow

**Buy 4  
get 1 FREE\***

# Maxcem™ Elite Universal

Self-Etch, Self-Adhesive Resin Dental Cement

- One-Peel™ clean-up, no hand mixing and no refrigeration required
- Dual-cure mechanism
- Compatible with all indirect restorations
- No colour shift
- Colour clean-up indicator

**Maxcem Elite Chroma**  
2x Automix Syringe 5g  
**36372** - Clear  
**36373** - White  
**36374** - Yellow

**Maxcem Elite**  
2x Automix Syringe 5g  
**34058** - Clear    **34060** - White Opaque  
**34059** - White   **34061** - Yellow  
**37003** - Dark



**Buy 3 get 1 FREE\***



# ZenFlex™

NiTi Rotary Shaping Files



## High Cutting Efficiency

Designed with triangular cross section with proprietary heat treatment, ZenFlex exhibits high cutting efficiency that lets you reach the full working length of even complex root geometries.



## Minimally Invasive

A maximum flute diameter of 1mm and a non-cutting tip work together to minimize transportation and preserve the tooth's structural integrity.



## Exceptional Strength

Resists torsional stress and cyclic fatigue resulting in minimal breakage and reduced transportation.



## Excellent Flexibility


Proprietary heat-treatment design delivers excellent controlled memory and flexibility for up to 90° curves.

**CONTACT YOUR KERR REP FOR SIZES**

**Buy 3 get 1 FREE\***

## Pro-Cup™

Splash-free Propphy Cups



**Buy 3  
get 1 FREE\***

## Coloured Propphy Cups



**Buy 3  
get 1 FREE\***

## Miniature Brushes and Pro-Brush™

Cleaning and polishing brushes



**Buy 1  
get 1 FREE\***

Not to be used in combination with any other offers or discount schemes. All free stock shipped directly from Dealer and claimed back quarterly along with transparency report. Dealer must purchase equivalent stock from Kerr UK during the promotional period to be eligible to claim promotional discount, or free of charge item credit. Credit will only be applied based on the quantity of promotional items purchased during the corresponding promotional period. Claims must be received within 30 days of promotion ending. Free item must be of same value within the product brand range. E & OE.

Kerr UK Ltd

c/o Omega Stockley Park | 4 Longwalk road | Stockley Park | Uxbridge | UB11 1FE | UK

Phone: 01733 892292 - [www.kerrdental.com](http://www.kerrdental.com)

MKT-25-0384 Rev0

**Kerr™**