



SAY HELLO

to the evolution in
ultrasonic and
powder therapy

NEW

Varios Combi Pro2
+ NSK's MIT guide
(Minimally Invasive Therapy)

FIND OUT MORE
ON PAGE 7



Special offers

Summer 2025

VALID UNTIL 31/08/2025

New NLZ E

Compact portable electric micromotor,
for general & endodontic applications.
Page 4 →



Restorative

New Ti-Max Z99L

New 1:5 Speed Increasing Contra-Angle with a 100° angled head
to enable easier access. Page 3 →



Implantology and Oral Surgery

New VarioSurg 4

Minimally invasive ultrasonic
bone surgery system
Plus option to connect to Surgic
Pro2 with one foot control.
Page 11 →



Decon and IPC

New iClave 24

+ 12 months additional
NSK handpiece warranty.
Page 12 →



Restorative - Turbines

Ti-Max Z HIGH PERFORMANCE TURBINES

2 + 1 **FREE*** + 1 FREE Z micro



4 years for iClave users



NEW

ULTIMATE POWER for quicker preps and cutting Zirconia

• Lightweight titanium body with Duragrip coating to reduce hand fatigue



- Optic
- **ULTIMATE Power 30-44w****
- Std head: Ø 12.1 x 12.7 mm

MODEL Z990L
CODE PA23730001
RRP £1,214

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®

**At 3 - 4.2 bar



- Optic
- **Power 28-41w****
- Mini head: Ø 10.6 x 10.8 mm

MODEL Z890L
CODE PA23720001
RRP £1,214

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®



- Optic
- **Power 20w**
- **Micro head: Ø 9 x 10.8 mm**

MODEL Z micro (ultra-mini head)
CODE PA23570001
RRP £1,214

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®

for regular or short shank burs

Scan/Click
For full info on models/order codes



S-Max M



3 years for iClave users



OPTIC 3 + 1 **FREE** + 1 **FREE Pico***



• Stainless steel • Enhanced power • 4 hole spray • DIY cartridge change • Clean head system



- Optic
- **Power 26w**
- **Std head: Ø 12.1 x 13.3 mm**

MODEL M900L
CODE P1254001
RRP £800

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®



- Optic
- **Power 23w**
- **Mini head: Ø 10.6x12.4mm**

MODEL M800L
CODE P1253001
RRP £800

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®



- Optic
- **Power 9w**
- **Ultra-mini head: Ø 8.6 x 9mm**

MODEL Pico
CODE P1140001
RRP £1,160

Also available in fittings for KaVo®, Dentsply Sirona®, W&H®, Bien-Air®

Used only with short shank burs.



- **Non-Optic**
- **Power 26w / 23w**
- **Std and mini heads**

MODEL M900
CODE P1256001
RRP £745

MODEL M800
CODE P1255001
RRP £745

M900 also available in fittings for KaVo®, W&H®

Pana-Max2



2 years for iClave users



1 + 1 **FREE***

• Non-Optic turbine with single hole spray • Clean head system
• Save time/money with DIY cartridge change • **12 month warranty**



- **Non-Optic**
- **Head: Ø 11.2 x 13.4 mm**
- **Coupling connection**

MODEL Pana-Max2 PTL
CODE P1211001
RRP £484
for NSK PTL coupling

MODEL Pana-Max2 KV
CODE P1212001
RRP £484
for KaVo® MULTIFlex® coupling



- **Non-Optic**
- **Head: Ø 11.2 x 13.4 mm**
- **Direct connection**

MODEL Pana-Max2 M4
CODE P1154001
RRP £439
For Midwest M4 hose

MODEL Pana-Max2 B2
CODE P1155001
RRP £439
for Borden hose

PANA SPRAY Plus

1 + 1 **FREE***



Specially formulated handpiece and air motor lubricant.

MODEL PANA SPRAY Plus
CODE Z182100
RRP £35.20

BUY NOW



*Terms and conditions apply. Offers only available through participating dealers.

Restorative - Contra Angles

Ti-Max Z



4 years for
iClave users



2 + 1 FREE + 1 FREE Z micro*



- Lightweight external titanium body • Optic with internal water • Small head/narrow neck
- Unique spray switch on Z99L, Z95L, Z45L, Z25L



- Optic
 - 1:5 Increasing
 - Head size
ø8.9 x 13.4mm,
at 100° angle
for **improved
access**
- MODEL Z99L
CODE CA13140001
RRP £1,666

NEW



- Optic
 - 1:5 Increasing
 - Head size
ø8.9 x 13.4mm
- MODEL Z95L
CODE C1034001-SW
RRP £1,666



- Optic
 - 1:5 Increasing
 - 10% shorter
& lighter than Z95L
 - Head size
ø8.9 x 13.4mm
- MODEL nano95LS
CODE C1099001
RRP £1,586



- Optic
 - 1:4.2 Increasing
+ Spray switch
 - 45° head
 - Head size
ø9.6 x 14.5mm
 - Ideal for minor oral
surgery too
 - For regular or extra
long (25mm) FG burs
- MODEL Z45L
CODE C1064001
RRP £1,666



- Optic
 - 1:1 Direct Drive
- MODEL Z25L
CODE C1038001-SW
RRP £1,320



Scan/Click

For a full range of
NSK Contra-angles

S-Max M



3 years for
iClave users



3 + 1 FREE + 1 FREE Pico*



- Stainless steel • All with internal water • Choice of optic & non-optic



- Optic
 - 1:5 Increasing
 - Head size 10 x H.14.6mm
- MODEL M95L
CODE C1023001
RRP £1,076

- Non-Optic
 - 1:5 Increasing
 - Head size 10 x H.14.6mm
- MODEL M95
CODE C1026001
RRP £893



- Optic
 - 1:1 Direct Drive
- MODEL M25L
CODE C1024001
RRP £871

- Non-Optic
 - 1:1 Direct Drive
- MODEL M25
CODE C1027001
RRP £572



- Optic
 - 4:1 Reducing
- MODEL M15L
CODE C1025001
RRP £1,019

- Non-Optic
 - 4:1 Reducing
- MODEL M15
CODE C1028001
RRP £681



- Non-Optic
 - 1:1 Straight
- MODEL M65
CODE H1008001
RRP £559

FXm



2 years for
iClave users



1 + 1 FREE*



- Non-Optic Straight
 - 1:1 Direct Drive
 - Push button chuck
 - No internal water
- MODEL FX65m
CODE H1014001
RRP £255



- Non-Optic Contra-angle
 - 1:1 Direct Drive
 - Push button chuck
 - No internal water
- MODEL FX25m
CODE C1052001
RRP £282



- Non-Optic Propphy
Contra-angle
 - 4:1 speed reducing
 - Max speed 1,000 rpm
 - Brushes/Cups:
Screw in
- MODEL FX57m
CODE C1054001
RRP £317

FX

3 FOR £201, 10 FOR £568*

New modern version of the classic NAC-EC



- Non-Optic Contra-angle
 - 1:1 Direct Drive
 - Latch-type chuck
 - No internal water
- MODEL FX22
CODE C1129001
RRP £165

NEW

Scan/Click

For a full range of
NSK Contra-angles



*Terms and conditions apply. Offers only available through participating dealers.

Restorative - Electric

NLZ E NEW

PORTABLE MICROMOTOR

Including reciprocating and rotary Endo features

NEW compact portable electric micromotor, for general & endodontic applications with contra angles (ie speed increasing/red band, 1:1 direct drive/blue band or speed reducing/green band).

- Lighter motor with enhanced power/torque
- Ideal for general (restorative, prophylaxis etc) and endo applications (rotary & reciprocating* modes)
- Enhanced patient safety function indicates when the handpiece needs repair to minimise it overheating
- Versatile installation options – portable or semi built-in

* The 'Reciprocating Endo Mode' can be operated only with the NLZ Endo Contra-angle and WAVE ONE® file. The NLZ requires access to a plug socket, and the use of an existing Midwest hose from the cart to gain access to coolant water/air and the existing foot control.

MODEL NLZ E
CODE Y1004004
RRP £3,363



FREE Ti-Max Z95L & NLZ Endo worth £2,031



NLX RETRO FITTED MICROMOTORS

Electric micromotors are the driving force behind speed increasing handpieces and can be fully integrated into the delivery system of most new and existing dental chairs. They can be retro fitted by an engineer or come as an option when buying new chairs. Controlling the specific speed and torque of micromotors within existing platforms is made easy thanks to built-in control pads.

Speed increasing contra-angles not only do everything a turbine does using the same FG burs, they provide a constant torque from the micromotor that does not decrease with resistance and higher load with no loss of power as they come in contact with the tooth.



To discuss the offers & options contact NSK:

OFFER 1 **BUY 3 x**  **AND GET OPTION A FREE** 

OFFER 2 **BUY 4 x**  **AND GET OPTION B FREE** 

OFFER 3 **BUY BUILT-IN OPTION A or B AND GET FREE Ti-Max Z95L** SPEED INCREASING CONTRA-ANGLE (1:5) 

- | | |
|-----------------|--|
| OPTION A | NLX nano micromotor & BUILT-IN KIT * (including MultiPad) RRP £2,387 |
| OPTION B | NLX plus micromotor (includes Rotary Endo functions) & BUILT-IN KIT * (including MultiPad) RRP £2,871 |
| OPTION C | Varios170 Lux ULTRASONIC HANDPIECE (tip options for scaling, perio, endo & restorative procedures) + BUILT-IN KIT * (including MultiPad) RRP £1,553 (Y1004145) |
| OPTION D | Option A + C RRP £3,781 |
| OPTION E | Option B + C RRP £4,238 |

Mobile dentistry

VIVA ace

MOBILE DENTISTRY AND SURGERY
BACK UP – ALL IN ONE SINGLE BOX

- Lightweight micromotor with additional endo functions
- Multipurpose ultrasonic handpiece
- Quick setup and storage
- 4-level adjustable height and angle



Width 462 mm, Height 284 mm, Depth 292 mm



MODEL VIVA ace
CODE Y1003773
RRP £13,408

Weight 8.6kg

FREE* 1 x



M95L



M25L



M65

Monthly RENTAL £433* Fully maintained for 3 years

Endodontics

Endo-Mate TC2

CORDLESS ENDO
MOTOR/HANDPIECE

The classic cordless endo motor/handpiece

- Use with rotary files


MODEL Endo-Mate TC2
CODE Y1001027
RRP £1,830



Varios Tips

ULTRASONIC TIPS FOR ENDO

Designed to aid agitation and irrigation. Compatible with stand alone
Varios 970 / Varios Combi Pro, Varios 170 (built-in).



- E5 & E6 for lateral condensation
- E4D (diamond coated) for cleaning/smoothing or widening canals
- E7 / E8 for removing filling material from the canal
- E15D (diamond coated) for locating root canals & decalcification

MODEL Varios Endo tip kit
(tip stand + 6 tips above)
CODE Z221086/6

OFFER £264

Equivalent of 3+3 free tips + free tip stand

Additional tips RRP £85.90 **Buy 1 get 1 free**
(diamond tips £92.20 Buy 1 get 1 free)

Scan/Click
For info article



Endo-Mate DT2

CORDED ENDO
MOTOR/HANDPIECE

- Use with rotary files

MODEL Endo-Mate DT2
CODE Y1004230
RRP £1,730



iPex II

APEX
LOCATOR

MODEL iPex II
CODE Y1002208
RRP £1,086



OFFER BUY Endo-Mate DT2 AND GET iPex II FREE WORTH £1,086*

Perio & oral hygiene

Scan/Click
For more info on
NSK Oral Hygiene



Varios 970 iPiezo engine® ULTRASONIC THERAPY



Ultrasonic scalers not only offer patients a more comfortable experience as less water is needed, but offer the user a wide range of 70 tips, not only for scaling/perio/implant maintenance, but for endo, extractions, restorative & MI applications.

Units come with:

- 1 handpiece
- G4, G8 & G16 tips
- 3 tip wrenches
- sterilisation case

• Optic

MODEL Varios 970 LUX
CODE Y1004161
RRP £3,103

• Non-Optic

MODEL Varios 970
CODE Y1004160
RRP £2,671

Scan/Click



For more info



OFFER

FREE

2 x



2 x



Varios2 handpiece

• Optic

MODEL VA2-LUX-HP
CODE E351050
RRP £568



Varios2 handpiece

• Non-Optic

MODEL VA2-HP
CODE E350050
RRP £317

1 + 1 FREE*

1 + 1 FREE*

Scaler Tips

1 + 1 FREE*



Scan/Click
To order



TIPS FOR VARIOS
ULTRASONIC SCALER




Scan/Click
Full range
of tips

TIPS FOR AIR SCALERS



Scan/Click
Tips for
Air Scaler



*Terms and conditions apply. Offers only available through participating dealers.

The Regeneration Course is back!

27th September 2025
Letchworth, London

Prof Luigi Nibali and
Dr Varkha Rattu present
a lecture and workshop
combo on MIST/MINST



**For more
details and
to book**
click or scan
the QR code



Perio & oral hygiene

Varios Combi Pro2 iPiezo engine®



NEW

THE EVOLUTION IN ULTRASONIC AND POWDER THERAPY

Gentle yet powerful, reliable and exceptionally durable, Varios Combi Pro2 has a host of advanced features designed to improve efficiency, control and enhance patient comfort. Engineered with meticulous attention to detail, Varios Combi Pro2 is an investment in clinical excellence and in the health of your patients.

- **Superbly comfortable water temperature, for both ultrasonic & powder**
- **Durable titanium ultrasonic handpiece**
- **Jet Nozzle for Perio Handpiece**
- **Auto Cleaning**
- **Removable Powder Chambers***
- **Multi-functional wireless foot control**

Elevate your practice and patient care with the Varios Combi Pro2 - the future of ultrasonic and powder therapy.

Ultrasonic Therapy

Powder Therapy

MODEL Varios Combi Pro2
CODE Y1500505
RRP £6,750

RENT Call for a price

Note: Connection to external 6mm air & water line sockets are needed to run this unit.

*Comes with detachable Prophy Chamber Set. Perio Chamber Set available as an optional extra.

iCart Prophy2

User-friendly and convenient trolley cart with the optional extra of a water bottle attachment, to create an independent water supply. This means that when the Varios Combi Pro2 is used with the iCart & bottle, you can treat patients without connecting to the external water mains supply.

MODEL iCart Prophy2
CODE S9098
RRP £1,077

MODEL Water Bottle Kit
CODE Y1500719
RRP £725

MODEL Perio Kit2
CODE Y1500470
RRP £1,274

Titanium Ultrasonic handpiece

- Durable & lightweight, with unrivaled ergonomics
- Non-slip & scratch-resistant Duragrip coating
- Slim LED handpiece enhances visibility



MODEL VA2-LUX-HP Ti
CODE EA1203X050
RRP £945

FREE 2ND TITANIUM ULTRASONIC HANDPIECE + FREE 2ND PROPHY HANDPIECE



Varios Combi Pro

ULTRASONIC & POWDER THERAPY UNIT

Comes with an ultrasonic handpiece and a detachable Prophy chamber.

Detachable Perio chamber available as an optional extra.

MODEL Varios Combi Pro
CODE Y1002843
RRP £5,331

Note: Connection to external 6mm air & water line sockets are needed to run this unit.

Ultrasonic Therapy

Powder Therapy

FREE Perio Kit
WORTH £1,274*



*Terms and conditions apply.
Offers only available through participating dealers.

MIT

Minimally Invasive Therapy

With Varios Combi Pro2 you can deliver NSK's Minimally Invasive Therapy (MIT).



Scan the QR code for full details on all the therapy steps

Perio & oral hygiene

Prophy-Mate neo

SUPRAGINGIVAL POWDER THERAPY

- Compact, lightweight & portable
- Includes 60° and 80° nozzles
- 360° rotation
- Includes 6 sachets of FLASH pearl® powder

CODE Y135029 for NSK PTL coupling*
 CODE Y135030 for KaVo® MULTiflex® LUX coupling
 CODE Y135031 for Sirona® quick coupling
 CODE Y135032 for W&H® Roto Quick® coupling
 CODE Y135033 for Bien-Air® Unifix® coupling

RRP £1,203



CHOICE OF OFFERS

Scan/Click
For more info



1. FREE SECOND HANDPIECE WORTH £456*



2. BUY WITH A Perio-Mate FOR £1,998* SAVE £576



3. BUY WITH A iProphy CORDLESS PROPHY HANDPIECE FOR £1,425* SAVE £939



Prophy-Mate neo nozzles



CODE T797012 60° nozzle
RRP £128.50



CODE T797013 80° nozzle
RRP £128.50

1 + 1 FREE*

IKIGAI

ORAL HYGIENE COMMUNITY

For all the latest in Perio & Prevention visit IKIGAI on



FACEBOOK



INSTAGRAM
@ikigaioralhygiene



Perio-Mate

SUBGINGIVAL POWDER THERAPY / BIOFILM REMOVAL

- Subgingival version of ProphyMate
- Eradicates biofilm in just 20 seconds
- Powder flow adjuster for delicate pockets
- Includes flexible plastic nozzle tips for deeper pockets
- Packed in a handy case with 160g powder
- Also used supragingivally on moderate staining and around composite restorations/orthodontic brackets



CODE Y1002654 for NSK coupling
 CODE Y1002655 for KaVo® MULTiflex® LUX coupling
 CODE Y1002656 for Sirona® quick coupling
 CODE Y1002657 for W&H® Roto Quick® coupling
 CODE Y1002658 for Bien-Air® Unifix® coupling

RRP £1,371

Scan/Click
For more info



+ FREE SECOND HANDPIECE WORTH £456*



Perio & oral hygiene

Powder Therapy

FLASH pearl

Supragingival powder for tough stains. Used with ProphyMate or the Varios Combi Pro (Prophy Chamber).



- Calcium carbonate powder (54µm)



4 x 300 g bottles
CODE Y900698
RRP £147.70
100 x 15 g sachets
CODE Y900693
RRP £201

SOFT pearl

Supragingival powder for moderate stains. Used with Prophy-Mate or the Varios Combi Pro (Prophy Chamber).



- Glycine based powder (65µm). Water soluble, making clean up easier.

NEW

4 x 250 g bottles
CODE Y1500721
RRP £142.40

Perio Mate Powder

For subgingival biofilm removal or mild supragingival stains. Used with Perio-Mate or the Varios Combi Pro (Perio Chamber).



- Glycine based powder (25µm)
2x 160 g bottles
CODE Y900938
RRP £102.70

BUY ANY 2 PACKS OF FLASH OR SOFT pearl
GET A FREE Prophy-Mate nozzle*



Air Scalers



Direct connection
air scaler

- Oscillation frequency 5,600-5,850Hz
- Kit includes 3 scaler tips

MODEL AS2000
CODE T698001 (Midwest M4 hose)
CODE T697001 (Borden B2/B3 hose)
RRP £510 each



Optic air scaler with
coupling connection

- Oscillation frequency 6,200-5,850Hz
- Scratch resistant with titanium body
- Kit includes 3 scaler tips

MODEL Ti-Max S970L
CODE T1005001**
RRP £1,433

**also available in KaVo® fitting T1013001
and Sirona® fitting T1015001

FREE FX22
WORTH £165



FREE 2 x FX22
WORTH £330



POWDER THERAPY PATIENT LITERATURE

Available on request with Prophy-Mate neo, Perio-Mate and Varios Combi Pro purchases.




Scan/
Click



FREE

Prophy Contra Angles



Paste defence system



FX57m



FX65m

- Max speed 1,000 rpm
- Brushes/Cups: **Screw in**
MODEL FX57m
CODE C1054001
RRP £317

- Max speed 40,000 rpm
- Brushes/Cups: **Disposable Prophy heads**
MODEL FX65m
CODE H1014001
RRP £255



AR-EC(k)



AR-EC(s)



FX22

- Max speed 5,000 rpm
- Brushes/Cups: **Snap on cups**

MODEL AR-EC(k)
CODE Y110165
RRP £160

- Max speed 5,000 rpm
- Brushes/Cups: **Screw in**

MODEL AR-EC(s)
CODE Y110164
RRP £157

NEW

- New modern version of the classic NAC-EC
- Max speed 30,000 rpm
- Brushes/Cups: **Latch grip chuck**

MODEL FX22
CODE C1129001
RRP £165

1 + 1 FREE*

3 FOR £201, 10 FOR £568*

*Terms and conditions apply. Offers only available through participating dealers.

Implantology & oral surgery

Surgic Pro and Surgic Pro2

MICROMOTORS FOR IMPLANTS / ORAL SURGERY

Scan/Click
For more about
Surgic Pro2



MODEL Surgic Pro (optic)
CODE Y1003586
RRP £6,429

MODEL Surgic Pro (non-optic)
CODE Y1003587
RRP £4,799

Scan/Click
For more about
Surgic Pro



Complete Kits contain:

- Control unit
- Motor/Cord
- 20:1 Implant contra angle
- Foot control
- 2 year warranty



MODEL Surgic Pro2 (optic)
CODE Y1004195
RRP £6,429

MODEL Surgic Pro2 (non-optic)
CODE Y1004196
RRP £4,799

FREE



Osseo 100+

OR



iSD900

OR 2 x



20:1 implant
contra-angles

ADDITIONAL MOTOR/CORD



• For Surgic Pro

MODEL SGL70M (optic)
CODE E1023001
RRP £2,854

MODEL SG70M (non-optic)
CODE E1025001
RRP £1,984



• For Surgic Pro2

MODEL SGL80M (optic)
CODE Y1004211
RRP £2,854

MODEL SG80M (non-optic)
CODE Y1004212
RRP £1,984

IRRIGATION TUBING

• Used with Surgic Pro / Pro2 / XT
and VarioSurg 3 / 4



NSK Eco Aware!



Irrigation tubing x5 (disposable)
CODE Y900079
RRP £66.50



Autoclavable Irrigation Tubing x1
CODE S900055
RRP £55.60

FREE 20:1 IMPLANT CONTRA-ANGLE OR SURGICAL STRAIGHT

WHEN YOU BUY AN ADDITIONAL MOTOR/CORD



SURGICAL STRAIGHTS 1:1

1 + 1 FREE*



- Used with standard HP burs
- **2 year warranty**
(3yrs for iClave users)

• Optic

MODEL Ti-Max X-SG65L
CODE H1009001
RRP £998

• Non-Optic

MODEL Ti-Max X-SG65
CODE H1038001
RRP £681

IMPLANT CONTRA-ANGLES 20:1

1 + 1 FREE*



- Max speed: 2,000min⁻¹
- **2 year warranty**
(3yrs for iClave users)

• Optic

MODEL Ti-Max X-SG20L
CODE C1003001
RRP £1,331

• Non-Optic

MODEL S-Max SG20
CODE C1010001
RRP £1,122

ZYGOMATIC CONTRA-ANGLE


1 + 1 FREE*

- **20:1 Reduction**
- Twist chuck
- **12 Month warranty**
(2yrs for iClave users)

Handpiece for Zygomatic & Pterygoid implants. Its unique angulated shape offers optimal access and visibility.

• Non-Optic

MODEL SGX-E20R
CODE HA1200001
RRP £1,637



*Terms and conditions apply. Offers only available through participating dealers.

Implantology & oral surgery

VarioSurg 4

ULTRASONIC BONE SURGERY SYSTEM

NEW

FREE 12 EXTRA TIPS OR AN OSSEO 100+

With its stylish design and wireless foot control, VarioSurg 4 builds on the reliability and control of previous models to deliver a modern ultrasonic bone surgery system, with Piezo technology to minimise damage to surrounding tissue and bone.

Contains

- Slimline lightweight LED handpiece
- Stylishly designed control box
- Wireless foot control**
- Tip Stand with 6 tips - SG11, SG6D and SG7D for sinus lifts. H-SG1 for cutting bone with SG3 and SG5 for scraping.



MODEL VarioSurg 4
CODES Y1500706 Full Kit (including cordless foot control)
Y1500712 (without foot control)
Y1500712 + Z1102003 (+ corded foot control)
RRP All £8,661 each

NEW

SG12 Tip

The ultra-slim profile of the SG12 has blades on 3 sides, to offer superior access from any angle to the mandibular ramus.

CODE ZT02010001
RRP £196



VarioSurg tips 1 + 1 FREE*



Scan/Click
Extensive range of tips



Scan/Click
For the core range of tips

LINK OPTIONS

Already have a Surgic Pro2?

Link it to a VarioSurg 4 and control both with 1 foot control**



VarioSurg 4 WITHOUT FOOT CONTROL - Y1500712

50% OFF FOR Surgic Pro2 USERS

+ FREE LinkStand 3 OR iCart Duo

Link Stand for VarioSurg 4 / Surgic Pro2

MODEL Link Stand3
CODE ZA16230001
RRP £321



iCart

MODEL iCart Duo
CODE S9090
RRP £1,077



Scan/Click
Watch VarioSurg 4 and
Surgic Pro2 link together



Osseo 100+

OSSEOINTEGRATION MONITORING DEVICE

Osseo 100+ measures implant stability and osseointegration to enhance decisions about when to load the implant. This is invaluable when using protocols with shorter treatment times and treating higher risk patients.

The simple steps to measure the ISQ value allow the implant loading period to be planned in advance.



30% OFF - ONLY £1,926*

MODEL Osseo 100+
CODE Y1004176
RRP £2,752

NSK Eco Aware!
Reusable pegs!

Model: Multipeg
RRP £44.70

Scan/Click



More info &
Multipeg List

“The Osseo device allows us to easily measure implant stability by assessing the increasing ISQ readings during the healing period. We use it routinely on sinus graft cases, when we place simultaneously, and also in cases with low stability at placement.”
The Campbell Clinic, Nottingham

*Terms and conditions apply. Offers only via participating dealers. **Only the Surgic Pro2 foot control (FC-81) is designed to link and control both units.

Decontamination & IPC

AFFORDABLE, RELIABLE & RAPID STERILISATION



iClave 24

NEW

Large B-type vacuum autoclave 24 litre



CODE DX92412B0100
RRP £8,274

iClave plus

B-type vacuum autoclave 18 litre



CODE DXY1003077
RRP £6,637

Sponsored
by NSK



Enter Your
Practice Now

- Perfectly combines performance, user experience, & technology, ensuring an effective, comfortable, and worry-free experience.
- 24L chamber with 5 trays of 2 different dimensions maximises capacity.
- High-efficiency steam generator allows for superior performance
- B-cycle 40min, S-cycle 31min
- User-friendly colour touchscreen
- Footprint: W451 x H492 x D669mm

- Easy to use + USB data storage
- Fast 25min cycle for unwrapped loads
- 35min for pouched loads (+econ cycle 20min – one large tray; 15min– two small trays)
- 18L cylinder-shaped chamber offers 20% more space than traditional chambers & is ideal for implant kits
- Includes Helix Test + Strips
- Footprint: W445mm x H428 x D532mm

FROM ONLY **£4,995**

+ install & delivery

RENT Call for a price

FROM ONLY **£3,995**

+ install & delivery

RENT for just £155pm[†]

12 MONTHS ADDITIONAL WARRANTY on new NSK handpieces or repairs for iClave users**

PROTECT IT

*Visit mynsk.co.uk/rental for more information. **Additional warranty applies to NSK handpieces (new or repaired by NSK), that are still within their warranty period.

All-Inclusive Maintenance and Breakdown Contracts for NSK autoclaves from just **£64pm**.

One annual maintenance visit or every 1500 cycles	Additional maintenance visit for units exceeding 1500 cycles	All service parts and service sundries included	All breakdown calls*	All breakdown parts including labour*	All engineer travel in the UK*	Free on-site operator training	Free technical software upgrades



More info on the
NSK decon range and the
Smart3 System of working

iCare

Handpiece lubrication unit



Precise dosing - less waste, better care

MODEL C2
CODE Y1002787
RRP £2,752

MODEL C3
CODE Y1002788
RRP £2,892

Connects to 6mm airline, 5-6 bar

NSK Eco Aware!

**£2,137 OF FREE
NSK HANDPIECES**

2 × Optic Turbine 1 × Straight 1 × 1:1 Contra-angle

OR HALF PRICE £1,376

RENT for just £73.95pm[†]

CONTACT NSK

National Sales Manager
Mark Beckwith
mbeckwith@nsk-uk.com 07900 246529

Scotland - Sales Manager
Angela Glasgow
aglasgow@nsk-uk.com 07525 911006

N. Ireland & ROI
info@nsk-uk.com +44 7435 840097

North West / North East / North Wales
Clare Morris
cmorris@nsk-uk.com 07901 117669

Derbys/Notts/Yorkshire/Lincs
Jo Tomlinson
jtomlinson@nsk-uk.com 07545 786619

West Midlands
Andy Hargreaves
ahargreaves@nsk-uk.com 07872 824175

Leics / Northants / CV & PE postcodes
Graham Brown
gbrown@nsk-uk.com 07541 864641

South Wales
Robert Musgrove
rmusgrove@nsk-uk.com 07798 918741

South West
Ben Guest
bguest@nsk-uk.com 07350 361543

South Central
Steve Barnes
sbarnes@nsk-uk.com 07889 181099

Central London
Teruyuki Nakanishi
tnakanishi@nsk-uk.com 07570 293317

South East & East Anglia
Francois Faro
ffaro@nsk-uk.com 07900 246822

Field Service Engineer - North
Jamie Duncombe
jduncombe@nsk-uk.com 07780 105457

Field Service Engineer - South
Lee Morris
lmorris@nsk-uk.com 07435 783707

NSK United Kingdom LTD www.mynsk.co.uk

Hertford House, Rutherford Close, Stevenage, Hertfordshire SG1 2EF, UK - info@nsk-uk.com - 0800 6341909



*Offers cannot be combined. NSK reserve the right to change prices and offers without notice. Published prices are RRP's and exclude VAT and installation. Offers only available through participating dealers. Additional warranties available from NSK UK on new products and NSK UK repairs. To qualify for specific offers please send a copy of invoice to NSK UK at Hertford House, Rutherford Close, Stevenage, Hertfordshire SG1 2EF, UK or email sales@nsk-uk.com. Offers end 31/08/25.