

Trycare

May/June Specials

Perfection Plus

Your partner in practice

Buy any 5 ErgoPLUS Instruments, Get Another 1 FREE*



* Cheapest item free

Mirror Handle

Code	Each
USUR003	£6.20

Probes

Single Ended

9.5mm Handles approx. 18g weight

No.	Code	Each
6	USUR004	£13.20
14W	USUR006	£13.20

Perio Screening Probes

9.5mm Handles approx. 18g weight

No.	Code	Each
CPITN-C	USUR005	£13.20
UNC15	USUR040	£13.20

Gracey Curettes

Double Ended Hollow

9.5mm Handles approx. 18g weight

No.	Code	Each
Gracey 1/2	USUR007	£17.90
3/4	USUR008	£17.90
5/6	USUR009	£17.90
7/8	USUR010	£17.90
9/10	USUR011	£17.90
11/12	USUR012	£17.90
13/14	USUR013	£17.90

Scalers

Double Ended Hollow

9.5mm Handles approx. 18g weight

No.	Code	Each
H6/H7	USUR014	£17.90
204S	USUR015	£17.90

Protect+ Ultrasonic Cleaner

- Specially formulated low foaming cleaner with proteolytic enzymes
- Triple enzyme for efficacy
- Non-Hazardous, unperfumed
- Suitable for use with Ultrasonics and Automatic washers in accordance with manufacturer's instructions
- Dilute 5ml per 2 litres warm water



Product	Code	Each
1 Litre Bottle	UPPUC01	£27.60

Protect+ Instrument Cleaning Sachets

- Individual pre-dosed soluble sachets, specially formulated surface-active cleaner
- Rinseable for rapid residue removal
- Removes dirt, dried blood and other organic matter from Instruments, glassware and hard surfaces
- Suitable for use in Ultrasonic
- Unperfumed
- 1 x 17.5g sachet makes up to 3 litre solution



Box 50

Product	Code	Each
17.5g Sachets	UPPIC55	£21.20



Non Latex Dam


Light Blue Synthetic Polyisoprene

- Non-sterile, single use
- High Tensile strength
- 6 x 6 inches (152mm x 152mm)
- Medium 0.18mm - 0.23mm
- Heavy 0.23mm - 0.28mm

Product	Quantity	Code	Each
Medium	20 sheets	UPPNLRM	£14.30
Heavy	20 sheets	UPPNLRH	£14.40
Medium	80 sheets	URDNLMD	£48.95
Heavy	80 sheets	URDNLHV	£51.45

Dental Dam Frames

Product	Code	Each
Plastic	URDF002	£5.30
Metal	URDF001	£8.95



Dental Dam Punch

(Ainsworth Type)

Code	Each
URDF004	£23.35





Dental Dam Clamp Forceps

(Brewer Type)

Code	Each
URDF003	£13.95




Dental Dam Clamps



Other sizes available. Use Dental Dam Forceps for placement of clamps.

Size	Code	Size	Code	Size	Code	Each
00	UPPRD00	7 Wingless	UPPRDW7	13A	UPPRD13	£3.50
1	UPPRD01	8	UPPRD08	14	UPPRD14	
2	UPPRD02	8A	UPPRD8A	14A	UPPRD14A	
2A	UPPRD2A	W8A	UPPRDW8A	26	UPPRD026	
3	UPPRD03	9	UPPRD09	29	UPPRD029	
4	UPPRD04	12A	UPPRD12	202	UPPRD202	
7	UPPRD07					



Dental Dam Template

Stainless steel

Code	Each
UPPDDT1	£16.75



Dental Dam Clamp Sterilisation Tray

(clamps sold separately)

Code	Each
UPPSCST	£12.92

EndoPLUS Alternate Path Files

- Single glide path file with Reciprocating motion enabling easier, safer faster pre-shaping of the root canal
- They are optimised for use with EndoPLUS Alternate files




Size	Pack Size	Code	Each
21mm	Pack of 4	UEPP210	£19.95
25mm	Pack of 4	UEPP250	£19.95
31mm	Pack of 4	UEPP310	£19.95

Sold in non-sterile blister pack



EndoPLUS Alternate Files

- Easy to use Reciprocating system for effective shaping of the root canal
- The premium NiTi alloy and advanced Gold technology enhance the file flexibility and increase the strength
- Superb cutting efficiency means significantly reduced shaping time
- Sold individually in packs of 4 or as an assorted pack of all taper sizes



20	25	35	45
<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>

Each
£19.95

ISO 35
Medium 0.06 Taper

Size	Code
21mm	UEPF212
25mm	UEPF252
31mm	UEPF312

Sold in non-sterile blister pack

ISO 20
Small 0.07 Taper

Size	Code
21mm	UEPF210
25mm	UEPF250
31mm	UEPF310

ISO 25
Primary 0.07 Taper

Size	Code
21mm	UEPF211
25mm	UEPF251
31mm	UEPF311

ISO 20-45

Size	Code
21mm	UEPF214
25mm	UEPF254
31mm	UEPF314

ISO 45
Large 0.05 Taper

Size	Code
21mm	UEPF213
25mm	UEPF253
31mm	UEPF313

EndoPLUS Cold Spray


- 10cm tube ensures precise application onto cotton wool rolls or pellets
- Temperature of -45°C
- Evaporates immediately without residue
- Mint
- 200ml



Code	Each
UEPCS00	£5.95


EndoPLUS Taper Gold

- Advanced gold metallurgy offers higher cutting efficiency and optimum flexibility
- Complete system with shaping files designed specifically for pre-enlarging the canals and finishing files with increased flexibility allowing easier shaping and finishing of the canal
- Variable Taper
- Handle length 13mm
- Sold individually in packs of 6 or assorted.



Shaping

Size	Length	Apical Diameter	Torque	Code	Each
VO	19mm	0.19mm	5.00N/cm	UEPT9V0	£24.50
V1	21mm	0.18mm	5.00N/cm	UEPT1V1	£24.50
V1	25mm	0.18mm	5.00N/cm	UEPT5V1	£24.50
V1	31mm	0.18mm	5.00N/cm	UEPT3V1	£24.50
V2	21mm	0.20mm	1.50N/cm	UEPT1V2	£24.50
V2	25mm	0.20mm	1.50N/cm	UEPT5V2	£24.50
V2	31mm	0.20mm	1.50N/cm	UEPT3V2	£24.50



Finishing

Size	Length	Apical Diameter	Torque	Code	Each
T1	21mm	0.20mm	1.50N/cm	UEPT1T1	£24.50
T1	25mm	0.20mm	1.50N/cm	UEPT5T1	£24.50
T1	31mm	0.20mm	1.50N/cm	UEPT3T1	£24.50
T2	21mm	0.25mm	3.0N/cm	UEPT1T2	£24.50
T2	25mm	0.25mm	3.0N/cm	UEPT5T2	£24.50
T2	31mm	0.25mm	3.0N/cm	UEPT3T2	£24.50
T3	21mm	0.30mm	3.0N/cm	UEPT1T3	£24.50
T3	25mm	0.30mm	3.0N/cm	UEPT5T3	£24.50
T3	31mm	0.30mm	3.0N/cm	UEPT3T3	£24.50
T4	21mm	0.40mm	2.5N/cm	UEPT1T4	£24.50
T4	25mm	0.40mm	2.5N/cm	UEPT5T4	£24.50
T4	31mm	0.40mm	2.5N/cm	UEPT3T4	£24.50
T5	21mm	0.50mm	2.5N/cm	UEPT1T5	£24.50
T5	25mm	0.50mm	2.5N/cm	UEPT5T5	£24.50
T5	31mm	0.50mm	2.5N/cm	UEPT3T5	£24.50

Assorted Pack of 6

Assortment	Code	Each
VO 19mm, V1, V2 & T1, T2, T3 21mm	UEPT921	£24.50
VO 19mm, V1, V2 & T1, T2, T3 31mm	UEPT925	£24.50
VO 19mm, V1, V2 & T1, T2, T3 25mm	UEPT931	£24.50

Protect+ Self-Seal Sterilisation Pouches

Easy to open pouches with triple seal end for extra protection

- Manufactured from 60g medical weight sterile packaging paper and transparent blue PET/PP-plastic composite laminate film
- Both paper and plastic are fully recyclable, refer to local/national recycling regulations
- Colour change confirms sterilisation

Packs of 200

Use	Size	Code	Each	10+
Burs & Endo Instruments	57 x 127mm	UPAF000	£5.30	£4.99
Burs & Endo Holders	85 x 165mm	UPAG000	£5.95	£5.65
Hand Instruments & HP's	85 x 253mm	UPAC000	£5.00	£4.75
Hand Instruments	90 x 230mm	UPSS007	£4.95	N/A
Multiple Hand Instruments	135 x 278mm	CNSM0026	£7.99	£7.50
Multiple Hand Instruments	140 x 260mm	UPSS008	£8.45	N/A
Push Bar Instrument Tray	190 x 330mm	CNSM0027	£17.00	£15.95
Full Instrument Tray	200 x 405mm	UPSS010	£21.25	N/A
Full Instrument Tray	300 x 450mm	UPSS009	£29.20	N/A
Full Instrument Tray	305 x 420mm	USBD000	£29.90	£27.95



For every box of Sterilisation Pouches sold, Perfection Plus will make a donation to SmileStar!



Scan the QR code and find out more about SmileStar



EasiJect Needles



Buy any 5 Packs of EasiJect Needles and get **CARTRIDGE SYRINGE* FREE**



Each	10+
£5.95	£4.95

* Code: USYC002
2.2ml Imperial Thread

Box 100

Imperial Thread	Code
27G x 25mm	UNEJ25S
27G x 35mm	UNEJ25L
30G x 12mm	UNEJ30X
30G x 21mm	UNEJ30S
30G X 25mm	UNEJ30G

- Single Use Sterile Dental Needles
- Silicone coating substantially reduces pain and tissue trauma
- Tri-bevel point for easy tissue penetration without coring
- High tech steel tube, non-toxic, pyrogenic free
- Round cutting edge to protect the gum

Protect+ 3-in-1 Air/Water Disposable Syringe Tips

- Effective replacement for standard metal syringe tips with no conversion adapter required
- Ensures instant dry air and prevents cross contamination



Pack 250

Product	Code	Each	3+
White Interior/Clear Exterior	UPER31W	£16.95	£14.95
Multi-Colour Interior/Clear Exterior	UPER31C	£21.95	£18.95

Scan the QR code to see the full range of Perfection Plus products



Follow Trycare on Social Media



trycareltd



trycaredental