

Surgical Essentials

& much more.

Bone Materials
Draping & Gowning
Surgical Sundries
Surgical Instruments
Irrigation Lines
Sutures
Sterile Gloves

Innovations | Exclusives | Specialities

Discover the Magic of Trycare

CONTENTS

4	BONE GRAFT MATERIALS		46	DRAPING & GOWNING	
	Bone Substitutes	5		Custom Kits	47
	Bone Blocks	8		Off The Shelf Kits	48
	Membranes & Barriers	10		Separates	49
	Allograft (Human)	12		Draping Accessories	51
	Synthetic Bone	13			
	•		52	SUTURES	
14	BONE TACS & SCREWS			Ethicon	52
				Omnia	54
	Bone Screws and Tacs	14			
	Bone Tacs & Trap	15	55	SUNDRIES	
16	TISSUE HARVESTING			Peri-Implantitis	55
10	HOOGE HARVEOTHO			Surgical Essentials	56
	PRF Equipment & Accessories	16		Anaesthetic	58
17	BONE HARVESTING		61	IMPLANTS -	
				Kerator Overdenture Attachments	61
	Bone Handling Instruments	17		DLC Implant Drills	63
	Chisels & Mallets Bone Curettes & Scrapers	. 19 . 20		Drill Stop Kits	64
	bone Corettes & Scrapers	. 20		Implant & Screw Removal Kits	65
21	SINUS GRAFTING		66	ROTARY	
	Kit and Individuals	21		Dana Duna Taankin aa 9 Duilla	//
				Bone Burs, Trephines & Drills Bur Butler & Total Bur Kit	
23	SURGICAL LIGHTING				
	Dr.Kim Hoodlight	23	68	IRRIGATION	
	Dr Kim Headlight	23			
0E	ODAL LIVEIENE			Giving Sets and Suction	68
25	ORAL HYGIENE		69	SUBCICAL FOLIDMENT	
	BioMin Toothpaste	25	07	SURGICAL EQUIPMENT	
	Clinisept+ Mouthwash	26		Sterilizers	69
				Decontamination	71
27	INICTOLIMENTO			Piezo Surgery Units	73
27	INSTRUMENTS			Surgical Motors & Handpieces	75
	Scalpel Blades and Handles	. 27		Osstell	77
	Retractors	. 27		W&H Premium Care	78
	Laschal Instruments	31			
	Needle Holders	35		OTHER CRECIALITIES	
	Tweezers and Pickups	37		OTHER SPECIALITIES	
	Scissors	39		True Light Smartphone Photography	22
	Periotomes	41		Surgical Courses	45
	Luxators	42		3Dalign Clear Aligners	58
	Forceps	43		BlancOne Whitening	59
	Periosteals	44		EMS	60

Osteo Biol by Tecnoss

OsteoBiol® represents a valid alternative to Autografts in most regenerative protocols(1-4)

OsteoBiol® xenografts cleverly retain most of natural bone Collagen/Mineral matrix. They reflect the massive progress in biomaterials technology. Collagen has been shown to stimulate bone regeneration processes⁽⁵⁾. Two experimental studies showed about 50% newly formed bone after one and two months in rabbit femurs⁽⁵⁾ and maxillary sinuses⁽⁶⁾, respectively. Generally, gradual graft resorption and adequate new bone formation^(1,4) can be expected. Scientific literature and clinicians report adequate bone regeneration and high implant success rates⁽¹⁻⁴⁾.

Truly, OsteoBiol® represents the next generation of dual-phase bone regenerative materials.

The OsteoBiol® Catalogue combines product details with technique advice(3). For clinical cases, visit the Osteobiol® website.



Scan to view or download the OsteoBiol® Catalogue

- Correia F et al. Materials (Basel), 2021 Jun;14(12):3439
 Felice P et al. Int J Periodontics Restorative Dent. 2017 Jul/Aug;37(4):469-4803.
- 3. Romasco T et al. J Funct Biomater, 2022 Aug 18;13(3):121
- 4. Correia F et al. Materials (Basel), 2023 Jan 31;16(3):1220
- 5. Falacho RI et al.Molecules, 2021 Mar 2;26(5):1339
- 6. Nannmark U et al. Clin Implant Dent Relat Res. 2008 Dec;10(4):264-70

Bone Substitutes




Fig. 2 Source: Nobil Bio Ricerche, Villafranca d'Asti, Italy | Fig. 3, 4, 6, 7 Author: Prof Ulf Nannmark, University of Goteborg, Sweden | Fig 10 Source: Politecnico di Torino, Italy

Gen-Os®

Mineral & Collagen granules

250-1000 μm **PARTICLE SIZE**

INDICATIONS Universal use

RE-ENTRY 4-5 months (approx.)

Gen-Os® bone granules are an appropriate mix of cortical and cancellous bone and are porcine or equine in origin. The unique production processes of Tecnoss® retain only naturally integrated mineral bone and collagen. The natural collagen stimulates bone regeneration. Gen-Os® is gradually resorbable, replicating autologous bone. Gen-Os® granules absorb blood and saline becoming even more like autologous bone. Pack them tightly. This will squeeze out any excess moisture from around the granules but not the useful blood within them. As well as ensuring the site is packed tight with bloody granules, the site integrity and contour will be maintained. Contraction will be minimal.



Porcine	Equine		
UGOP025	UG0E025	0.25G	£46.00
UGOP050	UG0E050	0.5G	£65.00
UGOP100	UGOE200	1.0G	£102.00
UGOP200	UG0E200	2.0G	£193.00

Apatos[®]

Microcrystalline hydroxyapatite cortico-cancellous & cortical granules

PARTICLE SIZE 600-1000 μm

INDICATIONS Alveolar regeneration, horizontal augmentation, peri-implant defect & sinus grafts

RE-ENTRY 4-5 months (approx.)

Apatos® is slower resorbing than Gen-Os®. There are two versions; one has more cortical bone than Gen-Os® but the other is all cortical bone and is slower still. Both have the same characteristics as Gen-Os® and human bone, making them an excellent alternative. The natural collagen content retained in Apatos® after processing stimulates bone regeneration. It is a universal bone material suitable for use any where you would use granulate bone.

Porcine	Equine		
UAPMP50	UAPME50	0.5G	£61.00
UAPMP10	UAPME10	1.0G	£95.00
UAPMP20	UAPME20	2.0G	£182.00
UAPMP01	UAPME01	1.0G <i>1000-2000</i> µm	£95.00



Cortical				
Porcine				
UAPCP50	0.5G	£58.00		
UAPCP10	1.0G	£91.00		

mp3

Ready to use

90% pre-hydrated dual-phase granules 10% collagen gel

PARTICLE SIZE

600-1000 μm

INDICATIONS

Lateral access sinus lift, alveolar regeneration & horizontal augmentation

RE-ENTRY

5 months (approx.)

mp3® contains pre-hydrated granules and requires no mixing: it can be simply syringed into the defect. Pack in site.

				Equine	Porcine
Buy 3 or more, £82 each	£86.00	Each	0.5CC	UTEC001	UTEC002
Buy 3 or more, £108 eac	£113.00	Each	1.0CC	UMP3E10	UMP3P10
0. (0	£180.00	Each	2.0CC	UMP3E62	UMP3P62
	£180.00	<i>1000-2000</i> μm	2.0CC	UMP3E12	UMP3P12

GTO[®]

The sticky xenograft.

80% pre-hydrated dual-phase granules 20% TSV gel

80% fine micronized dual-phase granules

PARTICLE SIZE

600-1000 μm

INDICATIONS

Peri-implant defects, alveolar regeneration, horizontal augmentation

RE-ENTRY 5 months (approx.)

With all the benefits of dual-phase biomaterials and being thermo sensitive, it is sticky bone which can be shaped and molded. Pack in site.

Porcine	Equine		
UGTOP50	UGTOE50	0.5CC	£105.00
UGTOP20	UGT0E20	2.00CC	£235.00



20% collagen gel

Putty

mp3

Ready to use

Up to 300 μ m PARTICLE SIZE

INDICATIONS

Perimplant defects & horizontal augmentation

RE-ENTRY

4 months (approx.)

Putty is designed to be used in peri implant defects with walls and any small cavities to regenerate bone. It has exceptional malleability and plasticity properties to facilitate this. Putty can be injected around the neck of implants in immediate extraction placement sites to achieve successful long term stability.

Porcine	Equine				
UPTP250	-	0.25CC	Each	£66.00	Buy 3 or more, £63 each
UPTP050	UPTE050	0.5CC	Each	£90.00	Buy 3 or more, £86 each
UPTP100	UPTE100	1.0CC*	Each	£148.00	

60% fine micronized dual-phase granules

ininia mi

Gel 40 Ready to use

PARTICLE SIZE Up to 300 μ m

INDICATIONS

periodontal regeneration & crestal access sinus lift

RE-ENTRY 4 months (approx.)

The addition of 40% type I and III collagen gel quickly stimulates the formation of primary blood clots and quick healing process. The collagen gel component is anti inflammatory and rapidly resorbed. The small syringe tip allows placement in crestal sinus procedures and deep peri implant defects.

Porcine UBGP050 0.5CC *Each* £78.00 Buy 3 or more, **£74 each**




40% collagen gel





Bone Blocks




























Fig. 2, 3, 9 Author: Prof Ulf Nannmark, University of Goteborg, Sweden | Fig 5,7 Source: Politecnico di Torino, Italy.

Sp-Block

Cancellous block for the inlay technique in the mandible





SEM image of OsteoBiol® Sp-Block

Courtesy of Prof Ulf Nannmark, University of Göteborg, Sweden

Sp-Block is a block of heterologous cancellous bone and preserved collagen. The rigid consistency maintains the original graft volume so it is suitable for use in large augmentation procedures and is particularly useful in the sandwich technique to gain ridge height.

RE-ENTRY

8 months (approx.)

Equine

USPE111	10x10x10mm	Each	£113.00
USPE112	10x10x20mm	Each	£143.00
USPE315	35x10x5mm	Each	£126.00

See pg14 for Bone Tacs & Screws...



Membranes & Barriers




Fig. 2, 4, 6 Source: Nobil Bio Ricerche, Villafranca d'Asti, Italy | Fig 7, 9 Source: Politecnico di Torino, Italy.

Evolution Heterologous mesenchymal tissue

HANDLING

It is completely resorbable. Shape with sterile scissors. The membrane should be rehydrated with lukewarm saline; unless the grafting site is already bleeding. Apply to the

INDICATIONS

- site when the desired plasticity is reached. • Peri-implant defect
- · Horizontal augmentation (2 wall defect) Intrabony defects
- Socket preservation
- Lateral sinus lift
- · Vertical augmentation

RE-ENTRY

- 3 months (approx.) fine
- 4 months (approx.) standard
- 2 months (approx.) extra fine



Porcine Standard 0.4mm	Equine Standard 0.4mm	Equine Fine 0.3mm		
UEVP020	UEVE020	-	20x20mm	£96.00
-	-	UEVFE20	20x20mm	£96.00
UEVP030	UEVE030	UEVFE30	30x30mm	£117.00
-	-	UEVFE40	40x40mm	£162.00



Special Extra Fine Membrane

HANDLING

It is completely resorbable. Shape with sterile scissors. The membrane should be rehydrated with lukewarm saline; unless the grafting site is already bleeding. Apply to the site when the desired plasticity is reached.

INDICATIONS

Protects hard & soft tissue grafts & sinus membrane perforations

RESORPTION

40 days (approx.)

Equine

USME222	20x20x0.2mm	£53.00
USME332	30x30x0.2mm	£65.00





Derma Soft Tissue Graft

HANDLING

Derma around the neck of an implant and with good closure will help produce thicker/stronger soft tissue which is important for good maintenance and the life of the implant. Partially exposed mp3® or Gen-Os® can be covered with Derma and the dehiscence will heal over it in approximately 3-4 weeks. During healing, while exposed, it is possible that the membrane thickness be reduced slightly by enzyme activity.

INDICATIONS

For soft tissue augmentation, covering dehiscences and implant sites with thin soft tissue!

RE-ENTRY

3 months (approx.) - fine

5 months (approx.) - standard

		())	,
Porcine Standard 2.0mm	Porcine Fine 0.9mm		
-	UDMF025	25x25mm	£92.00
UDMS050	UDMF050	50x50mm	£167.00
UDMS007	-	7x5mm	£87.00
UDMS015	-	15x5mm	£72.00
UDMS030	_	30x30mm	£92.00



Lamina® Cortical Bone

Resorbable & flexible bone barrier

Flexible cortical bone membranes. Margins are soft in order not to cause micro traumas to the surrounding tissues. Lamina® Soft (Curved) has a semi-rigid consistency and can be grafted without hydration, provided that it is previously shaped to fit the defect morphology.

Can be shaped with scissors. Apart from Lamina® Soft*, it must be hydrated for 5-10 minutes to the desired plasticity. Must be fixated either with tacks or screws. Lamina® Fine can be sutured directly to adjacent tissue.

Lamina® should only be removed if there is clear suprainfection, because its consistency is such as to allow it to achieve a complete second intention healing of the wound.

*Lamina® Soft (Curved) should not be hydrated.

INDICATIONS Lateral access sinus lift & horizontal augmentation (curved wall) in association with mp3[®].

RE-ENTRY 6 months (approx.) - Preformed curved 1.0mm

> 6 months (approx.) - Flat 0.5mm 8 months (approx.) - Flat 1.0mm

Porcine	Equine
ULMCP35	ULMCE35
ULMMP24	-
ULMFP25	-

HANDLING

1.0mm CURVED	35x35mm	£197.00
1.0mm MEDIUM	20x40mm	£197.00
0.5mm FINE	25x25mm	£155.00





Allograft Rocky Mountain

Why Use Rocky Mountain?

The allografts are from human bone donors. The bone is put through vigorous testing and sterilisation techniques to ensure that it is safe to use and is more effective than using synthetic materials. Allograft bone has clinically proven effectiveness and a great track record of safety.

- Cost effective No VAT.
- Re-entry in 6-9 months.
- Clinicians report faster more anticipated outcomes
- Ideal for sinus grafts & ridge expansion.



Cortical, Cancellous Granules

INDICATIONS

Use Rocky Mountain granules for any augmentation procedures

RE-ENTRY

4-6 months (approx.)

Rocky Mountain allograft granules can be used in all areas for augmentation and produce reliable results similar to that of patients own bone.



Cortical

CNCB10	0.5g	£94.00
CNCB11	1.0g	£188.00

Cancellous Fine Grind

Particle size 420-1000 microns

CNCBF01	0.25g	£47.97	
CNCBF02	0.5g	£94.00	

Periogrind Cancellous Particle size 1ml

Particle size Im

CNCBP02	0.25g	£48.05
CNCBP05	0.5g	£95.50
CNCBP1	1.0g	£189.15

Cancellous

Granule size 1-3mm

CNCB02	0.25g	£50.10
CNCB05	0.5g	£95.50
CNCB1	1.0g	£189.15
CNCB2	2.0g	£377.30
CNCB3	3.0g	£573.50

Block Graft



Use for any horizontal augmentation procedures to create more ridge height. Secure in place with bone screws.

6-9 months (approx.)

Rocky Mountain allograft blocks consist of cancellous bone with a cortical edge and producing excellent results in large augmentation procedures.



Block Graft

CNCBB2	5x10x10mm	£154.00
CNCBB3	5x10x20mm	£310.00
CNCBB4	5x10x30mm	£475.00
CNCBB5	10x10x10mm	£230.00
CNCBB6	10x10x20mm	£365.00
CNCBB7	10X10X30mm	£480.00

Rocky Mountain Tissue Bank is a highly respected non-profit organisation providing human graft tissue for dental professionals. Regulated by:- The UK Human Tissue Authority & The American Food & Drug Agency.

Accredited by:- The American Association of Tissue Banks EC ISO 9001:2000 & ISO 13485:2003.

Bond Apatite® in 60 seconds





1. Activate and place. 2. Sculpt and shape. 3. Press for instant set.



Bond Apatite is a synthetic bone graft cement consisting of 2/3 Biphasic Calcium Sulphate and 1/3 slow resorbing Hydroxyapatite granules presented in a 1.0cc ready to use syringe.

Biphasic Calcium Sulphate completely turns over in 4-10 weeks while slow-resorbing Hydroxyapatite maintains the graft volume and shape and gradually resorbs. Bond Apatite is delivered in a unique all-in-one syringe containing the granulated powder and physiological saline. When pushing the plunger, the Bond Apatite mixes with the liquid in the syringe and is then suitable for direct injection into the graft site. Simply shape and mold the material with a spatula and use a sterile gauze swab to press the material to set hard immediately. If full soft tissue closure can be achieved over the Bond Apatite, there is no need for a membrane. The whole procedure can take as little as one minute. Re-Entry - approx 3-4 months.



Augma Shield

Augma Shield™ is a non resorbable intraoral wound dressing. It protects the Bone Graft Cement during the first days, and allows soft tissue to close safely and predictably.

Benefits Include:

- A Non-Resorbable Barrier Using Bone Graft Cement requires a temporary barrier until soft tissue closes and healing begins.
- Additional Layer of Protection Use Augma Shield™ with every augmentation procedure. Make sure to cross-suture Augma Shield™ when it directly protects the graft.

dental@trycare.co.uk

An excellent range of Bone Tacs & Screws...

Perfectly designed for the fixation of membranes and bone blocks.



GBR IV Membrane & Block Kit

UGBRKIT4

Kit Contains:

1x Pilot Drill 1.0mm

1x CA Driver Tip for Handpiece

1x Tacking Tool

1x Driver tip for Handle

1x Autoclaveable case to house all contents

Other sizes of bone screw are available.

5 each Titanium Screws sizes:

Diameter	Length
1.2mm	3mm
1.2mm	4mm
1.4mm	4mm
1.6mm	8mm
1.6mm	10mm
1.8mm	10mm

12 Bone Tacks:

Diameter	Length
2.5mm	4.5mm



GBR III Membrane Tac & Screw Kit

UGBRKIT3

Kit Contains:

1 x Pilot Drill 1.0mm

1x CA Driver Tip for Handpiece

1x Tacking Tool

1x Driver tip for Handle

1x Autoclaveable case to house all contents

Other sizes of bone screw are available.

10 each Titanium Screws sizes:

Diameter	Length
1.2mm	3mm
1.2mm	4mm
1.4mm	4mm

9 each Bone Tacks:

Diameter	Length		
2.5mm	3.5mm		
2.5mm	4.5mm		



Bone Screw Refills

Ø mm	L: mm								
1.2	3	-	-	1.6	4	1.8	4	2.0	4
1.2	4	-	-	1.6	6	1.8	6	2.0	6
1.2	6	1.4	8	1.6	8	1.8	8	2.0	8
1.2	8	1.4	10	1.6	10	1.8	10	2.0	10
-	-	1.4	12	1.6	12	1.8	12	2.0	12

Pack of 5





UTBM032

Bone Tac Gun with 3.5mm Tacs

The automatic bone tack gun enables the user to secure membranes with the touch of a button. It is quick and easy to load the dispenser gun and allows accurate placement directly to the site. Note only works with MCT Bone Tacks. The kit is supplied with 21 bone tacks.



Bone Tac Tool with Tacs

Contains:

21 x Tacs 1x Tacking Tool 3.5mm or 4.5mm



Bone Tac Refills



UBT3521 Length 3.5mm Pack of 5

UBT4521 Length

Osseous Coagulum Trap

An improved method of collecting material from the aspirator tip into a filter. Excellent for retrieving root tips & precious metals. Can be chemiclaved, autoclaved or dry heat sterilised. Connects into a high volume cut-off valve or 1/4" surgical tubing & any surgical aspirator that fits a high volume cut-off valve will fit into the cap.

CGTRP Trap & 3 Filters **CGTRPFLT** Replacement Filter Pk3



Complete PRF Kit

USSOCOOHealthco PRF Complete
System

£975

Our custom built complete PRF kit is supplied with all the items required to produce a membrane from your patients own blood. Using a membrane produced from the patients own blood, a GBR procedure is considered to induce quicker healing and produce overall better results.


- Red Vacutainer Tubes (x100)
- Blood Collection Set (x25)
- Test Tube Rack (24 Holes)
- Centrifuge
- PRF Box




Individual PRF items...












Test Tube Rack

Plastic, 24 holes

Graft Placing & Condensing Instrument

[4](44]

Circumnavigates implants to evenly pack graft material



IP-1 Implant Bone Packer

- Pays for itself by saving allograft
- · Reduces fiddle factor
- Easy retrieval from any receptacle
- Conical shape allows greater site access



Sonic Steel®

contains zero traces of carbon, eliminating potential risk of implant

Hu-Friedy®

Bone Grafting Syringe



Hu-Friedy

Bone Graft Scoop



Ideal for packing OsteoBiol materials

Bone Condenser

Useful instrument for bone compacting with depth markings.

UBCGR1X f47.50



Perfect for OsteoBiol **Evolution**° membrane

Membrane Placement Instrument

Combination spike and curved ends for placement of guided tissue regeneration membrane material. Spike end holds membrane in place.


HDFY0024 f49.50







Ridge Expansion Instruments

Used for cutting or preparing bone.





curved





Hu-Friedy[®]

Surgical Curettes



Miller 10 USCM106

Lucas 2.8mm UHFCL00

Bone Scrapers



Buser 1/3 3mm/5mm UBSB136

£78

Buser 1/2 3mm/4mm UBSB126

Bone Chisel



Buser HDFY0005


6 x 135mm



Sinus Grafting Kit

Complete kit for sinus lift procedures contained in a fully autoclavable tray for ease of use, convenient storage and organisation.






Hu-Friedy[®]

Sinus Instruments - Black Line

Reduces light reflection enhancing visibility



£100

£100

UHSC01X **£98**

UHSC1SX **£98**

UHSLC2X **£98**

UHSLC2SX **£98**





True Light

Professional photography with your smartphone



More accurate shades



Clearer images



Ideal for patient communication






True Light helps you photograph all your clinical cases with the same before & after brightness. With an adaptability to most smartphones, True Light provides an excellent alternative to a professional camera.

- 6000K Daylight Illumination
- Compatible with most Smartphones
- Three banks of LEDs
- Diffusors
- Adjustable Intensity







Autoclavable VisionButler Perfect With True Light

UVBRMRB

The World's Best Shadow Free Headlight

Works with any through-the-lens loupe!


The patented Dr Kim Headlight produces chromatically-correct white light via dual LEDs, extinguishing harsh shadows to create a clear line of vision for efficient working conditions. With an extremely comfortable wireless, lightweight design and adjustable intensity dial, Dr Kim simplifies surgical lighting to effectively reduce strain and increase comfort.





Magnification available

£720








Scientifically proven fluoride & mineral therapy... ...not just toothpaste.

BioMin is unique in it's ability to utilise bioglass technology to deliver superior fluoride protection and sensitivity relief. BioMin adheres to teeth and slowly dissolves releasing calcium, phosphate and fluoride ions over a period of 8-12 hours to form acid resistant fluorapatite. This makes BioMin much more effective than conventional toothpastes where the active ingredients, including soluble fluoride, are washed away and become ineffective less than two hours after brushing.

- BioMin is the only toothpaste to gain accreditation by the Oral Health Foundation for sensitivity relief.
 - **Increased acid resistance** and **reduced sensitivity** by blocking exposed dentine tubules.



UBIOMINF 8 x 75ml tubes

Create extra revenue supplying your patients with the best toothpaste!

> Buy 8 tubes for £43.85 & sell for £10 a tube!

fluoride free toothpaste formulated to adhere to teeth & slowly release mineral ions, only without the Fluoride & much faster to produce a hydroxyapatite coating.



Dental

Goodbye Chlorhexidine

Hello Clinisept+

Clinisept+ is a biomimetic stabilised hypochlorous formulation offering unparalleled antimicrobial action with zero toxicity, achieving safe and effective oral care. It is 100% safe for patients and the environment and 99.99% effective against bacteria, viruses, spores and fungal infections. It cleanses and provides high levels of antimicrobial protection, but is also completely tissue compatible and doesn't harm, irritate or sensitise oral mucosa or dental enamel. Its ability to deliver effective antimicrobial cleansing without using harsh chemicals provides the optimum conditions for gums and mucosa to heal. Clinisept+ is highly recommended for use by those who suffer from mouth ulcers, bleeding gums and periodontitis.

Effective

Clinisept+ Dental Mouthwash removes 99.99% of the bacteria that cause gum disease and tooth decay. Meets EN1040, EN1276 and EN14476. Not subject to antimicrobial resistance.

Gentle

Contains no alcohol, chlorhexidine, artificial colouring or flavour. Non-irritant, non-sensitising, non-cytotoxic and doesn't deliver an 'arctic blast'. Clinically proven safe for use on sensitive mouths and gums.

Aids Recovery

Clinisept+ offers amazing wound healing properties in the oral environment. It is the perfect treatment for any Pre-Op and Post-Op implant, endodontic and surgical procedures to minimise risk of peri-implantitis and aid recovery.

Clinisept+ Skin

Professional skin cleansing prior to surgical procedures and facial aesthetics.

Clinisept+ Dental

Mouthwash

Professional oral cleanser and hygienic rinse. Suitable for those with normal and sensitive teeth, and recommended for use after implants.



Effective but gentle antimicrobial cleansing



Advanced hypochlorous technology



Causes no staining



Maintains soft tissue around implants



Provides optimum conditions for natural healing

400mle


Clinisept+ Dental Mouthwash has an 18-month shelf life and should be used within 3 months of opening.

In cases of

			00000 01	
UCSPMW4	Mouthwash	400ml	x12	£55.80
UCSSK01	Skin	100ml	x12	£78.00
UCSSK00	Skin	490ml	x6	£51.00



Standard Scalpel Handles B3 & 3 suitable for most blades except 20, 22, 23 which fit handle #4.



Standard Carbon Steel Sterile Blades









29



Self-Expanding Photographic Retractor

Orringer Multifunctional Cheek & Lip Spreader

- · Easy access to all parts of the dentition
- Troublesome cheek soft tissue are kept out of the way
 - "Hands Free" retraction
- Easy adjustment for the "front view" dental photo status

TD0230	Small	£45.00
TD0231	Medium	£45.00
TD0232	LargeX	£45.00

VisionButler Retractor

Reusable & Autoclavable!

UVBRMRB

Applications:

- Photography: Hands-free retraction & contrast.
- Orthodontist: Imaging, bond up & debond.
- Restorative: Composites, bonding, veneer placement, implant component placement, lip protection.
- Whitening: Lip protection.
- Cad-Cam: Bite registration & imaging.
- Hygienist: Lip control, spray control.





OptraGate Retractor Disposable

OptraGate latex free lip and cheek retractor is comfortable to wear over longer periods of time, as it has three-dimensional flexibility. OptraGate offers dentists and their assistants easier access to a considerably enlarged treatment field.

UOGJ000	Junior	Pk80	£135.00
UOGS000	Small	Pk80	£126.00
UOGR000	Regular	Pk80	£126.00
UOGASST	Assorted	Pk80	£126.00

Set Price

Suggested Basic Implant Instrument Set



Ridge Mapper STMTRC £90

Buser Periosteal UHLW001 f23.40

Castroviejo Needle Holder **UPPNHCS** f22.35

Minnesota Retractor UCRM155 £5.95

Lucas Curette UHFCL00 £74


Tissue **Tweezers** ETWT12 £14

Round Scalpel Handle USHRDST £12

Boley Gauge BLYG £26.35

Hu-Friedy®

Cheek/Tissue Retractors





01274 88 55 44



Snagless Needle Holders

Snagless® Featherlite Needle Holders with patented Thumlock* strongly grasps any size needle without splaying. Lighter than titanium, sonic steel will not break, loosen or rust. Your Snagless® needle holder will outlast others over and over again! Will not snag when tying sutures.

Ideal for delicate muco-gingival surgery & all suturing for that matter.

- Tungsten Carbide tips
- Thumlock* engages & disengages without fail
- Snagless joint design enabling suture tying with no snagging

Sonic Steel[™] has memory. It allows for a superior thumb-lock that doesn't lose its shape.

ULSI
15.5cm Straight:

£485

18cm Straight: ULSL012 £485

with **Castroviejo Lock**

Featherlite needle holders with **traditional locking mechanism** and **fine tips**, will strongly grasp any size needle without splaying.

- Tungsten Carbide tips
- Snagless Joint for Easy Tying of Suture
- Flexible Sonic Steel Ensures Long Term Durability & Grip of Any Needle

Sonic Steel

is elastic. The tips flex outwardly to hold every size needle without damaging them, or permanently splaying at the tips



SofTouch cut sutures not patients!

- Blades cross 300% more than all other scissors for a cleaner cut
- NO SERRATIONS or the damage they cause
- First they bite, then they cut. Tissue never slides

Trim tissue with confidence and create perfect margins for suturing

Sonic Steel™

fine scissors blades by an additional 300%. But instead of biting into each other, the blades glide together like liquid, reliably cutting all the way to the point and never crushing tissue


> Curved 14.25cm

ULSL017

Straight Blades 14.25cm

ULSL040

£315



Hands-Free Lasner Tissue Retractor

with Patented Thumlock

Unlock with a flick, without losing grip. It never fails!

Hands-free retraction. Non-penetrative gripping **preserves tissue** integrity and reduces bleeding & trauma.

> 17cm 45° Angle ULSL029

£510

18cm Straight ULSL006





Laschal Scissors... **Cut Cleaner!**

Soft Tissue Scissors

You can confidently trim soft tissue, creating perfect margins for suturing. Laschal scissors are made of Flexible Sonic Steel which enables the fine blades to cut cleaner because they are set at angles at least 300% greater than other scissors. Instead of biting into each other the fine blades cross each other gliding together in a fluid action cutting all the way to the tips and never crushing tissue which reduces the patients post operative trauma.

Watch me in



See the Laschal catalogue for the full range.



Classic Needle Holders

Olsen - Hegar			
ENH010	15cm TC	£78.00	
ENH015	17cm TC	£83.00	

Crile Wood		
ENH004	15cm Steel	£12.00

Mathieu		
ENH002	14cm Steel	£19.00
ENH003	17cm Steel	£21.00

Mayo - Heg	ar	
ENH011	16cm TC	£74.00

Gillies		
UVET049	16cm 6" Right	£36.90
UVET004	16cm 6" Left	£9.75

Castroviejo - Flat Handle		
ENH012	14cm TC Straight	£76.00
ENH013	14cm TC Curved	£76.00
UPPNHCS	14cm Steel Straight	£22.35

Castroviejo Gomel - Round Handle			
ENH014	16cm TC Straight	£86.00	
ENH009	18cm TC Straight	£81.00	
ENH006	14cm TC Curved	£78.00	
ENH007	18cm TC Curved	£81.00	







Hu-Friedy[®]

Premium Needle Holders







Classic Tweezers & Pickups

Tissue		
ETWT05	14cm Toothed	£5.00
ETWT03	20cm Toothed	£11.00



London College		
ETWT15	15cm Serrated	£6.00
ETWT14	15cm Locking Serrated	£9.00



Semkin-Taylor		
ETWT08	12.5cm Toothed	£6.00
ETWT13	12.5cm Serrated	£6.00



Adson		
ETWT01	12cm Serrated	£10.00



Adson Micro		
ETWT12	12cm Ultra Fine	£14.00
	Toothed	



Suture (Corn)		
ETWZS01	15cm Curved	£11.00

1	

Debakey Atraumatic		
ETWT10	15cm Straight	£21.00
ETWT11	15cm Curved	£21.00



Gomel Micro Surgical Atraumatic		
ETWT16	15cm Interlocking Toothed	£39.00



Hu-Friedy®

Premium Tweezers & Pickups



Semken-Taylor #33

12.5cm 1x2 Straight TD0934

f 25.50

Adson-Brown #42

12cm 7x7 Multiple Teeth UTPB431





Adson Anatomical

12cm TD2363

£131


DeBakey

15cm **UTPDBPS**

£220

Haemostats

Mosquito #3 Curved **UHMH000**



Mosquito #4 Straight TD1897 £31.80



Classic Scissors

La Grange		
ESC016	11.5cm Steel	£13.00
ESC018	13.5cm Steel	£13.00
ESC017	11.5cm TC	£55.00

Iris		
ESC001	11.5cm Steel Straight	£9.65
ESC007	11.5cm Steel Curved	£12.00
ESC006	11.5cm TC Straight	£53.00
ESC003	11.5cm TC Curved	£54.00

Goldman Fox		
ESC013	12.5cm Steel Straight	£12.00
ESC014	12.5cm Steel Curved	£13.00

Surgical		
ESC004	14cm Steel Sharp/ Blunt	£12.00
ESC005	14cm Steel Sharp/ Sharp	£12.00

Joseph		
ESC002	14cm Steel Curved	£14.00

Castroviejo-Gomel			
ESCM01	Steel Curved	£44.00	
ESCM02	TC Curved	£78.00	

Metzenbaum				
ESC008	12cm Straight Blunt	£13.00		
ESC009	14.5cm Curved Blunt	£14.00		



Hu-Friedy®



Goldman Fox

13cm Straight **UHFSGFS**

£189

Kelly 16cm Curved USK5001

£189



Perma Sharp

11.5cm Curved non-serrated £189



INSTRUMENTS 41 **Periotomes**

Periotome with Interchangeable Tips



Designed for use in the atraumatic extraction of teeth with minimal damage to the surrounding alveolar site.

- Sharp and delicate blades, so the marginal gingiva and periodontal ligament can be carefully detached.
- DALILA handle designed to ensure fast cleaning, easier sterilisation and a more comfortable and secure grip.

Handle with 4 **Blades** MSDY000

£75

Individual **Blades** MSDY004, MSDY001, MSDY002, MSDY003

£9.90



Pointed MSDY004

Straight MSDY001

Angled MSDY002

Angled Short MSDY003

Hu-Friedy®

Periotomes



HDFY0025 £53







Anterior HDFY0026




Dramatically less trauma... ...for you and your patient

Luxator® Original

Luxators are specially designed periodontal ligament knives with fine tapering blades that compress the alveolar bone, cut the membrane and gently ease the tooth from its socket. The whole operation is performed with minimal tissue damage and patient trauma.

Each £88




Luxator Kit of 4 3S, 3C, 5S, 5C **LXTK**

£285

Luxator Kit of 7 2S, 3S, 3C, 5S, 5C, 3CA, 3IC

LXTK6

f460



dental@trycare.co.uk

Hu-Friedy

Atraumatic Forceps

Designed for optimal access and enhanced crown and root grip, each forcep is meticulously crafted to the highest quality standard, ensuring superior performance and longevity.





UFA01XS #1 Standard Upper Incisors

UFAF1XS #1 Apical Upper Incisors

F35XS

UHFFA35 #35 Upper Canines & Premolars



UFA35XS #35 Forceps Atraumair Apical



UFA51XS #51 Upper Roots

£191



UHFFA36 #36 Lower Premolars, Canines & Inicsors



UFA36XS #36 Apical Premolars, Canines & Inicsors



UFA13XS #13 Lower Premolars



UFA79XS #79 Lower Molars



UFA73XS #73 Lower Molars



UFA23XS #23 Cowhorn Lower Molars



UHFF67XS #67 Upper Molars



UFA45XS #45 Lower Roots Serrated




Classic Periosteals



Hu-Friedy[®]

Premium Periosteals





Surgical Courses

Trycare is actively involved in supporting the dental community by hosting and sponsoring a range of professional events, including dental courses, workshops, symposiums, and other educational programs.



Earn CPD Points



UK & International



Early Bird Discounts





Lunch and learn training sessions also available...



Draping & Gowning



BlancOne Whitening



Dr Kim Headlamp



EMS Guided Biofilm Therapy



Kerator Overdenture Attachments



Tokuyama Spherical Restoratives



3Dalign Clear Aligners



BioMin Armour For Teeth



Adin Dental Implants



Bone Family & Surgical Essentials



Clinisept+ Mouthwash




without compromise



Our gowns and drapes are made using the most appropriate, no-compromise materials.

SMS(spunbond, meltblown, spunbond)

SMS is a breathable & repellent material; there is no film layer as this would make it too warm to wear for longer periods of time.



Anciliary Drapes

Our surface drapes and covers have an absorbant layer and a waterproof layer.

2ply Laminate

2ply laminate has a layer of PE film underneath that is impervious to fluids with an absorbant non-woven layer on top that collects fluids.



Design Your Own Kit



Scan the QR code

Fill in the order form and send to dental@trycare.co.uk for a personalised quote!



(A) Nurse Cap



(B) Aperture Drape



(C) Gauze Squares



(D) Tegaderm Site Dressing



(E) Hand Towel



F Surgeon Cap



(G) Hose Cover



(H) Plastic Tweezers



(I) Surgeon Gown



(J) Hygovac Aspirator Tip



(K) Adhesive Strips



(L) Plain Patient Drape



M Surface Drape



N Hygosurge Suction Tip





(P) Motor Unit Cover



(Q) Incineration Bag



(R) Mayo Stand Cover



S Adhesive Handle Cover



(T) Split Patient Drape



(U) Utility Drape

Renowned Professional Kits

From the biggest supplier to implant practices in the UK, these kits derive from the long term needs and experiences of Practices and Surgeons.

Surgikit

f**43**.50



- 1 Patient Cap (T)
- Aperture Patient Drape 160 x 210cm (B)
- **4** Gauze Squares (C)
- Tegaderm Site Dressing 15 x 20cm (D)
- 6 Hand Towels (E)
- Surgeons Cap (F)
- 3 Hose Covers (G)
- 3 Surgeons Gowns Large (I) 2 Nurses Hats (A)
- Plastic Tweezers (H)
- Hygovac Aspirator Tip (J)
- Adhesive Strips (S)
- Surface Drapes Adhesive 75 x 75cm (M)
- Instrument Tray (0)
- Motor Unit Cover (P)
- Incineration Bag (Q)
- 4 Handle Covers Adhesive (S)

Surgitwo DP0002

f**42**.50

- Patient Cap (T)
- Nurses Hats (A)
- Aperture Patient Drape 160 x 210cm (B)
- 8 Gauze Squares (C)
- Tegaderm Site Dressing 1 15 x 20cm (D)
- Hand Towels (E)
- Surgeons Cap (F)
- 3 Hose Covers (G)
- Plastic Tweezers (H)

- Hygovac Aspirator Tip (J)
- Adhesive Strips (S)
- Surface Drapes Adhesive 90 x 75cm (M)
- Hygosurge Suction Tip (N)
- Instrument Tray (0)
- Motor Unit Cover (P)
- Incineration Bag (Q)
- 1 Mayo Stand Cover 50 x

60cm (R)

2 Surgeons Gowns Large (I) 4 Handle Covers - Adhesive (S)

Endokit

DP0003

£**32.**15

- 1 Aperture Drape 160 x 210cm (B)
- 4 Hand Towels (E)
- 4 Hose Covers (G)
- 2 Surgeons Gown Large (I)
- 1 Hygosurge Suction Tip (N)
- 3 Surface Drape Adhesive 95 x 80cm (M)



- Hygovac Aspirator Tip (J)
- Instrument Tray (0)
- Incineration Bag (Q)
- Handle Covers Adhesive (S)
- Adhesive Strips (K)

Minimal **Implant**

DRPKT1

£18.05



- Patient Aperture Drape 150 x 200cm (B)
- Hose Cover 8 x 120cm (G)
- Adhesive Surface Drape 75 x 90cm (M)
- Handle Covers (S)
- Utility Drapes 90 x 75cm (V)

Universal Drape Kit

DC0001

- 2 Adhesive Strips (K)
- **2** Hose Cover 8 x 120cm (G)
- 2 Surface Drapes 75 x 75cm (M)
- 2 Handle Covers (S)
- 1 Split Patient Drape 57.5 x 80cm (U)



Dental Apicectomy Kit

DRPDDS

Buy x14 for

f16.80 OR f13.60 each



- Surface Drapes 175 X 100cm (M)
- Adhesive Strips 9 x 24.5cm (K)
- Utility Drapes 45 x 75cm (V)
- Hose Covers 8 x 120cm (S)



Economy Drape Kit

£5.35

- **2** Hose Cover *10 x 120cm* (G)
- Surface Drapes 75 x 75cm (M)
- 2 Handle Covers (S)
- 1 Plain Patient Drape 58 x 80cm (L)
- 2 Adhesive Strips (K)



Premium Sterile Separates

First for comfort, first for price

Patient Aperture Drape

Adhesive

cm Each

DC0022 125 X 95 £4.80 DC0007 160 X 210 £7.60



Adhesive Strips

Each cm

DRPDAS 49 X 9 £1.15



Mayo Stand Cover

X25 Each **DRPDMY** 145 X 79 £4.40 £3.10



Patient Split Drape

Each X20 DC0020 190 X 260 £12.80 £9.45

Motor Unit Cover

	cm	Each	X20	
DC0026	76 X 79	£2.10	£1.90	
DC0003	60 X 60	£1.95	-	

Patient Drape

Non-adhesive X20 Each cm DC0017 150 X 95 £2.95 £2.35

Handle Cover

	cm	Each	X40
DC0005	4 PACK	£2.20	-
DRPDHNDC	15 X 20 2 PACK	£1.10	£0.90

Surface Drape

Adhesive

	cm	Each	X20	- Constant
DC0024	50 X 47.5	£1.25	£1.15	
DC0023	75 X 75	£1.75	£1.60	

Hose Cover

	cm	Each	X20	
DC0027	10 X 122	£1.80	£1.60	R. C. C.
DC0004	10 X 120 3 PACK	£1.90	-	

Surface Drape

Non-Adhesive PK2

X20 cm Each DC0002 £4.45 70 X 75 £3.70

Clinical Waste Bag

PK25 cm DB610X765 15 X 28 X 39 £5.90



Economy Sterile Separates

Good quality, great price

Patient Aperture Drape

Adhesive

GDN013

cm Each 120 X 150 £3.80



Utility Drape

Adhesive

cm Each **GDN005** 90 X 150 £2.95 GDN006 160 X 175 £4.60



Patient Split Drape

Adhesive

Each cm GDN002 150 X 190 £6.90



Non-adhesive

GDN008 75 X 75 £1.05 **GDN009** 75 X 90 £1.20 **GDN011** 160 X 175 £4.05



Surface Drape

Adhesive

GDN004

Each cm 75 X 90 £1.25



Mayo Stand Cover

Each X20 **DRPDMY** 145 X 79 £4.40 £3.10



Surgical Gowns



All gowns include two hand towels.

Economy

Sterile

	Size	Each
GDN018	MEDIUM	£8.60
GDN017	LARGE	£8.60
GDN018	X-LARGE	£8.60

Premium

Sterile

	Size	Each	X24
GWNDM	MEDIUM	£10.65	£7.60
GWNDL	LARGE	£10.95	£7.60
GWNDXL	X-LARGE	£11.40	£8.40
GWNDXXL	XXL	£11.70	£8.80

Surgery Caps

Non-Sterile



Each **CPDASP** ₽**Ո.40**

PK150 **CPDA** f52.60



Non-Sterile

Each **MSKSESP** £0.30

Surgery Masks

Pk50 **MSKSE** f12.15

Annie



Each **CPDKSP** £**0**.75

Pk100 **CPDK** £46.35

Tie Band

Ear Loop

Each **MSKSTSP £1.60**

Pk60 **MSKST**

Kosack

Gloves

Nitrile

Nitrile powder-free sterile gloves

Sizes

PK50

S, M, L

£31.90



Gallipots

Stainless Steel

GLP4	57ml	£4.95
GLP6	170ml	£6.50
GLP5	284ml	£7.50



Biogel

Coated powder-free latex sterile gloves

Sizes

PK10

5.5 - 9.0

£19.10

Sizes: 5.5, 6.0, 6.5, 7.0, 7.5, 8.0, 8.5, 9.0



Disposable Plastic

GLPD	60ml	Each	£0.65
GALLIPOT60	60ml	X20	£8.60
GLPD1	120ml	Each	£0.90
GALLIPOT120	120ml	X12	£10.80



Towel Clamps











Non-penetrating towel forcep;

especially suitable for disposable drapes.



Ethicon Labelling Guide





Cutting Needle

Triangular in cross section, with the apical cutting edge on the outside of the curve improves strength and increases resistance to bending.



Triangular in cross section with the apex on the inside of the curve. The effective cutting edges are at front section and merge into a triangulated body which continues for half the length of the needle.

Cross geometry reduces cutting edge angles to give improved penetration and control. A square shank increases strength and

improved grip in the needle holder. Available as conventional or reverse cut.



Prime

Needle



Combines the initial penetration of a cutting needle with minimised trauma of a round bodied needle.

Mersilk Braided, Natural Material

Non-absorbable suture. Will require removal. Provides prolonged tensile strength retention in tissue for up to 3 months.



3/8 Circle Conventional Cut		Needle Size	Suture Length	Colour	Gauge Size	Pack	Code	£
	FS-2	19mm	45cm	black	4.0	36	W329H	£82.00
	CFS-3	16mm	75cm	black	4.0	36	W501H	£96.00
1/2 Circle Conventional Cut								
	CX-1	22mm	45cm	black	3.0	36	W577H	£83.00

Coated Vicryl Rapide Faster absorbing synthetic material.

Identical composition and ease of handling as coated vicryl. Excellent knot snug down. Absorption complete by 42 days.



3/8 Circle Reverse Cut P &	PS	Needle Size	Suture Length	Colour	Gauge Size	Pack	Code	£
	P-1	12mm	45cm	undyed	6.0	12	W9913	£91.00
	P-1	11mm	45cm	undyed	5.0	12	W9915	£88.00
	PS-2	19mm	75cm	undyed	4.0	12	W9922	£82.00
	PS	26mm	75cm	undyed	3.0	12	W9932	£79.00
3/8 Circle Reverse Cut								
	FS-2	19mm	75cm	undyed	4.0	12	W9930	£63.00
	FS-2	19mm	75cm	undyed	3.0	12	W9931	£59.00
3/8 Circle Conventional Cu	t - PC							
	PC-3	16mm	75cm	undyed	4.0	12	W9918	£86.00
1/2 Circle Conventional Cut	t-PC							
\ J	PC-22	22mm	45cm	undyed	4.0	12	W9926	£83.00
	PC-22	22mm	45cm	undyed	3.0	12	W9927	£68.00
3/8 Circle Conventional Cu	t							
	CFS-3	16mm	75cm	undyed	4.0	12	W9924	£61.00
	CFS-3	16mm	75cm	undyed	3.0	12	W9925	£59.00
1/2 Circle Conventional Cut	t							
	CX-1	22mm	45cm	undyed	3.0	12	W9935	£57.00
5/8 Circle Conventional Cu	t	Needle Size	Suture Length	Colour	Gauge Size	Pack	Code	£
	UX-25	26mm	75cm	undyed	3.0	12	W9929	£75.00

Compound Curved Reverse								
	PS-4C	16mm	45cm	undyed	5.0	12	W9910	£83.00
	PS-4C	16mm	45cm	undyed	4.0	12	W9911	£83.00
	PS	26mm	75cm	undyed	2.0	12	W9933	£77.00

Coated Vicryl Braided, synthetic, absorbable suture.

Excellent handling and knotting. Unique coating ensures minimal drag and sawing. Minimal tissue reaction. Mass absorption complete between 56 and 70 days.



3/8 Circle Reverse Cut - PC	C	Needle Size	Suture Length	Colour	Gauge Size	Pack	Code	£
$\overline{}$	P-1	11mm	45cm	undyed	6.0	24	W9500T	£134.00
1/2 Circle Conventional Cu	t							
	CX-1	22mm	45cm	voilet	4.0	12	W9471	£56.00
Compound Curved Reverse	e Cut							
L 1	PS-4C	16mm	45cm	undyed	5.0	12	W9436	£63.00
	PS-4C	16mm	45cm	undyed	4.0	12	W9437	£61.00
3/8 Circle Tapercut								
	V4	17mm	75cm	undyed	5.0	12	W9824	£57.00

Omnia PTFE Monofilament Sutures

Omnia Surgical Sutures are manufactured with PTFE (polytetrafluoroethylene) monofilament material. Omnia PTFE sutures are ideal for any implant, periodontal and bone graft surgery where the usage of a monofilament suture with low bacterial adhesion is recommended.




3/8 Circle Reverse Cut	Needle Size	Suture Length	Gauge Size	Pack	Code	£
	19mm	75cm	3.0	12	U0P3163	£111.00
	18mm	75cm	5.0	12	U0P5758	£111.00
	16mm	75cm	4.0	12	U0P4163	£111.00
1/2 Circle Oval Flat						
	12mm	75cm	5.0	12	U0P5131	£111.00

Titanium Brushes for Peri-implantitis

NiTi Brush for treatment of peri-implantitis is intended for mechanical debridgement and cleaning of titanium dental implant surfaces contaminated by plaque, calculus and infection resulting from peri-implantitis.





Omega NiTi Brush

Best access to every area, even if the space is limited.

UNIBOLG UNIBOLS

Standard (40mm long) Short (35mm long)

Pack of 4



Pocket NiTi Brush

Regular brush for normal peri-implantitis treatment.

UNIBTRG UNIBRSH

Standard (40mm long) Short (35mm long)

Pack of 4 f113.00




Nano NiTi Brush

Narrow brush for very precise and accurate operation. To treat initial stage of peri-implantitis.

UNIBTNR Standard (40mm long)


Pack of 4 **13**.00



Standard 40mm

Essentials & Accessories

Sterile Saline 0.9%





Sterile Water

WTRFRNJ01	5ml Ampoules	Pk10	£8.90
WTRRGST01	1 Litre Bottle	Pk1	£6.55
CNSM0064	1 Litre Bag	Each	£7.10

Sterile Topper

Swabs

Toppor o	O.G.I.G
UJJTS01	5x5cm 4ply 5s
UJJTS02	7.5x7.5cm 4ply 5

Topper 8 -	Sterile		-	
UJJTS01	5x5cm 4ply 5s	Pk25	£13.75	
UJJTS02	7.5x7.5cm 4ply 5s	Pk25	£9.70	
SW8105	10x10cm 4ply 5s	Pk40	£23.65	

Topper 12 - Sterile

SWBS10	10x10cm 6ply 5s	Pk30	£19.30
SWBS54	5x5cm	Pk45	£9.10
SWABST750	7.5x7.5cm	Pk25	£4.65

Rooko Suraitin

NOCK	Joong	lup		The same of
URS1000	1.2mm	Pk20	£45.50	
URS2000	2.8mm	Pk20	£41.50	
URS4000	4.8mm	Pk20	£37.10	
URNT000	4.8mm	Pk100	£63.00	

Disposable Scrub Brush





dental@trycare.co.uk

See P26 for Clinisept+

Pre-op & Post-op Mouthwash



Instant Cold Pack

Easy to use disposable instant cold pack. Ideal to use after treatment.



Pk60

£45.00



Normasol Sodium Chloride

Sterile topical irrigation solution. Contains 0.9% sodium chloride. For use on eyes, burns and wounds.



25ml Sachet Pk25 £13.90



Videne

Contains Povidone iodine

VDN01	Antiseptic Solution 250m
VDN02	Surgical Scrub 500ml
VDN03	Alcohol Solution 500ml

£20.50 £20.50 £15.00



Irripod Sodium Chloride

Can be used to wash out wounds thoroughly, rinse the mouth and remove objects and harmful substances from the eye.



20ml

Pk25 £13.40



Hydrogen Peroxide 12% Solution




01274 88 55 44

Debridement flushing liquid.



Reso-Pac Adhesive Perio Dressing

Ideal for the protection of wounds and suture sites. It stavs soft and is odourless and neutral in taste. Keeps out bacteria and remains for upto 30 hours. Self dissolving. Apply directly to wounds.





Hemocollagene

Hemostatic Sponges

Individually blister packed sterile sponges with collagen of bovine origin to promote local haemostasis. Hemocollagene is totally resorbable to help natural tissue regeneration.



Pk24

£45.95



Kaltostat

Calcium Sodium Alginate Dressing

Soft, white, non-woven dressings consisting of calcium alginate fibre. For moderately to highly exuding chronic and acute wounds or wounds with minor bleeding. When applied, the dressings convert exudates to a gel-like, fibrous mat creating a micro-environment that facilitates healing.



Pk10 £20.30

Sterilisation Paper

Pk250 £89.00





Autoclave Indicator Tape

UAUT000

18mm

Roll £5.20

Dentotemp Temporary Cement

Advantages

- Non eugenol No interferences with permanent
- High bonding strength, but with easy removal
- Leaves perfectly healthy gingiva after removal of the temporary crown
- Rebasing and cementation in one stage
- Can be used on vital teeth
- Thin film Very aesthetic result
- Radiopaque

Indications

- For temporary cementation in case of:
 - A small retention area
 - A long-term trial
- Permanent cementation of implant-retained crowns

UDTTA01	Automix	5ml	£30.60
UDTTA02	Automix	2x5ml	£59.00
UDTTS00	Syringe	2x10ml	£50.00





Anestadent Anaesthetic

Articaine hydrochloride/epinephrine tartrate

Anestadent is indicated for local infiltrative or conductive anaesthesia in both simple and complex dental procedures. This injectable anaesthetic includes a vasoconstrictor and complements the Dentsply Sirona range of anaesthetics.

Features & Benefits:

- Free from natural rubber latex
- Onset: Less than 10 minutes
- Duration: Approx. 1 hour for infiltrations, 2 hours for blocks
- Available in boxes of 50 x 2.2ml cartridges
- Each set of 10 cartridges is blister packed to avoid breakage





Horizontal and vertical movement!



Treatment options from £546.25

- Receive ongoing support & training.
- Horizontal AND vertical movement.
- Sculpted around teeth to reduce irritation.
- Treatment possible in 10 weeks.
- Almost invisible!

SCAN ME FOR MORE INFO!







Gentle to the enamel and the gingival tissues, BlancOne offer the only cosmetic in-office treatments which do not require gingival barriers, desensitising agents or remineralising agents.


Offer your patients an affordable, comfortable and pleasant experience which can be repeated often or combined with home treatments.





O1 ASSESSMENT AND INFECTION CONTROL

ASSESS EVERY CLINICAL CASE AND







AIRFLOW

- Clinically proven technology
- **Comfortable** patient experience
- Gold standard implant maintenance
- Swiss precision

Find out how comfortable...


Book your in practice demonstration 01274 88 55 44



Overdenture Attachment Systems

Whether Endodontically Treated Teeth or Implants...

KERATOR straight &



angled abutments for all implants!



Nobel Biocare, 3i, Astra, Straumann, Ankylos, BioHorizon, Megagen, & many others!



Kerator overdenture abutments are available straight (six heights) or 15 degree angle (two heights) and can accommodate up to 40 degrees implant displacement, achieved using different retentive inserts. All abutment packs contain all items for placement and fitting.

Straight Kit

Angled Kit

£86.00 £93.00

Retentive Inserts - Locator Interchangeable

For Angled Abutments:

Red Extra Light 0.8 - 1.0 lbs

Orange Light 1.8 - 2.0 lbs

Green Angled 2.8 - 3.0 lbs

Yellow Zero Retention

For Straight Abutments:

Blue 1.2 lbs

Pink 2.4 lbs

White 4 lbs

Black Zero Retention

Pack of 4



Drivers - Locator Interchangeable

KERAIOR	Tip Shape $\ \ \ \ \ \ \ \ \ \ \ \ \ $	UKLDTIP	Handpiece Driver Tip for use with Straight Kerator Abutments	Each	£36.30
	Tip Shape $igwedge$	UKTTS49	Wrench Torque Tip 4.9 Square Head for use with Straight Kerator Abutments	Each	£36.30
	Tip Shape \triangle	UKTTH63	Wrench Torque Tip 6.3 Hex Head for use with Straight Kerator Abutments	Each	£36.30
	Tip Shape $igwedge$	UKTTRND	Wrench Torque Tip Round Head for use with Straight Kerator Abutments	Each	£36.30
	Tip Shape	RS6080	Hex 1.27mm Short 10mm for use with Angled Kerator Abutments	Each	£54.00
	Tip Shape	RS6082	Hex 1.27mm Long 15mm for use with Angled Kerator Abutments	Each	£54.00
=	-	UKMAGTL	Magic Tool for insert fitting	Each	£58.00

For 'Straight Kerator Abutments' use Triangular Tipped Drivers.

For 'Angled Kerator Abutments' use the Hex Tipped Driver.

Accessories - Locator Interchangeable



UKAMH00	Metal Housing	Each	£13.00
UKWHSPR	Spacer Ring	Pack of 10	£3.60
UKIMPCO	Kerator Impression Coping	Each	£9.60
UKANALG	Kerator Analogue	Each	£9.60



Endodontically Treated Roots Complete Kerator Kit

Kerator Root is as simple as cementing a root post. The kit includes four straight kerator root posts, drills for root placement, the metal housing and a selection of retentive inserts.

UKRKITO Root Kit

Refill Kit Straight or Angled 10° & 20°

includes metal housing & inserts



DLC Drills

Last a lot, lot longer!

Natural diamonds exist in two crystal structures: common cubic lattice and rare hexagonal (lonsdaleite). Combining both forms creates a flexible, amorphous surface coating on drills, appearing smooth but microscopically textured like cobblestones. Pure DLC is as hard as diamond and slicker than Teflon, making DLC drills last longer than standard ones.

Tri-Step for one stage Osteotomies

Avoid numerous drill insertions & minimise invasiveness. Combines pilot tip drill with conventional two-step drill.

Tri-Step Drills 15mm Length £82.00 22mm Length £102.00

	1D1511
Dia: 1	1.8mm
Dia: 2	2.35mm
Dia: 3	3.2mm

	TD2218
Dia: 1	1.8mm
Dia: 2	2.35mm
Dia: 3	3.2mm

Two Step Drills

Short Twist Step Drills (Two step)

13mm Length

	DL1320	DL1328	DL1332	DL1336	DL1342	DL1346	DL1352	DL1356
Dia: 1	Pilot 2mm	2.3mm	2.8mm	3.2mm	3.8mm	4.2mm	4.8mm	5.2mm
Dia: 2		2.8mm	3.2mm	3.6mm	4.2mm	4.6mm	5.2mm	5.6mm

Long Twist Step Drills (Two step)

18mm Length

	DL1820	DL1828	DL1832	DL1836	DL1842	DL1846	DL1852	DL1856
Dia: 1	Pilot 2mm	2.3mm	2.8mm	3.2mm	3.8mm	4.2mm	4.8mm	5.2mm
Dia: 2		2.8mm	3.2mm	3.6mm	4.2mm	4.6mm	5.2mm	5.6mm



DLC Drill & Stop Kits


Control of Osteotomy depth at fantastic prices

A sure, quick and easy aid to ensuring the required depth without the tedium of referring to drill markings. Have full and easy control of the Osteotomy depth by quickly and simply snapping on a depth stop.

NOTE: To have drill stops for all diameters of drills, you must have both kits in long or short. Any customers using both long and short drills will need all four kits.

18mm Long Drill Stops


Drill Stop Lengths Drill Diameter 8 10 11.5 13 15 18 mm Ø 2mm Ø 2.8mm Ø 3.2mm Ø 3.6mm




S1001 / S1002 18mm Long Drill Stops

£**550**.00

13mm Short Drill Stops





S1003 / S1004 13mm Long Drill Stops

£**340**.00

Drill

Diameter

Ø 4.2mm

Ø 5.2mm

Ø 5.6mm

Fixture Removal Kit

Remove Implants Quickly & Easily with No Bone Loss!

Simply screw the fixture removal tip into any implant until tight, then use the hand adaptor or remover wrench to unscrew the implant from site. If it is not possible to remove the implant by this method then choose the correct diameter trephine for the implant diameter.



URKI000

Fixture Removal Kit

£**295**.00

Contains:

7 x Assorted trephines 3 x Fixture remover (1 of each size)

1 x Hand wrench

1x Wrench



Small

for 3-4mm diameter implants URTI003



Medium or 4-5mm

for 4-5mm diameter implants URTI001



Replacement Tips

£48.00

Large

for 5-8mm diameter implants URTI002

Abutment Screw Removal Kit

Remove Broken Screws

Select and locate drill guide that fits the head of the implant. Using screw remover drill through the guide, drill a channel in the broken screw. Insert the screw remover tip into the channel and rotate anticlockwise to remove the offending screw.



USRK000

Screw Removal Kit

£280.00

Contains:

5 x Assorted abutment guides

2 x Screw remover drills

1 x Tap drill

1x Wrench adaptor

1x Guide adaptor

1 Remover wrench

Fixture & Screw Removal Kit

Both kits in one!

Contains:

7 x Assorted trephines

3 x Fixture remover (1 of each size)

1x Hand wrench

1x Wrench

5 x Assorted abutment guides

2 x Screw remover drills

1 x Tap drill

1x Wrench adaptor

1x Guide adaptor

URKIS00

Fixture & Screw Removal Kit

£530.00



Bone Burs, Trephines, Pilot & Implant Drills




Lindemann Bone Burs






C162A 016	shank		
EDNT0085	HP	Pk3	£72.00
EDNT0102	RAXL	Pk3	£72.00



TC Round Burs








TC Tapered Fissure Burs



HP C33.104 No.009 No.010 No.012 No.016 Pack of 5 **£17.80**

Implant Removal Trephine Burs



TRPB374	3.7/4.4	
TRPB396	3.9/4.6	
TRPB430	4.3/5.0	
TRPB474	4.7/5.4	£27.60
TRPB529	5.2/5.9	
TRPB574	5.7/6.4	
TRPB629	6.2/6.9	

Ever dropped burs?
Bur Butler hasn't.

The BurButler is a **autoclavable silicone bur block** developed by dentists. Little pressure is needed to secure your burs safely & take them out again. Burs aren't dislodged; even in the washer disinfector!

- Tapered "soft" holes grip burs.
- · Snug fitting lid goes on & off easily with one hand.
- Contents are clearly visible through the lid.



All sizes available in 5 colours: Orange, White, Pink, Purple, Blue



5 Hole with Standard Lid

£19.40



10 Hole with Standard Lid

£**29**.20



25 Hole with Standard Lid

£34.10



25 Hole with Tall Lid

£46.00



Bone Expansion & Sinus Lifts **Made Easy!**



Total Bur Kit


The HaeNaem Total Bur Kit contains burs for **bone expansion** and **sinus treatments** to **increase bone depth** with **zero bone loss** due to the design which compacts bone. It also includes a full set of depth stops for the burs, ensuring **safe osteotomy preparation** to the **perfect depth** every time. Compacted bone around implants at placement is proven to offer better long term success.

Expansion Burs

The specially designed bone expansion burs compact bone to the sides of the osteotomy instead of cutting and ejecting bone, making bone expansion safe and easy.

Sinus Burs

The Sinus burs push the bone to the apex of the osteotomy extending the depth of bone without risk of perforating the sinus membrane.




Total Bur Kit £1200

Giving Sets & Suction

WSH Implant Motor Complete Giving Set

MPLNMDT5	Implantmed	PK6	£77.00
LCMDT05	Elcomed	PK6	£98.00



PK10 £62.00 MPLNMT1



W&H) 'Y' Connector

Plastic autoclavable Y connector for splitting irrigation tubing for internal and external irrigation through the handpiece.

MPLNMDTY3 £47.50



Omnia Complete Giving Set

To Fit Implantmed

OMN001	PK10	£65.00				
То	To Fit Elcomed					
OMN002	PK10	£65.00				



Yankaur Suction Tube

Tubing

YNKSTSP	Each	£2.05
YNKST	PK25	£46.50



Suction Tips

•			
YNKSTFSP	Fine ID 3.0mm	Each	£1.30
YNKSTF	Fine ID 3.0mm	PK50	£53.00
YNKSTMSP	Med ID 5.0mm	Each	£1.05
YNKSTM	Med ID 5.0mm	PK50	£41.95
YNKSTS	Suction Set	Each	£3.70

Nouvag Universal Giving Set

NVGT01	2 Way	PK10	£75.75
NVGT01SP	2 Way	Each	£8.60



NSK Surgi XT & Pro Complete Giving Set

Also used for NSK Variosurg bone surgery kit.

UNSKIT PK5 £64.30



Satelec Piezotome

for Bone Surgery



SAT0027	Autoclavable Irrigation Line	Each	£104.00
FDS245	Disposable Sterile Irrigation Line	PK10	£121.00
SAT0026	Autoclavable Irrigation Line (with 10 perforators)	Each	£138.00
SAT0028	Disposable Sterile Perforators	PK10	£34.50

Bien Air Implant Motor

Irrigation Line

Chiropro-L Irrigation Line

UBAILCL	PK10	£49.00		
Omnia Alternative				

OMN012	PK10	£65.00		



MI SCILICI IV	Cyci	
UKMIL10	PK10	£86.00

WSH>'L'Connector

Irrigation tubing L connector fits into the back of the handpiece head for internal irrigation purposes.

MPLNMDTL	Each	£35.10
----------	------	--------



Kavo INTRAsurg Implant Motor Irrigation Line Giving Set

 Б.	 	 	Б.	



Sterilizers

All W&H Sterilizers offer:



Intuitive user interface and touchscreen display for a smoother user experience

Wi-fi connection and built in data logger







2 Year Warranty or 2,000 Cycles (whichever occurs earlier)





Pre-programmable vacuum & helix test* function to run automatically out of surgery hours



Automatic Water Fill function

*for B-Type only

Lyla

S-Type Sterilizer



UWHL17S

UWHL21S

£3595

The **new** sterilizer concept



Capacity - 17 or 22 litre available

- Cycle Time from 23 minutes for an average 2kg
- **Traceability** built-in USB port, option cycle report printer and barcode label printer enable complete documentation.
- Mobile App available via upgrade for remote data storage
- Future-proof upgradability to a type B sterilizer using the new Activation Code concept
- **Activation Codes** Unique upgrade ability
 - B Cycle
 - Remote data storage
 - Traceability
 - All in One

Lara

B-Type Sterilizer

Designed for all load types with S Type capabilities

Equipped with fast cycle times and a secure data system, Lara simplifies and speeds up the sterilisation process. Lara offers the new activation code system allowing users to easily upgrade additional features.

- Eco Dry technology automatically adapts the drying time to the mass of the load
- Activation Codes Unique upgrade ability

- Fast Cycle - Performance - Traceability - All in One

Mobile App - for Real Time Remote Monitoring



ULR17BS £4495

ULR22BS f4995

Lara XL available ULR28BS £5295



28 litre chamber able to reprocess up to 7 kg wrapped instruments or 9 kg containers

Lisa B-Type Sterilizer

- EliTrace Full extended traceability
- State of the Art Vacuum B Sterilizer Technology for efficient, reliable, time saving, intuitive and ease of use
- Eco Dry technology automatically adapts the drying time to the mass of the load
- Elisense user interface generates status, temperature and other cycle information
- Cycle times from 13 minutes for specific loads
- Mobile App for Real Time Remote Monitoring

Efficient, reliable, time-saving & easy to use



UWHL17B

UWHL22B



Decon

Assistina One



Internal cleaning & lubrication of one handpiece in only 18 seconds. With Assistina One the aerosols stay in the HEPA filter and outcoming air is safe!

Includes care set and Quick RM and ISO adaptors.



UWHAT02 with RM+ISO Quick Adaptor £1395

Assistina TwinDual Chamber System



Achieve consistent and efficient handpiece lubrication. The Assistina Twin speeds up your handpiece care and maintenance process as you can connect a second handpiece while the first is being maintained. The HEPA filter thereby releases clean air.

Customer Feedback..

I did my research when it came to choosing equipment for my new decontamination room and W&H came out on top.

Now I have certainty that my handpiece can be reprocessed between treatments without additional waiting time for me or my patients. It's a win win situation!

The dual chamber means
I only need a fraction of the
time for the maintenance of
the handpieces. It is extremely
quiet and looks great too!

Seal2

Universal Pouch Sealing System



2 Second Perfect Seal!

Produces a 12mm hermetic seal in just two seconds. Roll holder can be wall mounted. Electronic temperature regulator guarantees perfect sealing and peace of mind.



EliReel

Sealing Pouches

Peace of mind for clinics and patients due to integrated steam indicator and validated drying system. 200m rolls available in different sizes.

Less wastage as the bags can be cut to the required size.



Seal2 Starter Pack UWHS2SS £1150



Sealing System 3 x Elireel **Pouches**



Motors & Handpieces

Piezomed

Precision and innovation, all in one device

The innovative Piezomed from W&H facilitates precise and safe treatment, with the Piezomed Module offering the flexibility to be combined with the new Implantmed SI-1023.

For optimum accuracy, the Piezomed employs highfrequency micro-vibrations, ensuring that only bone substance is resected, reducing pain for patients and speeding up the healing process.

The Piezomed has a unique, individual performance setting; this detects the instrument during insertion and automatically sets the correct power class. Clinicians can also save their last selected setting, allowing the Piezomed to adapt to their personal way of working. Other impressive features include an LED ring, which offers shadow-free illumination of the operating site, as well as an efficient supply of coolant that is sprayed close to the instrument's working area.

So, why not choose a surgical solution that ensures gold-standard treatment.

Stand-alone Optic Piezomed

- Piezomed Classic





Shadow-free illumination of operating site

Handpiece with LED ring

UWHPMW0

f5599

Piezomed Instruments

The W&H Piezomed instrument range includes different instruments for bone surgery, periodontology and extraction, for retrograde endodontics and for lateral sinus lift.



Compatible with all Piezomed units



Periodontology and extraction

Instruments for thorough cleaning in periodontitis treatment & for atraumatic tooth extraction

Tips: EX1, EX2, P1, P2LD, P2RD



Bone Surgery

The fine saw tooth design more teeth per instrument - sets new standards in the field of bone preparation

Tips: B1, B2R, B2L, B3, B4, B5



Lateral sinus lift

The instruments permit safe lateral preparation of the maxillary sinus wall & gentle mobilization of the Schneiderian membrane

Tips: S1, S2, S3, S4, S5



Retrograde endodontics

Diamond-coated, delicate instruments for easy access during retrograde root canal preparation

Tips: R1D, R2RD, R2LD, R3D, R4RD, R4LD



Preparation of the implant site and crestal sinus lift

The preparation instruments are specially matched to the bone qualities predominant in the maxilla.

Tips: I1, I3P, I4P, Z25P, Z35P

Add on modules...

Piezomed Classic Module



UPMCMDL £POA

- Power 18 W
- Handpiece without light SA-40
- Instrument changer
- 6 Irrigation tubes

Piezomed Plus Module



UPMPMDL

- Power 24W
- Handpiece with light SA-40 L
- Automatic instrument detection
- Instrument changer
- Instrument library
- 6 Irrigation tubes



Implantmed Plus SI-1023

Easy to use, powerful and reliable

The SI-1023 Implantmed Surgical Motor offers a range of innovative features, including a modular system for ultimate upgradeability. It is an efficient, versatile partner with cuttingedge functionality to meet the needs of every modern oral surgeon and implantologist. The Implantmed SI-1023 offers exceptional power and ease of use amongst many other benefits, including:

- Ergonomic, ultra-short, lightweight motor with 6.2 Ncm torque
- Precise torque control which can be adjusted from 5 to 80 Ncm, facilitating high safety standards when inserting implants
- Flexible speed range between 200 40,000 rpm
- Automatic thread cutting function
- Wireless foot control (optional) enhances flexibility with the ability to control multiple W&H devices
- The intuitive touch screen is user-friendly with a sleek glass surface

Implantmed Plus SI-1023

Enjoy the benefits of utilising a synergistic approach with the modular Implantmed surgical motor and Piezomed piezo surgery module alongside the Osstell Beacon implant stability measurement device for a full surgical solution. The Osstell Beacon used in conjunction with the Osstell Connect consultation tool provides ideal patient communication information along with an effective monitoring system to ensure predictable assessment of ideal loading.

Workflow driven **solution**



Implantmed Plus SI-1023 Kit 01 WHIMP1 £POA

- Non-optic motor
- Cable foot control

Implantmed Plus SI-1023 Kit 02 WHIMP3 FPOA

- Optic motor
- Cable foot control

Implantmed Plus SI-1023 Kit 03 WHIMP2 FPOA

- Non-optic motor
- Wireless foot control

Implantmed Plus SI-1023 **Kit 04** WHIMP4

- Optic motor
- Wireless foot control

Implantmed Classic SI-923

MPLNMD0 £2470



Easy to use, powerful and reliable

- Ultra short motor with 5.5 Ncm torque
- Speed range between 300 40,000 rpm
- Precise torque between 5 70 Ncm
- Automatic thread cutting function

Surgical Handpieces

Robust, reliable and comfortable



USUS11S S-11 f595

- Straight handpiece
- External single spray
- Dismantleable
- Without light

USU0S11L S-11 L £886

Mini LED+



USUW56E £875

- Contra-angle 1:1
- External single spray
- Dismantleable
- Without light

UWHWS56L **WS-56 L** £**995**

Mini LED+



USUS160 £945 **S-16**

- Angled handpiece
- External single spray
- Without light



USUSZ75 SZ-75 £1190

- Zygoma handpiece 20:1
- External single spray
- Without light

(MBH)

Osstell Beacon



Reduce treatment times & improve patient outcomes

Measure and evaluate the stability of the implant and control its osseointegration quickly and easily thanks to the easy to read colour light indicator and ISQ reading. Reduce treatment times and improve patient outcomes, especially when treating patients with high risk factors.

- Gives a diagnosis within seconds & presents the user with an easy to understand ISQ
- Wireless operation with Bluetooth
- Helps to non-invasively measure implant stability and osseointegration without jeopordising the healing process
- Improves the predictability of outcomes especially when early loading
- Unique technology based on 1400 scientific articles



UWHOSBN £1895



Osstell SmartPegs

The SmartPeg used together with Beacon measures the resonance frequency when it is attached to an implant. For reliable diagnosis and peace of mind always request the Osstell SmartPeg.

Pack of 5 **£114**

- Single-use ensures accuracy and implant integrity
 - Soft aluminium ensures no damage to implant

Osstell Connect

Free App for **Osstell Beacon**

Osstell Connect is the free online platform that offers dentists a new way to use data to improve their implant treatments. It allows the clinician to make decisions based on detailed and relevant information.





PremiumCare

The W&H PremiumCare Plans have been developed to ensure optimization of the performance and lifespan of your equipment.

They offer you the **confidence** of knowing that your W&H product is fully compliant with all current guidelines and is operating at peak performance.



Enjoy peace of mind with W&H PremiumCare.



Quality Guarantee

We only provide fully qualified factory-trained W&H accredited engineers. In the event of a breakdown, you will receive priority call out.



Plan Options

With a choice of plan options you have the flexibility to choose the plan that best suits your needs.

Contact Us **Trycare**

01274 88 55 44 dental@trycare.co.uk www.trycare.co.uk Mon - Fri 8.30am - 5.30pm



Rialto Buildings Station Road Clayton, Bradford West Yorkshire BD14 6JA