



Clip · Clip F · Clip Flow

LIGHT-CURING FILLING MATERIAL FOR
TEMPORARY RESTORATIONS



Clip · Clip F · Clip Flow

SAVING TIME IN YOUR DENTAL CLINIC

Clip, Clip F and Clip Flow are the ideal composite materials for temporary fillings of all types. Their indications reach from application in inlay / onlay technique to repeated temporary obturations in endodontic treatment. Clip Flow can, in addition, be used for relining, as a block out material, for covering the gingival margin, for fixing resin matrices during filling placement, and for sealing root canal orifices.

Saving time in your dental clinic

Clip is a restorative material for clean, temporary obturations.

The procedure of applying the material saves time as it can be cured using a halogen or LED polymerisation light. More over Clip remains elastic after curing and can be removed in one piece with a suitable instrument before permanent treatment of the cavity.



Easy removal in one piece for all indications

Source: Dr. Peter Bastobbe, Mainleus / Germany

Furthermore, the excellent seal provided by Clip protects from undesirable food particles and bacteria entering the cavity. Thus, Clip combines many advantages in one material. The process of treatment becomes simpler for you, and less time-consuming. Points in favour of Clip are its uncomplicated handling, its highly economical properties, the saving of material, and the simple removal of the provisional filling.

Clip in inlay / onlay technique

One advantage of all products in the Clip range is their simple use for the purposes of the inlay / onlay technique. As Clip can be removed in one piece, there is no need to use a drill to remove the temporary filling from the cavity. Thus, the cavity remains unchanged by the removal of Clip, and there is no adverse effect on the fit of the definitive restoration.



Application of Clip Flow into the cavity Result

Endodontic treatments

Clip products are also utilised in endodontics: in prolonged root treatment requiring more than one sitting Clip provides a saliva-proof seal for the cavity and thus protects it from germinal invasion. Here, too, the simple and time-saving handling of the material is advantageous during repeated sealing and re-opening in endodontic treatment.

Clip F – Fluoride-releasing

Clip F additionally contains fluorides; hence Clip F combines the advantages of Clip and those of a continued release of fluorides inside the freshly prepared cavity.

Clip Flow

VERSATILE AND FLOWABLE

Clip Flow as a relining material

Clip Flow is excellent for relining temporary crowns and bridges, guaranteeing perfect fit until the definite treatment is carried out.

Clip Flow as a block out material

Before taking impressions, for example, Clip Flow is the ideal material for blocking out retentive areas in the dental arch. With Clip Flow taking impressions can be considerably simplified.

Fixing of resin matrices with Clip Flow

Clip Flow allows the simple fixing of resin matrices to optimally hold them in place during placement of the filling.



Application of Clip Flow



Polymerisation



Blocked out interdental spaces



Removal of Clip Flow



Fixing the matrix



Restoring form and function of the tooth

Source: Dr. Walter Denner, Fulda / Germany

Implant Abutments

The aperture for the connecting screw on the implant-abutment interface is frequently given a temporary seal before the definitive luting of the restoration onto the abutment. This prevents luting cement from entering the aperture and thus sealing it permanently. With Clip Flow, this can be prevented in a quick, comfortable and reversible manner.

Gingiva protection with Clip Flow

In order to protect the teeth during, for example, the application of in-office whitening products, Clip Flow can be used for gingiva protection. It acts as a barrier for the whitening product and protects the surrounding soft tissue and teeth which are not to be whitened. Covering the gingival margin with Clip Flow is quick and simple.



Abutment with open aperture



Aperture sealed with Clip Flow



Application of gingiva protection



Removal in one piece



Light-curing of the material



Finished seal

Source: Dr. Walter Denner, Fulda / Germany

Clip · Clip F · Clip Flow

LIGHT-CURING FILLING MATERIAL FOR TEMPORARY RESTORATIONS

Indications

All types of temporary fillings
Temporary protection in inlay / onlay technique

Additional indications Clip Flow

Sealing of openings for implant screws
Relining material for temporary crowns and bridges
Block out material for retentive areas in the dental arch, e.g. before taking impressions
Covering of the gingival margin
Fixing of resin matrices during filling placement
Temporary filling of cavities and sealing of root canal orifices

Advantages

- Easy to place – easy to remove in one piece
- No damage to preparation boundaries
- Durable and saliva-proof
- Non-sticky
- Time-saving through light-curing
- Elastic
- Flowable (Clip Flow)
- NDT® syringe (Clip Flow)

Presentation

Clip

REF 1284	Syringe 2 × 4 g
REF 1285	Syringe 3 × 4 g
REF 1287	Caps 25 × 0.25 g
REF 9301	Dispenser – Caps

Clip F

REF 1289	Syringe 2 × 4 g
REF 1283	Syringe 3 × 4 g

Clip Flow

REF 1282	Syringe 2 × 1.8 g, accessories
REF 2148	Application cannulae type 46, 100 pcs.



VOCO GmbH
Anton-Flettner-Straße 1-3
27472 Cuxhaven
Germany

Freecall: 00 800 44 444 555
Fax: +49 (0) 4721-719-140

info@voco.com
www.voco.dental

Available from: