

AUTOCLAVABLE ELECTRODES

All the PRECISION You Require in Surgery with the VERSATILITY You Need.


All Ellman® electrodes are designed for use in high frequency (4 MHz) RF Generators. The insulated shafts can be bent or contoured to facilitate use in various surgical procedures.



Cat. No. XA



Cat. No. XB




Cat. No. XC



Cat. No. XD



Cat. No. XT



Cat. No. XP

CHOICE OF WIRE GAUGE
 — Fine - Used when minimal tissue alteration is desired
 — Regular - Used to produce more coagulation

	<u>Cat. No.</u>
	TA1
	TA2
	TA9
	TA3
	TF1
	TA8
	TBB
	TB1, TB4
	TB9, TB10
	TB2, TB5
	TB3, TB6
	TD3
	TD8
	TD9
	TD10
	TD11
	TC1, TC4
	TC3, TC6
	TC17, TC18
	TC19, TC20
	TC21, TC22
	TC9, TC10
	TC7, TC8
	TC11, TC12
	TC13, TC14
	TC15, TC16
	TP1, TP2
	TP3, TP4
	TP5, TP6
	TP7, TP8
	TP9, TP10