



## TOUAREG CloseFit™ - NP Ø3.0

OsseoFix™ Surface  
Titanium Grade 5

Catalog No.	Length mm		
ISPNP1030	10		
ISPNP1130	11.5		
ISPNP1330	13		
ISPNP1630	16		
ISPNP1830	18		

## Drill Sequence

BONE TYPE D-IV	BONE TYPE D-II & III	BONE TYPE D-I
1. 2.0	1. 2.0 2. (2.8)	1. 2.0 2. 2.8

Note: (x.xx) drill to the depth of the cortex only. All measurements are in mm.

⚠ CAUTION: The drill preparation is up to 1mm longer than the implants.

