

Less splatter. More comfort.

VivaDent® Polishing Paste

Prophy pastes



Less splatter. Highly efficient.

Professional, effective and gentle polishing without splatter.

Polished tooth surfaces are not only easier to clean, but it also reduces the bacterial load in the mouth, promoting long-term oral health. Polishing can help reduce the risk of future dental problems, particularly when combined with regular brushing and dental check-ups. In summary, professional tooth polishing plays a crucial role in both maintaining dental esthetics and promoting good oral hygiene by removing stains, plaque and bacteria.

Unique formulation

Thanks to their smooth consistency, the polishing pastes adhere well to teeth and instruments. Because of their unique formulation, VivaDent® polishing pastes are easy to rinse away. They are gluten-free and do not contain any microplastics.

Discover the gentle power and effective cleaning of our VivaDent® polishing pastes:

- ✓ **Our go-to solution** for effective and enamel-friendly polishing.
- ✓ **Designed to meet your patients' needs**, VivaDent® polishing pastes **ensure optimal results** in professional tooth cleaning procedures, regardless of the surface of natural teeth, restorations, implants or prosthetic reconstructions.
- ✓ Even high-quality restorations **maintain their lustre**.





Efficient and comfortable polishing with VivaDent®.



Less splatter. More comfort.

Designed to minimize splatter, the VivaDent® polishing pastes ensure a cleaner and more comfortable treatment environment for both practitioners and patients. The treatment area is kept cleaner, and visibility is improved.



Less splatter. More efficiency.

Thanks to their smooth consistency, VivaDent® polishing pastes polish with less product. No mixing and no preparation of the pastes is needed. This saves time during each procedure. The pastes spread easily, making them ideal for achieving a thorough polish in less time.



Less splatter. More variation.

We offer a range of abrasiveness levels and flavours, allowing you to tailor the polishing treatment to each patient's needs: Sensitive areas can be polished with VivaDent® Polishing Paste Medium (a fluoride-free option is available), while more intensive cleaning can be done with VivaDent® Polishing Paste Coarse; or you can simply enhance your tooth cleaning workflow with the non-scratching VivaDent® Polishing Paste Fine.



Less splatter. Easy handling.

The pastes are available in both tubes and single dose units, giving you flexibility in use. The tubes offer a convenient, long-lasting option for frequent treatments, ensuring you always have the right amount of paste on hand. Single dose units, however, are perfect for individual patients, providing a hygienic, mess-free option – ideal for busy practices and mobile treatments.



Less splatter. More variation.

Colour coded for easy selection of the grit

For better and easier identification of the different abrasion levels, the tubes are coded by the following colours:

Pink:	Fine - containing 500ppm sodium fluoride
Light Green:	Medium - containing 500ppm sodium fluoride
Dark Green:	Fluoride-free Medium
Blue:	Coarse - containing 500ppm sodium fluoride

New flavours

We are growing our range with two new flavours, Apple and Bubblegum. These new options provide a pleasant taste experience for patients while maintaining the high standards of care expected in professional dental cleaning.

Splatter free

Experience precision and cleanliness with VivaDent polishing pastes. Engineered to minimize splatter, our polishing pastes ensure smooth, efficient operation even at high speeds of up to 3000 rpms.



Aerosol-free tooth cleaning - elevating dental standards:

OptraGate®, VivaDent® Polishing Paste Coarse and VivaDent® Aerosol Reduction Gel in combination represent a significant advancement in dental care. They combine effectiveness with a commitment to personal and patient safety as well as environmental responsibility and better visibility of the treatment area. Ideal for routine dental cleaning and patients with sensitivities, this combination of products promotes a very high standard of care.



VivaDent Polishing Paste Fine

VivaDent Polishing Paste Fine is specifically formulated for polishing in a gentle yet effective manner. Designed for use after professional dental procedures, VivaDent Polishing Paste Fine ensures plaque removal and thorough polishing results without scratching or discomfort, on any surface, even composite resin or titanium.

VivaDent Polishing Paste Medium

VivaDent Polishing Paste Medium is designed for efficient plaque and stain removal, as well as polishing. It offers the ideal balance between cleaning power and gentleness. The paste is suitable for routine use on natural teeth and restorative materials, making it an essential part of regular dental hygiene treatment.

VivaDent Polishing Paste Fluoride-free Medium option without flavouring

Fluoride-free polishing paste is suitable for patients with allergies or sensitivities and those who prefer to avoid fluoride for personal reasons. It is also compatible with treatments involving self-etch adhesive bonding agents, as it does not interfere with polymerization or contain problematic flavorants. This option allows dental practices to cater to a broader range of patient needs while maintaining high-quality care.

VivaDent Polishing Paste Coarse

This paste is specifically formulated for intensive cleaning and polishing in cases of heavy staining and hard plaque. It is ideal for use in situations where thorough, deep cleaning is required: for example, stains caused by coffee, tea, tobacco or certain medications. VivaDent Polishing Paste Fine provides maximum cleaning power while still being safe for the enamel and therefore delivers the perfect final touch.

Delivery forms

Proxyt Paste	
761376	VD Polishing Paste Fine Mint 80g
761377	VD Polishing Paste Medium Mint 80g
761378	VD Polishing Paste Medium ff 80g
761379	VD Polishing Paste Coarse Mint 80g
759502	VD Polishing Paste Fine Apple 80g
759515	VD Polishing Paste Fine BGum 80g
759516	VD Polishing Paste Medium Apple 80g
759517	VD Polishing Paste Medium BGum 80g
759518	VD Polishing Paste Coarse Apple 80g
759519	VD Polishing Paste Coarse BGum 80g
761380	VD Polishing Paste Fine Mint 200x1.8g
761381	VD Polishing Paste Coarse Mint 200x1.8g
760115	VD Polishing Paste Fine Apple 200x1.8g
760116	VD Polishing Paste Fine BGum 200x1.8g
760117	VD Polishing Paste Coarse BGum 200x1.8g
760120	VD Polishing Paste Medium Apple 200x1.8g
760119	VD Polishing Paste Medium ff 200x1.8g

Ivoclar Vivadent AG | Bendererstr. 2 | 9494 Schaan | Liechtenstein
Tel. +423 235 35 35 | Fax +423 235 33 60

[ivoclar.com](https://www.ivoclar.com)