Warwick Range



SP2313 Silicone Pad - For trays IT/MIT2715 or small DIN Tray

Long Life Re-usable Silicone Protection Pad to protect Instruments during Washing, Sterilization and Transportation

Product Benefits



Long silicone fingers hold instruments securely in place. Pliable pad with perforations for maximum water penetration and

Autoclavable silicone protection pads to protect delicate instruments during washing, sterilization and transportation.

- steam circulation gives improved cleaning and drying results. Lightweight, strong and durable.
- Can be sterilized by Autoclave, Ethylene Oxide or Hydrogen Plasma.

The Warwick Range of polypropylene products from Warwick SASCo was originally designed to offer an alternative to stainless steel products traditionally used in Hospitals.

This Range is now widely used in Sterile Services and Decontamination Units, Operating Rooms and Clinics due to its low cost, effective product design, and repeat performance results in hospital autoclaves.

Easily recognised by its distinct blue colour, the Warwick Range includes bowls, gallipots, kidney dishes, instrument trays, instrument protection, instrument boxes, traceability tags, measuring jugs, scope and instrument



Product Specification

- Silicone Pad for IT/MIT2715 Instrument Tray
- 230mm x 132mm
- Pack Size 10
- The following standards are applicable:
 - ISO 9001:2008

Quality Management System

- ISO 13485:2003
 - Medical Devices Quality Management System
 - ISO 17664:2004

Sterilization of Medical Devices - Guidelines for the processing of re-sterilizable medical devices

Maximum processing temperature 200° C

Guide to Using Warwick Products from Warwick SASCo Cleaning

- Product should be washed before first use maximum temperature 200°C.
- Warwick Products are designed for a lengthy service and guaranteed for 1000 washes and sterilizations. Scouring and abrasive cleaners will reduce the performance of products. Acid/Amine solutions are extremely corrosive and must be used in the correct dilution according to manufacturers instructions.
- Silicone is not suitable for extended immersion in acid based cleaning solutions.

Storage

- This product has a shelf life of 12 months from date of manufacture providing the recommended storage conditions are adhered to.
- Product should be stored in its primary packaging, off the floor on a pallet or shelf, and in a clean dry ventilated room at a temperature of 4 °C - 25 °C.