


# Hole in One

Filling in a single step with **Fill-Up!**<sup>™</sup>

Try this amazing  
new product!



**Fill-Up!**<sup>™</sup> - in a single  
step to give a perfect result

- **Ideal for Amalgam replacement**
- Optimal depth polymerisation with minimal shrinkage due to dual curing system restoration
- Guaranteed single-layer technique - even in very deep cavities of 10 mm
- Optimised sealing of margins - reduced postoperative sensitivity
- Universal shade in a convenient Automix syringe for efficient placement

**Deep. Fast. Perfect.**

For further information:


**Freephone: 0500 295454**

Ext. 223 Julie or 224 Chris

[info.uk@coltene.com](http://info.uk@coltene.com) | [www.coltene.com](http://www.coltene.com)



**COLTENE**

 [www.facebook.com/coltene.uk](http://www.facebook.com/coltene.uk)




[@COLTENEUKLTD](https://twitter.com/COLTENEUKLTD)


 Follow us on  
Dental tubules

# Fill-Up!<sup>™</sup>


Dual curing bulk composite




etch




rinse




air dry




dispense & mix




massage




air dry




fill



wait for setting



light cure



finish