

4 reasons to use Luxator® Dual Edge

#1 Dual Edge

Integrated periotome and Luxator tip.



#3 Ideal for fractured teeth

Unique Dual Edge tip enables to find the root and margin and then start the luxating procedure.

#2 Safer insertion

The Dual Edge tip offers safer stepped insertion and reduces the risk of slipping.

#4 Perfect starter in difficult cases

The thin periotome part is ideal to use as a starter for difficult extractions.



**DE3
Dual Edge**
3/1.5 mm
Straight
Standard



**DE5
Dual Edge**
5/3 mm
Straight
Standard

Technique



Luxator Dual Edge Periotome Innovative Dual Purpose Tip

Integrated probing Luxator Periotome. Safer stepped insertion reduces the risk of slipping. Ideal for deeply fractured teeth and decayed roots. Available in 2 blade sizes. Sterilizable Fixus 1 case included.

Luxator® Dual Edge Periotome with Fixus 1

Prod No	Product Description	Quantity
506356	DE3 3mm/1,5mm, dual edge, grey	1 pc
506357	DE5 5mm/3mm, dual edge, brown	1 pc

