



Directa likes!



Preserved bone integrity and dramatically less trauma for your patients and their dentition

Luxator Periotomes are the preferred instruments for performing extractions.

These instruments were invented by a Swedish dentist to make extractions as trauma free as possible. He developed subtleties in the design that only a practicing dentist would appreciate. With an acclaimed, ergonomic handle Luxator instruments are discernibly different. The instruments are manufactured of high quality Swedish stainless steel.

Luxator Periotomes are specially designed periodontal ligament knives with fine tapering blades that compress the alveolar bone, cut the membrane and gently ease the tooth from its socket. The whole operation is performed with a minimum of tissue damage.

Luxator® is a registered trademark owned by Directa.

4 reasons to use Luxator® Periotome

#1

Ergonomic design

World acclaimed ergonomic design of the handle. Minimizes fatigue.

#2

Atraumatic extractions

The fine tapered blade compresses the alveolar bone and cuts the periodontal ligaments.

#3

Maximum tactility

The use of highest quality materials gives the instruments outstanding tactility.



#4

Wide range of sizes and shapes

Unique variety of sizes and shapes enables users to find the ideal instrument for any extraction.



Technique





L3CA
Periotome
3 mm
Contra
Standard

L3IC
Periotome
3 mm
Inverted
Standard

L5C
Periotome
5 mm
Curved
Standard

L5S
Periotome
5 mm
Straight
Standard

L1S
Periotome
1 mm
Straight
Standard

L2S
Periotome
2 mm
Straight
Standard

L3S
Periotome
3 mm
Straight
Standard

L3C
Periotome
3 mm
Curved
Standard

Luxator® Periotome with Fixus1

Prod No	Product Description	Quantity
506353	L3CA 3mm, contra angle, dark blue	1 pc
506354	L3IC 3mm, inverted curved, light blue	1 pc
506343	L5C 5mm, curved blade, brown	1 pc
506342	L5S 5mm, straight blade, tan	1 pc

Luxator® Periotome with Fixus1

Prod No	Product Description	Quantity
506355	L1S 1mm, straight blade, green	1 pc
506352	L2S 2mm, straight blade, dark green	1 pc
506341	L3C 3mm, curved blade, black	1 pc
506340	L3S 3mm, straight blade, grey	1 pc

Kit

506330	Kit L-K4	L3C, L3S, L5C, L5S, L10	4 pcs
506331	Kit L-K7	L2S, L3C, L3S, L3CA, L3IC, L5C, L5S, L10	7 pcs
506332	L-K VET	L1S, L2S, L3C, L3S, L3CA, L10	5 pcs
506399	L10	Sharpening Tool	1 pc