

Citrisil™

THE ONLY ALL IN ONE TABLET.

Citrisil™ is the **ONLY ALL IN ONE**, Shock and daily maintenance dental water line tablet with a quantified EPA Claim. **Citrisil™** and **Citrisil™ Blue** tablets continuously disinfect and maintain dental water lines. Simply start with one **Citrisil™ Shock** tablet and follow with one **Citrisil™** (white) or **Citrisil™ Blue** tablet per refill of each self-contained water bottle.

Using distilled or purified tap water, the tablet effervesces and provides a safe, pH balanced treatment. **Citrisil™ Blue** produces a light blue tint for visual compliance. EPA registered to produce less than or equal to 10 CFU per ml HPC purity.



PRODUCT HIGHLIGHTS

- + Eliminates the need to empty bottles and purge lines at the end of each day
- + Citrisil™ Shock included with all Citrisil™ and Citrisil™ Blue Daily maintenance boxes
- + Safe for patients and clinicians
- + Will not harm dental equipment
- + Does not adversely affect composite bond strength
- + Maintains BioFree® tubing
- + Contains no oxidizers and does not dissolve amalgam
- + Contains no iodine or chlorine
- + Citrisil™ Blue produces a light blue tint for visual Confirmation of compliance
- + Available for .70 to 1 liter bottles & 2 liter bottles
- + EPA registered to produce less than or equal to 10 CFU per ml HPC purity
- + FDA Cleared, EST.#3005844469
- + EPA Reg. No 83315-1
- + EPA Est. No. 83315-CO-001

For more information about our products contact your **Sterisil® Sales Representative** or local dealer.

For additional information call 01274 88 55 44 or visit www.trycare.co.uk



SELECTION CRITERIA

Part #	Description
C20-W	Citrisil™ White, (20) Maintenance Tablets & (1) Citrisil™ Shock Tablet, for .70 - 1 liter bottles.
CS-W	Citrisil™ Case Of (50) C20-W Boxes, White Tablets,.70-1 liter B.
C20-B	Citrisil™ Blue, (20) Maintenance Tablets & (1) Citrisil™ Shock Tablet, for .70 - 1 liter bottles.
CS-B	Citrisil™ Case of (50) C20-B Boxes, Blue Tablets, .70-1 liter B
C20-2W	Citrisil™ White, (20) Maintenance Tablets & (1) Citrisil™ Shock Tablet, for 2 liter bottles.
CS-SW	Citrisil™ Case of (50) C20-2W Boxes, White Tablets, 2 liter B.
C20-2B	Citrisil™ Blue, (20) Maintenance Tablets & (1) Citrisil™ Shock Tablet, for 2 liter bottles.
CS-2B	Citrisil™ Case of (50) C20-2B Boxes, Blue Tablets, 2 liter B.
ECS-10	Citrisil™ Shock Tablet, (10) Tablets.
CS-ECS-10	Citrisil™ Shock Tablet, Case of (50) Boxes of ECS-10.

Sterisil®, it doesn't get any easier, and it certainly doesn't get any better!

For more information about our products contact your **Sterisil® Sales Representative** or local dealer.

For additional information call **01274 88 55 44** or visit **www.trycare.co.uk**