

# Terumo Neolus<sup>®</sup>

*Hypodermic needle*



Unsurpassed comfort for patients  
and healthcare workers

- *Easy penetration and less painful puncture*
- *Smooth needle movement*
- *Increased flow rates*

## Easy penetration and less painful puncture

assured by the ultra-sharp lancet point with precisely ground edges.

## Smooth needle movement


guaranteed by high-quality silicon coating.

## Increased flow rates

realised by the Terumo-thin wall, facilitating the choice of more comfortable gauge sizes.

## Additional features and benefits:

- Easy handling with colour-coded hub and packaging
- Meeting customer requirements with wide range of products
- Reduced risk of vascular perforation with short bevels



## Products and specifications



Product Code	Colour	Cannula ø	Cannula Length	Wall Thickness
NN-1838R	Pink	18 G (1,2 mm)	1 ½" (40 mm)	Thin
NN-1838S	Pink	18 G (1,2 mm)	1 ½" (40 mm)	Thin
NN-1850R	Pink	18 G (1,2 mm)	2" (50 mm)	Regular
NN-1925R	Ivory	19 G (1,1 mm)	1" (25 mm)	Thin
NN-1938R	Ivory	19 G (1,1 mm)	1 ½" (40 mm)	Thin
NN-1938S	Ivory	19 G (1,1 mm)	1 ½" (40 mm)	Thin
NN-1950R	Ivory	19 G (1,1 mm)	2" (50 mm)	Regular
NN-1030R	Ivory	– (1,0 mm)	1 ¼" (30 mm)	Thin
NN-2013R	Yellow	20 G (0,9 mm)	½" (12 mm)	Thin
NN-2025R	Yellow	20 G (0,9 mm)	1" (25 mm)	Terumo-Thin
NN-2038R	Yellow	20 G (0,9 mm)	1 ½" (40 mm)	Terumo-Thin
NN-2038S	Yellow	20 G (0,9 mm)	1 ½" (40 mm)	Terumo-Thin
NN-2050R	Yellow	20 G (0,9 mm)	2" (50 mm)	Regular
NN-2055S	Yellow	20 G (0,9 mm)	2 ⅜" (55 mm)	Regular
NN-2070S	Yellow	20 G (0,9 mm)	2 ¾" (70 mm)	Regular
NN-2116R	Green	21 G (0,8 mm)	⅝" (16 mm)	Regular
NN-2125R	Green	21 G (0,8 mm)	1" (25 mm)	Terumo-Thin
NN-2132R	Green	21 G (0,8 mm)	1 ¼" (30 mm)	Terumo-Thin
NN-2138R	Green	21 G (0,8 mm)	1 ½" (40 mm)	Terumo-Thin
NN-2138S	Green	21 G (0,8 mm)	1 ½" (40 mm)	Terumo-Thin
NN-2150R	Green	21 G (0,8 mm)	2" (50 mm)	Regular
NN-2232R	Black	22 G (0,7 mm)	1 ¼" (30 mm)	Terumo-Thin
NN-2232S	Black	22 G (0,7 mm)	1 ¼" (30 mm)	Terumo-Thin
NN-2238R	Black	22 G (0,7 mm)	1 ½" (40 mm)	Terumo-Thin
NN-2250R	Black	22 G (0,7 mm)	2" (50 mm)	Regular
NN-2316R	Blue	23 G (0,6 mm)	⅝" (16 mm)	Terumo-Thin
NN-2325R	Blue	23 G (0,6 mm)	1" (25 mm)	Terumo-Thin
NN-2332R	Blue	23 G (0,6 mm)	1 ¼" (30 mm)	Terumo-Thin
NN-2332S	Blue	23 G (0,6 mm)	1 ¼" (30 mm)	Terumo-Thin
NN-2425R	Violet	24 G (0,55 mm)	1" (25 mm)	Terumo-Thin
NN-2516R	Orange	25 G (0,5 mm)	⅝" (16 mm)	Terumo-Thin
NN-2525R	Orange	25 G (0,5 mm)	1" (25 mm)	Terumo-Thin
NN-2613R	Brown	26 G (0,45 mm)	½" (12 mm)	Regular
NN-2623R	Brown	26 G (0,45 mm)	⅞" (23 mm)	Regular
NN-2716R	Grey	27 G (0,4 mm)	⅝" (16 mm)	Regular
NN-2719R	Grey	27 G (0,4 mm)	¾" (20 mm)	Regular
<b>Also available: BLUNT NEEDLE (mixing needle)</b>				
BN-1838	Pink	18 G (1,2 mm)	1 ½" (40 mm)	Thin

Unit Box: 100 pieces  
Shipping Carton: 3.000 pieces



R = Regular bevel 12°  
S = Short bevel 18,5°



**Terumo Corporation**  
44-1, 2-chome, Hatagaya  
Shibuya-ku, Tokyo 151-0072,  
Japan  
Tel.: +81 3 3374 8111  
Fax: +81 3 3374 8196

### European Headquarters

**Terumo Europe N.V.**  
Interleuvenlaan 40  
3001 Leuven, Belgium  
Tel.: +32 16 38 12 11  
Fax: +32 16 40 02 49

### European Offices

**Terumo Europe N.V.  
Benelux Sales Division**  
Interleuvenlaan 40  
3001 Leuven, Belgium  
Tel.: +32 16 39 25 80  
Fax: +32 16 39 25 99  
The Netherlands:  
Tel.: 0800 0220396  
Fax: 0800 0220414

**Laboratoires Terumo France S.A.**  
Parc Ariane - Bât. Uranus  
1, rue Hélène Boucher  
78284 Guyancourt Cedex, France  
Tel.: +33 130 96 13 00  
Fax: +33 130 43 60 85

**Terumo Deutschland GmbH**  
Ludwig-Erhard-Strasse 6  
65760 Eschborn, Germany  
Tel.: +49 6196 80 230  
Fax: +49 6196 80 23 200

**Terumo Deutschland GmbH**  
Zweigniederlassung Spreitenbach  
Bodenackerstrasse 3  
8957 Spreitenbach, Switzerland  
Tel.: +41 56 419 10 10  
Fax: +41 56 419 10 11

**Terumo Europe España SL**  
Avda. Juan Carlos I, Nº 13 - 7ª Planta  
Edificio Torre La Garena  
28806 Alcalá de Henares (Madrid),  
Spain  
Tel.: +34 9021 01 298  
Fax: +34 9021 01 358

**Terumo Italia S.r.l.**  
Via Simone Martini, 143/145  
00142 Roma RM,  
Italy  
Tel.: +39 0651 96 14 20  
Fax: +39 0650 30 407

**Terumo Sweden AB**  
Nya Varvet, Byggnad 90  
Sven Källfelts Gata 18  
426 71 Västra Frölunda,  
Sweden  
Tel.: +46 3174 85 880  
Fax: +46 3174 85 890

**Terumo Denmark  
Filial of Terumo Sweden AB**  
Tel.: +45 7020 93 80  
Fax: +45 7020 94 80

**Terumo UK Ltd**  
Tamesis  
The Causeway  
Egham, Surrey  
TW20 9AW, United Kingdom  
Tel.: +44 1784 476 200  
Fax: +44 1784 476 201

**Terumo Russia LLC**  
BC Northern Tower  
13 Floor, 10 Testovskaya Street  
Moscow 123317, Russia  
Tel.: +7 495 988 4740  
Fax: +7 495 988 4739

**Terumo Corporation Dubai Branch  
(Middle East and Africa)**  
Al Masraf Tower, 22nd Floor  
P. O. Box 20291 Dubai  
U.A.E.  
Tel.: +971 4 221 2220  
Fax: +971 4 221 3330

Addresses subject to change,  
please consult our Website.



© Registered Trademark  
Published by Terumo Europe N.V.  
HO 02GB-0702TE-VIII(08.13)E