

Max-i-Probe<sup>®</sup> becomes pro•rinse<sup>™</sup>



New brand, New packaging, Same Product.

+  
WE  
KNOW  
ENDO.

For better dentistry

DENTSPLY

PRORINSE™ is an irrigating probe with a closed, rounded end and a side port dispersal to assist the removal of debris in the root canal and periodontal pocket

- Improved debris removal (compared to standard open-ended needles).\*
- Sterilised probe and rounded-end for patient protection.
- With Luer-Lock system for convenient use.



**Periodontal Pockets**

Side port produces flushing motion to irrigate pockets


Rounded tip protects tissue




**Root Canals**

Design produces upward flushing motion for canal irrigation

Side port dispersal reduces risks of debris from being expressed through the apex compared to open-ended probes\*

Closed, rounded end reduces risk of apex damage compared to open-ended probes



|   | Description        | Box of 40 Probes | Outer Diameter | Gauge |
|---|--------------------|------------------|----------------|-------|
|  | PRORINSE 23g Pk 40 | A082100023400    | 0.63mm         | 23    |
|  | PRORINSE 25g Pk 40 | A082100025400    | 0.51mm         | 25    |
|  | PRORINSE 30g Pk 40 | A082100030400    | 0.30mm         | 30    |

\*Data on file

DENTSPLY (Australia) Pty Ltd ABN 15 004 290 322  
11-21 Gilby Road, Mount Waverley VIC 3149 Tel: 1300 55 29 29 Fax: 1300 55 31 31  
Email: [clientservices@dentsply.com](mailto:clientservices@dentsply.com)  
[www.dentsply.com.au](http://www.dentsply.com.au)

DENTSPLY New Zealand  
PO Box 113083, Newmarket, Auckland 1149 Tel: 0800 DENTSPLY (33 68 77)  
Fax: 0800 DENTFAX (33 68 32) Email: [clientservicesnz@dentsply.com](mailto:clientservicesnz@dentsply.com)  
[www.dentsply.co.nz](http://www.dentsply.co.nz)

**+  
WE  
KNOW  
ENDO.**