

4 reasons to use Luxator® Forte Elevator

#1

Ergonomic Design

World acclaimed handle. Ergonomic design. Minimizes fatigue.



#3

Cornered tip design

The unique cornered design grasps the crown or root for safe and controlled elevations.

#2

Thin, yet strong

The use of highest quality materials allows the tip to be very thin without its strength thus offering better access.


#4

Maximum tactility

Highest quality materials for instruments with outstanding tactility.

*The perfect match to Luxator Periotome to end a successful extraction!
/Directa*

Technique



F25 Forte
2.5 mm
Straight
Standard

F32 Forte
3.25 mm
Straight
Standard

F40 Forte
4 mm
Straight
Standard

Luxator Forte Elevator Ergonomic & Strong

Incorporates the world acclaimed ergonomic handle and a cornered, rigid tip that enables elevation of roots without risk of damaging the instrument. The unique cornered tip design grasps dentine or root cementum for controlled and safer tooth elevation. Available in 3 blade sizes. Sterilizable Fixus 1 case included.

Luxator® Forte Elevator with Fixus 1

Prod No	Product Description	Quantity
506370	F25 2,5 mm, straight blade, Green	1 pc
506371	F32 3,2 mm, straight blade, Grey	1 pc
506372	F40 4,0 mm, straight blade, Light Brown	1 pc

